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OPERATING INSTRUCTIONS

DISPOSAL U / B pro Type 90 safety storage cabinet

Retain for later use.

The operating instructions are available for download at info@dueperthal.com.
The German operating instructions are the original version. All other versions are translations of the original operating instructions.

Table of contents

1	General information	5
1.1	Notes on handling the operating instructions	5
1.2	Notes for reading	5
1.3	Type plate	6
2	Safety	7
2.1	Function of safety notices	7
2.2	Correct use	8
2.3	Misuse	9
2.4	The operator's responsibility	9
2.5	Requirements placed on employees	10
2.6	Stored goods	10
2.7	Residual risks	11
2.8	Safety concept for explosion protection	12
2.8.1	Active storage	12
2.8.2	Passive storage	14
2.8.3	Passive storage without industrial ventilation	14
2.9	Hazardous areas and their labelling	16
3	Technical specifications	18
3.1	Outdoors	18
3.2	Interior	20
3.3	Pressure drop during ventilation	22
4	Structure and function	23
4.1	Construction	23
4.2	Doors	23
4.2.1	Door options	23
4.2.2	Locking cylinder	23
4.3	Safety technology	24
4.3.1	Door closure in case of fire	24
4.3.2	Closure of venting cut-off flaps in case of fire	24
4.4	Pull-out tray with insert for collection container	25
4.4.1	Version for 10 litre collection container	25
4.4.2	Version for 5 litre collection container	26
4.5	Disposal system	27
4.5.1	Connection head	28
4.5.2	Connection pipe of the disposal system	29
4.6	Analysing unit (monitoring)	30
4.7	Penetration	31
4.8	Accessories (optional)	32
4.8.1	Variable base designs	32
4.8.2	Roller set	33
4.8.3	Extending the analysing unit	33
5	Transport and packaging	34
6	Installation and commissioning	35
6.1	Requirements for the installation location	35
6.2	Safety storage cabinet installation and alignment	36

6.3	Installation and operation of the locking screw on a roller set.....	38
6.4	Check the alignment of the safety storage cabinet	39
6.5	Connecting the safety storage cabinet to an exhaust air system	40
6.5.1	Exhaust air connection and air supply opening	40
6.5.2	Connection to an exhaust air system	41
6.6	Earthing the safety storage cabinet.....	42
6.7	Connecting the connection pipe	44
6.8	Commissioning the analysing unit.....	46
6.9	Calibrating the USB-C Extension Flex for level monitoring	47
6.9.1	Collection container > 5 litres.....	47
6.9.2	Collection container = 5 litres.....	48
7	Operation	49
7.1	Opening the safety storage cabinet.....	50
7.1.1	Opening the safety storage cabinet with the door handle.....	50
7.1.2	Opening a safety storage cabinet with push-to-open door operation.....	50
7.2	Filling and transferring.....	51
7.3	Changing the collection container	52
8	Procedure after a fire	54
9	Maintenance.....	55
10	Faults	58
11	Shutting down the safety storage cabinet.....	60
12	Spare parts and accessories	60
13	Disposal	61
14	Declaration of Conformity	62

1 General information

1.1 Notes on handling the operating instructions

These operating instructions must always be stored with the associated product.

The operating instructions must be available to all users and technically qualified personnel.

If the product is resold, these operating instructions must be included with it.

Subsequent amendments and updates made by the manufacturer must be incorporated into the operating instructions.

1.2 Notes for reading

The following symbols designate specific types of information

Table 1: Explanation of symbol

Symbol	Type of information
	Information for easier and more effective working
	Procedural step
	Result of a procedural step
	Link to another part of the document
	List

1.3 Type plate

The type plate is attached to the outside of the safety storage cabinet door.


Fig. 1: Type plate

Model

Maximum volume of individual containers

S/N: Serial number and year of manufacture
R/N: Item number

Maximum load per storage shelf

The details of the model, serial number and year of manufacture, as well as the item number, can also be found on another plate attached to the inside of the door of the safety storage cabinet.

2 Safety

2.1 Function of safety notices

Safety notices warn of personal injury or material damage and provide information on how such damage can be avoided.

The following signal words identify the degree of danger and the extent of the risk.

	WARNING The signal word " <i>WARNING</i> " indicates a potentially imminent hazard that can result in death or severe injuries.
	CAUTION The signal word " <i>CAUTION</i> " indicates a potentially imminent hazard that can result in minor or slight injuries.
	NOTE The signal word " <i>NOTE</i> " indicates a situation that can lead to damage to the safety storage cabinet.

2.2 Correct use

The safety storage cabinet is type-tested and classified as Type 90 in compliance with “EN 14470-1” with a fire resistance of 90 minutes.

The safety storage cabinet is to be used for the storage of flammable liquids in working areas. It is specially equipped with collection containers and disposal lines for the disposal of flammable liquids.

Any other use shall be deemed incorrect.

WARNING

Dangerous explosive atmosphere!

If flammable liquids are transferred or filled in safety storage cabinets, a potentially explosive atmosphere must be expected to develop in the safety storage cabinet.

The consequences can be death or serious injury.

- Carry out risk assessments and create an explosion protection document.
- Connect the safety storage cabinet to a permanently effective and monitored exhaust air system.
- Observe explosion protection regulations, such as Technical Rules for Hazardous Substances 727, etc.

Observe the safety instructions in these operating instructions to reduce health risks and avoid dangerous situations.

Any use that is not correct use as defined in these operating instructions involves a risk of accidents and a lack of fire protection.

2.3 Misuse

Any use that goes beyond the specified correct use is considered to be misuse.

DÜPERTHAL accepts no liability for damage caused by misuse.

The following constitute misuse:

- Storage of food in the safety storage cabinet.
- Storage of pressurised gas cylinders in the safety storage cabinet.
- Storage of lithium-ion batteries in the safety storage cabinet.
- Overloading the safety storage cabinet or individual storage levels.
- Storage of hazardous substances together, in contravention of national regulations.
- Storage of flammable hazardous substances together with electrical equipment.
- Objects on the cabinet roof.
- Blocking the doors of the safety storage cabinet or using objects to hold these doors open.
- Alteration and modification of the safety storage cabinet without the prior knowledge and approval of DÜPERTHAL.
- Using low-quality spare parts.
- Exceeding the maintenance intervals.
- Failure to collect and properly dispose of spilled liquids immediately in accordance with accident prevention regulations.
- Transferring and filling substances manually without using the connections provided for this purpose.
- Using unsuitable collection containers.

2.4 The operator's responsibility

The operator must ensure that:

- A risk assessment is prepared.
- An explosion protection document is prepared.
- Operating instructions are issued.
- The user has read and understood these operating instructions before using the safety storage cabinet for the first time.
- These operating instructions are available to the user at all times.
- Only trained users have access to and operate the safety storage cabinet.
- The user is trained in handling the stored goods.
- All activities and personal protective equipment are specified for the user.
- The safety storage cabinet is only used correctly.
- The safety storage cabinet is maintained in line with the specifications set out in the operating instructions.
- The safety storage cabinet is operated in accordance with local guidelines and health and safety regulations.
- All precautions are taken to avoid hazards caused by the safety storage cabinet.

2.5 Requirements placed on employees

WARNING

Unauthorised employees!

Employees who do not meet the requirements may contribute to the development of hazards when handling the safety storage cabinet.

Death or serious injuries.

- Designate employees who meet the requirements to carry out activities.

Employees need to meet specific requirements to carry out certain activities. These instructions for use assign activities to the following employees:

- Users
- Technically qualified personnel
- DÜPERTHAL service technicians

Users

Only people who meet the following requirements are approved as users:

- Instructed in the function and operation of the safety storage cabinet.
- Have basic knowledge in the general handling of chemicals and hazardous substances.
- In-depth knowledge of how to correctly store and remove flammable and liquid hazardous substances.

Technically qualified personnel

Technically qualified personnel also meet the following requirements in addition to the general requirements for users:

- Activity-specific professional training in the relevant technical field, which is recognised at the installation location.
- Basic knowledge and skills in the field of maintenance, repair and servicing.

DÜPERTHAL service technicians

DÜPERTHAL employees are specifically trained by DÜPERTHAL to carry out activities in relation to the safety storage cabinet.

2.6 Stored goods

Storage, handling and use of the stored goods must comply with the applicable national standards and regulations, e.g. "Technical Rule for Hazardous Substances 510"(TRGS 510) in Germany.

Transfer, filling or emptying of non-stationary containers must comply with the applicable national standards and regulations, e.g. "Technical Rule for Hazardous Substances 509" in Germany.

2.7 Residual risks

The following residual risks exist for users of the safety storage cabinet.

Risks when working in a potentially explosive area (ex-zone):

WARNING

Dangerous explosive atmosphere!

A lack of or insufficient air exchange can lead to formation of an explosive atmosphere in the safety storage cabinet.

Death or serious injuries due to explosion or deflagration.

- In a safety cabinet with ventilation system, for active storage the minimum hourly air exchange must be at least 30 m³/h when the doors are closed.

WARNING

Spark formation in the safety storage cabinet

Spark formation in the safety storage cabinet can lead to an explosion in the cabinet.

The consequences can be death or serious injury.

- Connect the safety storage cabinet to the earthing system.
- Use only suitable tools or equipment.
- Ventilate the safety storage cabinet thoroughly before working on or in it.
- Observe applicable national standards and regulations, e.g., “TRBS1112, Part 1” in Germany

Risks caused by the media to be transferred:

WARNING

Escaping hazardous substances

Escaping hazardous substances can result in health risks.

The consequences can be death or serious injury.

- Use only suitable storage or transport containers.
- Verify the resistance of the cabinet components to the hazardous substances with which they are to be filled.
- Immediately repair or replace damaged hoses, shut-off mechanisms, etc.
- Collect and properly dispose of spilled liquids immediately in accordance with accident prevention regulations.

2.8 Safety concept for explosion protection

The operator has an obligation to prepare a safety concept based on the operator's risk assessment.

2.8.1 Active storage

When all measures in the protection concept are observed, the ratio between the extraction flow rate and the fill flow rate guarantees that the lower explosion protection limit is exceeded only briefly.

When the safety concept is observed, it may be assumed that an ex-zone 1 inside the safety storage cabinet and at the collection container filling point becomes an ex-zone 2.

The operator is responsible for the classification and compliance of the explosion protection zones.

Compliance with the protection concept:

- Keep doors closed at all times during transfer and filling work.
- Use of a suitable, continuous monitoring system for the exhaust air flow rate.
- Ensure that the monitoring system is continuously operational (24/7).
- Ensure that the safety storage cabinet is vented continuously with a volumetric flow rate of at least 30 m³/h
- After commissioning, after changing containers and before maintenance work, the volume in the cabinet interior must be replaced at least 16 times per hour with doors closed.
- Earth the safety storage cabinet.

The continuous volumetric flow rate of at least 30 m³/h must be increased as necessary in accordance with the operator's safety concept.


Fig. 2: Possible classification of ex-zones in the cabinet interior in compliance with the protection concept

WARNING

Formation of a dangerous explosive atmosphere!

Unexpected malfunctions of the exhaust air system can lead to the development of an explosive atmosphere in the safety storage cabinet.

This may result in death or serious injuries due to explosion or deflagration.

- Stop any ongoing filling or transferring processes immediately.
- Do not carry out any filling or transferring processes until the malfunction has been corrected.
- Keep doors of the safety storage cabinet closed.

NOTE

Increasing the exhaust air flow

As part of the risk assessment, the operator must establish whether the recommended exhaust air flow rate is sufficient. In addition to explosion protection, other factors play a role, such as the substance to be stored etc.

2.8.2 Passive storage

If the safety storage cabinet is connected to industrial ventilation and is used exclusively for passive storage, the following measures must be observed:

- The air must be exchanged inside the safety storage cabinet with the doors closed at a rate of at least 10 times the internal volume of the cabinet.
- Use suitable, tightly sealed storage containers that are not wetted with flammable liquid.
- Regularly check storage containers for leaks.

2.8.3 Passive storage without industrial ventilation

If the safety storage cabinet is not connected to industrial ventilation and is used exclusively for passive storage, an explosive atmosphere may form in and around the safety storage cabinet. The operator must prepare a safety concept for this in compliance with the explosion protection regulations.

WARNING

Potentially explosive area on the safety storage cabinet

Risk of fire and explosion due to ignition of explosive mixtures in the safety storage cabinet.

The consequences can be death or serious injury.

- Determine the ex-zone in compliance with the applicable national standards and regulations, e.g. TRGS 722 in Germany, and mark it clearly and permanently.
- Naked flames and smoking in potentially explosive areas are prohibited.
- Do not use any tools that cause mechanically generated sparks.
- Avoid electrostatic charges.
- Avoid equipment with surface temperatures above the ignition temperatures of the flammable liquids stored.
- Only operate electrical equipment in potentially explosive areas if it complies with the requirements of the applicable national standards and regulations, e.g. Technical Rules for Hazardous Substances: TRGS 722 in Germany.
- Do not use any ignition sources in the safety storage cabinet.
- Earth the safety storage cabinet.

If an explosive atmosphere is to be expected, measures must be taken in compliance with the applicable national standards and regulations, e.g. TRGS 722 in Germany, and an explosion protection document must be drawn up.

Safety storage cabinets that are operated without industrial ventilation must be indicated using an appropriate notice.

The potentially explosive area must be marked clearly with the following warning sign in compliance with the applicable national standards and regulations, e.g. the European ATEX Directive “2014/34/EU”:

The size of all symbols and notices should be appropriate for the size of the safety storage cabinet.

2.9 Hazardous areas and their labelling

The following must be attached to the front of the safety storage cabinet and must be clearly visible:

- The notice “Keep door closed at all times”
- Fire resistance in minutes (e.g. “Type 90”)
- Name or trademark of the manufacturer
- Serial number and year of manufacture
- Information on the largest individual container volume that can be stored
- Specification of the shelves’ maximum load capacity

Furthermore, the following signs must be attached to the front of the safety storage cabinet and must be clearly visible:

Table 2: Prohibited action signs

Symbol	Meaning	Standard
	P003: No naked flames; fire, open ignition source and smoking prohibited	DIN EN ISO 7010

Table 3: Warning signs

Symbol	Meaning	Standard
	W012: Warning: Electricity	DIN EN ISO 7010
	W021: Warning: Flammable materials	DIN EN ISO 7010
	D-W012: Warning: Explosive atmosphere (only in zone 2 in the safety storage cabinet)	DIN 4844-2
	Warning signs: Electrostatic discharge (only in zone 2 in safety storage cabinet)	ANSI/ESD S8.1

Table 4: Mandatory action signs

Symbol	Meaning	Standard
	M002: Read the instructions	DIN EN ISO 7010
	Keep doors closed at all times	Manufacturer's instructions
	Connect the safety storage cabinet to a ventilation system.	Manufacturer's instructions

3 Technical specifications

3.1 Outdoors


Fig. 3: General diagram DISPOSAL U pro outside

01	Analysing unit (monitoring)	07	Connection pipe for the disposal system
02	Type plate	08	Exhaust air connection
03	Earthing connection	B	Width
04	Door handle	H	Height
05	Base	HS	Base height
06	Penetration	T	Depth

Table 5: External dimensions

Model	DISPOSAL U pro		DISPOSAL B pro
Size	S	LD	S
B (mm)	601	1100	601
H (mm)	596	596	770
HS (mm)*	35	35	35
T (mm)	594	594	594
Doors	1	2	1
Min. volumetric flow rate	30	30	30
Electrical connection for the analysing unit	110 – 230 V + 50 / 60 Hz	110 – 230 V + 50 / 60 Hz	110 – 230 V + 50 / 60 Hz

*Standard base height

3.2 Interior


Fig. 4: General diagram DISPOSAL U pro inside

- | | | | |
|----|---------------------------------|---|--|
| 01 | Disposal system trigger element | L | Maximum volume of largest individual container |
| 02 | Trigger element | M | Load-bearing capacity per standing surface (uniformly distributed) |
| 03 | Venting cut-off flap | V | Collection volume |
| 04 | Pull-out tray with insert | | |
| 05 | Spare container | | |
| 06 | Safety technology | | |
| 07 | Collection container | | |

Table 6: Dimensions internal

Model	DISPOSAL U pro		DISPOSAL B pro
Size	S	LD	S
L (l)	10	10 10*	20
M (kg)	45	45 45*	45
V (l)	18.7	18.1 18.1*	33.6
Max. storage volume (l)	20	20 20*	40
Max. payload (kg)	45	90	45
Empty weight (kg)	133	225	145
Max. total weight (kg)	178	315	190

*Safety storage cabinet with two (separate) storage areas: left storage area | right storage area

3.3 Pressure drop during ventilation

Industrial ventilation of the safety storage cabinet with a volumetric flow rate Q creates a pressure drop Δp at the exhaust air connection.

A required volumetric flow rate Q of 30 m³/h produces the following pressure drop ΔP .

Table 7: Pressure drop Δp with a volumetric flow rate $Q = 30$ m³/h

Model	DISPOSAL U pro		DISPOSAL B pro
Size	S	LD	S
Δp (Pa) exhaust air	68	80	86
Δp (Pa) cabinet interior	<1	<1	<1

The pressure drop Δp becomes larger as the volumetric flow rate Q increases. The pressure drop at the exhaust air connection, dependent on the volumetric flow rate averaged over all cabinet sizes, is illustrated in the following diagram.


Fig. 5: Average pressure drop

4 Structure and function

4.1 Construction

- Cabinet in multi-layer construction
- Outer casing: powder-coated sheet steel (conductive)*
- Interior surfaces: light grey coated decor panel or powder-coated sheet steel (conductive)*
- Safety components for closure in case of fire: Brass, spring steel (1.410)
- Parts in contact with media: Polytetrafluoroethylene (PTFE), polypropylene (PP), polyethylene (PE), fluorocarbon rubber (FKM) and stainless steel (1.4571 and 1.4408)

*Alternative: Stainless steel sheet in different qualities

4.2 Doors

4.2.1 Door options

The safety storage cabinet has one of the following door operations:

- Door handle
- Push-to-open

4.2.2 Locking cylinder

The door can be locked with the integrated locking cylinder. The key numbers are embossed on the locking cylinder and on the keys supplied, e.g. A007.

Safety storage cabinets with locking systems supplied as an option do not contain a key number.

4.3 Safety technology

4.3.1 Door closure in case of fire

At a temperature of approx. 50°C, open doors are closed automatically by the safety technology.

Safety storage cabinets with pull-out trays are equipped with a closure control system. This prevents the door from being blocked by the pull-out tray in case of fire.

4.3.2 Closure of venting cut-off flaps in case of fire


Fig. 6: Left: Ventilation ducts open; Right: Ventilation ducts closed

- 01 Venting cut-off flap
- 02 70°C trigger element

The safety storage cabinet is equipped with two venting cut-off flaps.

At a temperature of 70°C, the ventilation ducts are closed by the venting cut-off flaps installed in the cabinet.

4.4 Pull-out tray with insert for collection container

The function of the pull-out tray at the bottom of the safety storage cabinet is to collect leaking substances. In the DISPOSAL U and B pro models, it is equipped as standard with special inserts for the collection containers. The pull-out tray and insert are conductively connected to the outer casing.

The centring plate and the positioning plate ensure the correct positioning of the connection head. Pull-out trays are permanently installed in the factory and subsequent modification is not possible.

4.4.1 Version for 10 litre collection container


Fig. 7: Safety storage cabinet with pull-out tray with insert

- | | | | |
|----|-------------------------------------|----|---------------------------------|
| 01 | Centring plate | 04 | Insert for collection container |
| 02 | Equipotential bonding | 05 | Pull-out tray |
| 03 | Positioning plate with spring hinge | | |

4.4.2 Version for 5 litre collection container


Fig. 8: Safety storage cabinet with pull-out tray with insert

- | | | | |
|----|--------------------------|----|---------------------------------|
| 01 | Centring plate | 04 | Insert for collection container |
| 02 | Bracket for level sensor | 05 | Pull-out tray |
| 03 | Equipotential bonding | | |

4.5 Disposal system


Fig. 9: Disposal system

- | | | | |
|----|----------------------|----|--------------------------------|
| 01 | Connection pipe | 03 | Earth line with terminal clamp |
| 02 | Collection container | 04 | Connection head |

The disposal system consists of a collection container which is connected to the connection pipe via a connection head.

This collection container is connected to the cabinet by an earthing cable with terminal clamp.

Depending on the version, there are two or three canisters in the safety storage cabinet. Canisters up to 5 litres made of non-conductive material (polyethylene), and canisters larger than 5 litres made of electrostatically dissipative material (polyethylene).

Spare containers are not connected to the disposal system and must always be stored tightly sealed.

4.5.1 Connection head


Fig. 10: Connection head

01	Ventilation opening	05	Seal
02	Stopcock	06	Plastic ball
03	Screw plug cutting ring fitting	07	Screw plug
04	Cover cap		

The cover cap on the connection head must be securely connected to the collection container. The connection head is equipped with a stopcock and a ventilation opening with protection against overfilling (plastic ball).

The stopcock prevents any residues remaining in the connection pipe from escaping when the collection container is changed.

The connection pipe is connected via the cutting ring fitting (in the factory).

4.5.2 Connection pipe of the disposal system


Fig. 11: Connection pipe

- 01 PTFE hose (electrically conductive)
- 02 Anti-kinking sleeve
- 03 Elbow screw-in joint

The connection pipe consists of a PTFE hose reinforced with an anti-kinking sleeve and an elbow screw-in joint.

The connection pipe is attached with the elbow screw-in joint to a disposal line, from an HPLC system or a funnel for example.

4.6 Analysing unit (monitoring)


Fig. 12: Monitoring

- 01 Analysing unit display
- 02 USB-C Extension Flex for level monitoring (connection box, connecting cable with USB-C connector and I/O level sensor)

The analysing unit as supplied and connected serves to monitor the fill level in the collection container provided for this purpose. This consists of a display and one (DISPOSAL U and B pro S) or two (DISPOSAL U pro LD) USB-C Extension Flex units with a connected capacitive I/O level sensor.

The capacitive level sensor is suitable for use in an ex-zone 2.

The connection box and the connecting cable of the USB-C Extension Flex must not be located in the ex-zone.

In the event of an alarm or fault, the user is warned visually.

Observe the separate operating instructions for the analysing unit.

Optionally, 2 (DISPOSAL U pro LD) or 3 (DISPOSAL U and B pro S) additional sensors can be connected to the analysing unit. ↪ chapter 4.8.3 “Extending the analysing unit”

4.7 Penetration


Fig. 13 Penetration for hose and cables

The DISPOSAL U pro and DISPOSAL B pro type safety storage cabinets are equipped on the rear wall with one or two permanently installed pipe penetrations. The sensor cable and connection pipe are routed through the Type 90 type-tested penetrations into the safety storage cabinet, so that the holes do not impair the fire resistance.

NOTE

Subsequent installation of additional penetrations

Incorrect installation in the wrong location can cause damage to the safety storage cabinet and impair the fire resistance.

- Do not make any holes or penetrations before consulting with DÜPERTHAL.

NOTE

Routing and use of pipes, cables and hoses

The routing and use of pipes, cables and hoses are the responsibility of the operator.

Incorrect handling can cause damage to the safety storage cabinet.

- Carry out a separate risk assessment of the overall setup.

4.8 Accessories (optional)

Optionally, the DISPOSAL U pro and DISPOSAL B pro safety storage cabinets can be equipped with various upgrades.

4.8.1 Variable base designs


Fig. 14: Base designs

- 01 Outer base
- 02 Classic base
- 03 Toe-kick base

The safety storage cabinet is equipped with a classic 35 mm high base as standard. Optionally, the height and design of the of the base can be modified.

The total height of the safety storage cabinet varies depending on the design of the base.

The following versions are possible:

Outer base

In this case, there is no base and the safety storage cabinet stands directly on an outer base with an additional height of 3 mm.

Classic base

The classic base is adapted to the cabinet in terms of width and depth, and is provided with ventilation grilles on the front.

Height may be varied between 20 mm and 135 mm.

Toe-kick base

The toe-kick base is adapted to the cabinet in terms of width and depth. The front is set back by about 50 mm and is provided with ventilation grilles. This enables an ergonomic standing position when working on the bench or in the fume hood.

Height may be varied between 60 mm and 135 mm.

4.8.2 Roller set

The safety storage cabinet can optionally be fitted with a roller set. This enables the safety storage cabinet to be moved easily to its location and also to be pulled out for occasional cleaning or maintenance work.

The rollers are not designed for continuous mobile use of the safety storage cabinet.

Installing the roller set increases the cabinet height by 8 mm.

NOTE

Damaged floors

The rollers can leave permanent grooves on soft floors.

- Ensure appropriate floor quality.

4.8.3 Extending the analysing unit


Fig. 15: USB-C Extension Flex

- 01 USB-C Extension Flex for door status monitoring
- 02 USB-C Extension Flex for exhaust air monitoring

The analysing unit can optionally be extended with a USB-C Extension Flex for monitoring door and exhaust air status.

The USB-C Extension Flex units required for this can either be factory-fitted or retrofitted.

5 Transport and packaging

The safety storage cabinet is packaged for transportation and is protected against damage by transport restraints made of coloured plastic (white, yellow, red). The transport restraints should be refitted before any transportation.

For transportation, the safety storage cabinet can be placed on a transport pallet or other suitable transport equipment. It must be ensured that the safety storage cabinet is loaded as smoothly as possible.

WARNING

Crush hazard due to safety storage cabinet tipping over.

If the safety storage cabinet tips over when not transported with due caution, this can cause potentially fatal crushing.

- Wear personal protective equipment (PPE).
- Transport with at least two people.
- Only transport the safety storage cabinet upright and unladen.
- Only transport the safety storage cabinet using suitable transport equipment.
- Only lift the safety storage cabinet once it has been secured and there is no risk of it tipping over or slipping.

NOTE

Handling the transport restraints

Incorrect handling causes damage to the safety storage cabinet.

- The transport pallet should only be removed at the installation location.
- Do not remove the transport restraints until the safety storage cabinet has been taken off the pallet and is in position at its installation location.

6 Installation and commissioning

NOTE Handling the transport restraints

Incorrect handling causes damage to the safety storage cabinet.

- The transport restraints should only be removed at the installation location.
- Always replace transport restraints before lifting or moving the safety storage cabinet.

Install the safety storage cabinet so that the annual maintenance activities can be carried out without restriction.

6.1 Requirements for the installation location

The safety storage cabinet is only approved for installation in a building.


Fig. 16: Requirements for the installation location

Consider the following in relation to the installation location:

- The surface must be able to bear the weight of the safety storage cabinet when fully loaded.
- The foundation must be horizontal to ensure flawless functioning of the safety storage cabinet.
- The load-bearing capacity and stability of the surface must be assured both in normal situations and in the event of a fire.
- Do not install the safety storage cabinet near sources of heat.
- Protect the safety storage cabinet against moisture.
- At a relative humidity of >70 %, use in closed and heated buildings is permitted for a few weeks each year.
- The operating temperature must be between -5 °C and +40 °C.

6.2 Safety storage cabinet installation and alignment

NOTE

Damage to the safety storage cabinet from water or cleaning agents

Liquids such as cleaning water that get into or under the cabinet can damage the safety storage cabinet.

- Protect the safety storage cabinet on site, e.g., with joint sealant.
- Clean floors in a way that prevents liquids from running under the cabinet.
- Remove larger quantities of liquid immediately.

The alignment procedure described below is used for precision alignment. Major floor irregularities of more than 15 mm must be eliminated on site.


Fig. 17: Safety storage cabinet installation and alignment.

The safety storage cabinet is supplied without alignment elements. An alignment may be necessary in individual cases.

Install and align the safety storage cabinet:

Employees:

- Technically qualified personnel (2 to 3 people)

Tool:

- Suitable tool for lifting the cabinet

- At the installation location, carefully lift the safety storage cabinet off the pallet and set it down smoothly (Fig. 1).
- If the safety storage cabinet is equipped with a roller set, fit the locking screws.
 - ↳ Chapter 6.3“Installation and operation of the locking screw on a roller set”.
- Carefully push the safety storage cabinet to the installation location.
- Check the alignment of the safety storage cabinet with the spirit level (Fig. 2).
- If necessary, raise the safety storage cabinet slightly and place metal spacers underneath it to align the cabinet (image 2).
- Set the safety storage cabinet down carefully.
- If the safety storage cabinet is equipped with a roller set, pretension the locking screw.
 - ↳ Chapter 6.3“Installation and operation of the locking screw on a roller set”.
- Remove the coloured plastic transport restraints from the door gaps.
- Check the alignment. ↳ Chapter 6.4“Checking the alignment of the safety storage cabinet”.

6.3 Installation and operation of the locking screw on a roller set


Fig. 18: Installation and operation of the locking screw on a roller set.

The rollers allow easy installation of the safety storage cabinet under workplaces.

Installation and operation of the locking screw:

Employees:

- Technically qualified personnel

Tool:

- Suitable tool for lifting the cabinet
- Open-ended spanner SW 13 mm

- Raise the safety storage cabinet slightly.
- Fully tighten the accompanying locking screws in front of the left and right rollers.
- Set the safety storage cabinet back down smoothly.
- Push the safety storage cabinet to the required position.
- If necessary, raise the safety storage cabinet slightly and place metal spacers underneath it to align the cabinet.
- After aligning the safety storage cabinet, use the open-ended spanner to turn the locking screws until they reach the floor and pretension them slightly.

The locking screws prevent unintentional slipping of the safety storage cabinet. They are not designed for use as adjusting screws for alignment of the safety storage cabinet.

6.4 Check the alignment of the safety storage cabinet

If the safety storage cabinet is not aligned properly, the open wing doors will automatically close themselves or open fully, ↪ Chapter 6.2“Installing and aligning the safety storage cabinet”.


Fig. 19: Checking alignment

Correct alignment of the safety storage cabinet:

- The door gaps have a uniform width when the doors are closed.
- With two doors, the central gap and ceiling gap form an even “T”.

6.5 Connecting the safety storage cabinet to an exhaust air system

6.5.1 Exhaust air connection and air supply opening


Fig. 20: Exhaust air connection and air supply opening

01 Exhaust air connection 02 Air supply opening

On the DISPOSAL U pro S with door hinge on the left side, the supply air opening and the exhaust air connection are arranged mirror-symmetrically.

The exhaust air connection and air supply opening are located on the rear wall of the safety storage cabinet. The external exhaust air system can be connected to the exhaust air connection; it is designed to ensure the exchange of air and discharge the air from the safety storage cabinet to a safe region in the open air.

The layout of the ventilation ducts in the cabinet interior means that ventilation is effective directly above the bottom tray and the filling point.

6.5.2 Connection to an exhaust air system

The DISPOSAL U pro and DISPOSAL B pro safety storage cabinets are designed for active storage, particularly for disposal, of flammable liquids. Therefore, the safety storage cabinet must be connected to a continuously operating exhaust air system. The safety storage cabinet must be ventilated at a minimum of 30 m³/h. The ventilation must be monitored continuously. ↪ Chapter 2.8 “Safety concept for explosion protection”

Installation of industrial ventilation and connection to an existing exhaust air system must be carried out by a qualified company and is not a service provided by DÜPERTHAL.

The power of the exhaust air system can be determined using the technical specifications. ↪ Chapter 3.3 “Pressure drop during ventilation”

Connection to the exhaust air system:

Employees:

- Technically qualified personnel
- Connect the exhaust air line to the exhaust air connection socket and secure with a collar.
- Set the volumetric flow rate to at least 30 m³/h.
- After installing the safety storage cabinet, check the connection to the exhaust air system with smoke tubes.

6.6 Earthing the safety storage cabinet


Fig. 21: Earthing the safety storage cabinet

01 Earthing connection

Earthing the safety storage cabinet prevents ignition hazards.

The interior fittings are conductively connected to one another and to the outer casing. The connection for earthing is located on both sides on the rear of the cabinet.

Correct earthing is described in the applicable national standards and regulations, e.g. TRGS 727 in Germany.

Earthing connection:

Employees:

- Technically qualified personnel

Tool:

- Wide slotted screwdriver

- ➔ Bend open the earthing jumper with a slotted screwdriver.
- ➔ Connect the earthing cable (not included in scope of delivery).

6.7 Connecting the connection pipe


Fig. 22: Connection pipe with elbow screw-in joint

01	Connection angle	04	Screw cap
02	Reinforcing sleeve	05	PTFE hose (el)
03	Cutting ring	06	Anti-kinking sleeve

The PTFE hose (05) is fed into the interior at the rear of the safety storage cabinet through a penetration installed on the left-hand or right-hand side of the cabinet.

When delivered, the connection pipe is wound up and attached to the rear of the cabinet.

The elbow screw-in joint is supplied with the safety storage cabinet. It consists of a connection angle (01), a cutting ring (03) and a screw cap (04).


Fig. 23: Connecting the connection pipe

Installing the connection pipe:

Employees:

- Technically qualified personnel

Tool:

- Open-ended spanner SW 27 mm
- Open-ended spanner SW 22 mm
- Cutting pliers for plastic pipes

⚠ WARNING

Formation of an explosive atmosphere due to escaping flammable liquids

Incorrect installation or damage may allow flammable liquids to escape from leaking connection pipes.

The consequences may be death or serious injury due to explosion.

- Do not twist PTFE hose.
- To prevent kinks in the PTFE hose, do not exceed the minimum bending radius (121 mm) during installation.
- To protect the PTFE hose, use the anti-kinking sleeve.
- Tighten the cutting ring fitting correctly (1 ¼ turns).
- Cutting edge should always face the connection cone.
- Before commissioning, check the system for leaks using a suitable, non-hazardous liquid, e.g. water.

- Seal the supplied connection angle (G1/2") (01), e.g., with PTFE sealing tape, and connect to the external connection point.
- Completely unwind the PTFE hose (05) and place at the connection point.
- Cut the PTFE hose to the required length.
- Fit the anti-kinking sleeve (06) over the PTFE hose.
- Slide the reinforcing sleeve (02) into the PTFE hose (Fig. 1).
- Pass the PTFE hose through locking nut (04) and cutting ring (03).
- Observe the correct installation position of the cutting ring (Fig. 2).
- Slide the PTFE hose as far as possible into the open end of the connection angle.
- Slide cutting ring and locking nut forward and tighten by hand.
- Make sure that the PTFE hose does not slip out again.
- Securely tighten the locking nut by 1 ¼ turns using two open-ended wrenches (Fig. 3).

6.8 Commissioning the analysing unit

Commissioning the analysing unit:

Employees:

- Technically qualified personnel
- Insert suitable backup batteries (not included in the scope of delivery) into the battery compartment of the analysing unit.
- Connect the power supply unit to the analysing unit.
- Connect the power supply unit to the mains power supply.
- ✓ LED at the bottom right lights up green.
- Connect the USB-C Extension Flex for level monitoring to one of the connection sockets for USB-C Extension (Entry 1-4) on the analysing unit.
- In the operating and display menu, switch from the Device menu to the Flex display menu via the menu ribbon.
- ✓ The fill level status is displayed.

Observe the separate operating instructions for the analysing unit.

6.9 Calibrating the USB-C Extension Flex for level monitoring

The USB-C Extension Flex for level monitoring has already been calibrated by DÜPERTHAL with a water container prior to delivery. For fine calibration or in order to adjust the fill level limit, it must be re-calibrated with the substances that are to be removed and the required fill capacity.

Observe the separate operating instructions for the analysing unit, the Vi Hub or the level sensor.

6.9.1 Collection container > 5 litres

Calibrating the USB-C Extension Flex for level monitoring:

Employees:

- Technically qualified personnel

Tool:

- Small slotted screwdriver
- Fill the collection container with the required filling quantity of the substance to be removed, and place it in the safety storage cabinet.
- Place the USB-C Extension Flex for level monitoring into the designated bracket on the collection container.
- Open the cover flap of the potentiometer screw on the top of the USB-C Extension Flex for level monitoring.
- Adjust the sensitivity of the sensor with the potentiometer screw.
- If the LED on the USB-C Extension Flex for level monitoring does not light up, turn the potentiometer screw to the right using a small slotted screwdriver until the sensor switches.
 - ✓ LED on the sensor lights up.
- Then continue turning the potentiometer screw another $\frac{1}{4}$ turn to the right to ensure that the sensor responds reliably.
- When the LED on the sensor lights up, turn the potentiometer screw to the left until the sensor switches again.
 - ✓ LED on the sensor is no longer lit.
- Then turn the potentiometer screw $\frac{1}{4}$ turn to the right to ensure that the sensor responds reliably.
 - ✓ LED on the sensor lights up.
- Close flap again.

6.9.2 Collection container = 5 litres


Fig. 24: Calibrating the capacitive level sensor for level monitoring

Calibrating the capacitive level sensor.

Employees:

- Technically qualified personnel

Tool:

- Small slotted screwdriver

- Fill the collection container with the required filling quantity of the substance to be removed, and place it in the safety storage cabinet.
- Position the capacitive level sensor with the two fixing nuts in the bracket so that the sensor can switch safely (approx. 3 mm from the wall of the collection container and centred on the container fill level).
- Open the cover flap of the potentiometer screw on the top of the capacitive level sensor.
- The sensitivity of the sensor is set via the potentiometer screw.

LED on sensor does not light up:

- Turn the potentiometer screw to the right using a small slotted screwdriver until the sensor switches.
 - ✓ LED on the sensor lights up.
- Then continue turning the potentiometer screw another $\frac{1}{4}$ turn to the right to ensure that the sensor responds reliably.
- Close the flap.

LED on sensor lights up:

- Turn the potentiometer screw to the left until the sensor switches again.
 - ✓ LED on the sensor is no longer lit.
- Then turn the potentiometer screw $\frac{1}{4}$ turn to the right to ensure that the sensor responds reliably.
 - ✓ LED on sensor lights up again.
- Close the flap.

7 Operation

WARNING

Formation of an explosive atmosphere

Failure or impairment of the extraction can result in an explosive vapour-air mixture in the cabinet interior.

The consequences can be death or serious injury.

- Stop filling activities
- Have the malfunction corrected by technically qualified personnel.

WARNING

Risk of injury due to spark formation

Sparks can form due to electrostatic charge.

The consequences can be death or serious injury.

- Take precautions to prevent build-up of static electricity when transferring hazardous substances.
- Use conductive containers wherever possible.
- Connect containers to earth terminals.
- Do not bring any charges into the ex-zone on people, e.g. charged/unsuitable clothing.

7.1 Opening the safety storage cabinet

⚠ WARNING

Blocked doors

Door blocked by objects cannot close by themselves in case of fire.

This may result in death or serious injuries as a result of inadequate fire protection.

- Close the doors after every work process.
- Do not block the doors using any objects.

7.1.1 Opening the safety storage cabinet with the door handle

Opening the safety storage cabinet with door handle:

Employees:

- Users

- ➔ Open the safety storage cabinet by pulling on the door handle of the wing door.
- ✓ The wing door remains open in any position

7.1.2 Opening a safety storage cabinet with push-to-open door operation


Fig. 25: Door operation with push-to-open mechanism

Opening the safety storage cabinet with push-to-open door operation:

Employees:

- Users

- ➔ Press lightly on the door surface by hand in the area of the green dot.
- ✓ The door disengages and is pushed open slightly.
- ➔ Pull on the green handle at the top edge of the door to open it fully.
- ➔ To close, gently push the door.
- ✓ The door engages and is fully closed.

7.2 Filling and transferring

Filling process:

Employees:

- Users
- Before beginning the filling process, check the fill level in the collection container, verify that the safety storage cabinet is properly ventilated and the connection head is correctly seated on the collection container.
- Leave the door closed during the filling process.
- Carry out the filling process while observing all safety regulations.
- Interrupt the filling process as soon as the level indicator on the analysing unit switches from black "OK" text to red "NOT OK" text.

7.3 Changing the collection container


Fig. 26: Changing the collection container

Replacing the collection container:

Employees:

- Users
- Provide an absorber to collect drips.
- Open the door of the safety storage cabinet.
- Close the stopcock on the connection head (Fig. 1) to prevent residues from escaping.
- Pull the collection container backwards (10 l collection container: Figure 2.1 and 5 l collection container: Figure 2.2).
 - ✓ The connection head of the collection container is outside the centring plate.
- Loosen the connection head fitting.
- Lift the connection head out of the canister opening (Fig. 3).
- Wipe off any adherent residues and dispose of them properly.
- Seal the collection container tight with the supplied lid.
- Remove earthing cable.
- Remove the level sensor (only for 10 l collection containers).
- Remove the collection container.
- Place the connection head to one side, without kinking the PTFE hose.
- Place the empty collection container with the opening facing the rear wall of the safety storage cabinet.
- Place the connection head on the collection container so that the positioning plate presses the connection head into the centring plate and the connection head can no longer rotate (Figure 4).
- Screw the connection head tightly onto the empty collection container.
- Check that the round seal on the connection head is in perfect condition and properly in place.
- Slide the connection head of the collection container into the centring plate (Figure 5.1 or 5.2).
 - ✓ The positioning plate with spring hinge or the ramp of the insert plate continuously holds the connection head in the correct position in the centring plate.
- Connect earthing cable to the canister.
- Reinsert the level sensor into the bracket of the empty collection container (only for 10 l collection containers).
- Open stopcock.
- Remove the full collection container, or store it in the safety storage cabinet until it can be picked up.
- Close cabinet door.

8 Procedure after a fire

WARNING

Explosive vapour / air mixture

The consequences can be death or serious injury.

- Before opening the safety storage cabinet, remove all ignition sources within a 10m radius.
- Only open the safety storage cabinet with tools that do not cause any mechanically generated sparks.

WARNING

Damaged safety storage cabinet due to fire or extinguishing agents

The consequences can be death or serious injury.

- Do not use safety storage cabinets that have been damaged by fire or extinguishing agents.

After a fire, the safety storage cabinet must not be opened for at least 24 hours, and only with the utmost caution and by specialist employees.

In case of fire, the fire and extinguishing agents reduce the fire resistance of the safety storage cabinet to such an extent that it can no longer be used. It must be replaced with a new cabinet.

9 Maintenance

Safety equipment must be tested annually by a qualified person in compliance with BetrSichV and the maintenance interval stipulated by the manufacturer as set out in TRBS 1203.

Check the safety storage cabinet for any externally visible damage or defects.

Always perform checks:

- After installation.
- Before commissioning.
- After changes.
- After maintenance.

Maintenance intervals for individual components may differ from the maintenance intervals for the safety storage cabinet.

Follow the manufacturer's information in the operating instructions for the individual components.

The safety storage cabinet should also be inspected at the intervals listed.

Table 8: Maintenance work daily.

Interval	Maintenance work	Employees
Daily	Pull-out tray <ul style="list-style-type: none"> ▪ Check in accordance with regulations governing water legislation ▪ Collect and properly dispose of leaked liquids immediately. 	Users
	Exhaust air <ul style="list-style-type: none"> ▪ Check the function of the exhaust air system 	Users
	Earthing <ul style="list-style-type: none"> ▪ Check contact points that are created by clamp connections on non-stationary parts such as containers for secure fitting. 	Users
	Hoses <ul style="list-style-type: none"> ▪ Visual inspection of hoses and hose connections for damage and leakage. In case of defects, bring in technically qualified personnel. 	Users

Table 9: Maintenance work monthly.

Interval	Maintenance work	Employees
Monthly	Door function <ul style="list-style-type: none"> Check the ease of movement of the doors by opening and closing them. Check the door lock by locking and unlocking it. 	Technically qualified personnel
	Ventilation <ul style="list-style-type: none"> Check the effectiveness of the ventilation in the cabinet, upstream of the exhaust air duct, at the ventilation openings, using a smoke tube or similar. Remove any contamination and dust from the ventilation opening. 	Technically qualified personnel
	Seals <ul style="list-style-type: none"> Check that the fire protection and fabric seals on the carcass and doors are seated correctly, and in the event of visible damage, contact customer service. ↳ Fig. 27: Seals 	Technically qualified personnel
	Labelling <ul style="list-style-type: none"> Inspect the safety labels on the safety storage cabinet to ensure they are complete. ↳ Chapter 2.9“Hazardous areas and their labelling” 	Technically qualified personnel
	Alignment (push-to-open technology) <ul style="list-style-type: none"> Check ease of movement and engagement/disengagement of the doors. 	Technically qualified personnel
	Connection pipe and connections <ul style="list-style-type: none"> Inspect connection pipe and connections for leaks. Inspect connection pipe for kinks. If the hoses and hose connections are damaged, contact customer service. 	Technically qualified personnel

Table 10: Maintenance work yearly.

Interval	Maintenance work	Employees
Annually	Safety storage cabinet <ul style="list-style-type: none"> Inspection of the entire safety storage cabinet. 	DÜPERTHAL service technicians

Your contact person – Technical Service
 Phone: +49 6188 9139-166
service@dueperthal.com | service.dueperthal.com

For easier cleaning of the pull-out tray / bottom tray
 The shelf can be removed.


Fig. 27: Seals

- 01 Fabric seal
- 02 Fire protection seal around the carcass
- 03 Fabric seal, door (top and bottom)*
- 04 Fire protection seal, door**

* also to the centre gap for safety storage cabinets with two doors

** only for safety storage cabinets with two doors

10 Faults

Follow the manufacturer's information in the operating instructions for the individual components.

Table 11: (Part 1): Overview of possible faults

Fault description	Cause	Remedy	Employees
Doors do not close.	Safety storage cabinet is not aligned correctly.	Install and align the safety storage cabinet in a horizontal position. ↳ Chapter 6.4 "Checking the alignment of the safety storage cabinet".	Technically qualified personnel
	Doors are held open by objects.	Do not wedge or hold doors open with any objects.	Technically qualified personnel
	Safety storage cabinet is not correctly filled.	Make sure that products stored in the safety storage cabinet are uniformly distributed.	Technically qualified personnel
No extraction system available.	Venting cut-off flaps closed, as locking mechanism has been triggered	Replace the locking mechanism.	DÜPERTHAL service technicians
Doors do not move easily.	Moving parts, such as hinges, are dirty or corroded.	<ul style="list-style-type: none"> ▪ Remove rust. ▪ Lubricate parts. ▪ Remove corrosive substances from the safety storage cabinet. ▪ Notify technical customer service. 	Technically qualified personnel
Doors open again after closing and close again after opening.	Safety storage cabinet is not aligned correctly.	Install and align the safety storage cabinet in a horizontal position. ↳ Chapter 6.4 "Checking the alignment of the safety storage cabinet"	Technically qualified personnel

Table 11: (Part 2): Overview of possible faults

Fault description	Cause	Remedy	Employees
Doors cannot be opened	The “Door closure in case of fire” safety technology has been triggered and holds the door closed.	Open the door carefully and repair or replace the safety technology.	DÜPERTHAL service technicians
Liquids do not drain correctly.	Connection pipe blocked	<ul style="list-style-type: none"> ▪ Clean connection pipe ▪ Notify technical customer service 	Technically qualified personnel
	Ventilation opening clogged	<ul style="list-style-type: none"> ▪ Clean the ventilation opening ▪ Allow ventilation opening to dry 	Technically qualified personnel

If faults occur, assist the technical customer service by providing the cabinet model, production and key number, along with a description of the fault.

11 Shutting down the safety storage cabinet

The following items must have been carried out before shutting down a safety storage cabinet:

- Remove the stored goods.
- Remove any remaining residues of the stored goods.
- De-energise the analysing unit.
- Clean the cabinet.
- Close the doors and install transport restraints in the door gaps.
- A sign must be placed on the door indicating that the safety storage cabinet has been shut down.
- Pack up the cabinet before storage.

The safety storage cabinet must undergo a full maintenance cycle before it is recommissioned.

12 Spare parts and accessories

Only original parts from DÜPERTHAL are to be used for the safety storage cabinets.

- Pull-out tray
- Anti-slip mat made from rubber
- Door handle
- Perforated sheet insert
- Base
- Ventilators
- Exhaust air monitoring units
- Roller set
- Sensors for door monitoring
- Penetrations
- VISICON disposal system

13 Disposal

CAUTION

Dismantling the safety storage cabinet

Risk of injury due to improper dismantling of the safety storage cabinet.

- Ensure that the safety storage cabinet is dismantled by technically qualified personnel only.
- Use appropriate tools and wear personal protective equipment when dismantling the cabinet.

The safety storage cabinet can be completely dismantled by technically qualified personnel.

Recycle the individual material components separately.

Comply with national and local disposal regulations.

To save resources, do not place parts of the safety storage cabinet or the whole cabinet in bulky or domestic waste.

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hereby declares that the following products

Product designation: Safety storage cabinet type 90
Model designation: DISPOSAL U pro, DISPOSAL U standard
Model size: S, LD, XL

comply with all relevant provisions of the directives applied. This declaration refers only to the product in the condition as supplied to the customer. Parts and/or modifications subsequently fitted by the end user are not considered.

The following directives were applied:

Machinery Directive 2006/42/EG

The following harmonised standards were applied:

DIN EN ISO 1210:2011 Safety of machinery - General principles for design - Risk assessment and risk reduction

The following national standards and technical specifications have been applied:

DIN EN 14470-1:2004 Fire safety storage cabinets - Part 1: Safety storage cabinets for flammable liquids

DIN EN 16121:2017 Non-domestic storage furniture - Requirements for safety, strength, durability and stability

DIN EN 16122:2012 Domestic and non-domestic storage furniture - Test methods for the determination of strength, durability and stability

The person authorised to compile the technical documentation is: Mr Frank Backhaus, address see above.

Karlstein, 19.07.2024
(Location, date)

Signature
Franz-Josef Hagen / Managing Director

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hereby declares that the following products

Product designation: Safety storage cabinet type 90
Model designation: DISPOSAL B pro, DISPOSAL B standard
Model size: S

comply with all relevant provisions of the directives applied. This declaration refers only to the product in the condition as supplied to the customer. Parts and/or modifications subsequently fitted by the end user are not considered.

The following directives were applied:

Machinery Directive 2006/42/EG

The following harmonised standards were applied:

DIN EN ISO 1210:2011 Safety of machinery - General principles for design - Risk assessment and risk reduction

The following national standards and technical specifications have been applied:

DIN EN 14470-1:2004 Fire safety storage cabinets - Part 1: Safety storage cabinets for flammable liquids

DIN EN 16121:2017 Non-domestic storage furniture - Requirements for safety, strength, durability and stability

DIN EN 16122:2012 Domestic and non-domestic storage furniture - Test methods for the determination of strength, durability and stability

The person authorised to compile the technical documentation is: Mr Frank Backhaus, address see above.

Karlstein, 19.07.2024
(Location, date)

Signature
Franz-Josef Hagen / Managing Director

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Phone/WhatsApp +49 6188 9139-0

info@dueperthal.com | dueperthal.com

Your contact person – Technical Service

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Reproduction and duplication, including extracts, is not permitted. Violations will be prosecuted under civil law. We reserve the right to make technical modifications as a result of improvements and ongoing development, and to implement updates due to amended standards. We accept no liability for printing errors.

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DÜPERTHAL Sicherheitstechnik GmbH & Co. KG
DÜPERTHAL Sicherheitstechnik GmbH & Co. KG | gf | 55339_Rev00_EN