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User's Manual

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Contents

1	Introduction	5
1.1	Explaining Icons Used.....	5
1.2	Intended Use	6
1.3	Occupational Safety	6
2	Installing the Unit and Understanding its Functions	7
2.1	Installing and Checking out	7
2.2	Monitoring Unit Components.....	8
2.3	Understanding Operation, Switching on and off.....	8
2.4	Setting the Pressure Difference	9
3	Operation	9
3.1	Normal Operation and Error Alert.....	10
3.2	Regularly Checking/Replacing the Filters and Documenting It	10
3.3	Potential-free ("Dry") Contact	10
3.4	Status and Warning Messages by the LEDs of the Monitor.....	12
3.5	Online Monitoring over the Internet with DÜPERTHAL connect	12
3.6	Understanding Filter Operation and the User's Obligations	12
3.6.1	Stage 1 Filter	13
3.6.2	Servicing/Replacing the Particulate Filter	13
3.6.3	Particulate Filter Disposal	13
3.6.4	Technical Specifications of the Particulate Filter	13
3.6.5	Stage 2 Filter	14
3.6.6	Servicing/Replacing the Activated Carbon Filter	14
3.6.7	Activated Carbon Filter Disposal	14
3.6.8	Technical Specifications of the Activated Carbon Filter	14
3.6.9	Checking the Filters (every 3 months)	15
3.6.10	Replacing the Filters (every year)	16
3.6.11	Alert to Filter Service/Replacement by the Exhaust Air Monitor	17
4	Exhaust Air Monitor with Filter Technical Specifications.....	17
5	Troubleshooting.....	17
6	Contact Details	18
	Appendix: Substances with which the Product is Compatible.....	21

1 Introduction

Dear Customer:

Thank you very much for purchasing our product. We are confident that it will meet your expectations and perform satisfactorily in your daily routines.

The development, manufacture, and check of the ATEX-compliant exhaust airflow monitoring unit with fan and activated carbon and particulate filters has been guided by our concern to offer a device characterized by superior operation safety and user-friendliness. The unit was manufactured and checked according to state-of-the-art technology and widely recognized safety regulations.

1.1 Explaining Icons Used

The following expressions and icons are used in this User's Manual to refer to safety-related issues:

Caution

This icon indicates that non-compliance may result in a hazardous situation that possibly causes bodily injury or material damage.

Important

This icon indicates important messages related to the correct use of the product. Non-compliance may cause problems of operation and damage to the product.

Info

This icon indicates tips and useful information for using the product more efficiently and more economically.

Additionally, the following warning, prohibition, and mandatory signs provided for by ISO 7010 are used on the product with the meaning defined therein. Some of these signs are applicable not so much to the product itself as to its operation in an environment involving hazardous substances.

Electricity Hazard

Access restricted, qualified electrotechnicians only

Flammable Material Hazard

No Open Flame;

fire, open sources of ignition, and smoking prohibited

Wear protective gloves

Wear eye protection

Refer to instruction manual

Where used, the icon as per ISO 7010.

Caution

shall be considered equivalent to

1.2 Intended Use

The air extraction monitor with fan and activated carbon filter is used to extract air from safe-storage cabinets and chemical storage cabinets equipped with an air extraction connector. Two-stage filtering decontaminates the extracted air and releases it to the ambient air thereafter.

The monitor can be used for the extraction of a large number of chemicals. The repository in the appendix of this manual lists the substances for which appropriate filtering performance can be warranted. If substances not included in the repository are meant to be filtered, the filtering performance has to be confirmed by DÜPERTHAL prior to the start of operation. Feel free to consult with us on any questions related to product use.

From the model version capable of radio communication (P/N. 2.00.398), a wireless connection can be established to DÜPERTHAL. This requires a Bridge to be present as a gateway in the radio-covered area.

The product has to be installed outside of areas at risk of formation of explosive atmospheres. It shall be used exclusively for the positive extraction of air from safety cabinets of which the interior has been defined as Ex zone 2 at the maximum. Flammable substances only of the temperature classes T1 - T4 and also of the gas or explosion groups IIA and IIB shall be extracted. Refer to the respective Material Safety Data Sheet for information on gas group and temperature class.

The product is suitable only for fluid flows not transporting any particles or foreign bodies. When the unit is used, no EX atmosphere must form in any event.

Only such substances are acceptable for extraction that do not compromise the PPS connector sleeve, the PVC suction hose, the filters, and the fan. The relevant list of substances can be provided by the manufacturer in a digital format.

Any and all uses other than the above-described is prohibited and may cause hazardous situations, potentially danger of death, damage to the system and material property, and will cancel any and all warranty and liability for the product on the part of its manufacturer and distributor.

In particular, the use of the unit is prohibited for substances of which extracted air must not be released to the work area, e.g., cancer-producing, germ cell-mutagenic, and reproduction-toxic substances.

The unit allows the use of various types of activated carbon filter. Be sure that the installed filter is the proper one for the intended purpose **per the substance list** in the appendix.

!!! Filter without special mark = filter for AX/A B E K HG

!!! Filter with yellow mark = special filter for specific acids

1.3 Occupational Safety

The installation and operation of the product must be performed in compliance with any and all nationally applicable rules and regulations for the storage of hazardous substances, the ISO 13857 standard, and any and all other applicable codes and standards.

Prior to operating the product, the user is required to perform a risk assessment based on the substances that are going to be stored. **If the extracted air is released to the work area, the applicable occupational exposure limits must never be exceeded.** Current information

on occupational exposure limits can be found in the relevant Material Safety Data Sheet, in the databases of substances GESTIS (Germany) and ECHA (European Chemicals Agency), and Germany's technical rules on hazardous substances (series TRGS 900). Applicable national and international exposure limits shall be observed.

If cancer-producing, germ cell-mutagenic or reproduction-toxic hazardous substances of Categories 1A or 1B are stored in the cabinet to have its air extracted, this air must not be released to the work area under para. 10 of Germany's ordinance on hazardous substances. Internationally, similar provisions are contained in national laws and requirements.

Make sure you wear appropriate personal protective gear when replacing the filter. The kind and extent of required personal protective gear depend on the substances being filtered. Refer to the applicable Material Safety Data Sheet for information on this.

During the filter replacement process, the filtered substance may be released partially. Therefore, perform the replacement at a well ventilated location, e. g., the laboratory's extraction flue.

Vorsicht

As the filtered substances may partially attach to the filter body and may be released during filter replacement, the filter shall be handled in all instances as are the filtered substances.

The design of the product is such that the fan blades are off-touch. Modifications to the enclosure or the suction hose that may make the fan blades accessible to the touch when the unit is operating are prohibited.

Caution

The product must never be used if not completely installed, with properly set suction hose.

Electrical work on the product must only be performed when the air extraction monitor was previously disconnected from power and / or locked against unwanted reapplication of power. This also holds for electrical work on the extraction system of which the monitor may be a part.

Info

Keep in mind that along with the quoted German and international requirements, other regulations related to the design of the workplace and to substances may also apply, of which, then, the provisions have to be respected.


2 Installing the Unit and Understanding its Functions

2.1 Installing and Checking out

1. If needed or desired, connect the potential-free contact (so-called "dry" contact) to the terminal strip (see Sect. 3.3).
2. Connect the unit to the storage cabinet to be monitored:
Connect the exhaust nipple of the cabinet to the unit, e.g., with a flexible sleeve fitting, making sure no strain is imposed on the connection, and make a tight connection using an appropriate hose clamp. This is indispensable to ensure proper operation of the extraction fan.
3. Connect the exhaust airflow monitor to power. Take note of the behavior of the monitor when it is being switched on (see Sect. 2.3).
4. In some cases, it may be necessary to adapt the monitor to the situation found on site by setting the pressure difference at the relevant meter (see Sect. 2.4).

5. Register, if you so desire, the unit with the DÜPERTHAL connect system for monitoring the operation online over the internet (possible only with the model version capable of radio communication, P/N. 2.00.398; see Sect. 3.5).

2.2 Monitoring Unit Components

2.3 Understanding Operation, Switching on and off

To switch the monitor on, simply press the On/Off switch. To switch it off, this switch has to be held down for at least 4 sec. After the monitor was switched on, the fan located inside the enclosure draws air from the monitored system through the suction hose connected to the inlet nipple and extracts the air from the system via a filter setup through the output of the enclosure (see Figures in Sect. 2.2).

If an error condition develops for any reason, airflow through the line is restricted to a greater or lesser extent, and this causes the fluid flow velocity to drop. Then, an audible alarm and a visual alarm are triggered.

If an error condition occurs, first the problem has to be cleared to re-enable proper airflow through the exhaust line. Then certain errors and their correction have to be acknowledged using the On/Off switch (see Sect. 3.4).

2.4 Setting the Pressure Difference

At a metering orifice, the suction line from the storage cabinet is equipped with a pitot probe. The pressure difference it senses is measured by a differential pressure meter via the two transparent tubes routed to it. It rises as the airflow through the suction-side extraction line on site increases. With restricted or blocked airflow, the pressure difference sensed by the probe drops. If it drops below the threshold set at the differential pressure meter, a signal that an error condition has developed is sent via the cable to the monitoring system for further processing. For technical reasons (hysteresis), the threshold from which the monitor cancels the error condition when the pressure difference rises again, is a little higher than the threshold that triggers an error. Therefore, make sure at all times that the settings for the pressure difference and the fan speed of the air extraction system on site are matched to each other in such a way that both an error is indicated when a restriction is detected, and the error is canceled by the unit when the restriction is cleared. There is no upper threshold for a maximum pressure difference.

The value of the pressure difference threshold can be set, after opening the maintenance lid of the compartment with the pressure meter (refer to Sect. 2.2) from 20 Pa to 300 Pa at this meter. This is done by simply removing the transparent cover and using a screwdriver in the provided slot to turn the selector. The factory-set default threshold of the pressure difference is 40 Pa.

Setting the pressure difference threshold for monitoring the airflow requires the following steps:

1. Connect the exhaust air monitor properly to the extraction system (outlet nipple) on site.
2. Set the fan speed to the desired extraction airflow (see Sect. 3.2.5).



Caution

In a cabinet for flammable liquids, the extraction airflow has to be sufficient to replace the interior volume of the cabinet 10 times every hour.

3. Loosen the plastic cover of the differential pressure meter on the suction side (refer to the photographs to the right).
4. Using an appropriate screwdriver in the turning selector at the center of the differential pressure meter, set the threshold (reduce pressure difference) to make the status LEDs indicate normal operation without errors (refer to Sect. 3.4). If this is the case at the start, the pressure difference setting may first be increased until an error is indicated, then slowly reduced again. The higher the pressure difference threshold is set, the quicker / more sensitive the meter will respond by throwing an error message.
5. Using appropriate means on site (e.g. by covering the exhaust outlet with your hand), interrupt the extraction airflow temporarily.
6. If set correctly, the unit will trigger an alarm indicated by the status LEDs. After the airflow interruption is cleared, the unit should return to normal operation.
7. Reinstall the plastic cover on the differential pressure meter.

Point for screw-fastening the cover lid of the compartment for the meter
 (one of four points in total)

Selector for setting the pressure difference threshold

3 Operation

Caution

Be sure that all provisions of Sect. 1.1 and 1.2 above are strictly complied with for both the installation and the operation. Failure to do so may expose you to hazardous situations.


Immediately after switching the monitor on, all warning and status LEDs remain lit for 30 sec for verification.

3.1 Normal Operation and Error Alert

When it works normally, the unit lights the green status LEDs of the fan and the filter (see on the right). This indicates that both components are in normal, flawless operation.

In case of an error condition, the reason for it is indicated by a defined combination of sounds and behavior of the various LEDs (see Sect. 3.4 for reference).

After an error was cleared, the unit resumes normal operation after the fan “start-up” delay of 2 minutes.

3.2 Regularly Checking/Replacing the Filters and Documenting It

Caution

The filters used in the exhaust air monitoring unit have to be checked for proper state of wear at least once every quarter and replaced at least once every year. The self-monitoring functionalities of the unit support the user in this task. Failure to service the filters properly and as scheduled puts material and human resources at risk of severe damage and harm.

Caution

The filters have to be replaced immediately

- if a filter is broken;
- if a large spill happened (e.g., leaking container);
- if a maximum threshold was exceeded.



Important

Whenever the filters are checked or replaced, the fan has to be checked too. If it is dirty, damaged or shows signs of corrosion, it has to be cleaned or the defects have to be corrected, if needed by returning the unit for repair to the seller or manufacturer.

The filters have to be checked and/or replaced at determined intervals (see the appropriate warning indications by the LEDs in Sect. 3.4). During filter service the unit has to remain switched on. After the filter was reinstalled, this has in general to be confirmed on on the Service and Replacement Report (see the appendix) along with the date, description, operator's short identifier, and signature. The exact procedure and documentation is defined in the operating organization's rules for staff procedures.

3.3 Potential-free (“Dry”) Contact

Certain error conditions (see Sect. 3.4) cause a normally open, potential-free contact to be closed; for individually locating errors, their closure can be recorded at the switching panel from the backside of the enclosure.

Acknowledging the error and its correction and resetting the contact is possible when the status and warning LEDs exhibit their acknowledgment behavior (see Sect. 3.4), by pressing the On/Off switch twice.

No.	Event	Pri- ority	Fan LED Stat.	Fan LED Warn.	Filter LED Stat.	Filter LED Warn.	"Dry" Contact Warn.	Contact Status	Alarm Sound	Pushbutton Fct	Fan Status	Remark
1	30 sec wait when switching on, LED check		green lit	red lit	green lit	red lit	red lit				off	After fan start, wait another 2 min before relying on any error indications
2	Normal operation		green lit		green lit						on	
3	Excessive gas concentration	1	green lit		red lit	red lit	red lit	closed	4 x short every 30 sec	No sound for 1 hr after 1 x keypress*	off	After correction of issue acknowledgment/reset of contact needed; see at no. 11
4	Internal clock defective	2		red flash'g		red flash'g	red flash'g	closed	2 x short every 1 min		off	
5	Filter replacement overdue	3	green flash'g			red lit	red lit	closed	2 x long every 1 min		off	Error to block operation of unit when filter replacement overdue (error no. 10) for more than 1 month
6	No activated carbon filter installed	4				red lit		closed			off	After correction of issue acknowledgment/reset of contact needed; see at no. 11
7	RFID chip invalid	5		red lit		red lit					off	After correction (valid chip) unit resumes normal operation (no. 2) automatically
8	Insufficient extraction airflow velocity	6		red lit	green lit		red lit	closed	1 x short every 1 min	No sound for 24 hrs after 1 x keypress*	on	After correction of issue acknowledgment/reset of contact needed; see at no. 11
9	Filter replacement due (soon)	7	green flash'g			red flash'g			2 x short every 1 hour	No sound for 24 hrs after 1 x keypress*	on	Flashing to indicate error from 6 weeks before re-placement due date; alarm beep from 4 weeks before
10	Filter check/inspection due	8	green lit		green flash'g				1 x long every 1 hour	No sound for 24 hrs after 1 x keypress*	on	Alert at the end of every quarter after the first switching-on
11	"Dry" contact closed after error, awaiting acknowledgment		green lit		green lit		red flash'g			LED for contact off after 2 x keypress*	on	After acknowledgment unit resumes normal operation (no. 2)

* The key meant to be pressed is the On/Off switch

Important

The potential-free circuit is designed for **max. 24 V, 2 A**. If these values were exceeded, the unit may be irrevocably destroyed. Additionally, any such damage would not be covered by the product warranty.

3.4 Status and Warning Messages by the LEDs of the Monitor

The green status and red warning LEDs on the front panel of the monitor inform the user of its operation and alert him or her to possible issues and actions that might be needed. The indication has to be interpreted as a combination of LED behavior, status of the fan, and user action that are possible at that moment.

Refer to the table on the previous page to learn what these combined indicators are about.

3.5 Online Monitoring over the Internet with DÜPERTHAL connect

Using DÜPERTHAL connect, you can follow online the operating status and any malfunctions or alerts of the monitoring and air extraction units. A connection to DÜPERTHAL connect is possible only with units that are capable of radio communication (P/N 2.00.398).

Furthermore, a Bridge is required as the communication gateway in order to connect the unit to DÜPERTHAL connect. When a Bridge is within reach, the unit can be registered with DÜPERTHAL connect immediately. To do so, perform the following steps:


1. On the DÜPERTHAL connect home page, access the "Monitoring" menu item.
2. Use the "Add monitoring" button to start the installation process.
3. As the first step, select the cabinet in question.
4. Then, select the appropriate hub in the selection listing or add a new one with the help of the "Create new hub" button.
5. A new hub can be added by scanning the QR code present on the name plate tag of the unit itself or by entering its SN (serial number) and MAC address indicated on the name plate manually.
6. Select the desired functionality of the installation site. The data pertaining to the selected functionality will then be transmitted to the DÜPERTHAL connect online monitoring and will be displayed at the installation site.

Info

If you enter the SN and the MAC address of the unit manually in DÜPERTHAL connect, separate both by a semicolon ";" (in the case of the example in the margin, the required input is: 25839999;D8:47:8F:03:41:DA).

Installation
QR Code

3.6 Understanding Filter Operation and the User's Obligations

The filter enclosure can be accessed through the lateral panel of the unit at the side to which the inlet tube connects. A two-staged filter system is installed in this enclosure.

For servicing the particulate and activated carbon filters, observe in all cases the safety messages in Sect. 1 above, in particular:

Caution

Wear appropriate personal protective gear when replacing the filters, e.g., gloves or safety goggles. The kind and extent of required gear depend on the substances being filtered; the relevant Material Safety Data Sheet contains this information.

Caution

Replacing the filters at a well ventilated location, e. g., the flue, is recommended, as small quantities of the filtered substances may be released.

3.6.1 Stage 1 Filter

Stage 1 of the filter system is a particulate filter that comes under Class F5 pursuant to DIN EN 779. It is able to filter particles of 1 µm - 10 µm. The particulate filter works as a pre-filter pursuant to DIN 12927 for a main filter (Stage 2 Filter). The effective dust-carrying capacity of the filter, pursuant to ISO 12103-1 A4 (filter surface x dust-carrying capacity), amounts to approx. 108 g.

3.6.2 Servicing/Replacing the Particulate Filter

If the filter is subject to a final pressure drop of 450 Pa maximum, it has to be replaced or the absorbed dust has to be removed; it has to be replaced at least every year pursuant to DIN 6022. The final pressure drop is indicated by the exhaust air monitor by the appropriate error message (see Sect. 3.4). Whenever this error occurs, check also if there are any other restrictions of the airflow.

The filter element of the Stage 1 filter may be vacuumed or patted off carefully. While it is also water-resistant and may be wiped with water, this should be done only very exceptionally. After wiping it with water, it has to be left to dry completely before it is reinstalled. Only for replacement filters purchased through the distributor of the product, the full technical compatibility can be endorsed.

3.6.3 Particulate Filter Disposal

The way to dispose of the particulate filter depends on the filtered substances. Comply with the requirements in force where it is used.

Caution

Hazardous substance residues may attach to the filter body.

3.6.4 Technical Specifications of the Particulate Filter

Filter Class (pursuant to DIN EN 779)	F5
Avg. Filtration Level (pursuant to DIN EN 779)	> 90%
Avg. Efficiency Level (pursuant to DIN EN 779)	> 40%
Filter medium and Spacer	100% PP
Frame	100% metal
Maximum Operating Temperature	< 80°C
Recommended Operating Temperature	< 60°C
Short-time Temperature Peaks	< 110°C
Maximum Relative Ambient Humidity	< 100%
Characteristics in Fire (pursuant to DIN 53 438, part 3)	Class F1 (for PP filter medium and metal frame)

Filter Element	
Depth x Height	240 x 220 mm
Width	28 mm
Pleat Height	22 mm
Filter Surface	approx. 0,46 m ²
Glueing	Hotmelt (polyolefin)
Computed Initial Pressure Loss	approx. 40 Pa @ 200 m ³ /h
Recommended Final Pressure Loss	max. 450 Pa
Final Pressure Loss at Replacement	max. 450 Pa
Dust-carrying capacity measured pursuant to ISO 12103-1 A4:	168 g/m ²

3.6.5 Stage 2 Filter

Stage 2 of the filter system is a activated carbon filter. The substances it is able to filter are listed in the appendix (pages 20 ff.).

The activated carbon filter works as a main filter protected by a particulate filter as a pre-filter (Stage 1) pursuant to DIN 12927.

3.6.6 Servicing/Replacing the Activated Carbon Filter

The activated carbon filter needs replacing when its final weight exceeds 130% of its initial weight or, pursuant to DIN 6022, at least every year. The final weight has to be checked every 3 months (see Sect. 3.2). The very day the check interval expires, the exhaust air monitor indicates that it has. Then the filter has to be weighed and its weight has to be recorded for documentation purposes. This procedure serves as the filter check mandatory under DIN 12927. Only for replacement filters purchased through the distributor of the product, the full technical compatibility can be endorsed.

Caution

The filters have to be replaced immediately

- if a filter is broken;
- if a large spill happened (e.g., leaking container);
- if a maximum threshold was exceeded in the extracted air.

3.6.7 Activated Carbon Filter Disposal

The activated carbon filter comes under special/hazardous waste. It has to be disposed of via a local disposal agent or via the distributor. If the distributor is involved, return the filter in its original wrapping. Before returning it, a prior, written confirmation by the distributor is required, so filter disposal and its documentation can be carefully planned.

Caution


When returning the used filters, observe the instructions by the transporter and the requirements applicable to the transport of hazardous goods (ADR). Return shipments that are not properly documented and labeled will not be received.

3.6.8 Technical Specifications of the Activated Carbon Filter

Designation	HSL AK 66 (standard filter)
Type of Activated Carbon	ABEK Hg and AX (data sheet available upon request)
Designation	HSL AK S (special filter isobutane)
Type of Activated Carbon	A K and AX (data sheet available upon request)
Frame	100% PS
Maximum Operating Temperature	< 60°C
Recommended Operating Temp.	< 40°C
Short-time Temperature Peaks	< 80°C
Maximum Relative Ambient Humidity	< 70%
Filter Element	
Depth x Height	240 x 220 mm
Width	97 mm
Fill Dimensions	2 x 44 mm
Overall Activated Carbon Weight	approx. 2000 g
ABEK Hg Activated Carbon Weight	approx. 1000 g
AX Activated Carbon Weight	approx. 1000 g
Glueing	Hotmelt (polyolefin)

Grounding Wires

(open panel carefully, disconnect wires, then remove panel completely; prior to closing the panel, reconnect the grounding wires of both filters to the lugs)


Computed Initial Pressure Loss	not measured
Carrying Capacity per Activated Carbon w/w	30%
Overall Weight (act. carbon + enclosure)	see name plate

Important

Whenever filters are checked or replaced, the four screws in the filter enclosure, which ensure that the filters are in the optimum position along the sealing tape, have to be removed before the filters are pulled out. Failure to do so exposes the filters to a substantial risk of damage. When the check or the replacement is completed, the last step is to reinstall these screws. The strap around the filters helps pull them out of the enclosure more easily. The strap has to be reinstalled around the checked filters, or around the new filters in the case of a replacement, prior to the insertion of the filters into the enclosure.

Screws to Position
 Filters Correctly

(loosen prior to pulling filters out, tighten after reinstalling them)

3.6.9 Checking the Filters (every 3 months)

Every 3 months the exhaust air monitor indicates that the installed filters need checking. To check the filters, the lid panel has to be removed from the filter enclosure. A sticker on this panel also provides information on the filter stages and the used filters.

Behind the panel, the filters are accessible. Information texts on them make the filters easily identifiable. Information on the date of shipment and, possibly, the weight are also provided on the filters.

Important

For both filters, when checking them, 1 min must pass between the filter's being pulled out and its being reinserted and fastened in the filter compartment. Accordingly, wait for 1 min before reinstalling it even in those situations when the filter is clean and in order and nothing has to be done with it.

Checking the Stage 1 Filter (Particulate Filter)

Carefully pull the particulate filter out of the enclosure. Visual inspection is sufficient to check the level of contamination of the filter. If the filter is very dirty, it should be vacuumed, patted off or, very exceptionally, wiped off with water. If water is used, wait for the filter to dry completely before reinstalling it. If not, the activated carbon may absorb this humidity, which would be a risk of false activated carbon weight readings.

After the filter was reinstalled, this should be filled into the Service and Replacement Report (see the appendix) along with the date, description, short identifier, and signature of the operator. Comply with the rules for documentation established by the organization operating the monitor.

If the filter is so dirty that it cannot be cleaned anymore, the filter has to be replaced immediately. The replacement has to be filled into the Service and Replacement Report along with the date, description, short identifier, and signature of the operator.

Checking the Stage 2 Filter (Activated Carbon Filter)

Carefully pull the activated carbon filter out of the enclosure. The level of contamination of the filter can only be determined by way of weighing it. Its initial weight is indicated on the filter for reference. The current weight of the filter has to be determined with a sufficiently accurate precision scale. If the weight does not exceed 130% of the initial weight, no further action is required. If the current weight is heavier than 130% of the initial weight, both filter stages have to be replaced. The weight of the non-replaced filter or the replacement has to be documented in

writing with the Report in the appendix if the rules established by the operating organization say so.

3.6.10 Replacing the Filters (every year)

Caution

When the filters are replaced, be sure to avoid any contamination of the ambient air or other extraction units possibly existing in the system on site after the monitoring and filtering unit.

Important

The new filters, unless they are installed right after they were purchased, have to be stored dry and clean, protected from potential contamination and deterioration.

Every 12 months, the exhaust air monitor indicates that it is mandatory to replace the installed filters. To replace the filters, remove the lid panel from the filter enclosure. A sticker on this panel also provides information on the filter stages and the used filters.

Behind the panel, the filters are accessible. Information texts on them make the filters easily identifiable. Information on the date of shipment and, possibly, the weight are also provided on the filters.

Important

The manufacturer and distributor assume warranty and liability only for (original) filters they have approved for use beforehand. For other filters, neither full technical compatibility nor occupational safety can be endorsed.

Info

Both filters have to be replaced at least once a year. It makes sense to replace both of them always the same day. For filter replacement as indicated by the unit, priority should be given to the replacement at the correct intervals of the activated carbon filter, since this is the main filter of the unit. If, exceptionally, the particulate filter is replaced earlier or later than the activated carbon filter because it is unrecoverably dirty earlier or later, then the replacement that does not match the standard interval has to be documented separately in writing. **It is under the user's responsibility to make sure that the further replacements of the not-standard replaced filter are performed at the correct intervals.** In all and any cases, when one of the filters is replaced, the state of the other one has to be checked.

Replacing the Stage 1 Filter (Particulate Filter)

Carefully pull the particulate filter out of the enclosure and install a new particulate filter. Be sure that the filter sits right next to the sealing tape, making contact with it. Possibly the sealing tape has to be replaced too.

The replacement, the 12-month interval having expired, must be documented in writing. This is done by filling in the Service and Replacement Report (see the appendix) along with date, description, identifier, and signature of the operator.

Replacing the Stage 2 Filter (Activated Carbon Filter)

Carefully pull the activated carbon filter out of the enclosure and install a new activated carbon filter. Be sure that the filter sits right next to the sealing tape, making contact with it. Possibly the sealing tape has to be replaced too. The RFID chip on the filter has to point towards the antenna of the enclosure, so the exhaust air monitor will recognize the filter replacement and reset the replacement interval to the beginning.

The replacement, the 12-month interval having expired, must be docu-

mented in writing. This is done by filling in the Service and Replacement Report (see the appendix) along with date, description, identifier, and signature of the operator.

3.6.11 Alert to Filter Service/Replacement by the Exhaust Air Monitor

When the three-month filter service interval or the twelve-month filter replacement interval expires, the corresponding alert is triggered by the unit in the form of the behavior of the warning and status LEDs (no. 10 and 11 in the table in Sect. 3.4).

The service/check alert is triggered every quarter. The replacement alert starts visually (LEDs) six weeks from time and audibly (alarm beep) four weeks from time.

When the scheduled filter replacement date is overdue for more than 1 month, the alert no. 10 changes to alert no. 5 (per the table in Sect. 3.4), and the airflow monitor can no longer be used, before the filters were actually replaced.

The unit detects the new filters by the RFID chip inserted into every filter.

4 Exhaust Air Monitor with Filter Technical Specifications

Power Supply	110 V - 230 V, 50 Hz/60 Hz
Rated Power	26 W
Fluid velocity *)	2.19 m/s
Volumetric flow rate *)	30.8 m³/h
Weight (unit without filters)	17.2 kg
Weight of Particulate Filter	0.5 kg
Weight of Activ'd Carbon Filter	2.5 kf
Dimensions	310 x 460 x 415 mm
Ingress Protection	IP20
ATEX Rating	EX II (3) G [n T4]

*) as applying to approx. 3 m actual measured length, inside diameter 70.5 mm

5 Troubleshooting

In the course of using the appliance, malfunctions may occur. The table below is meant to help troubleshoot and correct malfunctions. If a malfunction cannot be traced to one of the possible causes or the correction given does not solve the issue, then use of the appliance has to be discontinued. Consult with the distributor on troubleshooting an repair.

Malfunction	Possible Cause	Verification of Cause	Correction of Malfunction
Appliance does not run	On/Off Switch Power Supply Cord Overvoltage Switch	- Is the appliance switched on at the On/Off switch? - Is the power supply cord damaged? - Is the power supply cord properly connected? - Is the fuse inside the appliance o.k.?	- Switch the appliance on at the On/Off switch. - Return unit to manufacturer/distributor for check/repair. - Connect the supply cord to the intended socket. - Replace the fuse (1 x 2 A medium slow blow).
Fan runs, but LEDs do not work	LED Electrics/ Electronics	Malfunction cannot be corrected by user.	Consult with the manufacturer/distributor on inspection/repair.
LEDs work, but fan does not run	Fan Control	Malfunction cannot be corrected by user.	Consult with the manufacturer/distributor on inspection/repair.

Malfunction	Possible Cause	Verification of Cause	Correction of Malfunction
Exhaust Airflow Error shows on the LEDs	Inlet or Outlet Hose Airflow Sensor Particulate Filter Activ'd Carbon F. Filter Enclosure Panel	- Is the suction or the extraction hose of the appliance blocked (airflow restriction) or very dirty? - Is the airflow sensor properly set? - Is the spiral hose connected tightly to the sleeve? - Is the particulate filter dirty? - Is the activated carbon filter dirty? - Is the panel covering the filter enclosure open?	- Clean the hose and/or clear the restriction. - Open the lid below which the airflow sensor is located and use the potentiometer to set it (only by a briefed person). - Tightly set the spiral hose onto the sleeve using the clamp. - Clean and replace the filter as needed. - Replace the activated carbon filter. - Close the filter enclosure with the panel.
Despite irregular extraction airflow no airflow error indicated	Airflow Sensor	Is the airflow sensor properly set?	Open the lid below which the airflow sensor is located and use the potentiometer to set it (only by a briefed person).
Error indicating absence of filter shown	Activated Carbon Filter	Is the activated carbon filter installed?	- Install the filter into the enclosure. - Ensure that activated carbon filter is correctly positioned with RFID chip towards the antenna.

6 Contact Details


DÜPERTHAL Sicherheitstechnik GmbH & Co. KG
 Frankenstr. 3 Tel.: +49 6188 9139 0
 63791 Karlstein, Germany Fax: +49 6188 9139 121
 Web: www.dueperthal.com E-mail: info@dueperthal.com

We reserve the right to change technical specs of the unit without prior notice.

Pursuant to the directive 2012/19/EU on Waste Electrical and Electronic Equipment (so-called WEEE Directive), equipment that was manufactured or distributed by ourselves may be returned to us. To discuss the exact procedure we will follow, please contact us with the details above.

We also declare that the equipment manufacture conforms to the directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (so-called RoHS Directive).

Service and Replacement Report

Serial Number of Exhaust Monitor:

Service	Replacem	Date	Filter Stage		Stage 2 Filter Weight	Filter Stage Date of Shipment	Identifier	Signature
			1.	2.				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				
☐	☐		☐	☐				

Please fill in or check as appropriate.

Remark:

The Service and Replacement Report has to be consistent with the data saved to system memory of the monitoring unit with that serial number.

Appendix: Substances with which the Product is Compatible

Substances not included in the below repository may be eligible for filtration using the product, too. However, this is subject to a prior verification case by case. Feel free to consult with us on this.

Prior to putting the product into service, it is imperative for the entity operating the extraction system to double-check whether the extracted substances are eligible for release to the work area (see Sect. 1.3).

A = Type of Filter

B = Main Area of Application

C = Harmful Substances

A

B

C

Substance list for standard filter, with no special mark

AX/A Gases and vapors of organic compounds, boiling point $\leq 65^{\circ}\text{C}$

Acetaldehyde	Chloro ethane
Acetic acid anhydride	2-chloro ethanol
Acetic acid 3-methyl butyl ester	n-Chloroformyl morpholine
Acetone	3-chloro 2-methyl propene
Acetone cyanhydride	1-chloronaphthalene
Acetonitrile	1-chloro 1-nitro propane
Acrylic aldehyde	1-chloropentane
Acrylic acid	1-chloropropane
Allyl chloride	2-chloropropene
1-aminobutane	5-chloro-o-toluidine
2-aminobutane	Cresols (isomer mixture)
2-amino 2-methyl propane	Cumene
Iso-amyl alcohol	Cyclohexane
Aniline	Cyclohexanol
p-Anisidine	Cyclohexanone
Benzaldehyde	Cyclohexene
Benzol amine	Cyclohexyl amine
Bromo benzol	1,3-cyclo pentadiene
2-bromobutane	Cyclopentanone
1-bromo 2-chloro éthane	n-Decane
Bromochloro methane	n-Decanol
2-bromo 2-chloro 1,1,1-trifluoro ethane	Demeton
Bromoform	Demeton methyl
2-bromopentane	Diazinon
n-Butanal	Dibenzylamine
1-butanol	Dibenzyl ether
2-butanol	Dibutyl ether
tert-Butanol	1,2-dichloro benzol
2-buténal	1,3-dichloro benzol
2-butoxy ethanol	1,4-dichloro benzol
2-butoxy ethyl acetate	2,2'-dichloro diethyl ether
1-butyl acetate	bis (2-chloro 1-methyl ethyl) ether
2-butyl acetate	Dichloro dimethyl ether
tert-Butyl acetat	1,1-dichloro ethane
n-butyl acrylate	1,1-dichloro ethene
Butyl formiate	1,2-dichloro ethene (cis/trans)
n-butyl methacrylate	Dichloro methane
Butyl stearate	1,1-dichloro 1-nitro ethane
p-tert-Butyl toluol	N,N-dimethyl aniline
Camphor	2,2-dimethyl butane
Chlorine dioxide	2,3-dimethyl butane
Chloro acetaldehyde (at 45% in water)	1,3-dimethyl butyl acetate
Chloroacetic acid ethyl ester	N,N-dimethyl cyclohexyl amine
Chlorobenzol	N,N-dimethyl ethanol amine

A

B

C

AX/A Gases and vapors of
organic compounds, boiling
point $\leq 65^{\circ}\text{C}$

Substance list for standard filter, with no special mark

Dimethyl ether	2-hexene (cis and trans isomeres)
N,N-dimethyl formamide	4-hydroxy 4-methyl pentan-2-on
Dimethyl hydrogen phosphite	Iodomethane
2,2-dimethyl propane	Isobutyl acetate
1,5-dinitro naphthalene	Isobutyl amine
2,6-dinitro naphthalene	Isobutyl formiate
1,4-dioxane	Isobutyraldehyde
1,3-dioxolane	Isoflurane
Dipentene	Isooctane
Diphenyl ether	Isopropyl acetate
Diphenyl ether/Biphenyl mixture	Isopropyl chloride
Dipropyl amine	Isopropyl nitrate
Dipropylene glycol methyl ether	Kerosene
Dipropyl ether	Lead tetraethyl
1,2-divinyl benzol	Lindane
1,3-divinyl benzol	Malathione
Divinyl ether	Methacryl nitrile
Dodecyl benzol	Methanol
Endrin	Methoxychlor
Enflurane	1-methoxy 2-propanol
EPN	2-methoxy 1-propanol
Ethanol	1-methoxy 2-propyl acetate
Ethanol amine	2-methoxy 1-propyl acetate
2-ethoxy ethanol	Methyl acetate
Ethyl acetate	Methyl acrylate
Ethyl acrylate	N-Methyl aniline
Ethyl benzol	Methyl bromide
Ethylene glycol	2-methyl butane
Ethylene glycol monoisopropyl ether	Methyl cyclohexane
Ethyl formiate	Methyl cyclohexanol
2-ethyl 1-hexanol	2-methyl cyclohexanol (isomer mixture)
2-ethyl hexyl amine	Methyl diethanol amine
5-ethylidene 2-norbornene	Methyl ethyl ketone
Ethyl mercaptane	Methyl formiate
Aircraft fuel F 34	Methyl isobutyl ketone
Fenthion	Methyl methacrylate
Ferbam	2-methyl pentane
Fluoro benzol	3-methyl pentane
Formamide	4-methyl pentan-2-ol
Furfurol	4-methyl pent-3-en-2-one
Furfuryl alcohol	2-methyl 1-propanol
Glycerin	2-methyl 2-propan thiol
Glyoxal, aqueous solution	Methyl propionate
Heptachlor	Methyl propyl ketone
n-Heptane	N-Methyl 2-pyrrolidone
2-heptanone	Methyl styrol
3-heptanone	Methyl tert-butyl ether
4-heptanone	Mevinphos
Hexachloro ethane	Morpholine
Hexal amine	Naled
Hexalene glycol	Nicotine
Hexamethylene diamine	2-nitro 4-amino phenol
n-Hexane	4-nitro aniline
n-Hexanol	Nitroethane
2-hexanone	Nitromethane
1-hexene	1-nitro naphthalene

C

AX/A Gases and vapors of organic compounds, boiling point $\leq 65^{\circ}\text{C}$

EN 23

A
B
C

B Inorganic gases and vapors, e.g. chlorine, hydrogen sulfide,

Substance list for standard filter, with no special mark

Acetyl chloride	Propyl mercaptane
Allyl propyl disulfide	Sodium cyanide
Antimony pentachloride	Sodium fluoracetate
Arsine	Sodium hypochlorite solution
Boro trifluoride	Stibine
Bromine	Sulphur dichloride
Bromocyan	Sulphur trioxide
2-butanone peroxide	Sulphuryl chloride
tert-Butyl hydroperoxide	Tetrahydrothiophene
1-butyl mercaptane	Thiourea
tert-Butyl peracetate	Thionyl chloride
Carbon sulfide	Thiram
Chlorine	Tin(IV) chloride
Chlorocyan	Tributyl tin benzoate
Chloro trifluoride	Tributyl tin chloride
Cyanacrylic acid methyl ester	Tributyl tin fluoride
Cyanuric chloride	Tributyl tin linoleate
Cyanuric chloride (aqueous suspension)	Tributyl tin methacrylate
Decaborane	bis(Tributyl tin)oxide
Demeton	Trimethyl amine
Demeton methyl	
Diborane	
Di-tert-butyl peroxide	
Diethyl sulfide	
Diisopropyl amine	
Dilauroyl peroxide	
Cumene hydroperoxide	
Dimethyl disulfide	
Dipropyl amine	
Disulphur dichlorid	
Disulfirame	
Fluorine	
Formic acid	
Hexamethylene tetramine	
Hydrogen bromide	
Hydrogen chloride	
Hydrogen sulfide	
Hydrogen selenide	
Hydroxyl amine	
Iodine	
Iron chloride	
2-mercapto ethanol	
Methyl mercaptane	
2-methyl 1-propanol	
2-methyl 2-propan thiol	
Nitric acid	
Oxalic acid dinitrile	
Phosphorated hydrogen	
Phosphorus oxide chloride	
Phosphorus pentachloride	
Phosphorus trichloride	
Potassium cyanide	
2-propan thiol	
Propionic acid	
Propoxur	

A

B

C

Substance list for standard filter, with no special mark

E Sulfur dioxide, hydrogen chloride

Acetic acid	Hydrogen chloride
Calcium bisulfide solution	Phosphorus trichloride
Formic acid	Sulphur dioxide
Hydrogen bromide	

K Ammonia

1-amino propane	N,N-Dimethyl ethyl amine
2-amino propane	Ethyl amine
Ammonia	Hexamethylene tetramine
Ammonia solution 25%	Hydroxyl amine, aqueous solution
Dimethyl amine	Methyl amine

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