SUPREME standard, ECO standard

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DÜPERTHAL®

Typ G90/G30

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Dokument: D00101441 Rev.: 01

SUPREME standard, ECO standard
Safety storage cabinets

Typ G90/G30

Operating instructions

DÜPERTHAL Sicherheitstechnik GmbH & Co. KG
Frankenstrasse 3
D-63791 Karlstein
Germany
Telephone: +49 6188 9139-0
Fax: +49 6188 9139-121
Email: info@dueperthal.com
Internet: www.dueperthal.com

© 09/2025 DÜPERTHAL 55228, D00098273, Rev. 02, en_GB

Table of contents

1	General information.....	71
1.1	Notes for reading.....	71
1.2	Type plate.....	71
2	Safety.....	72
2.1	Function of safety notices.....	72
2.2	Correct use.....	72
2.3	Misuse.....	72
2.4	Obligations of the operator.....	74
2.5	Demands on employees.....	74
2.6	Stored goods.....	74
2.7	Hazardous areas and their labelling.....	75
3	Technical specifications.....	77
3.1	General data.....	77
3.2	Dimensions and equipment.....	78
3.3	Pressure drop during ventilation.....	82
4	Structure and function.....	84
4.1	Construction.....	84
4.2	Earthing options.....	85
4.3	Exhaust air connection and feed opening.....	85
4.4	Doors.....	86
4.5	Closure of valve flaps in case of fire.....	86
4.6	Interior fittings.....	87
4.7	Extra load on the cabinet roof of the SUPREME standard (optional).....	95
5	Transport.....	96
6	Installation and commissioning.....	98
6.1	Requirements on the installation location.....	98
6.2	Installing the safety storage cabinet.....	99
6.3	Align the safety storage cabinet.....	100
6.4	Check the alignment of the safety storage cabinet.....	101
6.5	Venting the safety storage cabinet.....	102
6.6	Earth the safety storage cabinet.....	105
7	Operation.....	106
7.1	Open the safety storage cabinet.....	106
7.2	Fitting holder.....	107
7.3	Cylinder holder.....	111
7.4	Operating the roll-in flap.....	113
7.5	Connecting the purging gas pipe.....	114
7.6	Installing the storage shelf.....	115
7.7	Attaching the storage insert for screw caps (optional).....	115
7.8	Attaching penetrations.....	116
8	Opening the safety storage cabinet after a fire.....	118
9	Maintenance.....	119
10	Faults.....	121
11	Spare parts and accessories.....	122

12	Disposal.....	123
13	Certificates.....	124

1 General information

1.1 Notes for reading


The following symbols designate specific types of information.

Tab. 1: Explanation of symbol

Symbol	Type of information
	Information for easier and more effective working
	Procedural step
	Result of a procedural step
	Link to another part of the document

1.2 Type plate

The type plate is attached to the outside of the safety storage cabinet door.

Model

Serial number and year of
manufacture
Reference number

Maximum volume of individual containers

Maximum load per storage shelf

Fig. 1: Type plate

2 Safety

2.1 Function of safety notices

Safety notices warn against physical or material damage and provide information on how such damage can be avoided.

The following signal words identify the degree of danger and the extent of the risk.

WARNING!

The signal word '*WARNING*' refers to a potential hazard which could result in death or serious injury.

CAUTION!

The signal word '*CAUTION*' refers to a potential hazard which could result in slight or minor injury.

NOTICE!

The signal word '*NOTE*' indicates a situation that can lead to damage to the safety storage cabinet.

2.2 Correct use

Observe the safety instructions in these operating instructions to reduce health risks and avoid dangerous situations.

Any use that is not correct use as defined in these operating instructions involves a risk of accidents and a lack of fire protection.

The SUPREME standard safety storage cabinets are type tested and classified as '*Type G90*' in compliance with '*EN 14470-2*' with a fire resistance of 90 minutes.

The ECO standard safety storage cabinets are type tested and classified as '*Type G30*' in compliance with '*EN 14470-2*' with a fire resistance of 30 minutes.

The safety storage cabinet is to be used for the storage of pressurised gas cylinders in working spaces.

Nationally applicable regulations must be observed.

NOTICE!

The total volume of all pressurised gas cylinders in the safety storage cabinet must not be more than 220 litres. This includes pressurised gas cylinders for purging gas with a maximum volume of 10 litres.

2.3 Misuse

Any use that goes beyond the specified correct use is considered to be misuse.

DÜPERTHAL accepts no liability for damage arising from misuse.

In addition, the following safety notices must be observed:

WARNING!

Storage of living organisms in the safety storage cabinet

Living organisms can come into contact with the stored hazardous substances.

This may lead to death or serious injury.

- Use the safety storage cabinet exclusively for storage of pressurised gas cylinders.

WARNING!

Storage of food in the safety storage cabinet

Food can come into contact with the stored hazardous substances.

This may lead to death or serious injury.

- Use the safety storage cabinet exclusively for storage of pressurised gas cylinders.

WARNING!

Storage of flammable liquids

Flammable liquids may not be stored with gas cylinders in the safety storage cabinet.

This may lead to death or serious injury.

- Use the safety storage cabinet exclusively for storage of gas cylinders.

WARNING!

Objects on the cabinet roof

In case of fire, objects on the cabinet roof can impair the function of the safety technology.

The consequences may be death or serious injury.

- Do not store any objects on the cabinet roof.

WARNING!

Pipe or cable entries

In case of fire, too many or excessively large pipe or cable entries can impair the function of the safety equipment.

Pipe or cable entries outside the approved areas can also impair the safety equipment.

This may lead to death or serious injury.

- Only attach pipe or cable entries within the approved areas and in the specified number.

! NOTICE!**Alteration and modifications**

Do not alter or modify the safety storage cabinet.

This can lead to a lack of fire protection.

- If alteration or modification of the safety storage cabinet is required, contact DÜPERTHAL.

2.4 Obligations of the operator

The operator is obliged to comply with applicable legal regulations. This includes:

- Issuing operating instructions.
- Carrying out risk assessments.
- Creating explosion protection documents.
- Specifying activities by designated employees.

2.5 Demands on employees

! WARNING!**Employees who do not meet these requirements**

This may result in death or serious injury.

- Designate employees who meet the requirements to carry out activities.

These operating instructions set out the following employee activities:

- Specialist technical employees
- DÜPERTHAL service technicians

Only people who have been trained by the operator in use of the safety storage cabinet and handling of the stored goods are approved as specialist technical employees.

DÜPERTHAL service technicians

DÜPERTHAL employees are specifically trained by DÜPERTHAL to carry out their activities.

2.6 Stored goods

! WARNING!**Leaking fittings and connections**

Gases can escape from leaking fittings or connections.

This may lead to death or serious injury.

- Fittings and connections must be tested for leaks by the designated technical specialist personnel.

Valid national regulations regarding the storage of pressurised gas cylinders, such as 'TRGS 510' in Germany, must be observed. Information regarding storage volumes, combined storage of different hazardous substances and 'GHS' classifications is listed and defined in these regulations.

2.7 Hazardous areas and their labelling

The following must be attached to the front of the safety storage cabinet and must be clearly visible:

- The instruction "Close the door"
- Fire resistance in minutes (e.g. 'Type G90')
- Name or trademark of the manufacturer
- Serial number and year of manufacture

Furthermore, the following signs must be attached to the front of the safety storage cabinet and must be clearly visible:

Tab. 2: Prohibited action signs

Symbol	Meaning	Standard
	P003: No naked flames; fire, open ignition source and smoking prohibited	DIN EN ISO 7010

Tab. 3: Warning signs

Symbol	Meaning	Standard
	W029: Warning for pressurised gas cylinders	DIN EN ISO 7010

Tab. 4: Mandatory action signs

Symbol	Meaning	Standard
	M002: Read the instructions	DIN EN ISO 7010
	M004: Wear eye protection	DIN EN ISO 7010

Symbol	Meaning	Standard
	M009: Wear hand protection	DIN EN ISO 7010

The relevant national rules and legislation may require additional hazardous area markings for storing particular hazardous substances, such as acutely toxic gases.

3 Technical specifications

3.1 General data


Fig. 2: General diagram of safety storage cabinet Type G90

- 1 Exhaust air connection
- 2 Type plate
- 3 Installation rails
- 4 Cylinder holder
- 5 Roll-in flap

3.2 Dimensions and equipment

3.2.1 SUPREME standard


Fig. 3: Technical equipment for the SUPREME standard

H	Height	T1	Cabinet depth with open right doors
B	Width	T2	Cabinet depth with open left doors
T	Depth	LB	Clear width
H1	Adjustable height	L	Max. total volume of the pressurised gas cylinders (according to EN 14470-2)
BS	Standing surface width		
TS	Standing surface depth		

	S	M	L	XL	XXL
H (mm)	1375	2030	2030	2030	2030
B (mm)	595	595	895	1195	1195
T (mm)	610				
H1 (mm)	1171	1862			
BS (mm)	430	430	730	1030	1230
TS (mm)	268	268	249	249	249
T1 (mm)	-	-	965	1262	1365
T2 (mm)	1182	1182	1182	1182	1282
LB (mm)	415	415	699	999	1199
Empty weight (kg)	260	350	480	595	680
Max. payload (kg)*	100	100	300	400	400
L (l)	220				
Max. extra load on the cabinet roof (kg)	10	10	15	20	25
Intermediate shelf (optional) load-bearing capacity with uniformly distributed load (kg)	50	50	50	50	50
Load-bearing capacity per shelf for side attachment (optional) with uniformly distributed load (kg)	60	60	60	60	60
Load-bearing capacity of storage shelf for pressurised gas cylinders (optional) with uniformly distributed load (kg)	50	50	50	50	50

* When using an extra load, the max. payload is reduced.

3.2.2 ECO standard


Fig. 4: Technical equipment for the ECO standard

H	Height	T1	Cabinet depth with open right doors
B	Width	T2	Cabinet depth with open left doors
T	Depth	LB	Clear width
H1	Adjustable height	L	Max. total volume of the pressurised gas cylinders (according to EN 14470-2)
BS	Standing surface width		
TS	Standing surface depth		

	M	L	XL	XXL
H (mm)	2030			
B (mm)	595	895	1195	1395
T (mm)	596			
H1 (mm)	1884			
BS (mm)	500	800	1100	1294
TS (mm)	272			265
T1 (mm)	-	936	1236	1336
T2 (mm)	1168	1168	1168	1268
LB (mm)	428	705	1005	1199
Empty weight (kg)	285	395	480	560
Max. payload (kg)	200	300	400	400
L (l)	220			
Intermediate shelf (optional) load-bearing capacity with uniformly distributed load (kg)	50	50	50	50
Load-bearing capacity per shelf for side attachment (optional) with uniformly distributed load (kg)	60	60	60	60
Load-bearing capacity of storage shelf for pressurised gas cylinders (optional) with uniformly distributed load (kg)	50	50	50	50

3.3 Pressure drop during ventilation

Industrial ventilation of the safety storage cabinet results in a pressure drop at the exhaust air connection.

3.3.1 SUPREME standard


Fig. 5: SUPREME standard

- 1 Pressure drop
- 2 Average pressure drop from cabinet sizes S and M
- 3 Average pressure drop from cabinet sizes L, XL and XXL
- 4 Volumetric flow rate
- 5 Q with ten-fold air exchange (see table)

Model size	Q [m³/h]		ΔP [Pa]	
	10-fold air exchange	120-fold air exchange	10-fold air exchange	120-fold air exchange
S	1.8	21.8	<1	13
M	2.9	35.1	<1	41
L	5.3	63	<1	34
XL	7.5	89.4	<1	85
XXL	8.9	107.1	<1	131

3.3.2 ECO standard


Fig. 6: ECO standard

- 1 Pressure drop
- 2 Average pressure drop from all cabinet sizes
- 3 Volumetric flow rate
- 4 Q with ten-fold air exchange (see table)

Model size	Q [m³/h]		ΔP [Pa]	
	10-fold air exchange	120-fold air exchange	10-fold air exchange	120-fold air exchange
M	4.1	49.5	<1	27
L	6.6	79.4	<1	39
XL	8.9	107	<1	147
XXL	10.5	126.5	<1	220

4 Structure and function

4.1 Construction

- Cabinet carcass and doors in multi-layer construction
- External panelling: Powder-coated sheet steel
- Wall construction: Multi-layer design
- Interior surfaces: Light grey-coated decor panels
- Safety technology elements for closure of venting cut-off flaps in case of fire: Brass, spring steel (1.410)

4.2 Earthing options


Fig. 7: Earthing facility

1 Equipotential bonding saddle on the cabinet carcass

The equipotential bonding saddle is pre-fitted on the cabinet roof.

An optional earthing facility is possible in the cabinet interior.

4.3 Exhaust air connection and feed opening

The safety storage cabinets can be connected to an industrial exhaust air system which discharges air outside in a safe location. For this purpose, the exhaust air connection and the feed openings are located on the roof of the safety storage cabinet.

In normal operation, industrial ventilation of safety storage cabinets prevents the occurrence of a potentially explosive and hazardous atmosphere inside the cabinet.

Exhaust air connection NW 110 mm with reducer adapter NW 75 mm for adaptation to a ventilation system is possible.

The layout of the ventilation ducts in the cabinet means that ventilation is effective over the entire inner height.


Fig. 8: Exhaust air system (view from above)

- 1 Exhaust air connection
- 2 Air supply
- 3 Thermocouple
- 4 123 mm for SUPREME standard and 106 mm for ECO standard

4.4 Doors

Safety storage cabinets of the SUPREME standard and ECO standard type have the following door options:

- One-piece wing door with door handle
- Two-door opening with door handle and interior locking

4.4.1 Locking cylinder

The door is lockable using the locking cylinder integrated into the door handle.

Door handles with Euro profile half-cylinder are available as an option. Integration into a building locking system is possible.

4.5 Closure of valve flaps in case of fire

A valve flap is located below the exhaust air and air supply openings.

The valve flaps are kept continuously open by a fusible link.

At a temperature of approx. 60°C, the valve flaps close off the exhaust air and air supply openings.

The valve flaps and the fusible link may not have their function impaired and must not be covered or blocked.

4.6 Interior fittings

4.6.1 Fitting holder

4.6.1.1 SUPREME standard


Fig. 9: SUPREME standard fitting holder

The SUPREME standard type safety storage cabinet contains installation rails for mounting pressurised gas fittings.

4.6.1.2 ECO standard


Fig. 10: ECO standard fitting holder

The ECO standard type safety storage cabinet contains installation rails on the rear panel for mounting pressurised gas fittings.

4.6.2 Cylinder holder


Fig. 11: Cylinder holder

- 1 Cylinder holder
- 2 Retaining belts

The SUPREME and ECO standard safety storage cabinets are equipped with cylinder holders and retaining belts throughout as standard.

The retaining belts are designed to ensure safe mounting of the pressurised gas cylinders. The number of retaining belts depends on the size of the safety storage cabinet and its cylinder capacity.

4.6.3 Folding roll-in flap


Fig. 12: Roll-in flap

- 1 Roll-in flap
- 2 Locking bar

The SUPREME and ECO standard safety storage cabinets are equipped with a folding and lockable roll-in flap as standard.

4.6.4 Purging gas pipe connection


Fig. 13: Purging gas pipe connection

1 Purging gas pipe connection (4 x Ø10)

Below the covers of the closing valves for the ventilation are four holes for connection of the purging gas pipes.

4.6.5 Storage shelf (optional)


Fig. 14: Storage shelf

- 1 Mounting screw
- 2 Pedestal
- 3 Storage shelf
- 4 Support strip
- 5 Fixing screw

Storage shelves attached to the rear panel can be installed in the SUPREME and ECO standard safety storage cabinets as an option.

The heights of the storage shelves can be changed. The highest shelf must not be more than 1.75 m above the floor.

4.6.6 Shelf for side installation (optional)


Fig. 15: Shelf for side installation

- 1 Shelf
- 2 Retaining belt
- 3 Mounting screws (x4)
- 4 Mounting bracket (x2)

Shelves for side installation can be mounted on the side panel in the SUPREME and ECO standard safety storage cabinets as an option.

The height of the shelf can be freely selected up to a maximum of 1.75 m above the floor.

4.6.7 Storage shelves, diagonal layout (optional)


Fig. 16: Storage shelves, diagonal layout

1 Storage shelves for gas cylinders

Storage shelves for pressurised gas cylinders can be mounted on the side panels in the SUPREME and ECO standard safety storage cabinets as an option.

The height of the storage shelf can be freely selected up to a maximum of 1.75 m above the floor.

4.6.8 Storage insert for screw caps (optional)


Fig. 17: Storage insert for screw caps

1 Storage insert

Storage inserts for pressurised gas cylinder screw caps can be installed on the cylinder holder in the SUPREME and ECO standard safety storage cabinets as an option.

4.7 Extra load on the cabinet roof of the SUPREME standard (optional)


Fig. 18: SUPREME standard with sheet metal cladding as a load

Optionally, the cabinet top of SUPREME standard cabinets can be used as a shelf for light loads, such as gas detection systems or sheet metal cladding.

The maximum permissible load (☞ Chapter 3.2.1 'SUPREME standard' on page 78) must not be exceeded. The weight must be distributed evenly across the cabinet top.

The valve flaps and triggering thermocouples installed on the cabinet (☞ Chapter 4.5 'Closing the valve flaps in the event of a fire' on page 86) must remain freely accessible. Installed weights must not interfere with annual maintenance.

The load must be secured against slipping. It can be fastened to the cabinet ceiling with a maximum of 6 sheet metal screws with a maximum diameter of 4.8 mm and a maximum length of 12 mm.

If you have any questions about loads, please contact DÜPERTHAL before installation.

The extra load can only be installed for the SUPREME, not the ECO.

5 Transport

The safety storage cabinet is packaged for transport and is protected against damage by transport restraints. The transport restraints should be refitted before any transport.

WARNING!

Crush hazard due to safety storage cabinet tipping over

If the safety storage cabinet tips over when not transported with due caution, this can cause potentially fatal crushing.

- Wear personal protective equipment (PPE).
- Only transport with two people.
- Only transport the safety storage cabinet upright and unladen.
- When driving underneath it, only pick up the safety storage cabinet in the centre from the front or from the side.
- Only drive under the safety storage cabinet using suitable transport equipment.


Fig. 19: Transporting the safety storage cabinet

- 1 Pick up centrally from the front
- 2 Pick up centrally from the side

NOTICE!

Handling the transport restraints

Incorrect handling can damage the safety transport skids and the safety storage cabinet.

- Transport restraints and safety transport skids should only be removed at the installation location.

! NOTICE!

Tipping the safety storage cabinet over during transport

Damage to the safety storage cabinet caused by incorrect handling.

- Only pick up the safety storage cabinet from the side.
- Only pick up the safety storage cabinet using special and suitable transport or lifting equipment.
- Pick up the safety storage cabinet once it is securely lashed and is not at risk of slipping.

6 Installation and commissioning

Install the safety storage cabinet so that the annual maintenance activities can be carried out without restriction.

6.1 Requirements on the installation location

The safety storage cabinet is approved for installation in a building.


Fig. 20: Requirements on the installation location

Consider the following in relation to the installation location:

- The surface must be able to bear the weight of the safety storage cabinet when fully loaded.
- The surface must be horizontal in order to guarantee problem-free functioning of the safety storage cabinet.
- The load-bearing capacity and stability of the surface must be assured both in normal situations and in the event of a fire.
- Do not install the safety storage cabinet near sources of heat.
- Protect the safety storage cabinet against moisture.
 - At a relative humidity of >70 % use in closed and heated buildings is permissible for a few weeks each year.
- The operating temperature must be between -5°C and +40°C.

6.2 Installing the safety storage cabinet


Fig. 21: Installing the safety storage cabinet

6.3 Align the safety storage cabinet

The alignment procedure described below is used for precision alignment. Remedy any major floor unevenness of more than 1.5 mm on site.

Short adjusting screws are fitted in the corners of the base as standard. These are used to align the safety storage cabinet.


Fig. 22: Aligning the safety storage cabinet

1 Hexagon key

Alignment using the adjusting screws

Personal:

■ Technical specialist employees

➔ Screw the adjusting screws in or out using the hexagon key.

! NOTICE!

Align the safety storage cabinet.

Incorrect handling causes damage to the floor at the installation location.

– Place packing plates 90x90x3 under the adjusting screws.

6.4 Check the alignment of the safety storage cabinet


Fig. 23: Check the alignment (figure may differ from product).

Correct alignment of the safety storage cabinet:

- when the doors are closed, the door gaps are of equal width.
- With two doors, the central gap and ceiling gap form an even "T".

6.5 Venting the safety storage cabinet

6.5.1 Operating the safety storage cabinet with industrial ventilation

WARNING!

Insufficient safety storage cabinet air circulation

A lack of or insufficient air exchange can lead to formation of a hazardous atmosphere in the safety storage cabinet.

This may lead to death or serious injury.

In a safety storage cabinet with ventilation system, an air exchange must take place with the doors closed.

- If flammable or oxidising gases are being stored, the air exchange must be at least 10 times the cabinet interior volume per hour.
- If category 1-3 acutely toxic gases are being stored, the air exchange must be at least 120 times the cabinet interior volume per hour.

NOTICE!

Storage of corrosive gases can have an impact on the function of the shut-off mechanism for the ventilation openings for supply and exhaust air.

It may be necessary to adjust the maintenance cycle for the shut-off mechanisms.

Installation of industrial ventilation and connection to an existing exhaust air system must be carried out by a qualified company and is not a service provided by DÜPERTHAL.


Fig. 24: Connection to exhaust air

- 1 Exhaust air connection socket
- 2 Air supply opening
- 3 Thermocouples

Connection to the exhaust air system:

Personal:

- Technical specialist employees

1. ➔ Connect the exhaust air line to the exhaust air connection socket and secure with a collar.
2. ➔ After installing the safety storage cabinet, check the connection to the exhaust air system with smoke tubes.

The power of the exhaust air system can be determined using the technical specifications, ↗ Chapter 3.3 'Pressure drop during ventilation' on page 82.

6.5.2 Operating the safety storage cabinet without industrial ventilation

Non-acutely toxic gases and non-flammable gases can be stored in safety storage cabinets without industrial ventilation.

Safety storage cabinets that are operated without industrial ventilation must be indicated using an appropriate notice.

WARNING!

Safety storage cabinets without industrial ventilation

Danger due to formation of hazardous atmosphere in cabinet interior

This may lead to death or serious injury.

- Regularly check fittings and gas cylinders for leaks.
- If necessary, install additional gas warning equipment.

The size of all symbols and notices should be appropriate for the size of the safety storage cabinet.

6.6 Earth the safety storage cabinet


Fig. 25: Earthing connection

Earthing connection

Personal:

- Technical specialist employees

1. ➤ Bend the equipotential bonding saddle upwards.
2. ➤ Connect the earthing cable (not included in scope of delivery).

Only the external panels are earthed by the equipotential bonding saddle.

7 Operation

7.1 Open the safety storage cabinet


Fig. 26: Door handle position

- 1 Door closed
- 2 Door open

WARNING!

Blocked doors

Doors that are held open by objects impair the function of the safety technology.
This may result in death or serious injuries as a result of inadequate fire protection.

- Close the doors after every work process.

Door operation with door handle

Personal:

- Technical specialist employees

Opening the wing doors

1. ➤ Turn the door handle up by 90°.
 - ⇒ The door is unlocked
2. ➤ Open the wing door by pulling the door handle.
3. ➤ On two-door versions, then push forward the locking system on the end face of the left-hand door by 90°.
 - ⇒ The door is unlocked
4. ➤ Pull open the left-hand wing door
 - ⇒ The wing doors remain open in any position.

Closing wing doors

1. ➤ Close the left-hand wing door
2. ➤ Push back the locking system on the end face of the door by 90°
⇒ The door is locked.
3. ➤ Then close the right-hand door by pushing the door handle.
4. ➤ Turn the door handle down by 90°
⇒ The door is locked

7.2 Fitting holder

7.2.1 Adjust the fitting holder

7.1.1.1 SUPREME standard


Fig. 27: SUPREME standard fitting holder

- 1 Socket head screw
- 2 Bracket
- 3 Installation rail

Personal:

- Technical specialist employees

1. ➤ Loosen the socket head screws on the left-hand and right-hand side.
⇒ This detaches the bracket and installation rail.
2. ➤ Adjust the height of the installation rail to the preferred position.
3. ➤ Tighten the socket head screws on the left-hand and right-hand side.
⇒ The bracket and installation rail are securely mounted and can no longer be adjusted.

7.1.1.2 ECO standard


Fig. 28: ECO standard fitting holder

- 1 Countersunk screw
- 2 Installation rail

Personal:

- Technical specialist employees

1. ➤ Loosen the countersunk screws on the left-hand and right-hand side.
⇒ This detaches the installation rail.
2. ➤ Adjust the height of the installation rail to the preferred position.
3. ➤ Tighten the countersunk screws on the left-hand and right-hand side.
⇒ The installation rail is securely mounted and can no longer be adjusted.

7.2.2 Modifying the fitting holder (SUPREME standard)


Fig. 29: Modifying the fitting holder (SUPREME standard)

- A Inclined position
- B Vertical position
- C Rear panel installation

Modification to vertical position

Personal:

- Technical specialist employees

! NOTICE!

Damage due to incorrect mounting

Incorrect mounting of the fitting holder can result in damage to the safety storage cabinet.

- Only use the marking holes provided for mounting.
- Only use the tapping screws supplied for mounting (max. diameter: 3.5 mm; max. length: 10 mm).

1. ➤ Unscrew the built-in installation unit on the side panels.
2. ➤ Screw the installation unit back on vertically.

Rear panel installation

Personal:

- Technical specialist employees

! NOTICE!

Damage due to incorrect mounting

Incorrect mounting of the fitting holder can result in damage to the safety storage cabinet.

Permitted installation rail mounting options:

- In side area of rear panel (right or left) between side panels and air ducts.
Shown in Figure C (hatched area)

1. ➤ Unscrew the installation unit on the side panels.
2. ➤ Remove the connecting brackets for the vertical and horizontal installation rails.
3. ➤ Position the vertical rails behind the horizontal rails.

4. ➤ Join together the rails with the hammer head mountings supplied.
5. ➤ Attach the installation unit to the rear panel with short wood screws.

The wood screws are not supplied.

Screw requirements: max. diameter: 4 mm; max. length: 20 mm.

7.3 Cylinder holder

7.3.1 Moving the cylinder holder


Fig. 30: Moving the cylinder holder

- 1 Cylinder holder
- 2 Fixing screw
- 3 Area for positioning the cylinder holder

Tab. 5: SUPREME standard

	S	M	L	XL	XXL
a (mm)	273				
b (mm)	358	1000			
c (mm)	27				
d (mm)	316				

Tab. 6: ECO standard

	M	L	XL	XXL
a (mm)	273			
b (mm)	1000			
c (mm)	53			
d (mm)	341			335

Personal:

- Technical specialist employees

1. ➤ Loosen the fixing screws.
2. ➤ Move the cylinder holder to the required position.
3. ➤ Tighten the fixing screws (only use the supplied screws).

7.3.2 Securing the pressurised gas cylinders

Personal:

- Technical specialist employees

- Place the pressurised gas cylinders on the cylinder holder and secure each of them with a retaining belt.

7.4 Operating the roll-in flap


Fig. 31: Folding roll-in flap

- 1 Locking system
- 2 Roll-in flap
- 3 Handle

Personal:

- Technical specialist employees

Operate the roll-in flap with two hands.

1. ➤ Hold the roll-in flap in place with one hand at the designated position (3).
2. ➤ Pull up the locking system on the roll-in flap and unlock.
3. ➤ Carefully fold down the roll-in flap onto the floor.
4. ➤ Load or unload the safety storage cabinet.
5. ➤ Carefully fold the roll-in flap back up.
⇒ The locking system engages automatically and holds the roll-in flap in the vertical position.
6. ➤ Check the locking system.

! NOTICE!

Roll-in flap folded up

If the roll-in flap is not folded back up, jamming and damage can occur when closing the doors.

- Fold the roll-in flap back up after every work process.

7.5 Connecting the purging gas pipe


Fig. 32: Purging gas pipe connection

1 Connection

! WARNING!

Discharging the purging gases inside the safety storage cabinet

Unauthorised discharge of purging gases can lead to death or serious injury.

- The operator must check whether discharging the purging gases in the cabinet interior is permitted under the applicable regulations for operation of gas stations.

Personal:

- Technical specialist employees

1. ➤ Remove the cap from the hole.
2. ➤ Insert the purging gas pipe.
 - ⇒ Smaller purging gas pipes must be sealed after inserting.

! NOTICE!

Enlarging the bore holes

Expanding the bore holes results in damage to the safety storage cabinet and can impair the safety equipment.

- The holes must not be enlarged.
- The insertion depth of the pipe in the rear panel may not exceed 19 mm.

7.6 Installing the storage shelf

Personal:

- Technical specialist employees

The storage shelves for attachment to the rear panel or the side panel are attached using the mounting components supplied.

The installation instructions are enclosed with the storage shelves.

7.7 Attaching the storage insert for screw caps (optional)

Personal:

- Technical specialist employees

The storage insert is fitted in the recesses provided in the cylinder holder using the clips.

The storage insert can then be optionally secured using screws in the holes provided.

The mounting screws are not supplied.

7.8 Attaching penetrations


Fig. 33: Areas for pipe and cable penetrations

1 L marking as limitation for pipe and cable penetrations

The marked areas on the cabinet roof (inside and outside) and the side panel (inside only) of the safety storage cabinet are to be used for the pipe and cable penetrations.

Penetrations on the rear panel and doors are not approved. The number of pipe and cable penetrations is limited to the necessary minimum.

! NOTICE!

Subsequent drilling of pipe penetrations

Incorrect holes in the wrong location can cause damage to the safety storage cabinet.

- Only use pipe penetrations at the designated positions. Contact DÜPERTHAL for further information.

! NOTICE!

Unused penetration and holes

Improperly sealed penetrations and holes that are no longer used can lead to damage to the safety storage cabinet.

- Seal unused penetrations and holes properly.
- Contact DÜPERTHAL to check whether proper repair is possible.

7.8.1 Pipe penetrations

The maximum hole diameter for a gas pipe is 10 mm and may not generally exceed the nominal diameter of the pipe that runs through.

Producing the gas pipes from stainless steel or a material with similar heat conducting properties is recommended. This maintains the fire resistance classification of the safety storage cabinets.

The number of gas pipes must comply with EN 14470-2 (chapter 5.5) and the necessary risk assessment must be observed in its entirety.

7.8.2 Cable penetrations

The maximum hole diameter for cables is approx. 8 mm (electrical cable with 3 wire leads with a wire cross-section of 1.5 mm) and may not generally exceed the nominal diameter of the cable that runs through.

Cable installations with a different diameter and the number of cables must comply with EN 14470-2 (chapter 5.6) and the necessary risk assessment must be observed in its entirety.

8 Opening the safety storage cabinet after a fire

After a fire, the safety storage cabinet must not be opened for at least 24 hours, and only with the utmost caution and by specialist employees.

WARNING!

Explosive vapour-air mixture

This may result in death or serious injury.

- Before opening the safety storage cabinet, remove all ignition sources within a 10 metre radius.
- Only open the safety storage cabinet with tools that do not cause any mechanically generated sparks.

WARNING!

Damaged safety storage cabinet due to fire or extinguishing agents

This may result in death or serious injury.

- Do not use safety storage cabinets that have been damaged by fire or extinguishing agents.

9 Maintenance

Check the safety storage cabinet for any externally visible damage or defects.

Always perform checks:

- After installation.
- Before commissioning.
- After changes.
- After maintenance.

The safety storage cabinet should also be inspected periodically at the following intervals.

Interval	Maintenance work	Personnel
Monthly	Closing of the doors ■ Open the door and inspect the closure.	Technical specialist employees
	Ventilation ■ Check the effectiveness of the ventilation with a woollen thread or with a smoke tube in the cabinet in front of the exhaust air duct at the ventilation grilles.	Technical specialist employees
	Seals ■ Check the sealing strips are seated properly in the door rebate and on the end faces of the doors. ■ If visible damage is found, replace the sealing strips immediately.	Technical specialist employees
	Labelling ■ Inspect the safety labels on the safety storage cabinet to ensure they are complete.	Technical specialist employees

Interval	Maintenance work	Personnel
Annually	Safety storage cabinet ■ Check of the entire safety storage cabinet	DÜPERTHAL service technicians

If faults occur, assist the technical customer service by providing the cabinet model, production and key number, along with a description of the fault.

To comply with national regulations, e.g. BetrSichV in Germany, all safety equipment must be inspected annually by a qualified person as defined in the valid national standards and regulations, e.g. the technical rules for industrial safety: TRBS 1203 in Germany.

10 Faults

Fault description	Cause	Remedy	Personnel
Doors do not close.	Safety storage cabinet is not aligned correctly.	Install the safety storage cabinet so it is horizontal. <i>See Chapter 6.4 'Check the alignment of the safety storage cabinet' on page 101</i>	Technical specialist employees
	Doors are held open by objects.	Do not wedge or hold doors open with any objects.	Technical specialist employees
	Safety storage cabinet is not correctly filled.	Make sure that containers in the safety storage cabinet are uniformly distributed.	Technical specialist employees
No extractor present.	Venting cut-off flaps closed, as closing mechanism has been triggered.	Replace the locking mechanism.	DÜPERTHAL service technicians
Doors do not move easily.	Moving parts, such as hinges, are dirty or corroded.	■ Remove rust. ■ Lubricate parts. ■ Remove corrosive substances from the safety storage cabinet. ■ Notify technical customer service.	Technical specialist employees

11 Spare parts and accessories

Only original parts from DÜPERTHAL are to be used for the safety storage cabinets.

- Intermediate shelves
- Shelf for side installation
- Storage shelves for gas cylinders
- Door handle
- Locking cylinder
- Venting connection socket
- Ventilators
- Exhaust air monitoring units

12 Disposal

CAUTION!

Dismantling the safety storage cabinet

Risk of injury due to improper dismantling of the safety storage cabinet.

- Ensure that the safety storage cabinet is only dismantled by specialist technical employees.

The safety storage cabinet can be completely dismantled by specialist technical employees.

Recycle the individual material components separately.

Comply with national and local disposal regulations.

To save resources, do not place parts of the safety storage cabinet or the whole cabinet in bulky or domestic waste.

13 Certificates

CE Declaration of Conformity

in accordance with Machinery Directive 2006/42/EC, Annex II A

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hereby declare that the following machine

Machine designation: Safety storage cabinet for the storage of flammable gases
Machine model: ECO standard
Machine size: M, L, XL, XXL
Machine type: 73-20BB60-xyz
Fire resistance class: type G90

Key:

Machine size vs. machine type		M	L	XL	XXL
Width BB		06	09	12	14
Colour x		0, 4, 5, x			
Definition of fittings and features yz		11, 21, yz			
Special marks (ss)		L (only if necessary)			

If -xyz is mentioned in the machine number, the order number will define the colour, fittings and features and special marks.

complies with all relevant requirements of Machinery Directive 2006/42/EC.

Institution responsible for review of QS system according to annex X:

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Additionally, the machine complies with the following harmonised and national standards and specifications:

Transposed harmonised standards: DIN EN ISO 12100:2011


Transposed national standards and technical specifications:

DIN EN 14470-2:2006
DIN EN 16121:2017
DIN EN 16122:2012

Authorised person for compilation of technical documents:

(signee)
Frank Backhaus / CE-authorised person

Oerlinghausen, 28/04/2021
(place, date)

(signee)
Franz-Josef Hagen / Managing director

Document: D00101443 Rev.: 01

CE Declaration of Conformity

in accordance with Machinery Directive 2006/42/EC, Annex II A

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hereby declare that the following machine

Machine designation: Safety storage cabinet for the storage of flammable gases

Machine model: SUPREME standard

Machine size: S, M, L, XL, XXL

Machine type: 79-HHBB60-xyz

Fire resistance class: type G90

Key:

Machine size vs. machine type		S	M	L	XL	XXL
Height	HH	13	20	20	20	20
Width	BB	06	06	09	12	14
Colour	x	0, 4, 5, x				
Definition of fittings and features	yz	21, yz				
Special marks	(ss)	L (only if necessary)				

If -xyz is mentioned in the machine number, the order number will define the colour, fittings and features and special marks.

complies with all relevant requirements of Machinery Directive 2006/42/EC.

Institution responsible for review of QS system according to annex X:

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TÜV SÜD Management Service GmbH