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Modelle

Gasflaschendepot

**70-200740-001, 70-200740-003, 70-201040-001
70-201040-003, 70-201340-001, 70-201340-003**

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Operating instructions

DÜPERTHAL - Gas cylinder depots according to TRG 280

Models: 70-200740-001, 70-200740-003, 70-201040-001
70-201040-003, 70-201340-001, 70-201340-003

Please enter the following details:

(this information is important in the event of enquiries)

Cabinet model: _____ Production number: _____

Key number: _____ Location: _____

Responsible person: _____

1. Manufacturer and technical details

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E-mail: info@dueperthal.com
www.dueperthal.com

1.2 Name / application

Gas cylinder depots according to TRG280 for discharging, initialisation and storage of gas cylinders..

1.3 Possible uses

The height of the cabinet provides the application of standard gas cylinders including necessary fittings. Gas cylinders of different height can be placed by the insert of shelves.

1.4 Technical data

Models	70-200740-001		70-200740-003		70-201040-001		70-201040-003		70-201340-001		70-201340-003	
Dimensions:	Outer	Inner	Outer	Inner	Outer	Inner	Outer	Inner	Outer	Inner	Outer	Inner
Width in approx. mm	650	600	650	600	950	900	950	900	1250	1200	1250	1200
Depth in approx. mm	400	300	400	300	400	300	400	300	400	300	400	300
Height in approx. mm	2050	1970	2050	1970	2050	1970	2050	1970	2050	1970	2050	1970
Total depth incl. overhang	500		500		500		500		500		500	
No. of doors	1		1		2		2		2		2	
Places 50 l bottles Bottle-Ø = 230 mm	2		2		3		3		4		4	
Weight in approx. kg	71		71		98		98		117		117	
Number of windows	0		1		0		2		0		2	

2. Requirements for the installation location

2.1 Floor

A level assembly is the assumption for a correct function of the doors. Floor unevenness can be adjusted by flakes or similar.

3. Installation

3.1 Venting

Natural venting by slashes in the door.

3.2 Earthing

Earthing connection on the right of the top.

4. Designation

Warning "gas cylinder" on each lockable wing door.

5. Combination of different types of gases

According to TRG 280 it is acceptable to set up different gases in one cabinet. We recommend to store gases, which react together, separately.

6. Transport

The gas cylinder depot has to be protected against upset during transport. Only use adequate means of transport. Please contact the manufacturer in case of doubts.

7. Interior fitting

7.1 Fittings for gas cylinder

There are 2 adjustable installation rails for the assembly of gas cylinder stations.

7.2 Cylinder holder

The cylinders are secured against falling over by retaining belts on the cylinder holder

7.3 Line connection

There are 2 boreholes, each Ø 22 mm, per bottle place in the top. Those are closed with plastic plugs.

7.4 Cabinet bottom

The bottom of the cabinet is made of 3,5 mm thick checkered plate.

8. Other accessories

8.1 Intermediate shelf

We recommend adjustable intermediate shelves for application of small gas bottles. The shelf can be hinged in the sidewise adjustment

gib by carrier. Load bearing capacity under uniformly distributed load is 150 kg.

8.2 Shelf for hinge

If 50 l bottles and smaller bottles shall be put together in one cabinet, we recommend shelves for side installation. The shelf can be hinged directly in the sidewise adjustment gib.

8.3 Supplementary equipment

- Intermediate shelf, bearing capacity max. 150 kg, under uniformly distributed load.
- Shelf for side installation for hinge.
- additional installation rails for different cylinder sizes

9. Disposal

We recommend adjustable intermediate shelves for application of small gas bottles. The shelf can be hinged in the sidewise adjustment gib by carrier. Load bearing capacity under uniformly distributed load is 150 kg.

10. Warranty

All damages, which can be ascribed to improper handling and assembly as well as corrosion damages, are excepted from warranty.

All dimension are approximate
Subject to technical amendments. Printing errors excepted.

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All rights reserved. Reprint prohibited, even as extracts. No part of this document may be reproduced or processed using electronic systems, copied or distributed in any form whatsoever (photocopy, microfilm or other method, not even for teaching purposes, without the prior consent of DÜPERTHAL.

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© Copyright 2001
DÜPERTHAL Sicherheitstechnik GmbH & Co. KG