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TÜV®

COMFORT line
Safety storage cabinets

Typ 90

Operating instructions

DÜPERTHAL Sicherheitstechnik GmbH & Co. KG
Frankenstrasse 3
D-63791 Karlstein
Germany
Telephone: +49 6188 9139-0
Fax: +49 6188 9139-121
Email: info@dueperthal.com
Internet: www.dueperthal.com

© 08/2018 DÜPERTHAL 71289, 0, en_GB

Table of contents

1	General information	59
1.1	Cabinet model information.....	59
1.2	Notes for reading.....	59
1.3	Type plate.....	59
1.4	Guarantee.....	60
2	Safety	61
2.1	Function of safety notices.....	61
2.2	Correct use.....	61
2.3	Reasonably foreseeable misuse.....	62
2.4	Obligations of the operator.....	62
2.5	Demands on employees.....	63
2.6	Residual risks.....	64
2.7	Handling stored goods.....	64
2.7.1	Storage of flammable liquids.....	64
2.7.2	Formation of explosive atmosphere.....	64
3	Technical specifications	65
3.1	General technical specifications.....	65
3.2	Technical specifications for equipment.....	65
3.3	Pressure drop/volumetric flow rate.....	68
3.4	Noise emissions.....	68
4	Structure and function	69
4.1	Construction.....	69
4.2	Earthing options.....	69
4.3	Ventilation.....	70
4.3.1	Air flow system in the cabinet interior.....	70
4.3.2	Exhaust air connection and feed opening.....	71
4.4	Doors.....	72
4.4.1	Door operation.....	72
4.4.2	Locking cylinder.....	72
4.5	Safety technology.....	73
4.5.1	Door closure in case of fire.....	73
4.5.2	Closure of venting cut-off flaps in case of fire.....	73
4.6	Interior fittings.....	74
4.6.1	Storage shelves.....	74
4.6.2	Bottom tray.....	75
4.6.3	Perforated sheet insert.....	75
5	Transport	76
5.1	Transport before commissioning.....	76
5.2	Transport after commissioning.....	77
6	Installation and commissioning	78
6.1	Requirements on the installation location.....	78
6.2	Attaching the adjustable feet.....	79
6.3	Align the safety storage cabinet.....	80
6.3.1	Alignment (with plinth panel).....	80
6.3.2	Check the alignment of the safety storage cabinet.....	81

7	Commissioning	82
7.1	Ventilating the safety storage cabinet.....	82
7.1.1	Connection to an exhaust air system.....	82
7.1.2	Requirements for ventilators.....	84
7.2	Earth the safety storage cabinet.....	85
7.3	Piston rod cylinder installation.....	86
7.4	Plinth panel and compressed air installation.....	87
8	Operation	89
8.1	Open the safety storage cabinet.....	89
8.2	Changing the height of the storage shelves.....	90
8.3	Fitting additional storage shelves.....	91
8.4	Checking and cleaning the bottom tray.....	91
8.5	Inspect the safety labels on the safety storage cabinet to ensure they are complete.....	92
8.6	Operating the safety storage cabinet without industrial ventilation.....	93
8.7	Pipe penetration (optional).....	94
9	Opening the safety storage cabinet after a fire	95
10	Maintenance	96
10.1	Inspect daily.....	97
10.2	Inspect monthly.....	97
10.3	Inspect annually.....	97
11	Faults	98
12	Spare parts and accessories	99
13	Disposal	100
14	Certificates	101

1 General information

1.1 Cabinet model information

Enter the following information and have it ready in the event of queries:

Cabinet model	Production number	Key number	Location	Responsible person

1.2 Notes for reading

The following symbols designate specific types of information.

This graphical system identifies information for easier and more effective working.

➔ Procedural step

⇒ Result of a procedural step

Link to another part of the document: 🔗 *Chapter 1.2 'Notes for reading' on page 59*

1.3 Type plate

The type plate is attached to the outside of the safety storage cabinet door.


Fig. 1: Type plate

- 1 Model
- 2 Serial number and year of manufacture
- 3 Maximum volume of individual containers
- 4 Maximum load per storage shelf

1.4 Guarantee

Liability for defects in the delivered products is defined in the Terms and Conditions of Delivery. We accept no liability whatsoever for damage caused by failure to comply with the operating instructions and conditions of use.

Only flawless products with the properties pledged by DÜPERTHAL leave our factory. If the information in these operating instructions is followed, correct operation is guaranteed.

Annual review requirement

As safety equipment, DÜPERTHAL safety storage cabinets are subject to an annual safety review recommended by the manufacturer. This inspection is to be carried out by competent personnel qualified by the manufacturer. Without the recommended annual inspection the customer loses the right to make guarantee claims against the manufacturer.

2 Safety

2.1 Function of safety notices

Safety notices warn against physical or material damage and provide information on how such damage can be avoided.

The following signal words identify the degree of danger and the extent of the risk.

DANGER!

The signal word refers to an immediately imminent threat which can cause death or serious injury.

WARNING!

The signal word refers to a potential threat which can cause death or serious injury.

CAUTION!

The signal word refers to an potential threat which can cause slight or mild injury.

The following safety signs are used in these instructions for use:

WARNING!

Warning of a hazard.

NOTICE!

This symbol identifies a situation that can lead to damage to the safety storage cabinet.

2.2 Correct use

The COMFORT line safety storage cabinet is type tested and classified as *Type 90* in compliance with 'EN14470-1' with a fire resistance of 90 minutes.

The safety storage cabinet is suitable for passive storage of flammable and liquid hazardous substances, in accordance with applicable national technical regulations, such as 'TRGS 510' in Germany.

Passive storage is defined as storage in tightly closed containers. The tightness of the container seal must be guaranteed.

Any other use shall be deemed improper.

2.3 Reasonably foreseeable misuse

DÜPERTHAL accepts no liability for damage arising from misuse.

The following constitute misuse:

- Storage of living organisms in the safety storage cabinet.
- Storage of food in the safety storage cabinet.
- Performance of laboratory work in conjunction with flammable liquids in the safety storage cabinet.
- Performance of transfer and filling work in the safety storage cabinet.
- Conversion measures and modifications to the safety storage cabinet without the knowledge and consent of DÜPERTHAL and the irrespective use of such safety storage cabinets.
- Any utilisation beyond the terms of correct use.

2.4 Obligations of the operator

The operator is obliged to comply with applicable legal regulations. These include:

- Issuing operating instructions.
- Carrying out risk assessments.
- Creating explosion protection documents.
- Complying with specified maintenance and inspection work as well as recurring checks.
- Carry out documented safety training of employees in the correct storage and removal of hazardous substances and preparations at least 1x per year.
- Carry out documented safety training of employees in the correct storage and removal of water-endangering substances and preparations at least 1x per year.
- Define and comply with the segregated storage of hazardous substances and preparations.
- Prohibit access by unauthorised persons to the hazardous substances stored in the safety storage cabinet.
- Draw up rules and regulations concerning smoking, handling naked flames, working with naked flames and spark-generating work in and at the safety storage cabinet.

Other obligations of the operator:

- Safeguard the function of safety-relevant equipment in and at the cabinet.
- Store hazardous substances so that a visual inspection of the bottom tray for extraneous substances is possible on working days.
- Install the safety storage cabinet so that the annual maintenance activities can be carried out without restriction.
- Immobilise, collect and correctly dispose of leaked liquids in accordance with accident prevention regulations.

2.5 Demands on employees

WARNING!

Employees who do not meet these requirements

This may result in death or serious injury.

- Activities are to be carried out by designated employees who meet the requirements.

These operating instructions set out the following employee activities:

- Technically qualified personnel
- Laboratory / warehouse employees
- DÜPERTHAL service technicians

Only people who meet the following requirements are approved as Technical specialist employees and Laboratory / warehouse employees:

- Of the minimum age legally required at the installation location.
- Meet applicable legal health requirements for the installation location with regard to the substances being stored.
- Carrying out documented safety training in correct storage and removal of hazardous substances and preparations at least 1x per year.
- Carrying out documented safety training in correct storage and removal of water-polluting substances and preparations at least 1x per year.
- Have been instructed in the function and operation of the safety storage cabinet.
- Safety training directly related to the activity in question.
- Basic knowledge in handling of flammable and liquid hazardous substances.

Technically qualified personnel

Technical specialist employees also fulfil other requirements in addition to the general requirements:

- In-depth knowledge and skills in relation to maintenance and repair.

Laboratory / warehouse employees

Laboratory and warehouse employees also fulfil other requirements in addition to the general requirements:

- Scientific professional training or equivalent professional experience which is recognised at the installation location.
- Training and instruction in the activities to be carried out.

DÜPERTHAL service technicians

DÜPERTHAL employees are specifically trained by DÜPERTHAL to carry out their activities.

2.6 Residual risks

Despite design-related measures, residual risks can arise even if the safety storage cabinet is used correctly.

In order to reduce health risks and to avoid dangerous situations, observe the safety instructions listed here and the safety information contained in later chapters of these instructions.

WARNING!

Employees who do not meet these requirements

The consequences are death or serious injury.

- Carry out activities by named employees who meet the requirements.

2.7 Handling stored goods

2.7.1 Storage of flammable liquids

Flammable gases are substances and preparations with a low flash point.

The flash point of a substance is the lowest temperature at which an ignitable vapour/air mixture can form above that substance.

Storage of these liquids is set out in *'Annex 3 of TRGS 510 - Storage of flammable liquids in safety storage cabinets'*.

2.7.2 Formation of explosive atmosphere

An explosive atmosphere is defined as a mixture of flammable substances with oxygen.

Explosiveness depends on the substance used, its flammability and mixture with oxygen. An effective ignition source is an ignition source which is able to trigger an ignition in the explosive atmosphere under consideration.

Safety measures for prevention of an explosive atmosphere and ignition hazards set out in *'TRBS 2152-3 Hazardous explosive atmosphere - Preventing ignition of hazardous explosive atmosphere'* must be implemented.

3 Technical specifications

3.1 General technical specifications


Fig. 2: General diagram of Type 90 safety storage cabinet

- | | |
|--------------------------|--------------------------------------------------------|
| 1 Exhaust air connection | 4 Bottom tray |
| 2 Labelling | 5 Plinth panel including pneumatic unit and foot pedal |
| 3 Standing surface | 6 Piston rod cylinder |

3.2 Technical specifications for equipment

Model size	Door wings	COMFORT line with classic folding door technology	
		Number of storage shelves*	Number of bottom trays
M	1	4 units	1 unit
XL	2		

* Other quantities available on request.





Fig. 3: Technical specifications for COMFORT line standard

- | | | | |
|-------|-------------------------------------------------------------------|----|----------------------------------|
| H1-H4 | Height 1 - Height 4 | T | Inside depth of standing surface |
| M | Load-bearing capacity of standing surface (uniformly distributed) | B | Inside width of standing surface |
| L | Maximum volume of largest individual container | LB | Clear width |
| | | V | Maximum collection volume |

	M	XL
B (mm)	481	1081
LB (mm)	346	866
T (mm)	445	445
M (kg)	75	75
L (l)	10	30
V (l)	11	33
Empty weight (kg)	300	500
Max. payload (kg)	360	360

Maximum container height (mm)			
			M, XL
2 storage shelves		Height 1	766
		Height 4	739
3 storage shelves		Height 1	510
		Height 2	510
		Height 4	483
4 storage shelves		Height 1	382
		Height 2	382
		Height 3	382
		Height 4	355

3.3 Pressure drop/volumetric flow rate


Fig. 4: COMFORT line Pressure drop and volumetric flow rate

- 1 Pressure drop ΔP (Pa)
- 2 Volumetric flow rate Q (m^3/h)
- 3 Q with ten-fold air exchange (see table).
- 4 Averaged pressure drop from all cabinet sizes.

Tab. 1: Q with ten-fold air exchange

Tenfold air exchange	Volumetric flow rate (m^3/h)
M	4.3
XL	9.4

3.4 Noise emissions

The noise emitted by the door closing is < 70 dB(A).

4 Structure and function

4.1 Construction

- Cabinet carcass and doors in multi-layer construction
- External panelling: Powder-coated sheet steel
- Wall construction: Multi-layer design
- Interior surfaces: Light grey-coated decor panels
- Safety technology elements for closure in case of fire: Brass, spring steel (1.410)

4.2 Earthing options


Fig. 5: Earthing options

1 Equipotential bonding saddle on the cabinet carcass

Earthing the safety storage cabinet prevents ignition hazards.

The interior fittings are conductively connected to one another by an equipotential bonding saddle or equipotential bonding screw on the cabinet carcass.

This allows earthing in accordance with national regulations, such as TRGS 727 in Germany.

4.3 Ventilation

The safety storage cabinets can be connected to a technical exhaust air system which ducts outside at a danger-free location.

In normal operation, the industrial ventilation of safety storage cabinets prevents the occurrence of a potentially explosive atmosphere inside the cabinet.

! NOTICE!

In a safety storage cabinet with industrial ventilation, a minimum hourly air exchange must take place that is at least 10 times the internal volume of the cabinet when the doors are closed.

If the safety storage cabinet is connected to a separate industrial ventilation system, the ventilator must be installed as close as possible to the venting pipe leading to the open air.

4.3.1 Air flow system in the cabinet interior


Fig. 6: Ventilation ducts open

- 1 Green inspection cut-outs in the open ventilation ducts

The layout of the ventilation ducts in the cabinet means that ventilation takes place directly above the bottom tray and is effective on every cabinet level.

The cabinet is also equipped with an inspection window for visual checking of the connection sockets.

Below the suspended ceiling, green inspection cut-outs identify the open ventilation ducts.

4.3.2 Exhaust air connection and feed opening


Fig. 7: Exhaust air system in cabinet sizes M-XL (view from above)

- 1 Exhaust air connection socket
- 2 Air supply opening
- 3 Distance: 100 mm
- 4 Nominal width (NW)
- 5 Distance: 150 mm
- 6 Distance identical on both sides.

The exhaust air and air supply openings are located on the safety storage cabinet's roof.

The scope of delivery includes an exhaust air connection with nominal diameter 110

Also included in the scope of delivery is an adapting reducer connection with a nominal diameter 75, which allows connection to a ventilation system.

4.4 Doors

4.4.1 Door operation

There is a pedal on the cabinet plinth. Lightly depressing the pedal opens the folding doors automatically.

The folding doors remain open for as long as you hold down the pedal with your foot.

As soon as your foot leaves the pedal, the folding doors close automatically after a maximum of 20 seconds (adjustable).

4.4.2 Locking cylinder

The door can be locked with the integrated locking cylinder. The key number is imprinted on the locking cylinder and on the keys supplied, e.g. A007. Locks be subsequently adjusted to the operator's requirements.

4.5 Safety technology

4.5.1 Door closure in case of fire

At a temperature of approx. 50°C, the open doors are closed by the safety technology.

4.5.2 Closure of venting cut-off flaps in case of fire


Fig. 8: Closed ventilation ducts

1 Red inspection cut-outs on the closed ventilation ducts

At a temperature of 70°C, the open venting cut-off flaps are closed by the safety technology.

Below the suspended ceiling, red inspection cut-outs identify the closed ventilation ducts.

4.6 Interior fittings

4.6.1 Storage shelves


Fig. 9: Storage shelf

- 1 Fixing screw
- 2 Earthing screw
- 3 Storage shelf
- 4 Adjustable support
- 5 Support strip

The COMFORT line standard safety storage cabinet contains storage shelves that are distributed uniformly over the cabinet interior height.

The heights of the storage shelves can be changed.

The highest shelf must not be more than 1.75 m above the floor.

4.6.2 Bottom tray

Fig. 10: Bottom tray

1 Bottom tray

The function of the bottom tray in the floor area of the safety storage cabinet is to collect leaking substances in the cabinet interior. As a standard feature it cannot be used as additional standing surface.

4.6.3 Perforated sheet insert

The bottom tray in the COMFORT line can be used as an additional standing surface in accordance with the regulations if an optional perforated sheet is inserted.

! NOTICE!

Perforated sheet insert

The marked collection volume for the largest individual container that can be placed in the cabinet is only available if the perforated sheet that is also integrated is used as a standing surface.

5 Transport

5.1 Transport before commissioning

The safety storage cabinet is packaged for transport and is protected by transport restraints.

The packaging, transport restraints and safety transport skids protect the cabinet from damage during transport.

WARNING!

Risk of crushing from toppling safety storage cabinet if the centre of gravity is unbalanced

The consequences are death or serious injury.

- The safety storage cabinet must be transported by at least two employees.
- Only transport the safety storage cabinet standing and unladen.
- Only drive under the safety storage cabinet using suitable transport means.
- When driving underneath it, only pick up the safety storage cabinet in the centre from the front, rear or from the side.


Fig. 11: Transporting the safety storage cabinet

- 1 Pick up centrally from the front
- 2 Pick up centrally from the side

NOTICE!

Handling the transport restraints

Incorrect handling can damage the safety transport skids and the safety storage cabinet.

- Transport restraints and safety transport skids should only be removed at the installation location.
- Replace the safety transport skids after transport to the installation location with the enclosed adjustable feet.

5.2 Transport after commissioning

! NOTICE!

Tilting the safety storage cabinet

Tilting causes damage to the installed safety technology.

- Secure the doors of the safety storage cabinet using the transport restraints used for shipping in the door gaps.

! NOTICE!

Tipping the safety storage cabinet over during transport

Damage to the safety storage cabinet caused by incorrect handling.

- Only pick up the safety storage cabinet from the side or back wall.
- Only pick up the safety storage cabinet using special and suitable transport or lifting equipment.
- Pick up the safety storage cabinet once it is securely lashed and is not at risk of slipping.
- Do not damage the adjustable feet during transport.

! NOTICE!

Driving under the safety storage cabinet

Damage to the safety storage cabinet caused by incorrect handling.

- Unscrew the piston rod cylinder including mounting bracket.
- Remove the plinth panel.
- Disconnect the pneumatics.

Please also note the safety information detailed in [Chapter 5.1 'Transport before commissioning'](#) on page 76.

6 Installation and commissioning

Install the safety storage cabinet so that the annual maintenance activities can be carried out without restriction.

6.1 Requirements on the installation location

The safety storage cabinet is approved for installation in a building.

Installation in outdoor areas with or without roof covering is not permitted.

The surface must be able to bear the weight of the safety storage cabinet when fully loaded.

The surface must be horizontal in order to guarantee problem-free functioning of the safety storage cabinet.

The load-bearing capacity and stability of the surface must be assured both in normal situations and in the event of a fire.

The safety storage cabinet must be supplied with compressed air to enable operation. For this purpose, a fixed compressed air line with upstream maintenance unit must be connected to the safety storage cabinet. The compressed air line with upstream maintenance unit has an operating pressure of 6.5 bar. Compressed air quality conforming to ISO 8573-1:2010.


Fig. 12: Requirements for the installation location

Consider the following in relation to the installation location:

- Protect the safety storage cabinet from vehicles, shocks, etc.
- Do not install the safety storage cabinet near sources of heat.
- Protect the safety storage cabinet against moisture (service class 1). Use at > 70 % relative humidity for a few weeks per year in structures closed on all sides and heated.
- The operating temperature must be between - 5 °C and + 40 °C.

6.2 Attaching the adjustable feet


Fig. 13: Attaching the adjustable feet

- 1 Transport equipment for the safety storage cabinet
- 2 Safety storage cabinet
- 3 Safety transport skids
- 4 Attach the safety transport skids (4x screw, size 19 mm)
- 5 Adjustable foot

Personal:

- Technical specialist employees

1. ➤ Transport the cabinet to its place of use.
2. ➤ Remove packaging.
3. ➤ Remove the adjustable feet from the cabinet roof.
4. ➤ Raise the cabinet and loosen the screws for the safety transport skids (A).
⇒ The safety transport skids can be removed (B).
5. ➤ Screw the adjustable feet completely into the bottom of the cabinet from underneath (C-D).
6. ➤ Position the cabinet and set it down carefully.

6.3 Align the safety storage cabinet

The alignment procedure described below is used for precision alignment. Remedy any major floor unevenness of more than 1.5 mm on site.

6.3.1 Alignment (with plinth panel)

The safety storage cabinet with plinth panel has long adjustable feet in the corners of the base area. The adjustable feet are covered and protected by the plinth panel.

Alignment with plinth panel

Personal:

- Technical specialist employees

1. ▶ Lift the cabinet using suitable lifting equipment.
2. ▶ Screw the adjustable feet in or out by hand.
3. ▶ Set the safety storage cabinet back down.

6.3.2 Check the alignment of the safety storage cabinet


Fig. 14: Checking alignment

Correct alignment of the safety storage cabinet:

- when the doors are closed, the door gaps are of equal width.
- With all door elements, the central gap and ceiling gap form an even "T".

7 Commissioning

7.1 Ventilating the safety storage cabinet

7.1.1 Connection to an exhaust air system

! NOTICE!

Keep the exhaust air connection socket and air supply opening on the cabinet roof accessible and open at all times for maintenance and servicing purposes.

! NOTICE!

In the event of unplanned use of the cabinet roof as storage area, analyse the function of the safety elements and the weight load and include these in the risk assessment by the responsible operator.

! NOTICE!

Use of the cabinet roof is not planned in the fire safety concept of the European standard and has also not been tested!

The installation of industrial ventilation and/or connection to an existing exhaust air system must be carried out by a ventilation company and is not included in the DÜPERTHAL range of products and services.


Fig. 15: Connection to the exhaust air system

- 1 Exhaust air connection socket
- 2 Air supply opening
- 3 Thermocouples

Connection to the exhaust air system:

Personal:

- Technical specialist employees

1. ➤ Connect the exhaust air line to the exhaust air connection socket and secure with a collar.
2. ➤ After installing the safety storage cabinet, check the connection to the exhaust air system with smoke tubes.

7.1.2 Requirements for ventilators

If the safety storage cabinet is connected to a separate industrial ventilation system, the ventilator must be installed as close as possible to the venting pipe leading to the open air.

Only use ventilators that comply with the defined technical specifications.

In areas not at risk of explosion, this means:

- Inside: Zone 2
- Outside: area not at risk of explosion.

WARNING!

Use of incorrect ventilators

The ventilation function is not guaranteed if the ventilators used do not meet the necessary requirements.

This can create an explosive atmosphere inside the safety storage cabinet.

This can lead to death or serious injury.

- Only use ventilators that comply with the defined technical specifications.

NOTICE!

A ventilator must meet the requirements in the applicable national regulations for Zone 2 equipment.

7.2 Earth the safety storage cabinet


Fig. 16: Earthing connection

Earthing connection

Personal:

- Technical specialist employees

1. ➤ Place the screwdriver under the equipotential bonding saddle (A).
2. ➤ Bend up the equipotential bonding saddle.

7.3 Piston rod cylinder installation


Fig. 17: Piston rod cylinder

- 1 Door mounting bracket
- 2 Hexagon nuts and washers
- 3 Fork head
- 4 Piston rod cylinder
- 5 Piston rod cylinder mounting bracket
- 6 Hexagon head screw and serrated washers


Fig. 18: Installing the piston rod cylinder

Personal:

- Technical specialist employees

1. ➤ Screw the mounting bracket for the piston rod cylinder to the outer base of the cabinet using the hexagon head screws (M6x10) and serrated washers (A).

The position is defined by the M6 threaded holes in the outer base of the cabinet.

2. ➤ Screw the fork head of the piston rod cylinder to the threaded rod of the door mounting bracket using the hexagon nuts (M6) and washers (B).

7.4 Plinth panel and compressed air installation


Fig. 19: Plinth panel and compressed air

- 1 Plinth panel side section
- 2 Coupling socket (not included in scope of delivery)
- 3 Coupling plug
- 4 Throttle valve
- 5 Plinth panel front section (including pneumatic unit)


Fig. 20: Plinth panel and compressed air installation

Personal:

- Technical specialist employees

1. ▶ Connect the side sections of the plinth panel to the front panel (A-B).
2. ▶ Connect the pneumatic unit on the front section of the plinth panel to the piston rod cylinder (C).

Observe the markings on the hoses and standard cylinders.

3. ▶ Connect the external compressed air supply to the coupling plug using a size 4 coupling socket (D)
↳ Chapter 6.1 'Requirements on the installation location' on page 78.
4. ▶ Use the throttle valves to adjust the speed of the door closing and opening system.
5. ▶ Slide the three-part plinth panel under the cabinet from the front (E-F).

8 Operation

8.1 Open the safety storage cabinet


Fig. 21: Door operation

Door operation with foot pedal

Personal:

- Technical specialist employees
- Laboratory and warehouse employees

➔ Open the safety storage cabinet with your foot by lightly depressing the pedal on the cabinet plinth.

⇒ The folding door is opened fully to the defined end position and is held open in that position for as long as the foot pedal is pressed.

! NOTICE!

Keep the doors of the safety storage cabinet closed and close them again after each work process. Do not hold the doors of the safety storage cabinet open with a wooden wedge or similar objects. This blocks the safety technology!

8.2 Changing the height of the storage shelves


Fig. 22: Changing the height of the storage shelf

- 1 Fixing screw
- 2 Earthing screw
- 3 storage shelf
- 4 Adjustable support
- 5 Support strip

Changing the height of the storage shelf

Personal:

- Technical specialist employees

1. ➤ Loosen the earthing screw.
2. ➤ Loosen the fixing screw.
3. ➤ Remove the storage shelf.
4. ➤ Move the adjustable support into the support strips.
5. ➤ Insert the storage shelf.
6. ➤ Tighten the fastening screws.
7. ➤ Tighten the earthing screw.

8.3 Fitting additional storage shelves

Montage zusätzlicher Lagerböden: ↗ *Chapter 8.2 'Changing the height of the storage shelves' on page 90*

Personal:

- Technical specialist employees

8.4 Checking and cleaning the bottom tray

Personal:

- Laboratory and warehouse employees

➔ Carry out a visual inspection for extraneous substances on working days.

Storing objects of any kind in the bottom tray reduces the marked collection volume defined for the safety storage cabinet.

! NOTICE!

Store hazardous substances so that a visual inspection of the bottom tray for extraneous substances is possible on working days.

! NOTICE!

Remove any leaked liquids from the bottom tray straight away and dispose of them correctly.

! NOTICE!

Following dismantling for cleaning purposes, re-connect the bottom tray to the potential connection.

8.5 Inspect the safety labels on the safety storage cabinet to ensure they are complete

*Only clean cabinet stickers with conventional glass cleaners.
Other cleaners may damage them.*

The following must be attached to the front of the safety storage cabinet and must be clearly visible:

- The instruction "Close the door"
- Prohibited action sign or plaque bearing the words "Unauthorised use prohibited"
- Fire resistance in minutes (e.g. Type 90)
- Name or trademark of the manufacturer
- Serial number and year of manufacture
- Information on the largest individual container volume that can be stored
- Information on the shelves' maximum load capacity

Furthermore, the following signs must be attached to the front of the safety storage cabinet and must be clearly visible:

Tab. 2: Prohibited action signs

Symbol	Meaning	Standard
	P003: No naked flames; fire, open ignition source and smoking prohibited	DIN EN ISO 7010

Tab. 3: Warning signs

Symbol	Meaning	Standard
	W021: Warning: Flammable substances	DIN EN ISO 7010

Tab. 4: Mandatory action signs

Symbol	Meaning	Standard
	M002: Follow the instructions for use	DIN EN ISO 7010
	M004: Wear eye protection	DIN EN ISO 7010
	M009: Wear hand protection	DIN EN ISO 7010

8.6 Operating the safety storage cabinet without industrial ventilation

The safety storage cabinet can be operated in accordance with national regulations without industrial ventilation. Labelling through notice!

If an explosive atmosphere is to be expected, measures must be taken in compliance with national regulations and an explosion protection document must be drawn up.

! NOTICE!

- Determine the hazardous area in accordance with the hazardous substance regulations under all possible operating conditions and label it clearly, visibly and permanently.
- It is prohibited to handle naked flames or to smoke in areas at risk of explosion.
- Do not use any spark-generating tools in areas at risk of explosion.
- Ensure that no electrostatic charges can form in areas at risk of explosion and that no equipment is operated that has surface temperatures above the ignition temperatures of the stored, flammable liquids.
- Only operate electrical equipment in areas at risk of explosion if they comply with the requirements of European standards, such as DIN EN VDE.

Also note the following points from national regulations:

- Safety storage cabinets without industrial ventilation should, in case of fire, protect the stored goods from impermissible heating and ignition of possibly occurring explosive mixtures.
- There must be no ignition sources inside the safety storage cabinet. If this cannot be fully prevented, depending on the risk assessment, measures must be taken to avoid ignition sources. As a minimum requirement, the measures taken must be in accordance with those designated for Zone 2 areas as defined in the applicable national technical regulations, such as TRGS 727 in Germany.
- Safety storage cabinets without industrial ventilation must be earthed with equipotential bonding.

Clearly label the Zone 2 potentially explosive area according to national accident prevention regulations with the following warning symbol:

The size of all symbols and notices should be appropriate for the size of the safety storage cabinet.

8.7 Pipe penetration (optional)

Tested penetrations are optionally available for the safety storage cabinet with Type 90 classification. The penetrations are attached to the safety storage cabinets from the outside and can be equipped with holes of up to Ø 28 mm.

According to type testing, these have no negative effect whatsoever on the cabinet's fire resistance if used correctly. Unused, open holes in the pipe penetrations must be sealed. Contact DÜPERTHAL.

! NOTICE!

Later attachment of pipe penetrations

Incorrect attachment in the wrong location can cause damage to the safety storage cabinet.

- Attachment only at the locations provided.
- Holes only based on defined penetration profiles (see separate mounting instructions for each cabinet type).

! NOTICE!

Routing and using pipe penetrations

The routing and use of pipes, cables and tubes are the responsibility of the operator.

Incorrect handling can cause damage to the safety storage cabinet.

- Carry out a separate risk assessment of the overall setup.

9 Opening the safety storage cabinet after a fire

Only open the safety storage cabinet after 24 hours have elapsed, and only with extreme caution.

 **WARNING!**

Explosive vapour-air mixture

The consequences are death or serious injury.

- Before opening the safety storage cabinet, remove all ignition sources within a 10 metre radius.
- Open the safety storage cabinet only using spark-free tools.

Do not use safety storage cabinets that have been damaged by fire.

10 Maintenance

Check the safety storage cabinet for any externally visible damage or defects.

Always perform checks:

- After installation.
- Before commissioning.
- After changes.
- After maintenance.

The safety storage cabinet should also be inspected periodically at the following intervals.

Inspect daily

10.1 Inspect daily

Interval	Maintenance work	Personnel
Daily	Bottom tray and pull-out trays <ul style="list-style-type: none"> ■ Check in accordance with regulations governing water legislation ■ Collect and properly dispose of leaked liquids immediately. 	Laboratory and warehouse employees

10.2 Inspect monthly

Interval	Maintenance work	Personnel
Monthly	Closing of the doors <ul style="list-style-type: none"> ■ Open the door and inspect the closure. 	Technical specialist employees
	Ventilation <ul style="list-style-type: none"> ■ Check the effectiveness of the ventilation with a woollen thread or with a smoke tube in the cabinet in front of the left-hand exhaust air duct at the ventilation grilles. ■ Remove contamination on the flame screen. 	Technical specialist employees
	Seals <ul style="list-style-type: none"> ■ Check the sealing strips are seated properly in the carcass frame and on the end faces of the doors. ■ If visible damage is found, replace the sealing strips immediately. 	Technical specialist employees

10.3 Inspect annually

Interval	Maintenance work	Personnel
Annually	Safety storage cabinet <ul style="list-style-type: none"> ■ Check of the entire safety storage cabinet 	DÜPERTHAL service technicians

If faults should occur, assist the technical customer service by providing the cabinet model, production and key number, along with a description of the fault.

11 Faults

Fault description	Cause	Remedy	Personnel
Doors do not open	Safety storage cabinet is not aligned correctly.	Install the safety storage cabinet so it is horizontal. ↪ <i>Chapter 6.3.2 'Check the alignment of the safety storage cabinet' on page 81</i>	Technical specialist employees
	Cabinet is not connected to the compressed air supply	Connect the cabinet to the compressed air supply	Technical specialist employees
	Throttle valve for pneumatic unit is closed	Open the throttle valve	Technical specialist employees
	Doors are locked	Unlock the doors	Technical specialist employees
No extractor present.	Venting cut-off flaps closed, as closing mechanism has been triggered.	Replace the locking mechanism.	DÜPERTHAL service technicians
Doors do not move easily.	Air exchange in cabinet is too high. Moving parts, such as hinges, are dirty or corroded.	<ul style="list-style-type: none"> ■ Remove rust. ■ Lubricate parts. ■ Remove corrosive substances from the safety storage cabinet. ■ Notify technical customer service. 	Technical specialist employees

12 Spare parts and accessories

- Storage shelves
- Pull-out shelves
- Bottom tray
- PP insert
- Anti-slip mat made from rubber
- Door handle
- Perforated sheet insert
- Plinth panels
- Venting connection socket
- Ventilators
- Exhaust air monitoring units

13 Disposal

The safety storage cabinet can be completely dismantled.

Recycle the individual fractions, e.g. metal, insulation panels, door closer, etc. separately.

Comply with national and local disposal regulations.

To save resources, do not place parts of the safety storage cabinet or the whole cabinet in bulky or domestic waste.

14 Certificates

CE Declaration of Conformity

in accordance with Machinery Directive 2006/42/EC, Annex II A

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hereby declare that the following machine:

Machine designation: Safety cabinet type 90 for passive storage of flammable liquids
 Machine model: COMFORT
 Machine size: M, XL
 Machine type: 29-20BB66-xyz(ss)
 Key:

Model size vs. Machine type		M	XL
Width	BB	06	12
Colour	x	0, 4, 5, x	
Definition of fittings and features	yz	30, yz	30, 32, yz
Special marks	(ss)	A, D, L, R, S, T, U, V (only if necessary)	

If -xyz is mentioned in the machine number, the order number will define the colour, fittings and features and special marks.

complies with all relevant requirements of Machinery Directive 2006/42/EC.

Institution responsible for review of QS system according to annex X:

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
Additionally, the machine complies with the following harmonised and national standards and specifications:

Transposed harmonised standards: DIN EN ISO 12100:2011

Transposed national standards and technical specifications:
 DIN EN 14470-1:2004
 DIN EN 14727:2006

Authorised person for compilation of technical documents:

(signee)
 Frank Backhaus / CE-authorised person

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COMFORT line, 29-20BBTu-xyz(ss)
PREMIUM line, 29-20BBTu-xyz(ss)

Parameters:

Key code Art.-No.
BB width, stands for numbers: 06, 12,
T depth, stands for numbers: 6, 7
u door technique, stands for numbers: 2, 4, 6
x colour, stands for numbers: 0, 4, 5, 6, 7
yz equipment, stands for numbers: 00-99
ss special sign, stands for: L, R, B, U, A

Class of fire resistance: FWF 90
The article number defines exactly the material and the size of the product. A full description of the specifications is provided in the test report.


Tested according to:

DIN EN 14470-1:2004
DIN EN 14727:2006
EK5/AK4 09-10:2009
TRGS 510:2015 Anlage 3
AfPS GS 2014:01 PAK 3.1

The product meets the safety and health requirements of the German Product Safety Act section 20 to 22 ProdSG. The certification marks shown above can be affixed on the product. It is not permitted to alter the certification marks in any way. In addition the certificate holder must not transfer the certificate to third parties. This certificate is valid until the listed date, unless it is cancelled earlier. See also notes overleaf.

Test report no.: 713086011

Valid until: 2022-05-22

Date, 2017-05-29

(Horst Kristen)

Page 1 of 1