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Dokument: D00123770 Rev.: 00

UTS cabinet series, BENCH cabinet series
Safety storage cabinets

Typ 90

Operating instructions

DÜPERTHAL Sicherheitstechnik GmbH & Co.
KG Frankenstrasse 3
D-63791 Karlstein
Germany
Telephone: +49 6188 9139-0
Fax: +49 6188 9139-121
Email: info@dueperthal.com
Internet: www.dueperthal.com

© 07/2023 DÜPERTHAL 55213, 5, en_GB

Table of contents

1	General information	55
1.1	Notes for reading.....	55
1.2	Type plate.....	55
2	Safety	56
2.1	Function of safety notices.....	56
2.2	Correct use.....	56
2.3	Misuse.....	56
2.4	Obligations of the operator.....	58
2.5	Demands on employees.....	58
2.6	Stored goods.....	58
2.7	Hazardous areas and their labelling.....	59
3	Technical specifications	61
3.1	General data.....	61
3.2	Dimensions and equipment.....	62
3.3	Pressure drop during ventilation.....	68
4	Structure and function	69
4.1	Construction.....	69
4.2	Earthing options.....	69
4.3	Exhaust air connection and feed opening.....	69
4.4	Doors.....	70
4.5	Safety technology.....	71
4.6	Interior fittings.....	72
4.7	Additional pull-out level (optional).....	74
4.8	ToeKick base (optional).....	75
4.9	Roller set (optional).....	75
4.10	Storage shelves (optional).....	75
4.11	Pipe penetration (optional).....	76
5	Transport	77
6	Installation and commissioning	78
6.1	Requirements on the installation location.....	78
6.2	Align the safety storage cabinet.....	78
6.3	Installing the safety storage cabinet with roller set.....	79
6.4	Check the alignment of the safety storage cabinet.....	80
6.5	Venting the safety storage cabinet.....	81
6.6	Earth the safety storage cabinet.....	83
7	Operation	84
7.1	Opening the safety storage cabinet with the door handle.....	84
7.2	Opening a safety storage cabinet with push-to-open door operation.....	85
7.3	Pull-out tray.....	85
7.4	Check the bottom tray or pull-out tray.....	86
8	Opening the safety storage cabinet after a fire	87
9	Maintenance	88
10	Faults	90
11	Spare parts and accessories	91

12	Disposal.....	92
13	Certificates.....	93

1 General information

1.1 Notes for reading

The following symbols designate specific types of information.

Tab. 1: Explanation of symbol

Symbol	Type of information
	Information for easier and more effective working
	Procedural step
	Result of a procedural step
	Link to another part of the document

1.2 Type plate

The type plate is attached to the outside of the safety storage cabinet door.

Model

Serial number and year of manufacture
Reference number

Maximum volume of individual containers

Maximum load per storage shelf

Fig. 1: Type plate

2 Safety

2.1 Function of safety notices

Safety notices warn against physical or material damage and provide information on how such damage can be avoided.

The following signal words identify the degree of danger and the extent of the risk.

WARNING!

The signal word 'WARNING' refers to a potential hazard which could result in death or serious injury.

CAUTION!

The signal word 'CAUTION' refers to a potential hazard which could result in slight or minor injury.

NOTICE!

The signal word 'NOTE' indicates a situation that can lead to damage to the safety storage cabinet.

2.2 Correct use

Observe the safety instructions in these operating instructions to reduce health risks and avoid dangerous situations.

Any use that is not correct use as defined in these operating instructions involves a risk of accidents and a lack of fire protection.

The UTS and BENCH safety storage cabinet series are type tested and classified as *Type 90* in compliance with 'EN14470-1' with a fire resistance of 90 minutes.

The safety storage cabinet is to be used for passive storage of flammable liquids in working spaces.

Passive storage is defined as exclusively storage without work such as filling, mixing or transfer in tightly sealed containers. The tightness of the container seal must be guaranteed.

2.3 Misuse

Any use that goes beyond the specified correct use is considered to be misuse.

DÜPERTHAL accepts no liability for damage arising from misuse.

In addition, the following safety notices must be observed:

WARNING!

Storage of living organisms in the safety storage cabinet

Living organisms can come into contact with the stored hazardous substances.

This may result in death or serious injury.

- Use the safety storage cabinet exclusively for storage of flammable liquids.

 **WARNING!**

Storage of food in the safety storage cabinet

Food can come into contact with the stored hazardous substances.

This may result in death or serious injury.

- Use the safety storage cabinet exclusively for storage of flammable liquids.

 **WARNING!**

Transfer, filling and laboratory work in the safety storage cabinet

Inhalation of hazardous vapours can lead to life-threatening injuries to the respiratory system.

This may result in death or serious injury.

- Store flammable liquids in the safety storage cabinet only in closed containers.
- Do not perform any laboratory work in the safety storage cabinet.
- Do not perform any transfer or filling work in the safety storage cabinet.

 **WARNING!**

Storage of hazardous substances together

Risk of uncontrolled chemical reactions.

This may result in death or serious injury.

- Only store substances and preparations that are permitted to be stored together in the safety storage cabinet.

 **WARNING!**

Objects on the cabinet roof

In case of fire, objects on the cabinet roof may impair the functioning of the safety storage cabinet or affect its stability.

This may lead to death or serious injury.

- Do not store any objects on the cabinet roof.

 **WARNING!**

Spilled liquids

Inhalation of hazardous vapours from spilled liquids can lead to life-threatening injuries to the respiratory system.

Spilled liquids can lead to painful skin reactions.

- Collect and properly dispose of spilled liquids immediately in accordance with accident prevention regulations.

! NOTICE!**Alteration and modifications**

Do not alter or modify the safety storage cabinet.

This can lead to a lack of fire protection.

- If alteration or modification of the safety storage cabinet is required, contact DÜPERTHAL.

2.4 Obligations of the operator

The operator is obliged to comply with applicable legal regulations. This includes:

- Issuing operating instructions.
- Carrying out risk assessments.
- Creating explosion protection documents.
- Specifying activities by designated employees.

2.5 Demands on employees

! WARNING!**Employees who do not meet these requirements**

This may result in death or serious injury.

- Designate employees who meet the requirements to carry out activities.

These operating instructions set out the following employee activities:

- Specialist technical employees
- DÜPERTHAL service technicians

Only people who have been trained by the operator in use of the safety storage cabinet and handling of the stored goods are approved as specialist technical employees.

DÜPERTHAL service technicians

DÜPERTHAL employees are specifically trained by DÜPERTHAL to carry out their activities.

2.6 Stored goods

Storage, handling and use of the stored goods must comply with the applicable national standards and regulations, e.g. 'TRGS 510' in Germany.

2.7 Hazardous areas and their labelling

The following must be attached to the front of the safety storage cabinet and must be clearly visible:

- The instruction "Close the door"
- Fire resistance in minutes (e.g. 'type 90')
- Name or trademark of the manufacturer
- Serial number and year of manufacture
- Information on the largest individual container volume that can be stored
- Information on the shelves' maximum load capacity

Furthermore, the following signs must be attached to the front of the safety storage cabinet and must be clearly visible:

Tab. 2: Prohibited action sign

Symbol	Meaning	Standard
	P003: No naked flames; fire, open ignition source and smoking prohibited	DIN EN ISO 7010

Tab. 3: Warning sign

Symbol	Meaning	Standard
	W021: Warning: Flammable materials	DIN EN ISO 7010

Tab. 4: Mandatory action signs

Symbol	Meaning	Standard
	M002: Read operator's manual	DIN EN ISO 7010
	M004: Wear eye protection	DIN EN ISO 7010

Symbol	Meaning	Standard
	M009: Wear hand protection	DIN EN ISO 7010

3 Technical specifications

3.1 General data


Fig. 2: General diagram of Type 90 safety storage cabinet

- 1 Type plate
- 2 Exhaust air connection
- 3 Bottom tray or pull-out tray with standing surface
- 4 Partition wall with split storage areas

3.2 Dimensions and equipment

3.2.1 UTS standard cabinet series (bottom tray)


Fig. 3: UTS standard cabinet series with bottom tray dimensions

- H Cabinet height (without base)
- HS Base height
- T Cabinet depth
- T1 Max. cabinet depth with open doors
- B Cabinet width
- LB Clear width
- M Standing surface load-bearing capacity
- H1 Height between perforated sheet insert and cabinet roof
- TS Standing surface depth
- BS Standing surface width
- L Maximum volume of the largest individual container
- V Max. collection volume

Tab. 5: UTS standard cabinet series S-XL dimensions

	S	M	L	XL
H (mm)	596	596	596	596
H1 (mm)	380	380	380	380
HS (mm)*	20-135	20-135	20-135	20-135
B (mm)	601	889	1100	1400
BS (mm)	480	767	980	767 471**
LB (mm)	476	755	967	708 472**
T (mm)	594	594	594	594
T1 (mm)	1170	1013	1119	1119
TS (mm)	424	424	424	424
L (l)	12	19	24	19 12**
V (l)	13.2	21.3	27.2	21.3 13.2**
M (kg)	60	60	60	60
Empty weight (kg)	110	150	180	225
Max. payload (kg)	60	60	60	120

Tab. 6: UTS standard cabinet series S-5-XL-5 dimensions

	ST	MT	LT	XLT
H (mm)	596	596	596	596
H1 (mm)	380	380	380	380
HS (mm)*	20-135	20-135	20-135	20-135
B (mm)	601	889	1100	1400
BS (mm)	480	767	980	767 471**
LB (mm)	476	755	967	708 472**
T (mm)	494	494	494	494
T1 (mm)	1070	913	1019	1019
TS (mm)	324	324	324	324
L (l)	9	14	19	14 9**
V (l)	10.4	16.4	20.9	16.4 10.4**
M (kg)	60	60	60	60
Empty weight (kg)	100	135	165	205
Max. payload (kg)	60	60	60	120

* Cabinets with no base are equipped with a 3 mm thick bottom plate.

** Safety storage cabinet with two (separate) storage areas: left storage area | right storage area.

3.2.2 UTS pro cabinet series (pull-out tray)


Fig. 4: UTS pro cabinet series with pull-out tray dimensions

- H Cabinet height (without base)
- HS Base height
- T Cabinet depth
- T1 Max. cabinet depth with open doors
- B Cabinet width
- LB Clear width
- M Standing surface load-bearing capacity
- H1 Height between perforated sheet insert and cabinet roof
- TS Standing surface depth
- BS Standing surface width
- L Maximum volume of the largest individual container
- V Max. collection volume

Tab. 7: UTS pro cabinet series S-XL dimensions

	S	M	L	LD	XL
H (mm)	596	596	596	596	596
H1 (mm)	340	340	340	340	340
HS (mm)*	20-135	20-135	20-135	20-135	20-135
B (mm)	601	889	1100	1100	1400
BS (mm)	403	653	868	392 392**	607 392**
LB (mm)	476	755	967	472 472**	708 472**
T (mm)	594	594	594	594	594
T1 (mm)	1170	1013	1119	1119	1119
TS (mm)	418	418	418	418	418
L (l)	10	16	21	10 10**	15 10**
V (l)	11.3	18.1	23.9	11 11**	16.8 11**
M (kg)	45	45	45	45	45
Empty weight (kg)	115	155	185	195	235
Max. payload (kg)	45	45	45	90	90

Tab. 8: UTS pro cabinet series S-5-XL-5 dimensions

	S-5	M-5	L-5	LD-5	XL-5
H (mm)	596	596	596	596	596
H1 (mm)	340	340	340	340	340
HS (mm)*	20-135	20-135	20-135	20-135	20-135
B (mm)	601	889	1100	1100	1400
BS (mm)	403	653	868	392 392**	607 392**
LB (mm)	476	755	967	472 472**	708 472**
T (mm)	494	494	494	494	494
T1 (mm)	1070	913	1019	1019	1019
TS (mm)	318	318	318	318	318
L (l)	8	12	16	7 7**	11 7**
V (l)	8.9	13.9	18.3	8.4 8.4**	12.9 8.4**
M (kg)	40	40	40	40	40
Empty weight (kg)	100	135	165	170	215
Max. payload (kg)	40	40	40	80	80

* Cabinets with no base are equipped with a 3 mm thick bottom plate.

** Safety storage cabinet with two (separate) storage areas: left storage area | right storage area.

3.2.3 BENCH cabinet series


Fig. 5: BENCH cabinet series dimensions

- H Cabinet height (without base)
- HS Base height
- T Cabinet depth
- T1 Max. cabinet depth with open doors
- B Cabinet width
- LB Clear width
- M Standing surface load-bearing capacity
- H1 Height between perforated sheet insert and cabinet roof
- TS Standing surface depth
- BS Standing surface width
- L Maximum volume of the largest individual container
- V Max. collection volume

Tab. 9: BENCH cabinet series dimensions

	BENCH standard S	BENCH pro S
H (mm)	770	770
H1 (mm)	554	515
HS (mm)*	20-135	20-135
B (mm)	601	601
BS (mm)	480	403
LB (mm)	503	503
T (mm)	594	594
T1 (mm)	1170	1170
TS (mm)	424	418
L (l)	12	10
V (l)	13.2	11.3
M (kg)	60	45
Empty weight (kg)	111	134
Max. payload (kg)	60	45

* Safety storage cabinets with no base are equipped with a 3 mm thick bottom plate.

3.3 Pressure drop during ventilation

Industrial ventilation of the safety storage cabinet results in a pressure drop at the exhaust air connection, as shown in the following diagram.


Fig. 6: Average pressure drop for all cabinet sizes

- 1 Average pressure drop from all cabinet sizes
- 2 Q with ten-fold air exchange (see table)

Model size	Q [m³/h]	Δp [Pa]	Model size	Q [m³/h]	Δp [Pa]
UTS standard and pro cabinet series					
S	1.1	<1	S-5	0.9	<1
M	1.7	<1	M-5	1.4	<1
L	2.2	<1	L-5	1.7	<1
LD	2.2	<1	LD-5	1.7	<1
XL	2.7	<1	XL-5	2.2	<1
BENCH standard and pro cabinet series					
S	1.4			<1	

4 Structure and function

4.1 Construction

- Cabinet carcass and doors in multi-layer construction
- External panelling: Powder-coated sheet steel
- Wall construction: Multi-layer design
- Interior surfaces: Light grey-coated decor panels
- Safety technology elements for closure of venting cut-off flaps in case of fire: Brass, spring steel (1.410)

4.2 Earthing options


Fig. 7: Earthing options

1 Equipotential bonding saddle on the cabinet carcass

Earthing the safety storage cabinet prevents ignition hazards.

The interior fittings are conductively connected to one another by an equipotential bonding saddle or equipotential bonding screw on the cabinet carcass.

Correct earthing is stipulated in the applicable national standards and regulations, e.g. 'TRGS 727' in Germany.

4.3 Exhaust air connection and feed opening

The safety storage cabinets can be connected to a technical exhaust air system which ducts outside at a danger-free location.

The exhaust air connection and the feed opening are located in the rear panel of the safety storage cabinet for this purpose.

In normal operation, industrial ventilation of safety storage cabinets prevents the occurrence of a potentially explosive atmosphere inside the cabinet. The layout of the ventilation ducts in the cabinet means that ventilation is effective directly above the bottom tray.


Fig. 8: Connection socket (view from rear)

- 1 Air supply opening
- 2 Exhaust air connection: NW Ø 75 mm
- 3 UTS distance: 210 mm ; BENCH: 385 mm

4.4 Doors

4.4.1 Door options

The UTS and BENCH safety storage cabinet series have the following door options:

Door operation with door handle

- The wing door can be opened by pulling on the door handle and remains open in any position.

Push-to-open door operation

- Gently pressing the marked corner of the door opens the door slightly and it can then be pulled fully open with the hand. It remains open in any position.

4.4.2 Locking cylinder

The door can be locked with the integrated locking cylinder. The key numbers are embossed on the locking cylinder and on the keys supplied, e.g. A007. Locks be subsequently adjusted to the operator's requirements.

4.5 Safety technology

4.5.1 Door closure in case of fire

At an ambient temperature of approx. 50°C, open doors are closed by the safety technology.

On safety storage cabinets with a pull-out tray, the pull-out tray is automatically retracted when closing the door.

4.5.2 Closure of venting cut-off flaps in case of fire


Fig. 9: Left: Ventilation ducts open; Right: Ventilation ducts closed

1 Venting cut-off flaps

At a room temperature of 70°C the ventilation openings for air supply and exhaust air are sealed by the venting cut-off flaps installed in the cabinet.

4.6 Interior fittings

4.6.1 Bottom tray


Fig. 10: Safety storage cabinet with bottom tray

- 1 Bottom tray
- 2 Perforated sheet insert

The function of the bottom tray in the floor area of the safety storage cabinet is to collect leaking substances in the cabinet interior. It is equipped with a perforated sheet insert as standard. The perforated sheet insert allows the bottom tray to be used as a standing surface in accordance with the regulations.

4.6.2 Pull-out tray


Fig. 11: Safety storage cabinet with pull-out tray

- 1 Pull-out tray
- 2 Perforated sheet insert
- 3 Note: Maximum container height

The function of the pull-out tray at the bottom of the safety storage cabinet is to collect leaking substances in the cabinet interior. It is equipped with a perforated sheet insert as standard. The perforated sheet insert allows the pull-out tray to be used as a standing surface in accordance with the regulations. Pull-out trays are permanently installed in the factory and subsequent remodelling is not possible.

4.7 Additional pull-out level (optional)

The UTS pro and BENCH safety storage cabinet series can be equipped with an additional pull-out level. This changes the technical specifications set out in [Chapter 3.2 'Dimensions and equipment'](#) on page 62 as follows:


Fig. 12: Additional pull-out level

- H1 Max. container height
- TS Standing surface depth
- M Standing surface load-bearing capacity, $M = M1 + M2$ (uniformly distributed over both levels)
- L Maximum volume of the largest individual container
- V Max. collection volume

Tab. 10: UTS pro cabinet series S-XL dimensions

	S	M	L	ID	XL
H1 (mm)	190	190	190	190	190
TS (mm)	410	410	410	410	410
L (l)	1.2	2	3	1.2 1.2**	2 1.2**
V (l)	1.5	2.5	3.3	1.5 1.5**	2.3 1.5**
M (kg)	45	45	45	45 45**	45 45**

Tab. 11: UTS pro cabinet series S-5- XL-5 dimensions

	S-5	M-5	L-5	ID-5	XL-5
H1 (mm)	190	190	190	190	190
TS (mm)	310	310	310	310	310
L (l)	1	1.6	2	1 1**	1.5 1**
V (l)	1.1	1.9	2.5	1.1 1.1**	1.7 1.1**
M (kg)	40	40	40	40 40**	40 40**

** Safety storage cabinet with two (separate) storage areas: left storage area | right storage area.

Tab. 12: BENCH pro S dimensions

	S
H1 (mm)	270
TS (mm)	410
L (l)	2.5
V (l)	2.9
M (kg)	35

4.8 ToeKick base (optional)

The safety storage cabinet can be fitted with an optional ToeKick base. The front of the base is set back by 50 mm. This enables an ergonomic standing position when working on the bench or in the fume hood.

4.9 Roller set (optional)

The safety storage cabinet can be fitted with an optional roller set, which is integrated into the base. This enables the safety storage cabinet to be easily moved in its location and also to be pulled out for occasional cleaning or maintenance work.

The rollers are not designed for continuous mobile use of the safety storage cabinet.

Installing the roller set increases the cabinet height by 8 mm.

! NOTICE!

Ensure appropriate floor quality

The rollers can leave permanent grooves on soft floors.

! NOTICE!

Moving the safety storage cabinet.

Before moving the safety storage cabinet, completely screw in or remove the locking screws ↪ *Chapter 6.3 'Installing the safety storage cabinet with roller set' on page 79.*

4.10 Storage shelves (optional)

Additional storage shelves for storage of smaller containers are available on request.

The storage shelves are height adjustable in a 16 mm pitch. The maximum container height is reduced accordingly.

! NOTICE!

Subsequent installation of storage shelves

Incorrect installation in the wrong location can cause damage to the safety storage cabinet.

- They should only be installed on the approved surfaces (see additional installation instructions for storage shelves).
- Use only the attachment screws supplied for installation.

4.11 Pipe penetration (optional)

Tested penetrations are optionally available for the safety storage cabinet with Type 90 classification. The penetrations are attached to the safety storage cabinets from the outside and can be provided with holes. Refer to the separate instructions for the pipe penetrations.

When used correctly, it has no negative effect whatsoever on fire resistance. Unused, open holes in the pipe penetrations must be sealed.

! NOTICE!

Later attachment of pipe penetrations

Incorrect attachment in the wrong location can cause damage to the safety storage cabinet.

- They should only be attached on the approved surfaces (see additional instructions for pipe penetrations).
- Holes only based on defined penetration profiles (see separate mounting instructions for each cabinet type).

! NOTICE!

Routing and using pipe penetrations

The routing and use of pipes, cables and hoses are the responsibility of the operator.

Incorrect handling can damage the safety technology in the safety storage cabinet and cause it to fail.

- Carry out a separate risk assessment of the overall setup.

5 Transport

The safety storage cabinet is packaged for transportation and is protected against damage by transport restraints. The transport restraints should be refitted before any transportation. For transportation, the safety storage cabinet can be placed on a transport pallet or other suitable transport equipment. It must be ensured that the safety storage cabinet is loaded as smoothly as possible.

WARNING!

Risk of crushing due to tipping or falling safety storage cabinet

If the safety storage cabinet tips over or falls when not transported with due caution, this can cause potentially fatal crushing.

- Wear personal protective equipment (PPE).
- Only transport with two people.
- Only transport the safety storage cabinet upright and unladen.
- Only drive under the safety storage cabinet using suitable transport equipment.

NOTICE!

Handling the transport restraints

Removal of the transport restraints before transportation results in damage to the safety storage cabinet.

- The transport restraints and transport pallet should only be removed at the installation location.

6 Installation and commissioning

Install the safety storage cabinet so that the annual maintenance activities can be carried out without restriction.

6.1 Requirements on the installation location

The safety storage cabinet is approved for installation in a building.


Fig. 1.3: Requirements on the installation location

Consider the following in relation to the installation location:

- The surface must be able to bear the weight of the safety storage cabinet when fully loaded.
- The surface must be horizontal in order to guarantee problem-free functioning of the safety storage cabinet.
- The load-bearing capacity and stability of the surface must be assured both in normal situations and in the event of a fire.
- Do not install the safety storage cabinet near sources of heat.
- Protect the safety storage cabinet against moisture.
 - At a relative humidity of >70 % use in closed and heated buildings is permissible for a few weeks each year.
- The operating temperature must be between -5°C and $+40^{\circ}\text{C}$.

6.2 Align the safety storage cabinet

The alignment procedure described below is used for precision alignment. Remedy any major floor unevenness of more than 1.5 mm on site.

The safety storage cabinet is supplied without alignment elements. An alignment may be necessary in individual cases.

Personal:

- Technical specialist employees

Werkzeug:

- Suitable tool

1. ➤ Raise the safety storage cabinet slightly.
2. ➤ Place steel or stainless steel spacers underneath the safety storage cabinet.
3. ➤ Check the alignment. ↪ *Chapter 6.4 'Check the alignment of the safety storage cabinet' on page 80*
4. ➤ Set the safety storage cabinet down carefully.

6.3 Installing the safety storage cabinet with roller set

Rollers enable easy installation under workbenches and laboratory fume hoods.


Fig. 14: Installation and operation of locking screws

- 1 Jaw spanner size 13 mm

Installation and operation of locking screws

Personal:

- Technical specialist employees

1. ➤ Raise the safety storage cabinet slightly.
2. ➤ Completely screw in the enclosed locking screws in front of the left and right roller.
3. ➤ Set the safety storage cabinet back down smoothly.
4. ➤ Push the safety storage cabinet to the required position.
5. ➤ Use a jaw spanner to turn the locking screws until they reach the floor and slightly pretension them.

The locking screws prevent the safety storage cabinet from slipping.

Locking screws do not serve as adjustment screws for aligning the safety storage cabinet.

6.4 Check the alignment of the safety storage cabinet

If the safety storage cabinet is not aligned properly, the open wing doors will automatically close themselves or open fully, → Chapter 6.2 'Align the safety storage cabinet' on page 78.


Fig. 15: Checking alignment

Correct alignment of the safety storage cabinet:

- when the doors are closed, the door gaps are of equal width.
- With two doors, the central gap and ceiling gap form an even "T".

6.5 Venting the safety storage cabinet

6.5.1 Connection to an exhaust air system

WARNING!

Insufficient safety storage cabinet air circulation

A lack of or insufficient air exchange can lead to formation of an explosive atmosphere in the safety storage cabinet.

This may result in death or serious injury.

- In a safety storage cabinet with ventilation system, a minimum hourly air exchange must take place that is at least 10 times the internal volume of the cabinet when the doors are closed.

Installation of industrial ventilation and connection to an existing exhaust air system must be carried out by a qualified company and is not a service provided by DÜPERTHAL.

Connection to the exhaust air system:

Personal:

- Technical specialist employees

1. ➔ Connect the exhaust air line to the exhaust air connection socket and secure with a collar.
2. ➔ After installing the safety storage cabinet, check the connection to the exhaust air system with smoke tubes.

The power of the exhaust air system can be determined using the technical specifications, ↗ Chapter 3.3 'Pressure drop during ventilation' on page 68.

6.5.2 Operating the safety storage cabinet without industrial ventilation

Safety storage cabinets for passive storage can be operated without industrial ventilation.

Safety storage cabinets that are operated without industrial ventilation must be indicated using an appropriate notice.

WARNING!

Safety storage cabinets without industrial ventilation

Risk of fire and explosion due to ignition of explosive mixtures in the safety storage cabinet.

This may lead to death or serious injury.

- Do not use any ignition sources in the safety storage cabinet.
- Earth safety storage cabinets using equipotential bonding.

If an explosive atmosphere is to be expected, measures must be taken in compliance with applicable national standards and regulations, e.g. 'TRGS 722' in Germany, and an explosion protection document must be drawn up.

WARNING!

Potentially explosive area on safety storage cabinet

This may result in death or serious injury.

- Determine the ex-zone in compliance with the applicable national standards and regulations, e.g. 'TRGS 722' in Germany, and mark it clearly and permanently.
- Naked flames and smoking in potentially explosive areas are prohibited.
- Do not use any tools that cause mechanically generated sparks.
- Avoid electrostatic charges.
- Avoid equipment with surface temperatures above the ignition temperatures of the flammable liquids stored.
- Only operate electrical equipment in potentially explosive areas if it complies with the requirements of the applicable national standards and regulations, e.g. 'TRGS 722' in Germany.

Clearly mark the Zone 2 potentially explosive area with the following warning sign in compliance with the applicable national standards and regulations, e.g. the European ATEX product directive '2014/34/EU':

The size of all symbols and notices should be appropriate for the size of the safety storage cabinet.

6.6 Earth the safety storage cabinet


Fig. 16: Earthing connection

Earthing connection

Personal:

- Technical specialist employees

1. ➤ Place the screwdriver under the equipotential bonding saddle (A).
2. ➤ Bend the equipotential bonding saddle upwards.
3. ➤ Connect the earthing cable (not included in scope of delivery).

7 Operation

7.1 Opening the safety storage cabinet with the door handle


Fig. 17: Door operation

Door operation with door handle

- ➔ Open the safety storage cabinet by pulling on the door handle of the wing door.
- ⇒ The wing door remains open in any position.

7.2 Opening a safety storage cabinet with push-to-open door operation


Fig. 18: Door operation with push-to-open mechanism

Push-to-open door operation

1. ➤ Gently press on the surface of the door in the area of the green dot (A) so that the door disengages from the detent mechanism.
2. ➤ When you release the surface of the door, the door is pushed open slightly.
3. ➤ Pull on the green handle at the top edge of the door (B) to fully open.
4. ➤ To close, gently press on the door to engage it with the detent mechanism.

WARNING!

Blocked doors

Doors that are held open by objects impair the function of the safety technology.

This may result in death or serious injuries as a result of inadequate fire protection.

- Close the doors after every work process.

7.3 Pull-out tray

The pull-out tray is attached directly to the door, and is pulled out when the door is opened. When the door is closed, the pull-out tray is pushed into the safety storage cabinet.

NOTICE!

Observe the maximum container height

Exceeding the maximum container height can lead to damage to the container or prevent automatic closing of the door in case of fire.

- Observe the maximum container height marking on the inside of the door ↪ *Chapter 4.6.2 'Pull-out tray' on page 73.*

7.4 Check the bottom tray or pull-out tray

! NOTICE!

Store hazardous substances so that a visual inspection of the bottom tray for extraneous substances is possible on working days.

Personal:

- Technical specialist employees
- ➔ Perform a daily visual inspection for extraneous substances.

The perforated sheet insert can be removed for easier cleaning of the bottom tray.

Cleaning the bottom tray

Personal:

- Technical specialist employees
- 1. ➔ Detach the plug connection on the earthing cable.
- 2. ➔ Remove the perforated sheet insert.
- 3. ➔ Clean the bottom tray.
- 4. ➔ Re-insert the perforated sheet insert.
- 5. ➔ Reconnect the bottom tray to the equipotential bonding.

8 Opening the safety storage cabinet after a fire

After a fire, the safety storage cabinet must not be opened for at least 24 hours, and only with the utmost caution and by specialist employees.

 **WARNING!**

Explosive vapour-air mixture

This may result in death or serious injury.

- Before opening the safety storage cabinet, remove all ignition sources within a 10 metre radius.
- Only open the safety storage cabinet with tools that do not cause any mechanically generated sparks.

 **WARNING!**

Damaged safety storage cabinet due to fire or extinguishing agents

This may result in death or serious injury.

- Do not use safety storage cabinets that have been damaged by fire or extinguishing agents.

9 Maintenance

Check the safety storage cabinet for any externally visible damage or defects.

Always perform checks:

- After installation.
- Before commissioning.
- After changes.
- After maintenance.

The safety storage cabinet should also be inspected periodically at the following intervals.

Interval	Maintenance work	Personnel
Daily	Bottom tray and storage levels <ul style="list-style-type: none"> ■ Check in accordance with regulations governing water legislation ■ Collect and properly dispose of leaked liquids immediately. 	Technical specialist employees

Interval	Maintenance work	Personnel
Monthly	Closing of the doors closing system <ul style="list-style-type: none"> ■ Open the door and inspect the closure. 	Technical specialist employees
	Ventilation <ul style="list-style-type: none"> ■ Check the effectiveness of the ventilation with a woollen thread or with a smoke tube in front of the exhaust air opening. ■ Remove contamination at the exhaust air opening. 	Technical specialist employees
	Seals <ul style="list-style-type: none"> ■ Check the sealing strips are seated properly in the carcass frame and on the end faces of the doors. ■ If visible damage is found, replace the sealing strips immediately. 	Technical specialist employees
	Labelling <ul style="list-style-type: none"> ■ Inspect the safety labels on the safety storage cabinet to ensure they are complete. 	Technical specialist employees
	Alignment (push-to-open technology) <ul style="list-style-type: none"> ■ Check ease of movement and engagement/disengagement of the doors. 	Technical specialist employees

Interval	Maintenance work	Personnel
Annually	Safety storage cabinet <ul style="list-style-type: none"> ■ Check of the entire safety storage cabinet 	DÜPERTHAL service technicians

If faults occur, assist the technical customer service by providing the cabinet model, production and key number, along with a description of the fault.

To comply with national regulations, e.g. BetrSichV in Germany, all safety equipment must be inspected annually by a qualified person as defined in the valid national standards and regulations, e.g. the technical rules for industrial safety: TRBS 1203 in Germany.

10 Faults

Fault description	Cause	Remedy	Personnel
Doors do not close.	Safety storage cabinet is not aligned correctly.	Install the safety storage cabinet so it is horizontal. ↪ <i>Chapter 6.4 'Check the alignment of the safety storage cabinet' on page 80</i>	Technical specialist employees
	Doors are held open by objects.	Do not wedge or hold doors open with any objects.	Technical specialist employees
	Safety storage cabinet is not correctly filled.	Make sure that containers in the safety storage cabinet are uniformly distributed.	Technical specialist employees
No extractor present.	Venting cut-off flaps closed, as closing mechanism has been triggered.	Replace the locking mechanism.	DÜPERTHAL service technicians
Doors do not move easily.	Moving parts, such as hinges, are dirty or corroded.	<ul style="list-style-type: none"> ■ Remove rust. ■ Lubricate parts. ■ Remove corrosive substances from the safety storage cabinet. ■ Notify technical customer service. 	Technical specialist employees
Doors open again after being closed.	Safety storage cabinet is not aligned correctly.	<ul style="list-style-type: none"> ■ Align the safety storage cabinet so it is horizontal. ↪ <i>Chapter 6.2 'Align the safety storage cabinet' on page 78</i> 	Technical specialist employees
Doors close again after being opened.	Safety storage cabinet is not aligned correctly.	<ul style="list-style-type: none"> ■ Align the safety storage cabinet so it is horizontal. ↪ <i>Chapter 6.2 'Align the safety storage cabinet' on page 78</i> 	Technical specialist employees

11 Spare parts and accessories

- Pull-out tray
- Bottom tray
- PP insert
- Anti-slip mat made from rubber
- Door handle
- Perforated sheet insert
- Base
- Ventilators
- Exhaust air monitoring units
- Roller set
- Additional storage shelves
- Pipe penetrations

12 Disposal

 CAUTION!

Dismantling the safety storage cabinet

Risk of injury due to improper dismantling of the safety storage cabinet.

- Ensure that the safety storage cabinet is only dismantled by specialist technical employees.

The safety storage cabinet can be completely dismantled by specialist technical employees.

Recycle the individual material components separately.

Comply with national and local disposal regulations.

To save resources, do not place parts of the safety storage cabinet or the whole cabinet in bulky or domestic waste.

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hereby declares that the following products

Product designation: Safety storage cabinet type 90
 Model designation: UTS standard, UTS pro
 Model size: S-5, M-5, L-5, XL-5, XXL-5, S, M, L, XL, XXL

comply with all relevant provisions of the directives applied. This declaration refers only to the product in the condition as supplied to the customer. Parts and/or modifications subsequently fitted by the end user are not considered.

The following directives were applied:

Machinery Directive 2006/42/EG

The following harmonised standards were applied:

DIN EN ISO 1210:2011 Safety of machinery - General principles for design - Risk assessment and risk reduction

The following national standards and technical specifications have been applied:

DIN EN 14470-1:2004 Fire safety storage cabinets - Part 1: Safety storage cabinets for flammable liquids

DIN EN 16121:2017 Non-domestic storage furniture - Requirements for safety, strength, durability and stability

DIN EN 16122:2012 Domestic and non-domestic storage furniture - Test methods for the determination of strength, durability and stability

The person authorised to compile the technical documentation is: Mr Frank Backhaus, address see above.

Karlstein, 13.12.2022
 (Location, date)

 Signature
 Franz-Josef Hagen / Managing Director

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hereby declares that the following products

Product designation: Safety storage cabinet type 90
Model designation: BENCH standard, BENCH pro
Model size: S

comply with all relevant provisions of the directives applied. This declaration refers only to the product in the condition as supplied to the customer. Parts and/or modifications subsequently fitted by the end user are not considered.

The following directives were applied:

Machinery Directive 2006/42/EG

The following harmonised standards were applied:

DIN EN ISO 1210:2011 Safety of machinery - General principles for design - Risk assessment and risk reduction

The following national standards and technical specifications have been applied:

DIN EN 14470-1:2004 Fire safety storage cabinets - Part 1: Safety storage cabinets for flammable liquids

DIN EN 16121:2017 Non-domestic storage furniture - Requirements for safety, strength, durability and stability

DIN EN 16122:2012 Domestic and non-domestic storage furniture - Test methods for the determination of strength, durability and stability

The person authorised to compile the technical documentation is: Mr Frank Backhaus, address see above.

Karlstein, 13.12.2022
(Location, date)

Signature
Franz-Josef Hagen / Managing Director

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DÜPERTHAL Sicherheitstechnik GmbH & Co. KG