ACID C pro, ACID B pro, ACID U pro  
Sicherheitsschränke

**ACID**

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ACID C pro, ACID B pro, ACID U pro  
Safety storage cabinets

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© 05/2021 DÜPERTHAL 55251, D00110630, 0, en\_GB

## Table of contents


<b>1</b>	<b>General information</b> .....	<b>40</b>
1.1	Notes for reading.....	40
1.2	Type plate.....	40
<b>2</b>	<b>Safety</b> .....	<b>41</b>
2.1	Function of safety notices.....	41
2.2	Correct use.....	41
2.3	Misuse.....	41
2.4	Obligations of the operator.....	42
2.5	Demands on employees.....	43
2.6	Stored goods.....	43
2.7	Hazardous areas and their labelling.....	44
<b>3</b>	<b>Technical specifications</b> .....	<b>46</b>
3.1	General data.....	46
3.2	Dimensions and equipment.....	47
3.3	Pressure drop during ventilation.....	51
<b>4</b>	<b>Structure and function</b> .....	<b>52</b>
4.1	Construction.....	52
4.2	Exhaust air connection and feed opening.....	52
4.3	Doors.....	52
4.4	ACID C pro, ACID B pro and ACID U pro interior fittings.....	53
<b>5</b>	<b>Transport</b> .....	<b>54</b>
<b>6</b>	<b>Installation and commissioning</b> .....	<b>55</b>
6.1	Requirements on the installation location.....	55
6.2	Installing the safety storage cabinet.....	56
6.3	Wall mounting.....	58
6.4	Connection to an exhaust air system.....	59
<b>7</b>	<b>Operation</b> .....	<b>61</b>
7.1	Open the safety storage cabinet.....	61
7.2	Pull-out tray.....	61
7.3	Checking the pull-out trays.....	61
<b>8</b>	<b>Maintenance</b> .....	<b>62</b>
<b>9</b>	<b>Faults</b> .....	<b>64</b>
<b>10</b>	<b>Spare parts and accessories</b> .....	<b>65</b>
<b>11</b>	<b>Disposal</b> .....	<b>66</b>

# 1 General information

## 1.1 Notes for reading

The following symbols designate specific types of information.

Tab. 1: Explanation of symbol

Symbol	Type of information
	Information for easier and more effective working
	Procedural step
	Result of a procedural step
	Link to another part of the document

## 1.2 Type plate

The type plate is attached to the outside of the safety storage cabinet door.

Model

Maximum volume of individual containers

Serial number and year of manufacture

Maximum load per storage shelf

Fig. 1: Type plate

## 2 Safety

### 2.1 Function of safety notices

Safety notices warn against physical or material damage and provide information on how such damage can be avoided.

The following signal words identify the degree of danger and the extent of the risk.

#### WARNING!

The signal word refers to a potential threat which can cause death or serious injury.

#### NOTICE!

The signal word 'NOTE' indicates a situation that can lead to damage to the safety storage cabinet.

### 2.2 Correct use

*Observe the safety instructions in these operating instructions to reduce health risks and avoid dangerous situations.*

*Any improper use as defined in these operating instructions involves a risk of accidents and damage.*

The safety storage cabinet should be used for storage of solid and liquid, non-flammable corrosive chemicals in working spaces. It is particularly suitable for storage of acid and lye.

### 2.3 Misuse

Any use that goes beyond the specified correct use is considered to be misuse.

DÜPERTHAL accepts no liability for damage arising from misuse.

In addition, the following safety notices must be observed:

#### WARNING!

##### **Storage of living organisms in the safety storage cabinet**

Living organisms can come into contact with the stored hazardous substances.

This can lead to death or serious injury.

- Use the safety storage cabinet exclusively for storage of chemicals.

#### WARNING!

##### **Storage of food in the safety storage cabinet**

Food can come into contact with the stored hazardous substances.

This can lead to death or serious injury.

- Use the safety storage cabinet exclusively for storage of chemicals.

 **WARNING!**

**Storage of flammable liquids**

Flammable liquids may explode in case of fire.

This can lead to death or serious injury.

- Store flammable liquids in a Type 90 or Type 30 certified safety storage cabinet.

 **WARNING!**

**Storage of hazardous substances together**

Risk of uncontrolled chemical reactions.

This may result in death or serious injury.

- Only store substances and preparations that are permitted to be stored together in the safety storage cabinet.

 **WARNING!**

**Spilled liquids**

Inhalation of hazardous vapours from spilled liquids can lead to life-threatening injuries to the respiratory system.

Spilled liquids can lead to life-threatening burns.

- Collect and properly dispose of spilled liquids immediately in accordance with accident prevention regulations.

 **NOTICE!**

**Alteration and modifications**

Do not alter or modify the safety storage cabinet.

- If alteration or modification of the safety storage cabinet is required, contact DÜPERTHAL.

## 2.4 Obligations of the operator

The operator is obliged to comply with applicable legal regulations. This includes:

- Issuing operating instructions or laboratory regulations and defining the scope of application for the safety storage cabinet.
- Carrying out risk assessments.
- Defining and documenting activities of designated employees.

## 2.5 Demands on employees

### WARNING!

#### Employees who do not meet these requirements

This may result in death or serious injury.

- Designate employees who meet the requirements to carry out activities.

These operating instructions set out the following employee activities:

- Specialist technical employees
- DÜPERTHAL service technicians

Only people who have been trained by the operator in use of the safety storage cabinet and handling of the stored goods are approved as specialist technical employees.

#### DÜPERTHAL service technicians

DÜPERTHAL employees are specifically trained by DÜPERTHAL to carry out their activities.

## 2.6 Stored goods

### NOTICE!

#### Material damage due to escaping liquids or vapours

Corrosive liquids or vapours can escape from containers that are not tightly sealed or are unsuitable and cause damage to the safety storage cabinet.

- Only store the stored goods in suitable, tightly sealed containers.

Storage, handling and use of the stored goods must comply with the applicable national standards and regulations, e.g. 'TRGS 510' or 'TRGS 526' in Germany.

Hazardous substances must be labelled at least with the name of the substance, the components of the preparation as well as the hazard symbols with the associated indications of danger. Furthermore, they must be marked with the safety measures to be observed (H and P statements) or the corresponding safety data sheet must be provided.

Hazardous substances must be stored in such a way that no dangerous reactions are possible if the containers are damaged. This must be observed as part of the risk assessment. Detailed information on the combined storage of different classes of substances can be found in the German 'TRGS 510', for example.

## 2.7 Hazardous areas and their labelling

The following must be attached to the front of the safety storage cabinet and must be clearly visible:

- Name or trademark of the manufacturer
- Serial number and year of manufacture
- Information on the largest individual container volume that can be stored
- Information on the shelves' maximum load capacity

Furthermore, the following signs must be attached to the front of the safety storage cabinet and must be clearly visible:

Tab. 2: Prohibited action signs

Symbol	Meaning	Standard
	P011: Do not extinguish with water	DIN EN ISO 7010

Tab. 3: Warning signs

Symbol	Meaning	Standard
	W023: Warning for corrosive substances	DIN EN ISO 7010

Tab. 4: Mandatory action signs

Symbol	Meaning	Standard
	M002: Read the instructions	DIN EN ISO 7010
	M004: Wear eye protection	DIN EN ISO 7010
	M009: Wear hand protection	DIN EN ISO 7010

### 3 Technical specifications

#### 3.1 General data


Fig. 2: General diagram of ACID C pro safety storage cabinet

- 1 Exhaust air connection (DN 75)
- 2 Type plate
- 3 Standing surface
- 4 Plinth panel with adjustable feet


Fig. 3: General diagram of safety storage cabinet ACID B pro and ACID U pro

- 1 Exhaust air connection (DN 75)
- 2 Type plate
- 3 Standing surface

## 3.2 Dimensions and equipment

### 3.2.1 ACID C pro


Fig. 4: ACID C pro dimensions

- H Cabinet height
- T Cabinet depth
- B Cabinet width
- T1 Cabinet depth with open doors
- V Maximum collection volume
- L Maximum volume of largest individual container
- LB Clear width
- TS Standing surface depth
- BS Standing surface width
- H1-H6 Adjustable heights 1 to 6
- M Pull-out tray load-bearing capacity

Tab. 5: ACID C pro S – XL dimensions

ACID C pro				
Model	S	XS	M	XL
H (mm)	1385	1385	2045	2045
B (mm)	600	1200	600	1200
T (mm)	609	609	609	609
T1 (mm)	1207	1209	1207	1209
LB (mm)	553	2 x 563	553	2 x 563
TS (mm)	404	404	404	404
BS (mm)	474	474	474	474
L (l)	17	17	17	17
V (l)	19	19	19	19
M (kg)	40	40	40	40
Empty weight (kg)	75	145	110-115 <sup>1</sup>	205-215 <sup>1</sup>
Max. payload (kg)	240	2 x 240	240	2 x 240
Total permitted weight (kg)	315	625	350-355 <sup>1</sup>	685-695 <sup>1</sup>

<sup>1</sup> Empty weight varies depending on number of pull-out trays

Tab. 6: ACID C pro S - XL maximum adjustable heights

ACID C pro						
Model	S	XS	M	XL	M	XL
H1 (mm)	305	305	460	460	310	310
H2 (mm)	305	305	459	459	310	310
H3 (mm)	305	305	460	460	294	294
H4 (mm)	300	300	459	459	310	310
H5 (mm)	–	–	–	–	310	310
H6 (mm)	–	–	–	–	294	294

## 3.2.2 ACID B pro and ACID U pro


Fig. 5: ACID B pro and ACID U pro dimensions

- H Cabinet height
- T Cabinet depth
- B Cabinet width
- T1 Cabinet depth with open doors
- V Maximum collection volume
- L Maximum volume of largest individual container
- LB Clear width
- TS Standing surface depth
- BS Standing surface width
- H1-H2 Adjustable heights 1 to 2
- BH Base height

Tab. 7: ACID B pro S - L dimensions

ACID B pro					
Model	S	L		S	L
H <sup>1</sup> (mm)	770	770	H1 (mm)	350	350
HS (mm)	20-135	20-135	H2 (mm)	350	350
B (mm)	600	1100	L (l)	17	17
T (mm)	609	609	V (l)	19	19
T1 (mm)	1207	1159	M (kg)	40	40
LB (mm)	553	2 x 513	Empty weight (kg)	40	75
TS (mm)	404	474	Max. payload (kg)	80	2 x 80
BS (mm)	474	404	Total permitted weight (kg)	120	235

Tab. 8: ACID U pro S - L dimensions

ACID U pro					
Model	S	L		S	L
H <sup>1</sup> (mm)	569	569	H1 (mm)	264	264
HS (mm)	20-135	20-135	H2 (mm)	264	264
B (mm)	600	1100	L (l)	17	17
T (mm)	609	609	V (l)	19	19
T1 (mm)	1207	159	M (kg)	40	40
LB (mm)	553	2 x 513	Empty weight (kg)	35	65
TS (mm)	404	474	Max. payload (kg)	80	2 x 80
BS (mm)	474	404	Total permitted weight (kg)	115	225

<sup>1</sup> Height without base

### 3.3 Pressure drop during ventilation

Technical ventilation of the safety storage cabinet results in a pressure drop at the exhaust air connection, as shown in the following diagram.


Fig. 6: Pressure drop at exhaust air connection socket, averaged over all cabinet sizes

◆ Pressure drop at 30 m<sup>3</sup>/h

Tab. 9: Pressure drop  $\Delta p$  at the exhaust air connection socket and in the cabinet interior

ACID C pro			
Model	Q [m <sup>3</sup> /h]	$\Delta p$ at the exhaust air connection socket [Pa]	$\Delta p$ in the interior [Pa]
S	30	10	<1
XS	2 x 30	2 x 10	<1
M	30	14	<1
XL	2 x 30	2 x 14	<1
ACID B pro			
Model	Q [m <sup>3</sup> /h]	$\Delta p$ at the exhaust air connection socket [Pa]	$\Delta p$ in the interior [Pa]
S	30	10	<1
L	30	2 x 10	<1
ACID U pro			
Model	Q [m <sup>3</sup> /h]	$\Delta p$ at the exhaust air connection socket [Pa]	$\Delta p$ in the interior [Pa]
S	30	9	<1
L	30	2 x 9	<1

## 4 Structure and function

### 4.1 Construction

- Cabinet carcass and doors made of decor panels with melamine resin coating
- Plastic (HDPE) pull-outs with tray insert

### 4.2 Exhaust air connection and feed opening

The safety storage cabinets can be connected to an industrial exhaust air system which discharges air outside in a safe location.

There are one or two exhaust air connection sockets with a diameter of DN 75 for this purpose on the cabinet roof of the ACID C pro series and on the rear of the safety storage cabinet on the ACID B pro and ACID U pro series.

#### ! NOTICE!

##### Material damage due to aggressive and corrosive media

- If aggressive and corrosive media are stored in ACID series cabinets, a connection to a continuous ventilation system is mandatory.
- The volumetric flow rate must be at least 30 m<sup>3</sup>/h for each side panel.

### 4.3 Doors

The cabinets in the ACID U pro and ACID B pro cabinet series are equipped with one or two wing doors. The safety storage cabinets in the ACID C pro series are equipped with one to four wing doors.


The doors are lockable with the integrated locking cylinders. The key numbers are embossed on the locking cylinder and on the keys supplied.

#### ! NOTICE!

- Keep the doors of the safety storage cabinet closed and re-close them manually after each work process.

### 4.4 ACID C pro, ACID B pro and ACID U pro interior fittings

#### 4.4.1 Pull-out trays

*Fig. 7: Pull-out with plastic insert*

The safety storage cabinets in the ACID C pro, ACID B pro and ACID U pro series contain pull-outs with plastic inserts made of HDPE.

#### ! NOTICE!

The maximum height of the top shelf must not exceed 1750 mm.

#### 4.4.2 Bottom tray (optional)

The safety storage cabinet can be fitted with plastic (PP) bottom trays as an option. The plastic grid insert enables these to be used as an additional standing surface.

#### 4.4.3 Side storage inserts (optional)

As an alternative to the pull-out trays, plastic (PP) shelves can be attached to the side panels in addition to the bottom tray.

## 5 Transport

The safety storage cabinet is packaged for transport.

### WARNING!

#### Risk of crushing due to tipping or falling safety storage cabinet

If the safety storage cabinet tips over when not transported with due caution, this can cause potentially fatal crushing.


- Wear personal protective equipment (PPE).
- Transport with at least two people.
- Only transport the safety storage cabinet upright and unladen.

### NOTICE!

#### Tipping the safety storage cabinet over during transport

Damage to the safety storage cabinet caused by incorrect handling.

- Only pick up the safety storage cabinet using special and suitable transport or lifting equipment.
- Preferably transport the cabinet upright and tipp it over only in exceptional cases.
- Do not prop the the cabinet on the adjustable feet while tipping it over.


## 6 Installation and commissioning

### 6.1 Requirements on the installation location

*Install the safety storage cabinet so that the annual inspection and maintenance activities can be carried out without restriction.*

The safety storage cabinet is approved for installation in a building.

*Fig. 8: Requirements for the installation location*

Consider the following in relation to the installation location:

- The surface must be able to bear the weight of the safety storage cabinet when fully loaded.
- The surface must be horizontal in order to guarantee problem-free functioning of the safety storage cabinet.
- Protect the safety storage cabinet against moisture.
  - At a relative humidity of  $>70\%$  use in closed and heated buildings is permissible for a few weeks each year.
- The operating temperature must be between  $-5^{\circ}\text{C}$  and  $+40^{\circ}\text{C}$ .

## 6.2 Installing the safety storage cabinet

### 6.2.1 ACID C pro

Safety storage cabinets delivered to the place of use are normally supplied without packaging.


Fig. 9: ACID C pro installation

- 1 Accessories (Plinth panel, adjustable feet and exhaust air connection socket)
- 2 Safety storage cabinet
- 3 Transport pallet
- 4 Transport equipment
- 5 Adjustable foot

#### Personal:

- Technical specialist employees

1. ➤ Remove the packaging and detach the accessories (1) from the cabinet (A).
2. ➤ Lift the safety storage cabinet smoothly from the transport pallet (3) and put it upright onto a suitable transport equipment (4), for example a pallet truck (B).
3. ➤ Screw the standing feet completely into the prepared holes from below (C).
4. ➤ Position the safety storage cabinet at the installation location.
5. ➤ Align the safety storage cabinet (see Fig. 10).

## CAUTION!

The safety storage cabinet must be completely emptied for alignment. Any containers still inside can tip over, spill or break. This can lead to serious burns to operating personnel.

- Before alignment, check that the safety storage cabinet is empty.

The safety storage cabinet can be aligned by screwing the adjustable feet in or out. To do this, the safety storage cabinet can be raised slightly.


Fig. 10: Alignment of an ACID C pro.

1 4 x per cabinet (ACID C pro S and M) or 4x per half cabinet (ACID C pro XS and XL)

### Personal:

- Technical specialist employees

1. ➤ Align the safety storage cabinet by screwing the adjustable feet in or out.

- The adjustable feet can be screwed in or out by hand from the outside.
- The adjustable feet can be screwed in or out from the inside using a hexagon key (A).

2. ➤ Check the alignment of the safety storage cabinet with a spirit level.

3. ➤ Seal the holes on the base of the cabinet using protective caps (A).

4. ➤ Attach the safety storage cabinet to the wall (☞ Chapter 6.3 'Wall mounting' on page 58).

5. ➤ Put the three-part plinth together and mount it to the adjustable feet (B).

## 6.2.2 ACID B pro and ACID U pro

Personal:

- Technical specialist employees

1. ➤ Remove the packaging and detach the accessories.
2. ➤ Lift the safety storage cabinet smoothly from the transport pallet.
3. ➤ Position the safety storage cabinet at the installation location.
4. ➤ Check the alignment of the safety storage cabinet using suitable equipment, e.g. a spirit level.
5. ➤ To align the safety storage cabinet, raise it slightly and slide suitable spacer plates under the base.
6. ➤ If required, attach the safety storage cabinet to the wall ( ↪ *Chapter 6.3 'Wall mounting' on page 58*).

## 6.3 Wall mounting

To guarantee stability, ACID C pro safety storage cabinets must be attached to a wall using the mounting brackets supplied. The other safety storage cabinets in the series can be attached to a wall if required.

### CAUTION!

Unmounted safety storage cabinets, particularly the ACID C pro M and XL models, can tip over. This can lead to serious injuries or death due to crushing and burns due to escaping stored materials.

- Before use, ensure that the safety storage cabinet is mounted on a stable wall by technically qualified personnel.


Fig. 11: Wall mounting

- 1 Mounting bracket
- 2 Mounting screws

## Wall mounting

### Personal:

- Technical specialist employees

1. ➤ Position the mounting bracket on the cabinet roof.
2. ➤ Attach the mounting bracket to the cabinet roof using the mounting screws supplied.
3. ➤ Attach the second arm of the mounting bracket to the wall with appropriate mounting materials provided on site.

*On the ACID U pro and ACID B pro safety storage cabinet lines, the mounting brackets can also be attached to the side panel.*

### ! NOTICE!

The condition of the wall must always be checked before installation. Only use mounting materials that are appropriate for the relevant wall condition and for the load.

### CAUTION!

There must not be any electrical installation zones or supply and waste disposal lines in the mounting area.

- Use a suitable detector to check the mounting area before drilling. Note that plastic pipes cannot be detected.

## 6.4 Connection to an exhaust air system

### ! NOTICE!

To prevent damage caused by corrosive atmospheres inside the cabinet, the safety storage cabinet must be vented with a volumetric flow rate of min. 30 m<sup>3</sup>/h per side panel.

If large storage quantities, very corrosive substances or substances in sub-optimal containers are stored in the safety storage cabinet, a higher exhaust air flow than that recommended should be used.

If metal parts installed in the cabinet, e.g. the hinges, exhibit any signs of corrosion, the set air exchange rate is too low.

*Installation of industrial ventilation and connection to an existing exhaust air system must be carried out by a qualified company and is not a service provided by DÜPERTHAL.*

The power of the exhaust air system can be determined using the technical specifications, ↪ Chapter 3.3 'Pressure drop during ventilation' on page 51.

DÜPERTHAL stipulates that acid and lye cabinets must be connected to a continuous industrial ventilation system. If the cabinets are operated without continuous venting, DÜPERTHAL will not accept any complaints relating to the exhaust air.


Fig. 12: Cabinet roof with exhaust air connection socket (ACID C pro)

1 Exhaust air connection socket

## Connection to the exhaust air system

### Personal:

■ Technical specialist employees

1. ➤ Secure the exhaust air connection socket with reducing fitting in the specified hole in the cabinet roof (ACID C pro models only).
2. ➤ Connect the waste air line to the exhaust air connection socket on the cabinet roof (ACID C pro) or on the rear of the cabinet (ACID B pro and ACID U pro) and if necessary secure with a collar.
3. ➤ After installing the safety storage cabinet, check the connection to the exhaust air system with smoke tubes.

## 7 Operation

### 7.1 Open the safety storage cabinet

Unlock the locking cylinder with the key and open the safety storage cabinet by pulling the door handle on the wing door.

The wing door remains open in any position.

### 7.2 Pull-out tray

Personal:

- Technical specialist employees

1. ➤ Remove the pull-out from the safety storage cabinet using the edge of the plastic insert.
2. ➤ When work on the safety storage cabinet is complete, slide the pull-out back in.

### 7.3 Checking the pull-out trays

Personal:

- Technical specialist employees

#### ! NOTICE!

Store hazardous substances so that a visual inspection of the pull-out tray for extraneous substances is possible every working day.

1. ➤ Perform a daily visual inspection for extraneous substances.
2. ➤ Collect and properly dispose of spilled hazardous substances immediately in accordance with accident prevention regulations.

## 8 Maintenance

Check the safety storage cabinet for any externally visible damage or defects.

Always perform checks:

- After installation.
- Before commissioning.
- After changes.
- After maintenance.

The safety storage cabinet should also be inspected periodically at the following intervals.

Interval	Maintenance work	Personnel
Daily	Pull-out trays: <ul style="list-style-type: none"> <li>■ Check in accordance with regulations governing water legislation</li> <li>■ Collect and properly dispose of leaked liquids immediately.</li> </ul>	Technical specialist employees

Interval	Maintenance work	Personnel
Monthly	Closing of the doors closing system <ul style="list-style-type: none"> <li>■ Open the door and inspect the closure.</li> </ul>	Technical specialist employees
	Ventilation <ul style="list-style-type: none"> <li>■ Check the effectiveness of the ventilation with a woollen thread or with a smoke tube in the cabinet in front of the exhaust air duct at the ventilation grilles.</li> </ul>	Technical specialist employees
	Labelling <ul style="list-style-type: none"> <li>■ Inspect the safety labels on the safety storage cabinet to ensure they are complete.</li> </ul>	Technical specialist employees

Interval	Maintenance work	Personnel
Annually	Safety storage cabinet <ul style="list-style-type: none"> <li>■ Check of the entire safety storage cabinet</li> </ul>	DÜPERTHAL service technicians

*If faults occur, assist the technical customer service by providing the cabinet model, production and key number, along with a description of the fault.*

*Technical storage facilities must be tested annually by a qualified person (as defined in TRBS 1203 for example) in compliance with BetrSichV and the maintenance interval stipulated by the manufacturer. The applicable national regulations must be observed.*

## 9 Faults

Fault description	Cause	Remedy	Personnel
Doors do not close or open automatically.	Safety storage cabinet is not aligned correctly.	Install the safety storage cabinet in a horizontal position ↪ <i>Chapter 6.2 'Installing the safety storage cabinet' on page 56 and following.</i>	Technical specialist employees
Doors and pull-outs are sluggish.	Moving parts, such as hinges or pull-out rails, are dirty or corroded.	<ul style="list-style-type: none"> <li>■ Clean the pull-out rails, hinges and closing mechanism.</li> <li>■ Notify technical customer service.</li> </ul>	Technical specialist employees
Corrosion to metal parts.	Insufficient exhaust air flow. /Container not tightly sealed	<ul style="list-style-type: none"> <li>■ Increase exhaust air flow ↪ <i>Chapter 6.4 'Connection to an exhaust air system' on page 59</i></li> <li>■ Replace containers that are not tightly sealed.</li> </ul>	Technical specialist employees

## 10 Spare parts and accessories

*Only original parts from DÜPERTHAL are to be used for the safety storage cabinets.*

- PE tray inserts
- Pull-out rails
- Exhaust air connection socket
- Mounting bracket
- Ventilators
- Exhaust air monitoring units
- Filter systems
- Plinth panel
- Door handle
- Hinges
- Door lock
- Replacement key

## 11 Disposal

 CAUTION!

**Dismantling the safety storage cabinet**

Risk of injury due to improper dismantling of the safety storage cabinet.

- Ensure that the safety storage cabinet is only dismantled by specialist technical employees.

The safety storage cabinet can be completely dismantled by specialist technical employees.

Comply with national and local disposal regulations.

To save resources, do not place parts of the safety storage cabinet or the whole cabinet in bulky or domestic waste.

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