# ACID C pro, ACID B pro, ACID U pro

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© 09/2020 DÜPERTHAL 55251, 0, EN

## Table of contents

1	General information	31
1.1	Notes for reading	31
1.2	Type plate	31
2	Safety	32
2.1	Function of safety notices	32
2.2	Correct use	32
2.3	Misuse	32
2.4	Obligations of the operator	33
2.5	Requirements for employees	34
2.6	Stored goods	34
2.7	Hazardous areas and their labelling	35
3	Technical specifications	37
3.1	General data	37
3.2	Dimensions and equipment	38
3.3	Pressure drop during ventilation	42
4	Structure and function	43
4.1	Construction	43
4.2	Exhaust air connection and air supply opening	43
4.3	Doors	43
4.4	ACID C pro, ACID B pro and ACID U pro interior fittings	44
5	Transport	45
6	Installation and commissioning	45
6.1	Requirements for the installation location	45
6.2	Installing the safety storage cabinet	46
6.3	Wall mounting	48
6.4	Connection to an exhaust air system	50
7	Operation	51
7.1	Opening the safety storage cabinet	51
7.2	Pull-out tray	51
7.3	Checking the pull-out trays	51
8	Maintenance	52
9	Faults	53
10	Spare parts and accessories	53
11	Disposal	54

# 1 General information

## 1.1 Notes for reading

The following symbols designate specific types of information.

Table 10: Explanation of symbol

Symbol	Type of information
<b>i</b>	Information for easier and more effective working
➔	Procedural step
✓	Result of a procedural step
🔗	Link to another part of the document

## 1.2 Type plate

The type plate is attached to the outside of the safety storage cabinet door.

Fig. 13: Type plate

- 1 Model
- 2 Serial number and year of manufacture
- 3 Maximum volume of individual containers
- 4 Maximum load per pull-out tray

## 2 Safety

### 2.1 Function of safety notices

Safety notices warn against physical or material damage and provide information on how such damage can be avoided.

The following signal words identify the degree of danger and the extent of the risk.

#### **WARNING!**

The signal word "WARNING" refers to a potential hazard which could result in death or serious injury.

#### **NOTE!**

*The signal word "NOTE" indicates a situation that can lead to damage to the safety storage cabinet.*

### 2.2 Correct use

*Observe the safety instructions in these operating instructions to reduce health risks and avoid dangerous situations.*

*Any improper use as defined in these operating instructions involves a risk of accidents and damage.*

The safety storage cabinet should be used for storage of solid and liquid, non-flammable corrosive chemicals in working spaces. It is particularly suitable for storage of acid and lye.

### 2.3 Misuse

Any use that goes beyond the specified correct use is considered to be misuse.

DÜPERTHAL accepts no liability for damage arising from misuse.

In addition, the following safety notices must be observed:

#### **WARNING!**

##### **Storage of living organisms in the safety storage cabinet**

Living organisms can come into contact with the stored hazardous substances. This can lead to death or serious injury.

- Use the safety storage cabinet exclusively for storage of chemicals.

#### **WARNING!**

##### **Storage of food in the safety storage cabinet**

Food can come into contact with the stored hazardous substances. This can lead to death or serious injury.

- Use the safety storage cabinet exclusively for storage of chemicals.

**WARNING!****Storage of flammable liquids in the safety storage cabinet**

Flammable liquids may explode in case of fire.

This can lead to death or serious injury.

- Store flammable liquids in a Type 90 safety storage cabinet.

**WARNING!****Storage of hazardous substances together**

Risk of uncontrolled chemical reactions.

This can lead to death or serious injury.

- Only store substances and preparations that are permitted to be stored together in the safety storage cabinet (e.g. see TRGS 510).

**WARNING!****Spilled liquids**

Inhalation of hazardous vapours from spilled liquids can lead to life-threatening injuries to the respiratory system.

Spilled liquids can lead to life-threatening burns.

- Collect and properly dispose of spilled liquids immediately in accordance with accident prevention regulations.

**NOTE!****Alteration and modifications**

*Do not alter or modify the safety storage cabinet.*

- *If alteration or modification of the safety storage cabinet is required, contact DÜPERTHAL.*

## 2.4 Obligations of the operator

The operator is obliged to comply with applicable legal regulations. This includes:

- Carrying out risk assessments.
- Issuing operating instructions or laboratory regulations and defining the scope of application for the safety storage cabinet.
- Defining and documenting activities of designated employees.

## 2.5 Requirements for employees

### WARNING!

#### Employees who do not meet these requirements

This can lead to death or serious injury.

- Designate employees who meet the requirements to carry out activities.

These operating instructions set out the following employee activities:

Technically qualified personnel:

Only people who have been trained by the operator in use of the safety storage cabinet and handling of the stored goods are approved as specialist technical employees.

DÜPERTHAL service technicians:

DÜPERTHAL employees are specifically trained by DÜPERTHAL to carry out their activities.

## 2.6 Stored goods

Storage, handling and use of the stored goods must comply with the applicable national standards and regulations, e.g. "TRGS 510" or "TRGS 526" in Germany.

Hazardous substances must be labelled at least with the name of the substance, the components of the preparation as well as the hazard symbols with the associated indications of danger. Furthermore, they must be marked with the safety measures to be observed (H and P statements) or the corresponding safety data sheet must be provided.

Hazardous substances must be stored in such a way that no dangerous reactions are possible if the containers are damaged. This must be observed as part of the risk assessment. Detailed information on the combined storage of different classes of substances can be found in the German "TRGS 510", for example.

### NOTE!

*Only store the stored goods in suitable, tightly sealed containers.*

- *Corrosive liquids or vapours can escape from containers that are not tightly sealed or are unsuitable and cause damage to the safety storage cabinet.*

## 2.7 Hazardous areas and their labelling

The following must be attached to the front of the safety storage cabinet and must be clearly visible:

- Name or trademark of the manufacturer
- Serial number and year of manufacture
- Information on the largest individual container volume that can be stored
- Information on the shelves' maximum load capacity

Furthermore, the following signs must be attached to the front of the safety storage cabinet and must be clearly visible:

Table 11: Prohibited action signs

Prohibited action signs	Meaning	Standard
	P011: Do not extinguish with water	DIN EN ISO 7010

Table 12: Warning signs

Warning signs	Meaning	Standard
	W023: Warning for corrosive substances	DIN EN ISO 7010

Table 13: Mandatory action signs

Mandatory action signs	Meaning	Standard
	M002: Read operator's manual	DIN EN ISO 7010
	M004: Wear eye protection	DIN EN ISO 7010
	M009: Wear hand protection	DIN EN ISO 7010

### 3 Technical specifications

#### 3.1 General data


Fig. 14: General diagram of ACID C pro safety storage cabinet

- 1 Exhaust air connection (DN 75)
- 2 Type plate
- 3 Standing surface
- 4 Plinth panel with adjustable feet


Fig. 15: General diagram of ACID B pro safety storage cabinet

- 1 Exhaust air connection (DN 75)
- 2 Type plate
- 3 Standing surface

### 3.2 Dimensions and equipment

#### 3.2.1 ACID C pro


Fig. 16: ACID C pro dimensions

H	Cabinet height	LB	Clear width
T	Cabinet depth	TS	Standing surface depth
B	Cabinet width	BS	Standing surface width
T1	Cabinet depth with open doors	H1-H6	Adjustable heights 1 to 6
V	Maximum collection volume	M	Pull-out tray load-bearing capacity

L Maximum volume of largest individual container

Table 14: ACID C pro S – XL dimensions

ACID C pro				
Model	S	XS	M	XL
H (mm)	1385	1385	2045	2045
B (mm)	600	1200	600	1200
T (mm)	609	609	609	609
T1 (mm)	1207	1207	1207	1207
LB (mm)	553	2x553	553	2x553
TS (mm)	404	404	404	404
BS (mm)	474	474	474	474
L (l)	17	17	17	17
V (l)	19	19	19	19
M (kg)	40	40	40	40
Empty weight (kg)	75	150	110-115 <sup>1</sup>	220-230 <sup>1</sup>
Max. payload (kg)	240	2x240	240	2x240
Total permitted weight (kg)	315	630	350-355 <sup>1</sup>	700-710 <sup>1</sup>

<sup>1</sup> Empty weight varies depending on number of pull-out trays

Table 15: ACID C pro S - XL maximum adjustable heights

ACID C pro						
Model	S	XS	M	XL	M	XL
H1 (mm)	305	305	460	460	310	310
H2 (mm)	305	305	459	459	310	310
H3 (mm)	305	305	460	460	294	294
H4 (mm)	300	300	459	459	310	310
H5 (mm)	---	---	---	---	310	310
H6 (mm)	---	---	---	---	294	294

### 3.2.2 ACID B pro and ACID U pro


Fig. 17: ACID B pro and ACID U pro dimensions

H	Cabinet height	LW	Clear width
T	Cabinet depth	TS	Standing surface depth
B	Cabinet width	BS	Standing surface width
T1	Cabinet depth with open doors	H1-H2	Adjustable heights 1 to 2
V	Maximum collection volume	M	Pull-out tray load-bearing capacity
L	Maximum volume of largest individual container	HS	Base height

Table 16: ACID B pro S - L dimensions

ACID B pro					
Model	S	L		S	L
H <sup>1</sup> (mm)	770	770	H1 (mm)	350	350
HS (mm)	20-135	20-135	H2 (mm)	350	350
B (mm)	600	1200	L (l)	17	17
T (mm)	609	609	V (l)	19	19
T1 (mm)	1207	1207	M (kg)	40	40
LB (mm)	553	2x553	Empty weight (kg)	40	80
TS (mm)	404	404	Max. payload (kg)	80	2x80
BS (mm)	474	474	Total permitted weight (kg)	120	240

<sup>1</sup> Height without base

Table 17: ACID U pro S - L dimensions

ACID U pro					
Model	S	L		S	L
H <sup>1</sup> (mm)	596	596	H1 (mm)	264	264
HS (mm)	20-135	20-135	H2 (mm)	264	264
B (mm)	600	1200	L (l)	17	17
T (mm)	609	609	V (l)	19	19
T1 (mm)	1207	1207	M (kg)	40	40
LB (mm)	553	2x553	Empty weight (kg)	35	70
TS (mm)	404	404	Max. payload (kg)	80	2x80
BS (mm)	474	474	Total permitted weight (kg)	115	230

<sup>1</sup> Height without base

### 3.3 Pressure drop during ventilation

Technical ventilation of the safety storage cabinet results in a pressure drop at the exhaust air connection, as shown in the following diagram. The maximum pressure drop during ventilation must not exceed 150 Pa.


Fig. 18: Pressure drop at exhaust air connection socket, averaged over all cabinet sizes

◆ Pressure drop at 30 m<sup>3</sup>/h

Table 18: Pressure drop at exhaust air flow rate of 30 m<sup>3</sup>/h

ACID C pro		
Model	Q [m <sup>3</sup> /h]	Δp [Pa]
S	30	8
XS	2 x 30	2x8
M	30	10
XL	2 x 30	2 x 10
ACID B pro		
Model	Q [m <sup>3</sup> /h]	Δp [Pa]
S	30	10
L	30	2 x 10
ACID U pro		
Model	Q [m <sup>3</sup> /h]	Δp [Pa]
S	30	12
L	30	2x12

## 4 Structure and function

### 4.1 Construction

- Cabinet carcass and doors made of decor panels with melamine resin coating
- Plastic (HDPE) pull-outs with handle and tray insert

### 4.2 Exhaust air connection and air supply opening

The safety storage cabinets can be connected to an industrial exhaust air system which discharges air outside in a safe location.

There are one or two exhaust air connection sockets with a diameter of DN 75 for this purpose on the cabinet roof of the ACID C pro series and on the rear of the safety storage cabinet on the ACID B pro and ACID U pro series.

#### **NOTE!**

*If aggressive and corrosive media are stored in ACID series cabinets, a connection to a continuous ventilation system is mandatory.*

*The volumetric flow rate must be at least **30 m<sup>3</sup>/h** for each side panel.*

### 4.3 Doors

The cabinets in the ACID U pro and ACID B pro cabinet series are equipped with one or two wing doors. The safety storage cabinets in the ACID C pro series are equipped with one to four wing doors.

The doors are lockable with the integrated locking cylinders. The key numbers are embossed on the locking cylinder and on the keys supplied.

#### **NOTE!**

*Keep the doors of the safety storage cabinet closed and re-close them manually after each work process.*

## 4.4 ACID C pro, ACID B pro and ACID U pro interior fittings

### 4.4.1 Pull-out trays


Fig. 19: Pull-out with plastic insert

#### 1 Pull-out with plastic insert

The safety storage cabinets in the ACID C pro, ACID B pro and ACID U pro series contain pull-outs with handles and plastic inserts made of HDPE.

#### **NOTE!**

*The maximum height of the top shelf must not exceed 1750 mm.*

### 4.4.2 Bottom tray (optional)

The safety storage cabinet can be fitted with plastic (PP) bottom trays as an option. The plastic grid insert enables these to be used as an additional standing surface.

### 4.4.3 Side storage inserts (optional)

As an alternative to the pull-out trays, plastic (PP) shelves can be attached to the side panels in addition to the bottom tray.

## 5 Transport

The safety storage cabinet is packaged for transport.

### WARNING!

**Crush hazard due to safety storage cabinet tipping over**

**If the safety storage cabinet tips over when not transported with due caution, this can cause potentially fatal crushing.**

- Wear personal protective equipment (PPE).
- Transport with at least two people.
- Only transport the safety storage cabinet upright and unladen.

### NOTE!

*Tipping the safety storage cabinet over during transport*

*Damage to the safety storage cabinet caused by incorrect handling.*

- Only pick up the safety storage cabinet using special and suitable transport or lifting equipment.
- Pick up the safety storage cabinet once it is securely lashed and is not at risk of slipping.

## 6 Installation and commissioning

### 6.1 Requirements for the installation location

*Install the safety storage cabinet so that inspection and maintenance activities can be carried out without restriction.*

The safety storage cabinet is approved for installation in a building.


Fig. 20: Requirements for the installation location

Consider the following in relation to the installation location:

- The surface must be able to bear the weight of the safety storage cabinet when fully loaded.
- The surface must be horizontal in order to guarantee problem-free functioning of the safety storage cabinet.
- Protect the safety storage cabinet against moisture.
  - Use at a relative humidity of >70 % in closed and heated buildings is permissible for a few weeks each year.
- The operating temperature must be between -5 °C and + 40 °C.

## 6.2 Installing the safety storage cabinet

### 6.2.1 ACID C pro

**i** Safety storage cabinets delivered to the place of use are normally supplied without packaging and transport skids.


Fig. 21: ACID C pro installation

- 1 Transport equipment for the safety storage cabinet
- 2 Safety storage cabinet
- 3 Transport skid
- 4 Attach the safety transport skids (4x screws, size 19 mm)
- 5 Adjustable foot

- ➔ Remove the packaging and detach the accessories.
- ➔ Lift the safety storage cabinet smoothly using suitable transport equipment (A).
- ➔ Loosen the screws for the transport skids and remove the transport skids (A and B).
- ➔ Screw the standing feet completely into the prepared holes from below (C and D).
- ➔ Position the safety storage cabinet at the installation location.
- ➔ Align the safety storage cabinet (↺ Figure 10 on page 47).

**WARNING!**

The safety storage cabinet must be completely emptied for alignment. Any containers still inside can tip over, spill or break. This can lead to serious burns to operating personnel.

- Before alignment, check that the safety storage cabinet is empty.

The safety storage cabinet can be aligned by screwing the adjustable feet in or out. To do this, the safety storage cabinet can be raised slightly.


Fig. 22: Alignment of a two-part ACID C pro XL

- Align the safety storage cabinet by screwing the adjustable feet in or out.
  - The adjustable feet can be screwed in or out by hand from the outside.
  - The adjustable feet can be screwed in or out from the inside using a hexagon key (A).
- Check the alignment of the safety storage cabinet with a spirit level.
- Seal the holes on the base of the cabinet using protective caps (A).
- Join together two-part cabinets using connecting screws (B).
- Attaching the safety storage cabinet to the wall (☞ Chapter 6.3: “Wall mounting” on page 48).

### 6.2.2 ACID B pro and ACID U pro

- Remove the packaging and detach the accessories.
- Lift the safety storage cabinet smoothly from the transport pallet.
- Position the safety storage cabinet at the installation location.
- Check the alignment of the safety storage cabinet using suitable equipment, e.g. a spirit level.
- To align the safety storage cabinet, raise it slightly and slide suitable spacer plates under the base.
- If required, attach the safety storage cabinet to the wall (☞ Chapter 6.3: “Wall mounting” on page 48).

### 6.3 Wall mounting

To guarantee stability, ACID C pro safety storage cabinets must be attached to a wall using the mounting brackets supplied. The other safety storage cabinets in the series can be attached to a wall if required.

**WARNING!**

**Unmounted safety storage cabinets, particularly the ACID C pro M and XL models, can tip over. This can lead to serious injuries or death due to crushing and burns due to escaping stored materials.**

- Before use, ensure that the safety storage cabinet is mounted on a stable wall by technically qualified personnel.


Fig. 23: Wall mounting

- 1 Mounting bracket
- 2 Mounting screws

Wall mounting:

Personnel:

- Technically qualified personnel
- ➔ Position the mounting bracket on the cabinet roof.
- ➔ Attach the mounting bracket to the cabinet roof using the mounting screws supplied.
- ➔ Attach the second arm of the mounting bracket to the wall with appropriate mounting materials provided on site.

*On the ACID U pro and ACID B pro safety storage cabinet lines, the mounting brackets can also be attached to the side panel.*

**NOTE!**

*The condition of the wall must always be checked before installation. Only use mounting materials that are appropriate for the relevant wall condition and for the load.*

**WARNING!**

**There must not be any electrical installation zones or supply and waste disposal lines in the mounting area.**

- Use a suitable detector to check the mounting area before drilling. Note that plastic pipes cannot be detected.

## 6.4 Connection to an exhaust air system

**i** Installation of industrial ventilation and connection to an existing exhaust air system must be carried out by a qualified company and is not a service provided by DÜPERTHAL.

**i** The power of the exhaust air system can be determined using the technical specifications, ↪ Chapter 3.3 “Druckabfall bei Entlüftung” on page 42.

### NOTE!

To prevent damage caused by corrosive atmospheres inside the cabinet, the safety storage cabinet must be vented with a volumetric flow rate of **min. 30 m<sup>3</sup>/h** per side panel.

**!** If large storage quantities, very corrosive substances or substances in sub-optimal containers are stored in the safety storage cabinet, a higher exhaust air flow than that recommended should be used.

If metal parts installed in the cabinet, e.g. the hinges, exhibit any signs of corrosion, the set air exchange rate is too low.


Fig. 24: Cabinet roof with exhaust air connection socket (ACID C pro)

1 Exhaust air connection socket

Connection to the exhaust air system:

Personnel:

- Technically qualified personnel
- ➔ Secure the exhaust air connection socket with reducing fitting in the specified hole in the cabinet roof (ACID C pro models only).
- ➔ Connect the waste air line to the exhaust air connection socket on the cabinet roof (ACID C pro) or on the rear of the cabinet (ACID B pro and ACID U pro) and if necessary secure with a collar.
- ➔ After installing the safety storage cabinet, check the connection to the exhaust air system with smoke tubes.

*DÜPERTHAL stipulates that acid and lye cabinets must be connected to a continuous industrial ventilation system. If the cabinets are operated without continuous venting, DÜPERTHAL will not accept any complaints relating to the exhaust air.*

## 7 Operation

### 7.1 Opening the safety storage cabinet

Unlock the locking cylinder with the key and open the safety storage cabinet by pulling the door handle on the wing door.

The wing door remains open in any position.

### 7.2 Pull-out tray

Personnel:

- Technically qualified personnel
- ➔ Extract the pull-out from the safety storage cabinet using the handle.
- ➔ When work on the safety storage cabinet is complete, slide the pull-out back in.

### 7.3 Checking the pull-out trays

Personnel:

- Technically qualified personnel
- ➔ Perform a daily visual inspection for extraneous substances.
- ➔ Collect and properly dispose of spilled hazardous substances immediately in accordance with accident prevention regulations.

#### **NOTE!**

*Store hazardous substances so that a visual inspection of the pull-out tray for extraneous substances is possible every working day.*

## 8 Maintenance

Check the safety storage cabinet for any externally visible damage or defects.

Always perform checks:

- After installation.
- Before commissioning.
- After modifications or maintenance work

The safety storage cabinet should also be inspected periodically at the following intervals.

Interval	Maintenance work	Personnel
Daily	Pull-out trays: <ul style="list-style-type: none"> <li>▪ Check in accordance with regulations governing water legislation</li> <li>▪ Collect and properly dispose of leaked liquids immediately.</li> </ul>	Technically qualified personnel
Interval	Maintenance work	Personnel
Monthly	Closing of the doors closing system <ul style="list-style-type: none"> <li>▪ Open the door and inspect the closure.</li> </ul>	Technically qualified personnel
	Ventilation <ul style="list-style-type: none"> <li>▪ Check the effectiveness of the ventilation with a woollen thread or with a smoke tube in the cabinet in front of the exhaust air duct at the ventilation grilles.</li> </ul>	Technically qualified personnel
	Labelling <ul style="list-style-type: none"> <li>▪ Inspect the labelling on the safety storage cabinet for completeness.</li> </ul>	Technically qualified personnel
Interval	Maintenance work	Personnel
Annually	Safety storage cabinet Test on the entire safety storage cabinet.	DÜPERTHAL service technicians

**i** *If faults occur, assist the technical customer service by providing the cabinet model, production and key number, along with a description of the fault.*

**i** *Technical storage facilities must be tested annually by a qualified person (as defined in TRBS 1203 for example) in compliance with BetrSichV and the maintenance interval stipulated by the manufacturer. The applicable national regulations must be observed.*

## 9 Faults

Fault description	Cause	Remedy	Personnel
Doors do not close or open automatically.	Safety storage cabinet is not aligned correctly.	Install the safety storage cabinet so it is horizontal. ↪ Chapter 6.2: "Sicherheitsschrank aufstellen" on page 46	Technically qualified personnel
Doors and pull-outs are sluggish.	Moving parts, such as hinges or pull-out rails, are dirty or corroded.	<ul style="list-style-type: none"> <li>• Clean the pull-out rails, hinges and closing mechanism.</li> <li>• Notify technical customer service.</li> </ul>	Technically qualified personnel
Corrosion to metal parts	Insufficient exhaust air flow	Increase exhaust air flow ↪ Chapter 6.4: "Anschluss an ein Abluftsystem" on page 50	Technically qualified personnel
	Container not tightly sealed	Replace containers that are not tightly sealed	Technically qualified personnel

## 10 Spare parts and accessories

**i** Only original parts from DÜPERTHAL are to be used for the safety storage cabinets.

- PP tray inserts
- Pull-out rails
- Exhaust air connection socket
- Mounting bracket
- Ventilators
- Exhaust air monitoring units
- Filter systems
- Plinth panel
- Door handle
- Hinges

## 11 Disposal

**WARNING!**

**Risk of injury due to improper dismantling of the safety storage cabinet.**

- Ensure that the safety storage cabinet is dismantled by technically qualified personnel only.

The safety storage cabinet can be completely dismantled by specialist technical personnel.

Comply with national and local disposal regulations.

To save resources, do not place parts of the safety storage cabinet or the whole cabinet in bulky or domestic waste.