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User's Manual

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
Contents

1	Introduction	5
1.1	Intended Use	5
1.2	Occupational Safety	5
2	Installing the Unit and Understanding its Functions	5
2.1	Installing and Checking out	5
2.2	Exhaust Air Extractor Components	6
2.3	Understanding Operation	6
3	Operation	6
3.1	Normal Operation	6
4	Technical Specifications	7
5	Troubleshooting	7
6	Contact Details	8

1 Introduction

Dear Customer:

Thank you very much for purchasing our product. We are confident that it will meet your expectations and perform satisfactorily in your daily routines.

The development, manufacture, and check of the ATEX-compliant exhaust airflow unit with fan has been guided by our concern to offer a device characterized by superior operation safety and user-friendliness. The blower was manufactured and checked according to state-of-the-art technology and widely recognized safety regulations. It is approved for use in explosible atmospheres.

1.1 Intended Use

The exhaust air extractor with fan has to be used in EX-free areas and is intended exclusively for the positive extraction of exhaust air from safety cabinets for storage installed indoors, zone 2 for gases belonging to temperature classes T1 - T4, gas type IIA or gas type IIB. The product is suitable only for fluid flows not transporting any particles or foreign bodies. The maximum temperature of extracted medium is 50°C (122°F).

Only such substances are acceptable for extraction that do not compromise the PP's connector sleeve, the PVC suction hose, and the fan. The relevant list of substances can be provided by the manufacturer in a digital format.

Caution

Any and all uses other than the above-described is prohibited and may cause hazardous situations, potentially danger of death, damage to the system and material property, and will cancel any and all warranty and liability for the product on the part of its manufacturer and distributor.

1.2 Occupational Safety

The installation and operation of the product must be performed in compliance with any and all nationally applicable rules and regulations for the storage of flammable liquids, the ISO 13857 standard, and any and all other applicable codes and standards. The design of the product is such that the fan blades are off-touch. Modifications to the enclosure or the suction and exhaust hoses that may make the fan blades accessible to the touch when the unit is in normal use are prohibited.

Caution

The product must never be used if not totally installed, with properly set suction and exhaust hoses. Electrical work on the product must only be performed when the exhaust air extractor was previously disconnected from power and/or locked against unwanted reapplication of power. This also holds for electrical work on the extraction system of which the extractor may be a part.

2 Installing the Unit and Understanding its Functions

2.1 Installing and Checking out

1. Connect the extractor to the air extraction system on site or the safety cabinet to be served:
Connect the external ventilation system or the exhaust outlet of the cabinet to the unit, e.g., with a flexible sleeve fitting, making sure

the connection is under no strain, and make a tight connection using an appropriate hose clamp. This is indispensable to ensure proper operation of the extraction fan. As a general rule, be sure that the system is set up with air lines that are as short as possible.


Caution

In no circumstances, it is acceptable to transfer the air extracted by the blower fan back into the work room.

2. Connect the exhaust airflow extractor to power.

2.2 Exhaust Air Extractor Components

2.3 Understanding Operation

After the On/Off switch was pressed to turn the extractor on, the fan located inside its enclosure draws air from the served system through the suction hose connected to the inlet nipple and extracts the air from the system through the blower nipple on the enclosure (see Figures in section 2.2).

3 Operation

Caution

Be sure that all provisions of Sect. 1.1 and 1.2 above are strictly complied with for both the installation and the operation. Failure to do so might expose you to hazardous situations.

3.1 Normal Operation

In normal operation, the fan extracts the air from the connected system to the outside after it was switched on.


It is under the operator's responsibility to ensure proper monitoring of the unit, to be sure it performs as intended.

The operation characteristics of the fan are given on the graph to the right.

Important


The fan has to be checked at regular intervals. If it is dirty, damaged or shows signs of corrosion, it has to be cleaned or the defects have to be corrected, if needed by returning the blower box for repair to the seller or manufacturer.

4 Technical Specifications

Power Supply Voltage	110 V
Frequency	50 Hz
Current Rating	1.4 A (H model version) 1.6 A (L model version)
Motor Power	100 W (H model version) 80 W (L model version)
Revolutions	2810 1/min (H model version) 1300 1/min (L model version)
Fluid Velocity	5.1 m/s (H model version) 4.3 m/s (L model version)
Volume Flow Rate	72 m³/h (H model version) 60 m³/h (L model version)
Weight	approx. 7.0 kg (H model version) approx. 7.5 kg (L model version)
Dimensions (enclosure edges)	260 x 290 x 275 mm
Ingress Protection	IP20
ATEX Rating	EX II (3) G [c II B T4]

see also the component name plate, as to the right

Important

Values given for the exhaust airflow were determined without subjecting the fan box to any load. When it is connected to a storage cabinet, the velocity and the volume flow rate may differ, potentially substantially.

5 Troubleshooting

In the course of using the appliance, malfunctions may occur. The table below is meant to help troubleshoot and correct malfunctions. If a malfunction cannot be traced to one of the possible causes or the correction given does not solve the issue, then use of the appliance has to be discontinued and it has to be returned for inspection and repair.

Malfunction	Possible Cause	Verification of Cause	Correction of Malfunction
Appliance does not run	On/Off Switch Power Supply Cord Overvoltage Switch	- Is the appliance switched on at the On/Off switch? - Is the power supply cord damaged? - Is the power supply cord properly connected? - Are the fuses inside the appliance o.k.?	- Switch the appliance on at the On/Off switch. - Replace the supply cord with a new one. - Connect the supply cord to the intended socket. - Replace the fuses (2 x 1.25 A slow blow).

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Web: www.dueperthal.com

E-mail: info@dueperthal.com

We reserve the right to change technical specs of the unit without prior notice.

Pursuant to the directive 2012/19/EU on Waste Electrical and Electronic Equipment (so-called WEEE Directive), equipment that was manufactured or distributed by ourselves may be returned to us. To discuss the exact procedure we will follow, please contact us with the details above.

We also declare that the equipment manufacture conforms to the directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (so-called RoHS Directive).

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