ACID pure
Sicherheitsschränke

**ACID
pure**

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Table of contents

1	General information	27
1.1	Notes for reading	27
1.2	Type plate	27
2	Safety	28
2.1	Function of safety notices	28
2.2	Correct use	28
2.3	Misuse	28
2.4	Obligations of the operator	29
2.5	Requirements for employees	30
2.6	Stored goods	30
2.7	Hazardous areas and their labelling	31
3	Technical specifications	33
3.1	General data	33
3.2	Dimensions and equipment	34
3.3	Pressure drop during ventilation	36
4	Structure and function	37
4.1	Construction	37
4.2	Exhaust air connection and air supply opening	37
4.3	Doors	37
4.4	ACID pure interior fittings	38
5	Transport	39
6	Installation and commissioning	40
6.1	Requirements for the installation location	40
6.2	Installing the safety storage cabinet	40
6.3	Connection to an exhaust air system	41
7	Operation	43
7.1	Opening the safety storage cabinet	43
7.2	Storage shelves	43
7.3	Checking the bottom tray and storage shelves	44
7.4	Cleaning the safety storage cabinet	44
8	Maintenance	45
9	Faults	46
10	Spare parts and accessories	47
11	Disposal	47

12 General information

12.1 Notes for reading

The following symbols designate specific types of information.

Table 8: Explanation of symbol

Symbol	Type of information
i	Information for easier and more effective working
➔	Procedural step
✓	Result of a procedural step
🔗	Link to another part of the document

12.2 Type plate

The type plate is attached to the outside of the safety storage cabinet door.

Fig. 9: Type plate

- 1 Model
- 2 Serial number and year of manufacture
- 3 Maximum volume of individual containers
- 4 Maximum load per standing surface

13 Safety

13.1 Function of safety notices

Safety notices warn against physical or material damage and provide information on how such damage can be avoided.

The following signal words identify the degree of danger and the extent of the risk.

WARNING!

The signal word "WARNING" refers to a potential hazard which could result in death or serious injury.

NOTE!

The signal word "NOTE" indicates a situation that can lead to damage to the safety storage cabinet.

13.2 Correct use

Observe the safety instructions in these operating instructions to reduce health risks and avoid dangerous situations.

Any improper use as defined in these operating instructions involves a risk of accidents and damage.

The safety storage cabinet should be used for storage of solid and liquid, non-flammable corrosive chemicals in working spaces. It is particularly suitable for storage of acid and lye.

13.3 Misuse

Any use that goes beyond the specified correct use is considered to be misuse.

DÜPERTHAL accepts no liability for damage arising from misuse.

In addition, the following safety notices must be observed:

WARNING!

Storage of living organisms in the safety storage cabinet

Living organisms can come into contact with the stored hazardous substances. This can lead to death or serious injury.

- Use the safety storage cabinet exclusively for storage of chemicals.

WARNING!

Storage of food in the safety storage cabinet

Food can come into contact with the stored hazardous substances. This can lead to death or serious injury.

- Use the safety storage cabinet exclusively for storage of chemicals.

WARNING!

Storage of flammable liquids in the safety storage cabinet

Flammable liquids may explode in case of fire.

This can lead to death or serious injury.

- Store flammable liquids in a Type 90 safety storage cabinet.

WARNING!

Storage of hazardous substances together

Risk of uncontrolled chemical reactions.

This can lead to death or serious injury.

- Only store substances and preparations that are permitted to be stored together in the safety storage cabinet (e.g. see TRGS 510).

WARNING!

Spilled liquids

Inhalation of hazardous vapours from spilled liquids can lead to life-threatening injuries to the respiratory system.

Spilled liquids can lead to life-threatening burns.

- Collect and properly dispose of spilled liquids immediately in accordance with accident prevention regulations.

NOTE!

Alteration and modifications

Do not alter or modify the safety storage cabinet.

- *If alteration or modification of the safety storage cabinet is required, contact DÜPERTHAL.*

13.4 Obligations of the operator

The operator is obliged to comply with applicable legal regulations. This includes:

- Carrying out risk assessments.
- Issuing operating instructions or laboratory regulations and defining the scope of application for the safety storage cabinet.
- Defining and documenting activities of designated employees.

13.5 Requirements for employees

WARNING!

Employees who do not meet these requirements

This can lead to death or serious injury.

- Designate employees who meet the requirements to carry out activities.

These operating instructions set out the following employee activities:

Technically qualified personnel:

Only people who have been trained by the operator in use of the safety storage cabinet and handling of the stored goods are approved as specialist technical employees.

DÜPERTHAL service technicians:

DÜPERTHAL employees are specifically trained by DÜPERTHAL to carry out their activities.

13.6 Stored goods

Storage, handling and use of the stored goods must comply with the applicable national standards and regulations, e.g. "TRGS 510" or "TRGS 526" in Germany.

Hazardous substances must be labelled at least with the name of the substance, the components of the preparation as well as the hazard symbols with the associated indications of danger. Furthermore, they must be marked with the safety measures to be observed (H and P statements) or the corresponding safety data sheet must be provided.

Hazardous substances must be stored in such a way that no dangerous reactions are possible if the containers are damaged. This must be observed as part of the risk assessment. Detailed information on the combined storage of different classes of substances can be found in the German "TRGS 510", for example.

NOTE!

Only store the stored goods in suitable, tightly sealed containers.

- *Corrosive liquids or vapours can escape from containers that are not tightly sealed or are unsuitable and cause damage to the safety storage cabinet.*

13.7 Hazardous areas and their labelling

The following must be attached to the front of the safety storage cabinet and must be clearly visible:

- Name or trademark of the manufacturer
- Serial number and year of manufacture
- Information on the largest individual container volume that can be stored
- Information on the shelves' maximum load capacity

Furthermore, the following signs must be attached to the front of the safety storage cabinet and must be clearly visible:

Table 9: Prohibited action signs

Prohibited action signs	Meaning	Standard
	P011: Do not extinguish with water	DIN EN ISO 7010

Table 10: Warning signs

Warning signs	Meaning	Standard
	W023: Warning for corrosive substances	DIN EN ISO 7010

Table 11: Mandatory action signs

Mandatory action signs	Meaning	Standard
	M002: Read the instructions	DIN EN ISO 7010
	M004: Wear eye protection	DIN EN ISO 7010
	M009: Wear hand protection	DIN EN ISO 7010

14 Technical specifications

14.1 General data


Fig. 10: General diagram of ACID pure safety storage cabinet

- 1 Exhaust air connection (adapter DN100)
- 2 Air supply connection (adapter DN100)
- 3 Type plate
- 4 Standing surfaces

14.2 Dimensions and equipment


Fig. 11: ACID pure dimensions

H	Cabinet height	LB	Clear width
T	Cabinet depth	TS	Standing surface depth
B	Cabinet width	BS	Standing surface width
T1	Cabinet depth with open doors	H1–H4	Adjustable heights 1 to 4
V	Maximum collection volume	M	Load-bearing capacity per standing surface
L	Maximum volume of largest individual container		

Table 12: ACID pure M and XL dimensions

ACID pure		
Model	M	XL
H (mm)	1955	1955
B (mm)	595	1195
T (mm)	601	601
T1 (mm)	1142	1176
LB (mm)	504	1073
TS (mm)	474	474
BS (mm)	529	1129
L (l)	12	25
V (l)	14	28
M (kg)	60	60
Empty weight (kg)	155	255
Max. payload (kg)	240	240
Total permitted weight (kg)	395	495

Table 13: ACID pure M and XL maximum adjustable heights

ACID pure		
Model	M	XL
H1 (mm)	375	375
H2 (mm)	479	479
H3 (mm)	479	479
H4 (mm)	459	459

14.3 Pressure drop during ventilation

Technical ventilation of the safety storage cabinet results in a pressure drop at the exhaust air connection, as shown in the following diagram. The maximum pressure drop during ventilation must not exceed 150 Pa.


Fig. 12: Pressure drop at exhaust air connection socket, averaged over all cabinet sizes

◆ Pressure drop at 30 m³/h

Table 14: Pressure drop Δp at the exhaust air connection socket and in the cabinet interior

ACID pure			
Model	Q [m³/h]	Δp at the exhaust air connection socket [Pa]	Δp in the interior [Pa]
M	30	10	<1
XL	30	10	<1

15 Structure and function

15.1 Construction

- Outer carcass and door trim made of powder-coated sheet steel
- Air ducts, door handles and fittings made of plastic
- Interior trim made of plastic (PP)
- Interior fittings made of plastic (PP, PA)

15.2 Exhaust air connection and air supply opening

The safety storage cabinets can be connected to an industrial exhaust air system which discharges air outside in a safe location.

For this reason, cabinets in the ACID pure series are equipped with a connection with a DN100 adapter for a feed line and another for an exhaust air line on the cabinet roof.

NOTE!

If aggressive and corrosive media are stored in ACID series cabinets, a connection to a continuous ventilation system is mandatory.

*The volumetric flow rate must be at least **30 m³/h**.*

15.3 Doors

The cabinets in the ACID pure cabinet series are equipped with one or two wing doors.

The doors are lockable with an integrated locking cylinder.

NOTE!

Keep the doors of the safety storage cabinet closed and re-close them manually after each work process.

15.4 ACID pure interior fittings

15.4.1 Storage shelf


Fig. 13: Storage shelves with mounting components

The ACID pure safety storage cabinets are equipped with storage shelves made of plastic (PP). The adjustable heights are variable in a grid of approx. 128 mm.

NOTE!

The maximum height of the top shelf must not exceed 1750 mm.

15.4.2 Bottom tray

The ACID pure safety storage cabinets are equipped with a bottom tray made of plastic (PP). This can also be used as an additional standing surface. If it is used as a standing surface, the maximum collection volume is reduced.

NOTE!

If the bottom tray is used as a standing surface, the maximum collection volume is reduced by the proportional volume of the containers placed on it.

16 Transport

The safety storage cabinet is packaged for transport on a pallet.

WARNING!

Crush hazard due to safety storage cabinet tipping over

If the safety storage cabinet tips over when not transported with due caution, this can cause potentially fatal crushing.

- *Wear personal protective equipment (PPE).*
- *Transport with at least two people.*
- *Only transport the safety storage cabinet upright and unladen.*

NOTE!

Damage to the safety storage cabinet caused by incorrect handling.

- *Only pick up the safety storage cabinet using special and suitable transport or lifting equipment.*
- *Pick up the safety storage cabinet once it is securely lashed and is not at risk of slipping.*

17 Installation and commissioning

17.1 Requirements for the installation location

i Install the safety storage cabinet so that inspection and maintenance activities can be carried out without restriction.

The safety storage cabinet is approved for installation in a building.


Fig. 14: Requirements for the installation location

Consider the following in relation to the installation location:

- The surface must be able to bear the weight of the safety storage cabinet when fully loaded.
- The surface must be horizontal in order to guarantee problem-free functioning of the safety storage cabinet.
- Protect the safety storage cabinet against moisture.
 - Use at a relative humidity of >70 % in closed and heated buildings is permissible for a few weeks each year.
- The operating temperature must be between -5 °C and + 40 °C.

17.2 Installing the safety storage cabinet

Personnel:

- Technically qualified personnel
- Remove the packaging and detach the accessories.
- Lift the safety storage cabinet smoothly from the transport pallet.
- Position the safety storage cabinet at the installation location.
- Check the alignment of the safety storage cabinet using suitable equipment, e.g. a spirit level.
- To align the safety storage cabinet, raise it slightly and slide suitable spacer plates under the base.

i The gap between the base of the safety storage cabinet and the floor of the clean room and between the rear or side panel and the wall of the clean room can be sealed using a suitable sealant, such as silicone.

17.3 Connection to an exhaust air system

i Installation of industrial ventilation and connection to an existing exhaust air system must be carried out by a qualified company and is not a service provided by DÜPERTHAL.

i The power of the exhaust air system can be determined using the technical specifications, ↪ Chapter 3.3 "Druckabfall bei Entlüftung" on page 14.

NOTE!

To prevent damage caused by corrosive atmospheres inside the cabinet, the safety storage cabinet must be vented with a volumetric flow rate of **min. 30 m³/h**.

! If large storage quantities, very corrosive substances or substances in sub-optimal containers are stored in the safety storage cabinet, a higher exhaust air flow than that recommended should be used.

If the cabinet's integral components show signs of corrosion, the air exchange rate is set too low.

i DÜPERTHAL stipulates that acid and lye cabinets must be connected to a continuous industrial ventilation system. If the cabinets are operated without continuous venting, DÜPERTHAL will not accept any complaints relating to the exhaust air.


Fig. 15: Cabinet roof with air supply and exhaust air connection sockets

- 1 Air supply connection socket
- 2 Exhaust air connection socket

Connection to the exhaust air system:

Personnel:

- Technically qualified personnel
- ➔ Secure the exhaust air connection socket and air supply connection socket with reducing fitting in the specified opening in the cabinet roof.
- ➔ Connect the exhaust air line to the exhaust air socket on the cabinet roof and secure it with a collar if necessary.
- ➔ If necessary, a feed line can be connected to the air supply connection socket on the cabinet roof. This may need to be secured with a collar.
- ➔ After installing the safety storage cabinet, check the connection to the exhaust air system and the effectiveness of the system.

18 Operation

18.1 Opening the safety storage cabinet

Unlock the locking cylinder with the key and open the safety storage cabinet by pulling the door handle on the wing door.

The wing door remains open in any position.

18.2 Storage shelves

The height of the storage shelves can be individually adjusted.


Fig. 16: Storage shelves with mounting components

- 1 Locking system
- 2 Retaining bolts

Personnel:

- Technically qualified personnel
- ➔ Open the door(s) of the safety storage cabinet fully (approx. 180°).
- ➔ Loosen the locking system on the rear of the cabinet (1).
- ➔ Remove the storage shelf and re-position it on the retaining bolts (2).
- ➔ Fasten the locking system on the rear of the cabinet (1).

18.3 Checking the bottom tray and storage shelves

Personnel:

- Technically qualified personnel
- ➔ Perform a daily visual inspection for extraneous substances.
- ➔ Collect and properly dispose of spilled hazardous substances immediately in accordance with accident prevention regulations.

NOTE!

Store hazardous substances in such a way that a visual inspection of the standing surfaces to check for extraneous substances can be performed every working day.

18.4 Cleaning the safety storage cabinet

All surfaces of the safety storage cabinet can be cleaned with a soft, non-abrasive cloth. It is recommended that a small amount of natural cleaning agent be diluted in water for cleaning. Alternatively, isopropyl alcohol (IPA) can be used to clean the surfaces and ethanol can be used to clean interior fittings made of PP.

The application time of the cleaning agent must not exceed 1 hour.

NOTE!

The surfaces of the safety storage cabinet may be damaged by unsuitable cleaning agents, such as acidic or alkaline substances, or agents containing esters, ketones, polyvalent alcohols, aromatic substances, glycol ethers or halogen hydrocarbons.

19 Maintenance

Check the safety storage cabinet for any externally visible damage or defects.

Always perform checks:

- After installation.
- Before commissioning.
- After modifications or maintenance work

The safety storage cabinet should also be inspected periodically at the following intervals.

Interval	Maintenance work	Personnel
Daily	Storage shelves and bottom tray: ▪ Check in accordance with regulations governing water legislation ▪ Collect and properly dispose of leaked liquids immediately.	Technically qualified personnel
Interval	Maintenance work	Personnel
Monthly	Closing of the doors closing system ▪ Open the door and inspect the closure.	Technically qualified personnel
	Ventilation ▪ Check the effectiveness of the ventilation in the cabinet in front of the exhaust air opening.	Technically qualified personnel
	Labelling ▪ Inspect the labelling on the safety storage cabinet for completeness.	Technically qualified personnel
Interval	Maintenance work	Personnel
Annually	Safety storage cabinet Test on the entire safety storage cabinet.	DÜPERTHAL service technicians

i If faults occur, assist the technical customer service by providing the cabinet model, production and serial number, along with a description of the fault.

i Technical storage facilities must be tested annually by a qualified person (as defined in TRBS 1203 for example) in compliance with BetrSichV and the maintenance interval stipulated by the manufacturer. The applicable national regulations must be observed.

20 Faults

Fault description	Cause	Remedy	Personnel
Doors do not close or open automatically.	Safety storage cabinet is not aligned correctly.	Install the safety storage cabinet so it is horizontal. ↪ Chapter 6.2: "Sicherheitsschrank aufstellen" on page 18	Technically qualified personnel
Doors are sluggish.	Moving parts, such as hinges, are dirty or corroded.	• Clean the hinges and closing mechanism. • Notify technical customer service.	Technically qualified personnel
Corrosion on components	Insufficient exhaust air flow	Increase exhaust air flow ↪ Chapter 6.3: "Anschluss an ein Abluftsystem" on page 19f.	Technically qualified personnel
	Container not tightly sealed	Replace containers that are not tightly sealed	Technically qualified personnel

21 Spare parts and accessories

Only original parts from DÜPERTHAL are to be used for the safety storage cabinets.

- Exhaust air connection socket
- Exhaust air monitoring units
- Door handle
- Hinges

22 Disposal

WARNING!

Risk of injury due to improper dismantling of the safety storage cabinet.

- Ensure that the safety storage cabinet is dismantled by technically qualified personnel only.

The safety storage cabinet can be completely dismantled by specialist technical personnel. Comply with national and local disposal regulations.

To save resources, do not place parts of the safety storage cabinet or the whole cabinet in bulky or domestic waste.

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© 05/2021 DÜPERTHAL 55280, D00112786, 0, DE