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Dokument: D00121082 Rev.: 00

Table of contents

1	General information	40
1.1	Notes for reading.....	40
1.2	Type plate.....	41
2	Safety	42
2.1	Function of safety notices	42
2.2	Correct use.....	42
2.3	Misuse	43
2.4	Obligations of the operator	43
2.5	Requirements for employees	43
2.6	Stored goods	44
2.7	Hazardous areas and their labelling.....	44
3	Technical specifications	46
3.1	General data.....	46
3.2	Dimensions and equipment.....	48
3.3	Pressure drop during ventilation	50
4	Structure and function	51
4.1	Construction	51
4.2	Doors	51
4.3	Safety technology.....	52
4.3.1	Door closure in case of fire	52
4.3.2	Closure of venting cut-off flaps in case of fire	52
4.4	Interior fittings	52
5	Transport and packaging	53
6	Installation and commissioning.....	54
6.1	Requirements for the installation location	54
6.2	Attaching the adjustable feet.....	55
6.3	Aligning the safety storage cabinet	56
6.4	Checking the alignment of the safety storage cabinet	57
6.5	Fitting the plinth panel	58
6.6	Venting the safety storage cabinet.....	59
6.6.1	Exhaust air connection and air supply opening	59
6.6.2	Connection to an exhaust air system.....	60
6.6.3	Operating the safety storage cabinet without industrial ventilation.....	61
6.7	Earthing the safety storage cabinet.....	62
6.8	Changing the height of the storage shelves.....	63

Table of contents

7	Operation.....	64
7.1	Storage and removal	64
7.2	Checking and cleaning the bottom tray.....	65
8	Procedure after a fire	66
9	Maintenance.....	67
10	Faults	68
11	Spare parts and accessories	69
12	Disposal	70
13	Declaration of Conformity.....	71

1 General information

1.1 Notes for reading

The following symbols designate specific types of information

Table 7: Explanation of symbol

Symbol	Type of information
	Information for easier and more effective working
	Procedural step
	Result of a procedural step
	Link to another part of the document
	List

1.2 Type plate

The type plate is attached to the outside of the safety storage cabinet door.


Fig. 16: Type plate

- 1 Model
- 2 Item number (R/N), serial number (S/N) and year of manufacture
- 3 Maximum volume of individual containers
- 4 Maximum load per storage shelf

2 Safety

2.1 Function of safety notices

Safety notices provide warnings about physical or material damage and provide information on how such damage can be avoided.

The following signal words identify the degree of danger and the extent of the risk.

WARNING

The signal word “*WARNING*” refers to a potential hazard which could result in death or serious injury.

CAUTION

The signal word “*CAUTION*” refers to a potential hazard which could result in slight or minor injury.

NOTE

The signal word “*NOTE*” indicates a situation that can lead to damage to the safety storage cabinet.

2.2 Correct use

Observe the safety instructions in these operating instructions to reduce health risks and avoid dangerous situations.

Any use that is not correct use as defined in these operating instructions involves a risk of accidents and a lack of fire protection.

The BASIC standard line safety storage cabinets are type-tested and classified as Type 30 according to EN 144701 with a fire resistance of 30 minutes.

The safety storage cabinet is to be used for passive storage of flammable liquids in working spaces.

Passive storage is defined as storage exclusively in tightly sealed containers in a safety storage cabinet with no work such as filling, mixing or transfer.

2.3 Misuse

Any use that goes beyond the specified correct use is considered to be misuse.

DÜPERTHAL accepts no liability for damage arising from misuse.

The following constitute misuse:

- Storage of food in the safety storage cabinet.
- Transfer, filling and laboratory work in the safety storage cabinet.
- Incorrect storage of hazardous substances together.
- Objects on the cabinet roof.
- Blocking the doors of the safety storage cabinet or using objects to hold these doors open.
- Alteration and modification of the safety storage cabinet without the prior knowledge and approval of DÜPERTHAL.
- Using low-quality spare parts.
- Non-compliance with specified maintenance intervals.
- Collect and properly dispose of spilled liquids immediately in accordance with accident prevention regulations.

2.4 Obligations of the operator

The operator is obliged to comply with applicable legal regulations. This includes:

- Issuing operating instructions.
- Carrying out risk assessments.
- Creating explosion protection documents.
- Specifying activities by designated employees.
- Specifying personal protective equipment for employees.
- Ensuring that all employees have read and understood these operating instructions before using the safety storage cabinet for the first time.
- Ensuring that these operating instructions are available to employees at all times.

2.5 Requirements for employees

WARNING

Employees who do not meet these requirements!

This can lead to death or serious injury.

- Designate employees who meet the requirements to carry out activities.

Employees need to meet specific requirements in order to carry out certain activities. These instructions for use assign activities to the following employees:

- User
- Technically qualified personnel
- DÜPERTHAL service technicians

User

Only people who meet the following requirements are approved as employees:

- Instructed in the function and operation of the safety storage cabinet
- In-depth knowledge of how to correctly store and remove flammable and liquid hazardous substances.

Technically qualified personnel

Technical specialist employees also fulfil the following requirements in addition to the general requirements:

- Activity-related professional training in the technical area which is recognized at the installation location.
- In-depth knowledge and skills in relation to maintenance and repair.

DÜPERTHAL service technicians

DÜPERTHAL employees are specifically trained by DÜPERTHAL to carry out activities in relation to the safety storage cabinet.

2.6 Stored goods

Storage, handling and use of the stored goods must comply with the applicable national standards and regulations, e.g. "Technical Rule for Hazardous Substances 510" (TRGS 510) in Germany.

2.7 Hazardous areas and their labelling

The following must be attached to the front of the safety storage cabinet and must be clearly visible:

- The instruction "Close the door"
- Fire resistance in minutes (e.g. "Type 30")
- Name or trademark of the manufacturer
- Serial number and year of manufacture
- Information on the largest individual container volume that can be stored
- Specification of the shelves' maximum load capacity

Furthermore, the following signs must be attached to the front of the safety storage cabinet and must be clearly visible:

Table 8: Prohibited action signs

Symbol	Meaning	Standard
	P003: No naked flames; fire, open ignition source and smoking prohibited	DIN EN ISO 7010

Table 9: Warning signs

Symbol	Meaning	Standard
	W021: Warning: Flammable materials	DIN EN ISO 7010

Table 10: Mandatory action signs

Symbol	Meaning	Standard
	M002: Read the instructions	DIN EN ISO 7010
	Keep doors closed	Manufacturer's instructions
	Ventilation system connection: It is recommended to connect the safety storage cabinet to a ventilation system.	Manufacturer's instructions

3 Technical specifications

3.1 General data


Fig. 17: General view of BASIC standard safety storage cabinet

- 1 Exhaust air connection
- 2 Type plate
- 3 Storage shelf
- 4 Bottom tray

3.2 Dimensions and equipment


Fig. 18: BASIC standard line dimensions

H	Cabinet height	V	Maximum collection volume
L	Maximum volume of largest individual container	H1-H3	Adjustable heights 1 to 3, adjustable in 16mm grid
M	Standing surface load-bearing capacity	LB	Clear width
BS	Standing surface width	TS	Standing surface depth
B	Cabinet width	T	Cabinet depth
T1	Cabinet depth with open doors		

Technical specifications

Table 11: BASIC standard line dimensions

	M	XL
H (mm)	2045	2045
B (mm)	594	1194
BS (mm)	481	1081
LB (mm)	374	974
T (mm)	612	612
T1 (mm)	1100	1140
TS (mm)	444	444
H1 (mm)	526	526
H2 (mm)	526	526
H3 (mm)	538	538
L (l)	10	30
V (l)	11	33
M (kg)	75	75
Empty weight (kg)	270	450
Max. payload (kg)	360	360
Max. total weight (kg)	630	810

3.3 Pressure drop during ventilation

Technical ventilation of the safety storage cabinet results in a pressure drop at the exhaust air connection, as shown in the following diagram.


Fig. 19: Average pressure drop for all cabinet sizes

- 1 Pressure drop
- 2 Average pressure drop from all cabinet sizes
- 3 Volumetric flow rate
- 4 Q with ten-fold air exchange (see table)

Table 12: Volumetric flow rate Q and pressure drop ΔP with 10-fold air exchange

Model size	Q [m³/h]	ΔP [Pa]
M	4.3	<1
XL	9.5	<1

4 Structure and function

4.1 Construction

- Cabinet carcass and doors in multi-layer construction
- External panelling: Powder-coated sheet steel
- Wall construction: Multi-layer design
- Interior surfaces: Light grey-coated decor panels
- Safety technology elements for closure in case of fire: Brass, spring steel (1.410)

4.2 Doors

The BASIC standard line safety storage cabinets are equipped with one (model size M) or two (model size XL) wing doors, which can be opened by pulling the door handle. The doors remain open in any position.

Each door is lockable using the integrated locking cylinder. The key numbers are embossed on the locking cylinder and on the keys supplied, e.g. A007. Locks can subsequently be adjusted to the operator's requirements.

4.3 Safety technology

4.3.1 Door closure in case of fire

At an ambient temperature of approx. 50°C, open doors are closed automatically by the safety technology.

WARNING

Propping open the safety storage cabinet doors

Blocked doors mean that the safety storage cabinet loses its protection function. This can lead to death or serious injury.

- Do not block the doors.

4.3.2 Closure of venting cut-off flaps in case of fire

The safety storage cabinet is equipped with an inspection window for visual inspection of the ventilation openings for air supply and exhaust air. Above the intermediate roof, inspection cut-outs identify the position of the closing mechanism for the ventilation openings.


Fig. 20: Opened air ducts

- 1 Green inspection cut-outs in the open ventilation ducts

At an ambient temperature of 70°C, the closing mechanism is closed by the safety technology. The inspection cut-outs turn red.

4.4 Interior fittings

The BASIC standard line safety storage cabinets are equipped with several uniformly distributed and fixed storage shelves. The storage shelves are used as a standing surface for storage of closed containers.

The heights of the storage shelves can be changed. ↪ Chapter 6.8. “Changing the height of the storage shelves” on page **Fehler! Textmarke nicht definiert.**

The function of the bottom tray in the floor area of the safety storage cabinet is to collect leaking substances in the cabinet interior.

5 Transport and packaging

The safety storage cabinet is packaged for transportation and is protected against damage by transport restraints. The transport restraints should be refitted before any transportation.

WARNING

Crush hazard due to safety storage cabinet tipping over.

If the safety storage cabinet tips over when not transported with due caution, this can cause potentially fatal crushing.

- Wear personal protective equipment (PPE).
- Transport with at least two people.
- Only transport the safety storage cabinet upright and unladen.
- When driving underneath it, only pick up the safety storage cabinet in the centre from the front or from the side.
- Only drive under and transport the safety storage cabinet using suitable transport equipment.


Fig. 21: Transporting the safety storage cabinet

1 Pick up centrally from the front

2 Pick up centrally from the side

NOTE

Handling the transport restraints

Incorrect handling can damage the safety transport skids or the safety storage cabinet.

- Transport restraints and safety transport skids should only be removed at the installation location.
- Replace the safety transport skids after transport to the installation location with the enclosed adjustable feet.

6 Installation and commissioning

Install the safety storage cabinet so that the annual maintenance activities can be carried out without restriction.

6.1 Requirements for the installation location

The safety storage cabinet is approved for installation in a building.


Fig. 22: Requirements for the installation location

Consider the following in relation to the installation location:

- The surface must be able to bear the weight of the safety storage cabinet when fully loaded.
- The surface must be horizontal in order to guarantee problem-free functioning of the safety storage cabinet.
- The load-bearing capacity and stability of the surface must be assured both in normal situations and in the event of a fire.
- Do not install the safety storage cabinet near sources of heat.
- Protect the safety storage cabinet against moisture.
 - At a relative humidity of >70 % use in closed and heated buildings is permissible for a few weeks each year.
- The operating temperature must be between -5°C and +40°C.

6.2 Attaching the adjustable feet


Fig. 23: Attaching the adjustable feet

- 1 Transport equipment for the safety storage cabinet
- 2 Safety storage cabinet
- 3 Safety transport skids
- 4 Attach the safety transport skids (4x screw, size 19 mm)
- 5 Adjustable foot

Personnel:

- Technically qualified personnel
- ➔ Transport the cabinet to its place of use.
- ➔ Remove the packaging.
- ➔ Take the adjustable feet out of the cabinet.
- ➔ Raise the cabinet and loosen the screws for the safety transport skids. (A)
 - ✓ The safety transport skids can be removed (B).
- ➔ Screw the adjustable feet completely into the bottom of the cabinet from underneath (C-D).
- ➔ Position the cabinet and set it down carefully.

6.3 Aligning the safety storage cabinet

The alignment procedure described below is used for precision alignment. Remedy any major floor unevenness of more than 15 mm on site.

Long adjustable feet are fitted in the corners of the base as standard. These are used to align the safety storage cabinet.

Aligning with adjustable feet

Personnel:

- Technically qualified personnel
- ➔ Lift the cabinet using suitable lifting equipment.
- ➔ Screw the adjustable feet in or out by hand.
- ➔ Set the safety storage cabinet back down.

Aligning without adjustable feet

Personnel:

- Technically qualified personnel

Special tool:

- Suitable tool

Optionally, the safety storage cabinet can be supplied without adjustable feet.

The safety storage cabinet is supplied without alignment elements. An alignment may be necessary in individual cases.

- ➔ Raise the safety storage cabinet slightly.
- ➔ Place steel or stainless-steel spacers underneath the safety storage cabinet.
- ➔ Set the safety storage cabinet down carefully.
- ➔ Check the alignment of the safety storage cabinet and repeat the steps if necessary.
(🔗 Chapter 6.4 “Checking the alignment of the safety storage cabinet” on page 57)

6.4 Checking the alignment of the safety storage cabinet

In the event of the safety storage cabinet not being aligned properly, the open wing doors will automatically close themselves or open fully.

🔗 Chapter 6.3 “Aligning the safety storage cabinet” on page 56.


Fig. 24: Checking alignment

Correct alignment of the safety storage cabinet:

- when the doors are closed, the door gaps are of equal width.
- With two doors, the central gap and ceiling gap form an even “T”.

6.5 Fitting the plinth panel

The adjustable feet are covered and protected by the plinth panel.


Fig. 25: Fitting the plinth panel

- 1 Three-part plinth panel

Fitting the plinth panel

Personnel:

- Technically qualified personnel
- ➔ Connect the side sections of the plinth panel to the front panel (A-B).
- ➔ Slide the three-part plinth panel under the cabinet from the front and fit the spring clips (C-F) on the front adjustable feet.

6.6 Venting the safety storage cabinet

6.6.1 Exhaust air connection and air supply opening

The safety storage cabinets can be connected to an industrial exhaust air system which discharges air outside in a safe location. For this purpose, the exhaust air connection and the feed opening are located on the safety storage cabinet roof.

In normal operation, industrial ventilation of safety storage cabinets prevents the formation of a potentially explosive atmosphere inside the cabinet.

The exhaust air connection NW 110 mm has a reducer connection NW 75 mm for adaptation to a ventilation system.

The layout of the ventilation ducts in the cabinet means that ventilation takes place directly above the bottom tray and is effective on every cabinet level.


Fig. 26: Exhaust air system (connection from top)

- 1 Exhaust air connection
- 2 Air supply opening
- 3 Thermocouple

6.6.2 Connection to an exhaust air system

WARNING

Insufficient air circulation in the safety storage cabinet

A lack of or insufficient air exchange can lead to formation of an explosive atmosphere in the safety storage cabinet.

This can lead to death or serious injury.

- In a safety storage cabinet with ventilation system, a minimum hourly air exchange must take place that is at least 10 times the internal volume of the cabinet when the doors are closed.

Installation of industrial ventilation and connection to an existing exhaust air system must be carried out by a qualified company and is not a service provided by DÜPERTHAL.

Connection to the exhaust air system:

Personnel:

- Technically qualified personnel
- ➔ Connect the exhaust air line to the exhaust air connection socket and secure with a collar.
- ➔ After installing the safety storage cabinet, check the connection to the exhaust air system with smoke tubes.

The power of the exhaust air system can be determined from the technical specifications,  Chapter 3.3 “Pressure drop during ventilation” on page 50.

6.6.3 Operating the safety storage cabinet without industrial ventilation

Safety storage cabinets for passive storage can be operated without industrial ventilation.

Safety storage cabinets that are operated without industrial ventilation must be indicated using an appropriate notice.

WARNING

Safety storage cabinets without industrial ventilation

Risk of fire and explosion due to ignition of explosive mixtures in the safety storage cabinet.

This can lead to death or serious injury.

- Do not use any ignition sources in the safety storage cabinet.
- Earth safety storage cabinets using equipotential bonding.

If an explosive atmosphere is to be expected, measures must be taken in compliance with applicable national standards and regulations, e.g. TRGS 722 in Germany, and an explosion protection document must be drawn up.

WARNING

Potentially explosive area in safety storage cabinet

The consequences are death or serious injury.

- Determine the ex-zone in compliance with the applicable national standards and regulations, e.g. “TRGS 722” in Germany, and mark it clearly and permanently.
- Naked flames and smoking in potentially explosive areas are prohibited.
- Do not use any tools that cause mechanically generated sparks.
- Avoid electrostatic charges.
- Avoid equipment with surface temperatures above the ignition temperatures of the flammable liquids stored.
- Only operate electrical equipment in potentially explosive areas if it complies with the requirements of the applicable national standards and regulations, e.g. “TRGS 722” in Germany.

Clearly mark the Zone 2 potentially explosive area with the following warning sign in compliance with the applicable national standards and regulations, e.g. the European ATEX Directive “2014/34/EU”:

The size of all symbols and notices should be appropriate for the size of the safety storage cabinet.

6.7 Earthing the safety storage cabinet

Earthing the safety storage cabinet prevents ignition hazards.

The interior fittings are conductively connected to one another and to the carcass of the cabinet.

The equipotential bonding saddles for attaching an earthing cable are located on the cabinet roof on BASIC standard safety storage cabinets.

Correct earthing is stipulated in the applicable national standards and regulations, e.g. TRGS 727 in Germany.


Fig. 27: Earthing connection

Earthing connection

Personnel:

- Technically qualified personnel

Tool:

- Slotted screwdriver

- ➔ Bend the equipotential bonding saddle upwards.
- ➔ Connect the earthing cable (not included in scope of delivery).

6.8 Changing the height of the storage shelves


Fig. 28: Open safety storage cabinet with storage shelf

- | | |
|------------------|----------------------|
| 1 Earthing screw | 4 Adjustable support |
| 2 Fixing screw | 5 Storage shelf |
| 3 Support strip | |

The storage shelves in the BASIC standard safety storage cabinets are height adjustable. Additional storage shelves can be installed or removed.

The highest shelf must not be more than 1.75 m above the floor.

Personnel:

- Technically qualified personnel

Tool:

- Cordless screwdriver with crossed slot bit Z2 and AW20

Removing a storage shelf

- ➔ Remove the earthing screw including serrated washer and the fixing screws.
- ➔ Remove the storage shelf.
- ➔ Remove the shelf support.

Installing storage shelves

- ➔ Fit the shelf supports at the desired positions in the shelf support rails.
- ➔ Place the storage shelf on the shelf supports.
- ➔ Screw in the three front and rear left and front right fixing screws.
- ➔ Screw the earthing screw with serrated washer through the support strip at the rear right.

The storage shelf is conductively connected due to the attachment to the interior fittings and the equipotential bonding saddle on the cabinet roof.

7 Operation

7.1 Storage and removal


Fig. 29: Door operation

WARNING

Blocked doors

Doors that are held open by objects impair the function of the safety technology. This may result in death or serious injuries as a result of inadequate fire protection.

- Close the doors after every work process.

Personnel:

- User

Door

- ➔ Open the door by pulling the door handle as far as it will go.
 - ✓ The revolving door remains open in any position.
- ➔ Store or remove the solvents in tightly sealed containers.
- ➔ Press the door handle to close the door.
 - ✓ The safety technology holds the door in its original position.

7.2 Checking and cleaning the bottom tray

Personnel:

- User
- ➔ Store hazardous substances in such a way that a visual inspection of the bottom tray is possible on working days.
- ➔ Perform a daily visual inspection for extraneous substances.
- ➔ Collect and properly dispose of spilled liquids immediately in accordance with accident prevention regulations.

The bottom tray can be removed for easier cleaning.


Fig. 30: Removing the bottom tray

Personnel:

- Technically qualified personnel

Tool:

- Wrench size 5.5 mm

Removing the bottom tray

- ➔ Detach and set aside the lip seals (2).
- ➔ Remove the fixing screws (4).
- ➔ Remove the bottom tray (1).

Installing the bottom tray

- ➔ Adjust the bottom tray (1).
- ➔ Securely screw the fixing screws (4) through the shelf support rails (3).
- ➔ Replace the lip seals (2).

- ✓ The bottom tray is conductively connected to the interior fittings and the equipotential bonding saddle on the outside of the cabinet carcass.

NOTE

Using the bottom tray as a standing surface

Using the bottom tray as a standing surface reduces the collection volume. Escaping liquids can no longer be completely collected.

- Do not use the bottom tray as a standing surface.

8 Procedure after a fire

After a fire, the safety storage cabinet must not be opened for at least 24 hours, and only with the utmost caution and by specialist employees.

WARNING

Explosive vapour / air mixture

Risk of fire and explosion due to ignition of explosive mixtures in the safety storage cabinet.

This can lead to death or serious injury.

- Before opening the safety storage cabinet, remove all ignition sources within a 10-metre radius.
- Only open the safety storage cabinet with tools that do not cause any mechanically generated sparks.

WARNING

Damaged safety storage cabinet due to fire or extinguishing agents

A damaged safety storage cabinet can no longer provide its protection function.

This can lead to death or serious injury.

- Do not use safety storage cabinets that have been damaged by fire or extinguishing agents.

9 Maintenance

Check the safety storage cabinet for any externally visible damage or defects.

Always perform checks:

- After installation.
- Before commissioning.
- After changes
- After maintenance.

The safety storage cabinet should also be inspected periodically at the following intervals.

Interval	Maintenance work	Personnel
Daily	Bottom tray and storage levels ▪ Check in accordance with regulations governing water legislation ▪ Collect leaked liquids and hazardous substances immediately and dispose of them properly.	User

Interval	Maintenance work	Personnel
Monthly	Closing of the doors closing system ▪ Open the door and inspect the closure.	Technically qualified personnel
	Ventilation ▪ Check the effectiveness of the ventilation with a woollen thread or with a smoke tube in the cabinet in front of the exhaust air duct at the ventilation grilles. ▪ Remove contamination at the exhaust air opening.	Technically qualified personnel
	Seals ▪ Check that the sealing strips are positioned correctly on the carcass and around the doors. ▪ Check that the fabric seals on the doors are positioned correctly. ▪ If there is any visible damage, replace the seals immediately.	Technically qualified personnel
	Labelling ▪ Inspect the safety labels on the safety storage cabinet to ensure they are complete. ↪ <i>Chapter 2.7. "Hazardous areas and their labelling" on page 44</i>	Technically qualified personnel

Interval	Maintenance work	Personnel
Annually	Safety storage cabinet ▪ Test on the entire safety storage cabinet.	DÜPERTHAL service technicians

If faults occur, assist the technical customer service by providing the cabinet model, production, and key number, along with a description of the fault.

Safety equipment must be tested annually by a qualified person in compliance with BetrSichV and the maintenance interval stipulated by the manufacturer as set out in TRBS 1203.

10 Faults

Description of fault	Cause	Remedy	Personnel
Doors do not close.	Safety storage cabinet is not aligned correctly.	Install the safety storage cabinet in a level, horizontal position. ↳ <i>Chapter 6.4 "Check the alignment of the safety storage cabinet" on page 57.</i>	Technically qualified personnel
	Doors are held open by objects.	Do not wedge or hold doors open with any objects.	Technically qualified personnel
	Safety storage cabinet is not correctly filled.	Make sure that containers in the safety storage cabinet are uniformly distributed.	Technically qualified personnel
No extractor present.	Venting cut-off flaps closed, as locking mechanism has been triggered.	Replace the locking mechanism.	DÜPERTHAL service technicians
Doors do not move easily.	Moving parts, such as hinges, are dirty or corroded.	- Remove rust. - Lubricate parts. - Remove corrosive substances from the safety storage cabinet. - Notify technical customer service.	Technically qualified personnel

11 Spare parts and accessories

Only original parts from DÜPERTHAL are to be used for the safety storage cabinets.

- Storage shelf
- Perforated sheet insert
- PP insert
- Anti-slip mat made from rubber
- Plinth panel
- Venting connection socket
- Ventilators
- Exhaust air monitoring units

Additional spare parts and accessories are available on request.

12 Disposal

CAUTION

Dismantling the safety storage cabinet

Risk of injury due to improper dismantling of the safety storage cabinet.

- Ensure that the safety storage cabinet is dismantled by technically qualified personnel only.

The safety storage cabinet can be completely dismantled by specialist technical employees.

Recycle the individual material components separately.

Comply with national and local disposal regulations.

To save resources, do not place parts of the safety storage cabinet or the whole cabinet in bulky or domestic waste.

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hereby declares that the following products

Product designation: Safety storage cabinet type 30
Model designation: BASIC standard
Model size: M, XL

comply with all relevant provisions of the directives applied. This declaration refers only to the product in the condition as supplied to the customer. Parts and/or modifications subsequently fitted by the end user are not considered.

The following directives were applied:

Machinery Directive 2006/42/EG

The following harmonised standards were applied:

DIN EN ISO 1210:2011 Safety of machinery - General principles for design - Risk assessment and risk reduction

The following national standards and technical specifications have been applied:


DIN EN 14470-1:2004 Fire safety storage cabinets - Part 1: Safety storage cabinets for flammable liquids

DIN EN 16121:2017 Non-domestic storage furniture - Requirements for safety, strength, durability and stability

DIN EN 16122:2012 Domestic and non-domestic storage furniture - Test methods for the determination of strength, durability and stability

The person authorised to compile the technical documentation is: Mr Frank Backhaus, address see above.

Karlstein, 20.10.2022
(Location, date)

Signature
Franz-Josef Hagen / Managing Director

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DÜPERTHAL Sicherheitstechnik GmbH & Co. KG