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
OPERATING INSTRUCTIONS

COOL standard/dual LL
Type 90 safety storage cabinet

Retain for later use.

The operating instructions are available for download at info@dueperthal.com.
The German operating instructions are the original version. All other versions are translations of the original operating instructions.

Table of contents

1	General information	5
1.1	Notes on handling the operating instructions.....	5
1.2	Notes for reading.....	5
1.3	Type plate	6
2	Safety	7
2.1	Function of safety notices.....	7
2.2	Correct use	7
2.3	Misuse	8
2.4	The operator's responsibility	8
2.5	Requirements placed on employees	9
2.6	Stored goods.....	9
2.7	Residual risks.....	10
2.8	Safety concept for explosion protection.....	11
2.9	Hazardous areas and their labelling	13
2.10	Type testing and product labeling	14
3	Technical specifications	15
3.1	Exterior	15
3.2	Interior.....	16
3.3	Cooling.....	18
3.4	Pressure drop during ventilation.....	18
4	Structure and function	19
4.1	Construction.....	19
4.2	Doors	19
4.2.1	Door options	19
4.2.2	Locking cylinder	19
4.3	Safety technology.....	20
4.3.1	Door closure in case of fire	20
4.3.2	Closure of venting cut-off flaps in case of fire.....	20
4.4	Standing surfaces	21
4.5	Penetration.....	22
4.6	Extra load adapter	22
5	Transport and packaging	23
6	Installation and commissioning	24
6.1	Requirements for the installation location	24
6.2	Fitting the standing skids and feet	25
6.3	Installing and aligning the safety storage cabinet	26
6.4	Checking the alignment of the safety storage cabinet	27
6.5	Fitting the plinth panel	28
6.6	Connecting the safety storage cabinet to an exhaust air system	29
6.7	Earthing the safety storage cabinet	31
6.8	Connecting the safety storage cabinet to the power supply	31

7	Operation	32
7.1	Cooling unit operation	33
7.2	Opening the safety storage cabinet.....	33
7.3	Storing containers	33
8	Procedure after a fire	34
9	Maintenance.....	35
10	Faults	37
11	Shutting down the safety storage cabinet.....	39
12	Spare parts and accessories	39
13	Disposal	40
14	Declaration of Conformity	41

1 General information

1.1 Notes on handling the operating instructions

These operating instructions must always be stored with the associated product.

The operating instructions must be available to all users and technically qualified personnel.

If the product is resold, these operating instructions must be included with it.

Subsequent amendments and updates made by the manufacturer must be incorporated into the operating instructions.

1.2 Notes for reading

The following symbols designate specific types of information.

Table 1: Explanation of symbol

Symbol	Type of information
	Information for easier and more effective working
	Procedural step
	Result of a procedural step
	Link to another part of the document
	List

1.3 Type plate

The type plate is attached to the outside of the safety storage cabinet door.


Fig. 1: Type plate

Model
max. total load and volume

Maximum volume of individual containers

S/N: Serial number and year of
manufacture
R/N: Item number

Maximum load per storage shelf

The details of the model, serial number and year of manufacture, as well as the item number, can also be found on another plate attached to the inside of the door of the safety storage cabinet.

2 Safety

2.1 Function of safety notices

Safety notices warn of physical injury or material damage and provide information on how such injury and damage can be avoided.

The following signal words identify the degree of danger and the extent of the risk.

	WARNING The signal word " <i>WARNING</i> " refers to a potential hazard that could result in death or serious injuries.
	CAUTION The signal word " <i>CAUTION</i> " refers to a potential hazard that could result in slight or minor injuries.
	NOTE The signal word " <i>NOTE</i> " indicates a situation that could lead to damage to the safety storage cabinet.

2.2 Correct use

COOL standard/dual LL safety storage cabinets are type tested and classified as Type 90 in compliance with EN14470-1 with a fire resistance of 90 minutes.

The safety storage cabinet is to be used for cooled passive storage of flammable liquids in working spaces.

Any other use shall be deemed incorrect.

	Observe the safety instructions in these operating instructions to reduce health risks and avoid dangerous situations. Any use that is not correct use as defined in these operating instructions involves a risk of accidents and a lack of fire protection.
---	--

2.3 Misuse

Any use that goes beyond the specified correct use is considered to be misuse.

DÜPERTHAL accepts no liability for damage arising from misuse.

The following constitute misuse:

- Storage of food in the safety storage cabinet.
- Storage of pressurised gas cylinders in the safety storage cabinet.
- Operating the safety storage cabinet with insufficient ventilation.
- Overloading the safety storage cabinet or individual storage levels.
- Transfer, filling and laboratory work in the safety storage cabinet.
- Storage of hazardous substances together, in contravention of national regulations.
- Use of the safety storage cabinet in potentially explosive areas.
- Blocking the doors of the safety storage cabinet or using objects to hold these doors open.
- Alteration and modification of the safety storage cabinet without the prior knowledge and approval of DÜPERTHAL.
- Using low-quality spare parts.
- Failure to comply with specified maintenance intervals.
- Failure to collect and properly dispose of spilled liquids immediately in accordance with accident prevention regulations.

2.4 The operator's responsibility

The operator must ensure that:

- A risk assessment is prepared.
- An explosion protection document is prepared.
- Operating instructions are issued.
- The user has read and understood these operating instructions before using the safety storage cabinet for the first time.
- These operating instructions are available to the user at all times.
- Only trained users have access to and operate the safety storage cabinet.
- The user is trained in handling the stored goods.
- All activities and personal protective equipment are specified for the user.
- The safety storage cabinet is only used correctly.
- The safety storage cabinet is maintained in line with the specifications set out in the operating instructions.
- The safety storage cabinet is operated in accordance with local guidelines and health and safety regulations.
- All precautions are taken to avoid hazards caused by the safety storage cabinet.

2.5 Requirements placed on employees

WARNING

Unauthorised employees!

Employees who do not meet the requirements can contribute to the development of hazards when handling the safety storage cabinet.

Death or serious injuries.

- Designate employees who meet the requirements to carry out activities.

Employees need to meet specific requirements to carry out certain activities. These instructions for use assign activities to the following employees:

- Users
- Technically qualified personnel
- DÜPERTHAL service technicians

Users

Only people who meet the following requirements are approved as users:

- Instructed in the function and operation of the safety storage cabinet.
- Have basic knowledge in the general handling of chemicals and hazardous substances.
- In-depth knowledge of how to correctly store and remove flammable and liquid hazardous substances.

Technically qualified personnel

Technically qualified personnel also meet the following requirements in addition to the general requirements for users:

- Activity-specific professional training in the relevant technical field, which is recognised at the installation location.
- In-depth knowledge and skills in maintenance, repair and servicing.

DÜPERTHAL service technicians

DÜPERTHAL employees are specifically trained by DÜPERTHAL to carry out activities in relation to the safety storage cabinet.

2.6 Stored goods

Storage, handling and use of the stored goods must comply with the applicable national standards and regulations, e.g. "Technical Rule for Hazardous Substances 510" (TRGS 510) in Germany.

Furthermore, the respective operating instructions of the manufacturer of the integral refrigerators must be observed.

2.7 Residual risks

The following residual risks exist for users of the safety storage cabinet.

WARNING

Dangerous explosive atmosphere!

A lack of air exchange in the cooled storage area can lead to formation of an explosive atmosphere.

This may result in death or serious injuries due to explosion or deflagration.

- Ventilate the safety storage cabinet thoroughly before working on or in it.
- Do not perform any filling or transferring in the safety storage cabinet
- Check the stored materials regularly for leaks and damage.
- Only store tightly sealed containers.
- Remove any adherent flammable liquids on the exterior of the containers before storage.
- Mark the ex-zone.
- Create explosion protection documents.

WARNING

Spark formation in the safety storage cabinet

Spark formation in the safety storage cabinet can lead to an explosion in the cabinet.

This can lead to death or serious injuries.

- Do not use any naked flames or ignition sources in or near the safety storage cabinet.
- Do not wear statically charged clothing.
- Connect the safety storage cabinet to the earthing system.
- Use only suitable tools or equipment.
- Observe applicable national standards and regulations, e.g., "TRBS1112, Part 1" in Germany

WARNING

Electrical voltage

Touching damaged cables, plugs or sockets can cause a dangerous electric shock.

This can lead to death or serious injuries.

- Discontinue use of damaged electrical components, and have them replaced immediately by qualified personnel.
- Switch off damaged components by removing the relevant fuse or pulling out the mains plug.

2.8 Safety concept for explosion protection

As the cooling areas of the COOL standard/dual LL safety storage cabinet cannot be connected to a continuously effective exhaust air system, they should be viewed as potentially explosive areas.

The operator of the cabinet must divide them into zones – **taking into account the results of the hazard assessment and the applicable national regulations.**

Based on the applicable national regulations in Germany, e.g TRGS 722, the cooling compartment can be defined as zone 2. The following requirements must be met for this:

- Only store tightly sealed containers.
- There must be no adherent liquids on the exterior of the containers.
- Check the containers regularly for leaks and damage.
- Do not open the containers in the safety storage cabinet.
- No sampling, filling or transferring in the safety storage cabinet.


Fig. 2: Possible classification of ex-zones in compliance with the protection concept

⚠ WARNING**Potentially explosive area on the safety storage cabinet**

Risk of fire and explosion due to ignition of explosive mixtures in the safety storage cabinet.

This can lead to death or serious injuries.

- Determine the ex-zone in compliance with the applicable national standards and regulations, e.g. TRGS 722 in Germany, and mark it clearly and permanently.
- Naked flames and smoking in potentially explosive areas are prohibited.
- Do not use any tools that cause mechanically generated sparks.
- Avoid electrostatic charges.
- Only use tools, devices and equipment that are suitable for use in an ex-zone and marked accordingly.
- Do not use any ignition sources in the safety storage cabinet.
- Earth safety storage cabinets.

Because an explosive atmosphere is to be expected, measures must be taken in compliance with applicable national standards and regulations, e.g. "TRGS 722" in Germany, and an explosion protection document must be drawn up.

As the cooling areas of the COOL standard/dual LL cannot be operated with an industrial exhaust system, they must be labelled using an appropriate notice.

The Zone 2 potentially explosive area must be marked clearly with the following warning sign in compliance with the applicable national standards and regulations, e.g. the European ATEX Directive "2014/34/EU":

The size of all symbols and notices should be appropriate for the size of the safety storage cabinet.

2.9 Hazardous areas and their labelling

The following must be attached to the front of the safety storage cabinet and must be clearly visible:

- The notice “Keep door closed at all times”
- Fire resistance in minutes (e.g. “Type 90”)
- Name or trademark of the manufacturer
- Serial number and year of manufacture
- Information on the largest individual container volume that can be stored
- Information on the shelves’ maximum load capacity
- Information on the maximum total load and maximum total volume

Furthermore, the following signs must be attached to the front of the safety storage cabinet and must be clearly visible:

Table 2: Prohibited action signs

Symbol	Meaning	Standard
	P003: No naked flames; fire, open ignition source and smoking prohibited	DIN EN ISO 7010

Table 3: Warning signs

Symbol	Meaning	Standard
	W012: Warning: Electricity	DIN EN ISO 7010
	W021: Warning: Flammable materials	DIN EN ISO 7010
	D-W021: Warning: Explosive atmosphere	DIN 4844-2

Table 4: Mandatory action signs


Symbol	Meaning	Standard
	M002: Read the instructions	DIN EN ISO 7010
	Keep doors closed at all times	Manufacturer's instructions
	Connect the safety storage cabinet to a ventilation system.	Manufacturer's instructions

2.10 Type testing and product labeling

The COOL standard LL / dual LL safety storage cabinet is type-tested in accordance with DIN EN 14470-1, DIN EN 16121 and DIN EN 16122 and labelled with the GS mark ('tested safety') and the high-quality-mark of TÜV SÜD.

An explosion-proof refrigerator (CE labelled) is installed in the safety storage cabinet. In addition, the safety storage cabinet is equipped with an exhaust air monitoring unit (CE labelled).

Table. 5: Labelling on the type plate

Symbol	Meaning
	Safety storage Typ 90 (GS mark („tested safety”) und high-quality-mark)
	Exhaust air monitoring (CE declaration of conformity)
	Explosion-proof refrigerator (CE declaration of conformity)

3 Technical specifications

3.1 Exterior


Fig. 3: General view of the COOL standard/dual LL

01	Exhaust air monitor (with or without ventilator)	07	Penetration
02	Extra load adapter	08	Trigger element
03	Type plate	09	Earthing connection
04	Door handle	B	Width
05	Plinth panel	H	Height
06	Foot & standing skid	T	Depth

Table 6: External dimensions of the COOL standard/dual LL

Designation	COOL standard/dual LL
B (mm)	894
H (mm)	2050
T (mm)	791
Doors	1

*Height without extra load adapter

3.2 Interior


Fig. 4: Interior view of the COOL standard LL


Fig. 5: Interior view of the COOL dual LL

01	Cooling compartment	BS2	Shelf standing surface width
02	Glass shelf	TS1	Glass shelf standing surface depth
03	Bottom tray with storage insert	TS2	Shelf standing surface depth
H1	Adjustable height above shelf	M	Load-bearing capacity per standing surface (uniformly distributed)
H2	Adjustment range for glass shelves – adjustable in a 32 mm grid	L	Maximum volume of largest individual container
BS1	Glass shelf standing surface width	V	Max. collection volume

Table 7: Technical specifications of the COOL standard/dual LL

Designation	COOL standard LL	COOL dual LL
H1 (mm)	1373	628
H2 (mm)	896	320
BS1 (mm)	477	460
TS1 (mm)	440	294
BS2 (mm)	470	450
TS2 (mm)	282	410
M (kg)	40	40
L (l)	10	5
V (l)	11	5.5
Max. total volume (l)	110	2x 55
Number of glass shelves	4	2x 3
Number of shelves	1	2x 1
Empty weight (kg)	485	490
Max. payload of refrigerator (kg)	180	2x 90
Max. payload of extra load adapter (kg)	20*	20*
Max. total weight (kg)	685	690
Emission sound pressure level L_p (dB)	< 70	< 70

*The weight of the exhaust air monitor must be taken into account.

3.3 Cooling

Further cooling technical and electrical specifications can be found in the manufacturer's operating instructions.

Table 8: Cooling technical and electrical specifications

Designation	COOL standard LL	COOL dual LL
Cooling area volume (m ³)	0.261	2x 0.130
Refrigerant	R 600a	R 600a
Temperature range (°C)	+3 to +16	+3 to +16
Climate category	7	7
Connection voltage (V)	220-240	220-240
Frequency (Hz)	50	50
Power consumption (W)	120	2x 100

3.4 Pressure drop during ventilation

The safety storage cabinet must be operated with a continuously effective exhaust air system so that the heat produced during cooling can be discharged.

Industrial ventilation of the safety storage cabinet with a volumetric flow rate Q of 30 m³/h creates the following pressure drop Δp at the exhaust air connection socket and in the cabinet.

The pressure drop Δp becomes larger as the volumetric flow rate Q increases.

Table 9: Exhaust air values of the COOL standard/dual LL

Model	ΔP in cabinet (Pa)	Δp at exhaust air socket (Pa)
COOL standard LL	4	55
COOL dual LL	4	55

*Inner diameter of exhaust air duct: 70.5 mm

4 Structure and function

4.1 Construction

- Cabinet in multi-layer construction
- Outer casing: powder-coated sheet steel
- Refrigerator interior surfaces: White plastic
- Safety components for closure in case of fire: Brass, spring steel (1.410)
- Extra load adapter: Powder-coated sheet steel, connection to the safety storage cabinet by means of drilling screws

4.2 Doors

4.2.1 Door options

The COOL standard/dual LL has the following doors:

- External door
- Refrigerator door

4.2.2 Locking cylinder

The external door can be locked using the integrated locking cylinder. The key numbers are embossed on the locking cylinder and on the keys supplied, e.g. A007. Locks can subsequently be adjusted to the operator's requirements.

4.3 Safety technology

4.3.1 Door closure in case of fire

At an ambient room or cabinet interior temperature of approx. 50 °C, open doors are automatically closed by the safety technology.

4.3.2 Closure of venting cut-off flaps in case of fire


Fig. 6: Open ventilation openings for air supply and exhaust air (symbolic image)

01 Inspection cut-outs of ventilation openings for air supply and exhaust air

At a temperature of 70 °C, the ventilation openings for air supply and exhaust air are sealed by the venting cut-off flaps installed in the cabinet.

The safety storage cabinet is also equipped with an inspection window for visual inspection of the ventilation openings for air supply and exhaust air. Above the intermediate roof, inspection cut-outs identify the position of the locking mechanism. Under normal circumstances, these are green.

When the ventilation openings for air supply and exhaust air are closed by the locking mechanism, the inspection cut-outs change to red.

4.4 Standing surfaces


Fig. 7: Standing surfaces of the COOL standard/dual LL

- 01 Glass shelf
- 02 Bottom tray with storage insert

The COOL standard LL is equipped with one and the COOL dual LL with two integral refrigerators.

The refrigerator in the COOL standard LL has four glass shelves and the two refrigerators in the COOL dual LL each have three glass shelves, which serve as standing surface. The height of the glass shelves in the refrigerators can be adjusted in a 32 mm grid. For a secure hold, they must snap into place in the rails.

Each refrigerator is also equipped with a bottom tray with storage insert. The bottom tray catches any substances escaping in the cabinet interior. The integrated storage insert serves as additional standing surface.

The bottom tray and the storage insert can be removed from the respective refrigerator for easier cleaning.

4.5 Penetration

The roof of the COOL standard/dual LL safety storage cabinet is fitted with a tested penetration that does not impair the fire resistance. The connection cables of the refrigerators and additional cable sensors are fed into the safety storage cabinet through the prepared penetration.

NOTE

Using the penetration

Improper use of the penetration can cause damage to the safety storage cabinet or the electrical installations.

- Carry out a separate risk assessment of the overall setup.

NOTE

Extending the penetration

Attaching additional penetrations or drilling additional holes in existing penetrations is not deemed to be correct use and leads to a lack of/restricted fire protection.

- Do not make any holes or penetrations before consulting with DÜPERTHAL.

4.6 Extra load adapter

The COOL standard/dual LL is equipped with a tested extra load adapter. The extra load adapter is used for supporting loads on the cabinet roof. Proper use does not impair the fire resistance of the safety storage cabinet.

5 Transport and packaging

The safety storage cabinet is packaged for transportation and is protected against damage by transport restraints. The transport restraints should be refitted before any transportation.

The safety storage cabinets can be fitted with safety transport skids for transport. They must be removed at the installation location and replaced by the feet and standing skids supplied.

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⚠ WARNING

Crush hazard due to safety storage cabinet tipping over.

If the safety storage cabinet tips over when not transported with due caution, this can cause potentially fatal crushing.

- Wear personal protective equipment (PPE).
- Transport with at least two people.
- Only transport the safety storage cabinet upright and unladen.
- Only transport the safety storage cabinet using suitable transport equipment.
- Only lift the safety storage cabinet once it has been secured and there is no risk of it tipping over or slipping.

NOTE

Handling the transport restraints

Incorrect handling causes damage to the safety storage cabinet.

- The transport restraints should only be removed at the installation location
- After transporting the safety storage cabinet to the installation location, remove the safety transport skids (if fitted) or replace them with the feet and standing skids supplied.


Fig. 8: Transport: Image 1: Pick up centrally from the front / Image 2: Pick up centrally from the side

6 Installation and commissioning

Install the safety storage cabinet so that the annual maintenance activities can be carried out without restriction.

6.1 Requirements for the installation location

The safety storage cabinet is approved for installation in a building.


Fig. 9: Requirements for the installation location

Consider the following in relation to the installation location:

- The surface must be able to bear the weight of the safety storage cabinet when fully loaded.
- The surface must be horizontal to guarantee problem-free functioning of the safety storage cabinet.
- The load-bearing capacity and stability of the surface must be assured both in normal situations and in the event of a fire.
- Do not install the safety storage cabinet near sources of heat.
- Protect the safety storage cabinet against moisture.
- At a relative humidity of >70 %, use in closed and heated buildings is permitted for a few weeks each year.
- The operating temperature must be between +10 °C and +35 °C.

6.2 Fitting the standing skids and feet

Safety storage cabinets delivered to the place of use by a haulage firm are already fitted with standing skids and feet.


Fig. 10: Fitting the standing skids and feet

Personnel:

- Technically qualified personnel

Tool:

- Suitable tool for lifting the cabinet
- Open-ended wrench SW19

- ➔ Transport the cabinet to the place of use using a suitable means of transport, e.g. a hand pallet truck, and set it down carefully.
- ➔ Remove the packaging from the cabinet.
- ➔ Remove the feet and standing skids from the cabinet roof or from the cabinet (to do this, remove the transport restraints from the doors before opening them and replace them after closing the doors).
- ➔ Raise the cabinet and loosen the screws for the safety transport skids (image 1).
- ➔ Remove the safety transport skids (image 2).
- ➔ Fit the standing skids to the existing attachment points on the base of the cabinet using the feet (image 3). In doing so, ensure that the magnets are positioned at the front and outside and that the standing skids are centred.

6.3 Installing and aligning the safety storage cabinet

The alignment procedure described below is used for precision alignment. Remedy any major floor unevenness on site.

Aligning with feet and standing skids:


Fig. 11: Aligning the safety storage cabinet

Personnel:

- Technically qualified personnel

Tool and material:

- Suitable tool for lifting the cabinet
- steel spacers
 - 6x 1,5mm steel spacers are included
 - additional steel spacers can be reordered

- ➔ Position the cabinet and set it down carefully.
- ➔ Remove the transport restraints from the doors.
- ➔ Check the alignment of the safety storage cabinet with the spirit level (image 1).
- ➔ Replace the transport restraints if necessary and lift the safety storage cabinet using a suitable lifting tool, e.g. hand pallet truck.
- ➔ Position the steel spacers supplied underneath the foot and standing skid (image 1).
- ➔ Set the safety storage cabinet back down carefully.
- ✓ The foot and standing skid rest on spacers after setting the cabinet down.
- ➔ Remove the transport restraints from the doors.
- ➔ Check the alignment. ↪ Chapter 6.4 "Checking the alignment of the safety storage cabinet".

6.4 Checking the alignment of the safety storage cabinet

In the event of the safety storage cabinet not being aligned properly, the open wing doors will automatically close or open fully. ↪ Chapter 6.3 “Installing and aligning the safety storage cabinet”


Fig. 12: Checking alignment

Correct alignment of the safety storage cabinet:

- When the doors are closed, the door gaps are of equal width.
- With two doors, the central gap and ceiling gap form an even “T”.

6.5 Fitting the plinth panel


Fig. 13: Fitting the plinth panel

Personnel:

- Technically qualified personnel

- ➔ Hook the side plinth panel onto the rear of the respective standing skid (image 1).
- ➔ Attach the side plinth panel at the front to the respective magnet (image 2).
- ➔ Attach the front plinth panel to the front magnets of the standing skids (image 3).
- ✓ The plinth panels are fully assembled (image 4).

6.6 Connecting the safety storage cabinet to an exhaust air system

⚠ WARNING

Potentially explosive area

Ventilation of the safety storage cabinet **does not** prevent escape of a potentially explosive atmosphere in the cooling area of the safety storage cabinet.

This can lead to death or serious injuries.

- Do not use any ignition sources in the cooling area of the safety storage cabinet.
- Check the stored materials regularly for leaks and damage.
- Only use tightly sealed containers.
- Remove any adherent hazardous substances before storage.

NOTE

Operating the safety storage cabinet without an industrial exhaust system

If the safety storage cabinet is operated without a continuously effective exhaust air system, the heat produced during cooling cannot be discharged.

This can damage the installed electronics or the cooling unit and impair safe operation of the cabinet.

- Connect the safety storage cabinet to a continuous effective industrial exhaust air system and vent with an air exchange of at least 30 m³/h.

Installation of industrial ventilation and connection to an existing exhaust air system must be carried out by a qualified company and is not a service provided by DÜPERTHAL.


Fig. 14: Connection to an exhaust air system with ventilator (image 1) or without ventilator (image 2)

- 01 Exhaust air connection
- 02 Differential pressure sensor with connecting sleeve

Connection to an exhaust air system with ventilator (image 1) or without ventilator (image 2)

Personnel:

- Technically qualified personnel
- ➔ Attach a differential pressure sensor with connecting sleeve to the exhaust air connection socket. An adapter is required for the exhaust air system with ventilator.
- ➔ Connect the exhaust air line to the differential pressure sensor with connecting sleeve and secure it.
- ➔ After installation of the safety storage cabinet, check the effectiveness of the exhaust air system in the area of the supply air opening using a smoke tube or similar.
- ➔ Start up the differential pressure sensor and exhaust air monitoring unit as described in the operating instructions.

For the exact range of functions and operation, refer to the operating instructions supplied with the respective exhaust air monitoring system.

The power of the exhaust air system can be determined using the technical specifications. (See ➔ Chapter 3.4: “Pressure drop during ventilation”)

6.7 Earthing the safety storage cabinet


Fig. 15: Earthing connection

Earthing the safety storage cabinet prevents ignition hazards.

The external panels are conductively connected to one another. The connection for earthing is located on the extra load adapter.

Correct earthing is described in the applicable national standards and regulations, e.g. TRGS 727 in Germany.

Earthing connection:

Personnel:

- Technically qualified personnel
- ➔ Connect the earthing cable (not included in the scope of delivery) to the equipotential bonding saddle.

6.8 Connecting the safety storage cabinet to the power supply

The COOL standard/dual LL safety storage cabinet is supplied with voltage by the exhaust air monitoring unit supplied.

To connect the power supply, connect the connecting cable on the cabinet roof to the power connection on the rear of the exhaust air monitoring unit housing.

For the exact range of functions and operation, refer to the operating instructions supplied with the respective exhaust air monitoring system.

7 Operation

WARNING

Dangerous explosive atmosphere!

A lack of air exchange in the cooled storage area can lead to formation of an explosive atmosphere.

This may result in death or serious injuries due to explosion or deflagration.

- Ventilate the safety storage cabinet thoroughly before working on or in it.
- Do not perform any filling or transferring in the safety storage cabinet
- Check the stored materials regularly for leaks and damage.
- Only store tightly sealed containers.
- Remove any adherent flammable liquids on the exterior of the containers before storage.
- Mark the ex-zone.
- Create explosion protection documents.

WARNING

Spark formation in the safety storage cabinet

Spark formation in the safety storage cabinet can lead to an explosion in the cabinet.

This can lead to death or serious injuries.

- Do not use any naked flames or ignition sources near the safety storage cabinet.
- Do not wear statically charged clothing.
- Connect the safety storage cabinet to the earthing system.
- Use only suitable tools or equipment.
- Observe applicable national standards and regulations, e.g., "TRBS1112, Part 1" in Germany.

7.1 Cooling unit operation

For the exact range of functions and operation, refer to the manufacturer's operating instructions supplied with the respective cooling unit.

7.2 Opening the safety storage cabinet

WARNING

Blocked doors

Door blocked by objects cannot close by themselves in case of fire.

This may result in death or serious injuries as a result of inadequate fire protection.

- Close the doors after every work process.
- Do not block the doors using any objects.

If the safety storage cabinet is not aligned properly, the wing doors will automatically close themselves, open fully or do not close securely.

→ Chapter 6.3 "Installing and aligning the safety storage cabinet".

Opening the doors:

Personnel:

- Users

- Open the door of the safety storage cabinet by pulling the door handle.
 - ✓ The wing door remains open in any position.
- After opening the safety storage cabinet, open the cooling area by pulling the respective refrigerator door.

7.3 Storing containers

Personnel:

- Users

- Open the doors.
- Check the containers for leaks, damage and adherent liquids.
- Store or remove containers.
- Close the doors again.

8 Procedure after a fire

WARNING

Explosive vapour / air mixture

This can lead to death or serious injuries.

- Before opening the safety storage cabinet, remove all ignition sources within a 10m radius.
- Only open the safety storage cabinet with tools that do not cause any mechanically generated sparks.

WARNING

Damaged safety storage cabinet due to fire or extinguishing agents

This can lead to death or serious injuries.

- Do not use safety storage cabinets that have been damaged by fire or extinguishing agents.

After a fire, the safety storage cabinet must not be opened for at least 24 hours, and only with the utmost caution and by specialist employees.

In case of fire, the fire and extinguishing agents reduce the fire resistance of the safety storage cabinet to such an extent that it can no longer be used. It must be replaced with a new cabinet.

9 Maintenance

Safety equipment must be tested annually by a qualified person in compliance with BetrSichV and the maintenance interval stipulated by the manufacturer as set out in TRBS 1203.

Check the safety storage cabinet for any externally visible damage or defects.

Always perform checks:

- After installation.
- Before commissioning.
- After changes
- After maintenance.

Maintenance intervals for individual components may differ from the maintenance intervals for the safety storage cabinet.

Follow the manufacturer's information in the operating instructions for the individual components.

The safety storage cabinet should also be inspected at the intervals listed.

Table 10: Maintenance work daily

Interval	Maintenance work	Personnel
Daily	Bottom tray ▪ Check in accordance with regulations governing water legislation ▪ Collect and properly dispose of leaked liquids immediately.	Users
	Cooling areas ▪ Check for ice formation. ▪ Defrost if necessary.	Users
	Exhaust air monitor ▪ Check error indication on the exhaust air monitor.	Users

Table 11: Maintenance work monthly

Interval	Maintenance work	Personnel
Monthly	Closing of the doors closing system Open the door and inspect the closure.	Technically qualified personnel
	Ventilation - Check the effectiveness of the ventilation with a smoke tube or similar in the cabinet in front of the exhaust air duct at the ventilation grilles. - Remove any contamination and dust from the ventilation opening.	Technically qualified personnel
	Seals Check that the fire protection and fabric seals on the carcass and doors are seated correctly, and in the event of visible damage replace immediately -> Contact customer service	Technically qualified personnel
	Refrigerator seal Check that the refrigerator seal is in the correct position on the inside of the door and replace immediately if there is any visible damage.	Technically qualified personnel
	Labelling Inspect the safety labels on the safety storage cabinet to ensure they are complete. ↗ Chapter 2.9 “Hazardous areas and their labelling”	Technically qualified personnel

Table 12: Maintenance work yearly

Interval	Maintenance work	Personnel
Annually	Safety storage cabinet Test of the entire safety storage cabinet.	DÜPERTHAL service technicians

Your contact person – Technical Service

Phone: +49 6188 9139-166

service@dueperthal.com | service.dueperthal.com

10 Faults

Table 13: Overview of possible faults

Fault description	Cause	Remedy	Personnel
Doors do not close.	Safety storage cabinet is not aligned correctly.	Install and align the safety storage cabinet in a horizontal position. ↳ Chapter 6.4 “Checking the alignment of the safety storage cabinet”.	Technically qualified personnel
	Doors are held open by objects.	Do not wedge or hold doors open with any objects.	Technically qualified personnel
	Safety storage cabinet is not correctly filled.	Make sure that containers do not obstruct closing of the doors.	Technically qualified personnel
No extractor available.	Venting cut-off flaps are closed because locking mechanism has been triggered.	Replace the locking mechanism.	DÜPERTHAL service technicians
Doors do not move easily.	Moving parts, such as hinges, are dirty or corroded.	- Remove rust. - Lubricate parts. - Remove corrosive substances from the safety storage cabinet. - Notify technical customer service.	Technically qualified personnel
Doors open again after closing and close again after opening	Safety storage cabinet is not aligned correctly.	Install and align the safety storage cabinet in a horizontal position. ↳ Chapter 6.4 “Checking the alignment of the safety storage cabinet”	Technically qualified personnel
Doors cannot be opened	The “Door closure in case of fire” safety technology has been triggered and holds the doors closed.	Open the doors carefully and repair or replace the safety technology.	DÜPERTHAL service technicians

Fault description	Cause	Remedy	Personnel
Cooling unit does not cool	Temperature controller not switched on	Check the temperature controller	Technically qualified personnel
	Power supply interrupted	Check the power supply and mains plug	Technically qualified personnel
Temperature too high	Temperature controller set incorrectly	Adjust the temperature controller	Technically qualified personnel
	Cooling area or cold generator iced up	Defrost and if necessary set the temperature controller to a lower value	Technically qualified personnel
	Heat source too close to installation location	Change the installation location	Technically qualified personnel

If faults occur, assist the technical customer service by providing the cabinet model, production and key number, along with a description of the fault.

11 Shutting down the safety storage cabinet

If the unit is to be shut down for a lengthy period, the following points must be observed:

Personnel:

- Technically qualified personnel
- ➔ Remove the stored goods.
- ➔ De-energise the COOL standard/dual LL.
- ➔ Defrost the unit.
- ➔ Clean the cooling area and bottom tray and allow to dry thoroughly.
- ➔ Leave the doors open to prevent odours and mould formation.
- ➔ A sign must be placed on the door with the information that the safety storage cabinet has been shut down.

The safety storage cabinet must undergo a full maintenance cycle before it is recommissioned.

12 Spare parts and accessories

Only original parts from DÜPERTHAL are to be used for the safety storage cabinets.

- Bottom tray
- Storage insert
- Anti-slip mat made from rubber
- Door handle, lock
- Plinth panels
- Venting connection socket
- Monitoring systems

13 Disposal

⚠ CAUTION

Damage to cooling circuit during disposal

If the cooling circuit is damaged, the pressurised refrigerant can escape and cause serious injuries.

- Do not damage the cooling circuit when transporting an old unit for disposal.

⚠ CAUTION

Dismantling the safety storage cabinet

Risk of injury due to improper dismantling of the safety storage cabinet.


- Ensure that the safety storage cabinet is dismantled by technically qualified personnel only.
- Use appropriate tools and wear personal protective equipment when dismantling the cabinet.

The safety storage cabinet can be completely dismantled by technically qualified personnel.

Recycle the individual material components separately.

Comply with national and local disposal regulations.

To save resources, do not place parts of the safety storage cabinet or the whole cabinet in bulky or domestic waste.

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DÜPERTHAL Sicherheitstechnik GmbH & Co. KG
DÜPERTHAL Sicherheitstechnik GmbH & Co. KG | mf | 55332_Rev00_FR