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
Dokument: D00091030 Rev.: 01

CLASSIC line pure  
Safety storage cabinets

**Typ 90**

Operating instructions

DÜPERTHAL Sicherheitstechnik GmbH & Co. KG  
Frankenstraße 3  
63791 Karlstein  
Germany  
Fon: +49 6188 9139-0  
Fax: +49 6188 9139-121  
Email: [info@dueperthal.com](mailto:info@dueperthal.com)  
Internet: [www.dueperthal.com](http://www.dueperthal.com)

© 01/2024 DÜPERTHAL 71276, 2, en\_GB

## Table of contents

|          |                                                                        |           |
|----------|------------------------------------------------------------------------|-----------|
| <b>1</b> | <b>General information.....</b>                                        | <b>61</b> |
| 1.1      | Cabinet model information.....                                         | 61        |
| 1.2      | Notes for reading.....                                                 | 61        |
| 1.3      | Labelling.....                                                         | 61        |
| 1.4      | Guarantee.....                                                         | 62        |
| 1.5      | Requirements on the installation location.....                         | 62        |
| <b>2</b> | <b>Safety.....</b>                                                     | <b>63</b> |
| 2.1      | Function of safety notices.....                                        | 63        |
| 2.2      | Correct use.....                                                       | 63        |
| 2.3      | Reasonably foreseeable misuse.....                                     | 63        |
| 2.4      | Obligations of the operator.....                                       | 64        |
| 2.5      | Demands on employees.....                                              | 64        |
| 2.6      | Residual risks.....                                                    | 65        |
| 2.7      | Handling stored goods.....                                             | 65        |
| 2.7.1    | Storage of flammable liquids.....                                      | 65        |
| 2.7.2    | Presence of explosive atmosphere.....                                  | 66        |
| <b>3</b> | <b>Technical specifications.....</b>                                   | <b>67</b> |
| 3.1      | General technical specifications.....                                  | 67        |
| 3.2      | Technical specifications CLASSIC line pure equipment.....              | 67        |
| 3.3      | Pressure drop / volumetric flow rate with cabinet depth of 610 mm..... | 70        |
| 3.4      | Noise emissions.....                                                   | 70        |
| <b>4</b> | <b>Structure and function.....</b>                                     | <b>71</b> |
| 4.1      | Application.....                                                       | 71        |
| 4.2      | Construction.....                                                      | 71        |
| 4.3      | Pipe penetration.....                                                  | 71        |
| 4.4      | Earthing options.....                                                  | 72        |
| 4.5      | Ventilation.....                                                       | 72        |
| 4.5.1    | Air flow system in the cabinet interior.....                           | 73        |
| 4.5.2    | Choice of ventilator.....                                              | 73        |
| 4.5.3    | Exhaust air and air supply opening.....                                | 74        |
| 4.6      | Automatic closure of the ventilation openings.....                     | 75        |
| 4.7      | Doors.....                                                             | 75        |
| 4.7.1    | Door operation without arrest system.....                              | 75        |
| 4.7.2    | Locking cylinder.....                                                  | 75        |
| 4.7.3    | Door closure in case of fire.....                                      | 75        |
| 4.8      | CLASSIC line pure interior fittings.....                               | 76        |
| 4.8.1    | Storage shelves.....                                                   | 76        |
| 4.8.2    | Bottom tray.....                                                       | 77        |
| 4.8.3    | Perforated sheet insert.....                                           | 77        |
| <b>5</b> | <b>Transport.....</b>                                                  | <b>79</b> |
| 5.1      | Transport measures.....                                                | 79        |
| 5.2      | Measures for in-house transport.....                                   | 80        |
| <b>6</b> | <b>Install.....</b>                                                    | <b>81</b> |
| 6.1      | Attaching the adjustable feet.....                                     | 81        |
| 6.2      | Mount the plinth panel.....                                            | 82        |

|           |                                                                                          |            |
|-----------|------------------------------------------------------------------------------------------|------------|
| 6.3       | Align the safety storage cabinet.....                                                    | 82         |
| 6.3.1     | Alignment (with plinth panel).....                                                       | 82         |
| 6.3.2     | Alignment (without plinth panel).....                                                    | 83         |
| 6.3.3     | Check the alignment of the safety storage cabinet.....                                   | 84         |
| <b>7</b>  | <b>Commissioning.....</b>                                                                | <b>85</b>  |
| 7.1       | Connection to an exhaust air system.....                                                 | 85         |
| 7.2       | Earth the safety storage cabinet.....                                                    | 86         |
| <b>8</b>  | <b>Operation.....</b>                                                                    | <b>87</b>  |
| 8.1       | Open the safety storage cabinet.....                                                     | 87         |
| 8.2       | Changing the height of the storage shelves.....                                          | 88         |
| 8.3       | Fitting additional storage shelves.....                                                  | 88         |
| 8.4       | Checking the bottom tray.....                                                            | 89         |
| 8.5       | Inspect the safety labels on the safety storage cabinet to ensure they are complete..... | 89         |
| 8.6       | Operating the safety storage cabinet without industrial ventilation.....                 | 90         |
| <b>9</b>  | <b>Opening the safety storage cabinet after a fire.....</b>                              | <b>93</b>  |
| <b>10</b> | <b>Maintenance.....</b>                                                                  | <b>95</b>  |
| 10.1      | Inspect daily.....                                                                       | 95         |
| 10.2      | Inspect monthly.....                                                                     | 95         |
| 10.3      | Inspect annually.....                                                                    | 95         |
| <b>11</b> | <b>Faults.....</b>                                                                       | <b>97</b>  |
| <b>12</b> | <b>Spare parts and accessories.....</b>                                                  | <b>99</b>  |
| <b>13</b> | <b>Disposal.....</b>                                                                     | <b>101</b> |
| <b>14</b> | <b>Certificates.....</b>                                                                 | <b>103</b> |

# 1 General information

## 1.1 Cabinet model information

Enter the following information and have it ready in the event of queries:

| Cabinet model: | Production number: | Key number: | Location: | Responsible person: |
|----------------|--------------------|-------------|-----------|---------------------|
|                |                    |             |           |                     |

## 1.2 Notes for reading

The following symbols designate specific types of information.

*This graphical system identifies information for easier and more effective working.*

Procedural step

⇒ Result of a procedural step

Link to another part of the document:  Chapter 1.2 „Notes for reading“ on page 61

## 1.3 Labelling

|  |                                                                                                                                                                                                       | TRbF                                                                                                                    | DIN EN ISO 9001                                                                                                                                                                                                                                                                             |  | CE |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|----|
|                                                                                     |                                                                                                                                                                                                       | DIN EN 14470-1                                                                                                          | DIN EN 14727                                                                                                                                                                                                                                                                                |                                                                                       |    |
|                                                                                     |                                                                                                                                                                                                       | max. Liter*                                                                                                             |                                                                                                                                                                                                                                                                                             | max. kg*                                                                              |    |
| ①                                                                                   | ②                                                                                                                                                                                                     | ③                                                                                                                       | ④                                                                                                                                                                                                                                                                                           |                                                                                       |    |
| Modell<br>Model<br>Modèle<br>Modelo<br>Modello<br>Model                             | Seriennr./Baujahr<br>Serial number/year of model<br>Numéro de série/année de fabrication<br>Número de serie/año de fabricación<br>Numero di serie/anno di fabbricazione<br>Seriennummer/Jaar van bouw | Einzelgebinde<br>Individual container<br>Récipient individuel<br>Recipiente individual<br>Pacco singolo<br>Enige bundel | pro Boden gleichmäßig verteilte Last<br>loading capacity per shelf, uniformly distributed<br>charge par étagère, uniformément répartie<br>carga por balda, peso distribuido uniformemente<br>carico dei ripiani, peso distribuito simmetrico<br>per legbord gelijkmatig verdeelde belasting |                                                                                       |    |

Fig. 1: Type plate

- 1 Model
- 2 Serial number and year of manufacture
- 3 Maximum volume of individual containers
- 4 Uniformly distributed load for each storage shelf

The type plate is attached outside the safety storage cabinet's door.

## 1.4 Guarantee

Liability for defects in the delivered products is defined in the Terms and Conditions of Delivery. We accept no liability whatsoever for damage caused by failure to comply with the operating instructions and conditions of use.

Only flawless products with the properties pledged by DÜPERTHAL leave our factory. If the information in these operating instructions is followed, correct operation is guaranteed.

### *Annual review requirement*

*As safety equipment, DÜPERTHAL safety storage cabinets are subject to an annual safety review recommended by the manufacturer. This inspection is to be carried out by competent personnel qualified by the manufacturer. Without the recommended annual inspection the customer loses the right to make guarantee claims against the manufacturer.*

## 1.5 Requirements on the installation location

The safety storage cabinet is approved for installation in a building.

Installation in outdoor areas with or without roof covering is not permitted.

The surface must be able to bear the weight of the safety storage cabinet when fully loaded.

The surface must be horizontal in order to guarantee the problem-free functioning of the safety storage cabinet.

The load-bearing capacity and stability of the surface must be assured both in normal situations and in the event of a fire.


Fig. 2: Requirements on the installation location

Consider the following in relation to the installation location:

- Protect the safety storage cabinet from vehicles, shocks, etc.
- Do not install the safety storage cabinet near sources of heat.
- Protect the safety storage cabinet against moisture (service class 1). Use at >70 % relative humidity for a few weeks per year in structures closed on all sides and heated.
- The operating temperature must be between -5 °C and + 40 °C.

## 2 Safety

### 2.1 Function of safety notices

Safety notices warn against physical or material damage and provide information on how such damage can be avoided.

The following signal words identify the degree of danger and the extent of the risk.

#### **DANGER!**

The signal word refers to an immediately imminent threat which can cause death or serious injury.

#### **WARNING!**

The signal word refers to a potential threat which can cause death or serious injury.

#### **CAUTION!**

The signal word refers to an potential threat which can cause slight or mild injury.

The following safety signs are used in these instructions for use:

#### **WARNING!**

Warning of a hazard.

#### **NOTICE!**

This symbol identifies a situation that can lead to damage to the safety storage cabinet.

### 2.2 Correct use

The safety storage cabinet is to be used for the passive storage of flammable liquids in working areas.

Passive storage is defined as storage in tightly closed containers. The tightness of the container seal must be guaranteed. Nationally applicable regulations must be observed.

Any other use shall be deemed incorrect.

### 2.3 Reasonably foreseeable misuse

DÜPERTHAL accepts no liability for damage arising from misuse.

The following constitute misuse:

- Storage of living organisms in the safety storage cabinet.
- Storage of food in the safety storage cabinet.

- Storage of radioactive substances in the safety storage cabinet.
- Storage of liquid, combustible and non-combustible hazardous substances in the safety storage cabinet.
- Performance of work in connection with repairs, maintenance and electrical charge in the safety storage cabinet.
- Performance of work on lithium batteries in the safety storage cabinet.
- Conversion measures and modifications to the safety storage cabinet without the knowledge and consent of DÜPERTHAL and the irrespective use of such safety storage cabinets.
- Any utilisation beyond the terms of correct use.

## 2.4 Obligations of the operator

The operator is obliged to comply with applicable legal regulations. These include:

- Issuing operating instructions.
- Carrying out risk assessments.
- Creating explosion protection documents.
- Complying with specified maintenance and inspection work as well as recurring checks.
- Conducting documented safety instruction for employees involved in operating and working with the safety storage cabinet at least once a year.
- Completing documented assessment of chemical and thermal reactions between the lithium batteries to be stored and the effects thereof.
- Defining and complying with the segregated storage of hazardous substances and hazardous goods.
- Prohibiting access by unauthorised persons to the lithium batteries stored in the safety storage cabinet.
- Drawing up rules and regulations concerning smoking, handling naked flames, working with naked flames and spark-generating work in and at the safety storage cabinet.

Other obligations of the operator:

- Safeguard the function of safety-relevant equipment in and at the cabinet.
- Store hazardous substances so that a visual inspection of the bottom tray for extra-neous substances is possible on working days.
- Install the safety storage cabinet so that the annual maintenance activities can be carried out without restriction.
- Immobilise, collect and correctly dispose of leaked liquids in accordance with accident prevention regulations.

## 2.5 Demands on employees

Activities make demands on employees. These instructions for use assign activities to the following employees:

- Technical specialist employees
- Laboratory and warehouse employees
- DÜPERTHAL service technicians

Only people who meet the following requirements are approved as employees:

- Of the minimum age legally required at the installation location.
- Meet applicable legal health requirements for the installation location with regard to the substances being stored.
- Have been instructed in the function and operation of the safety storage cabinet.
- In-depth knowledge of what to do when storing and removing flammable and liquid hazardous substances.

## Technical specialist employees

Technical specialist employees also fulfil other requirements in addition to the general requirements:

- Activity-related professional training in the technical area which is recognized at the installation location.
- In-depth knowledge and skills in relation to servicing and repair.

## Laboratory and warehouse employees

Laboratory and warehouse employees also fulfil other requirements in addition to the general requirements:

- Scientific professional training or equivalent professional experience which is recognised at the installation location.
- Training and instruction in the activities to be carried out.

## DÜPERTHAL service technicians

DÜPERTHAL employees are specifically trained by DÜPERTHAL to carry out their activities.

## 2.6 Residual risks

Despite design-related measures, residual risks can arise even if the safety storage cabinet is used correctly.

In order to reduce health risks and to avoid dangerous situations, observe the safety instructions listed here and the safety information contained in later chapters of these instructions.

### WARNING!

Employees who do not meet these requirements

The consequences are death or serious injury.

- Carry out activities by named employees who meet the requirements.

## 2.7 Handling stored goods

### 2.7.1 Storage of flammable liquids

Valid national regulations, such as TRGS 510 in Germany, regarding the storage of flammable liquids, must be observed.

Information regarding storage volumes, the common storage of different hazardous substances and GHS classifications is listed and defined in these regulations.

### **2.7.2 Presence of explosive atmosphere**

The applicable national technical regulations stipulate safety measures regarding the prevention of explosive atmospheres and ignition hazards, as well as ventilation requirements.

## 3 Technical specifications

### 3.1 General technical specifications


Fig. 3: General diagram of safety storage cabinet Type 90

- |                          |                    |
|--------------------------|--------------------|
| 1 Type class             | 4 Standing surface |
| 2 Exhaust air connection | 5 Bottom tray      |
| 3 Labelling              | 6 Plinth panel     |

### 3.2 Technical specifications CLASSIC line pure equipment

| Model size | Door wings | Number of storage shelves* | Number of bottom trays |
|------------|------------|----------------------------|------------------------|
| M          | 1          | 4                          | 1                      |
| XL         | 2          |                            |                        |

\* Other quantities available on request.





Fig. 4: Technical specifications CLASSIC line pure

- |   |                                                                       |    |                                                        |
|---|-----------------------------------------------------------------------|----|--------------------------------------------------------|
| 1 | Height 1                                                              | 6  | Maximum volume of the largest individual container [L] |
| 2 | Height 2                                                              | 7  | Inside depth of the standing surface [T]               |
| 3 | Height 3                                                              | 8  | Inside width of the standing surface [B]               |
| 4 | Height 4                                                              | 9  | Clear width [LB]                                       |
| 5 | Load-bearing capacity of standing surface (uniformly distributed) [M] | 10 | Maximum collection volume [V]                          |

|                   | M   | XL   |
|-------------------|-----|------|
| B (mm)            | 479 | 1079 |
| LB (mm)           | 372 | 972  |
| T (mm)            | 435 | 435  |
| M (kg)            | 75  | 75   |
| L (l)             | 10  | 30   |
| V (l)             | 11  | 33   |
| Empty weight (kg) | 315 | 530  |
| Max. payload (kg) | 360 | 360  |

| Maximum container height (mm) |                                                                                   |          |       |
|-------------------------------|-----------------------------------------------------------------------------------|----------|-------|
|                               |                                                                                   |          | M, XL |
| 2 storage shelves             |  | Height 1 | 782   |
|                               |                                                                                   | Height 4 | 790   |
| 3 storage shelves             |  | Height 1 | 526   |
|                               |                                                                                   | Height 2 | 526   |
|                               |                                                                                   | Height 4 | 518   |
| 4 storage shelves             |  | Height 1 | 382   |
|                               |                                                                                   | Height 2 | 398   |
|                               |                                                                                   | Height 3 | 398   |
|                               |                                                                                   | Height 4 | 390   |

### 3.3 Pressure drop / volumetric flow rate with cabinet depth of 610 mm


Fig. 5: Cabinet depth 610 mm

- 1 Pressure drop  $\Delta P$  (PA)
- 2 Volumetric flow rate  $Q$  (M<sup>3</sup>/H)
- 3  $Q$  with ten-fold air exchange (see table).
- 4 Averaged pressure drop from all cabinet sizes.

Tab. 1:  $Q$  with ten-fold air exchange

| Ten-fold air exchange | Volumetric flow rate |
|-----------------------|----------------------|
| S                     | 2.8                  |
| M                     | 4.3                  |
| L                     | 6.9                  |
| XL                    | 9.5                  |
| XS                    | 6.1                  |

### 3.4 Noise emissions

The noise emitted by the door closing is < 70 dB(A).

## 4 Structure and function

### 4.1 Application

The CLASSIC line pure safety storage cabinet is tested and classified as Type 90 according to EN 14470-1 with a fire resistance of 90 minutes.

The CLASSIC line pure is suitable for the passive storage of flammable and liquid hazardous substances, in accordance with applicable national technical regulations, such as TRGS 510 in Germany.

### 4.2 Construction

- Cabinet carcass and doors in multi-layer construction
- External panelling: Powder-coated sheet steel
- Wall construction: Sandwich construction
- Interior surfaces: Stainless steel
- Safety technology elements for closure in case of fire: Brass, spring steel (1.410)

### 4.3 Pipe penetration

Tested penetrations are optionally available for the safety storage cabinet with Type 90 classification. The penetrations are attached to the safety storage cabinets from the outside and can be equipped with holes of up to Ø 28 mm.

According to type testing, these have no negative effect whatsoever on the cabinet's fire resistance if used correctly. Unused, open holes in the pipe penetrations must be sealed. Contact DÜPERTHAL.

#### NOTICE!

##### Later attachment of pipe penetrations

Incorrect attachment in the wrong location can cause damage to the safety storage cabinet.

- Attachment only at the locations provided.
- Holes only based on defined penetration profiles (see separate mounting instructions for each cabinet type).

#### NOTICE!

##### Routing and using pipe penetrations

The routing and use of pipes, cables and tubes are the responsibility of the operator.

Incorrect handling can cause damage to the safety storage cabinet.

- Carry out a separate risk assessment of the overall setup.

## 4.4 Earthing options


Fig. 6: Earthing options

1 Equipotential bonding saddle on the cabinet carcass

Earthing the safety storage cabinet prevents ignition hazards.

The interior fittings are conductively connected to one another by an equipotential bonding saddle or equipotential bonding screw on the cabinet carcass.

This allows earthing in accordance with national regulations, such as TRGS 727 in Germany.

## 4.5 Ventilation

The safety storage cabinets can be connected to a technical exhaust air system which ducts outside at a danger-free location.

In normal operation, the industrial ventilation of safety storage cabinets prevents the occurrence of a potentially explosive atmosphere inside the cabinet.

### NOTICE!

In a safety storage cabinet with industrial ventilation, a minimum hourly air exchange must take place that is at least 10 times the internal volume of the cabinet when the doors are closed.

*If the safety storage cabinet is connected to a separate industrial ventilation system, the ventilator must be installed as close as possible to the venting pipe leading to the open air.*

### 4.5.1 Air flow system in the cabinet interior


Fig. 7: Ventilation ducts open

1 Green inspection cut-outs in the open ventilation ducts

The layout of the ventilation ducts in the cabinet means that ventilation takes place directly above the bottom tray and is effective on every cabinet level.

The cabinet is also equipped with an inspection window for visual checking of the connection sockets.

Above the suspended ceiling, green inspection cut-outs identify the open ventilation ducts.

### 4.5.2 Choice of ventilator

Only use ventilators that comply with the defined technical specifications.

In areas not at risk of explosion, this means:


- Inner: Zone 2
- Outer: area not at risk of explosion.

#### NOTICE!

As Zone 2 equipment, a ventilator must fulfil the requirements of Directive 2014/34/EU (ATEX) equipment group II, category 3 G.

## 4.5.3 Exhaust air and air supply opening

*Fig. 8: Exhaust air system (view from above)*

- 1 Exhaust air connection socket for cabinet size S-XL
- 2 Supply air opening for cabinet size S-XL
- 3 Thermocouple for cabinet size S-XL
- 4 Nominal width (NW)
- 5 Distance: 100 mm
- 6 Distance: 150 mm

The exhaust air and air supply openings are located on the safety storage cabinet's roof.

The scope of supply includes one exhaust air connection socket with a nominal diameter of 110.

Also included in the scope of supply is an adapter reducer with a nominal diameter of 75, which enables connection to a ventilation system.

### 4.6 Automatic closure of the ventilation openings


Fig. 9: Closed ventilation ducts

1 Red inspection cut-outs on the closed ventilation ducts

At a temperature of 70°C, the open venting cut-off flaps are closed by the safety technology.

Above the suspended ceiling, red inspection cut-outs identify the closed ventilation ducts.

### 4.7 Doors

#### 4.7.1 Door operation without arrest system

The wing door can be opened by pulling on the door handle and remains open in any position.

#### 4.7.2 Locking cylinder

The door can be locked with the integrated locking cylinder. The key number is imprinted on the locking cylinder and on the keys supplied, e.g. A007. Locks be subsequently adjusted to the operator's requirements.

#### 4.7.3 Door closure in case of fire


At a temperature of approx. 50°C, the open doors are closed by the safety technology.

## 4.8 CLASSIC line pure interior fittings

CLASSIC line pure models are equipped with multiple, uniformly distributed standing surfaces for the storage of closed containers.

The standing surfaces are permanently installed storage shelves.

### 4.8.1 Storage shelves

*Fig. 10: Storage shelf*

- 1 Fixing screw / Earthing screw
- 2 Storage shelf
- 3 Adjustable support
- 4 Support strip

The CLASSIC line pure safety storage cabinet contains storage shelves that are distributed uniformly over the cabinet's interior height.

The heights of the storage shelves can be changed.

The highest shelf must not be more than 1.75 m above the floor.

### 4.8.2 Bottom tray

Fig. 11: Bottom tray

#### 1 Bottom tray

The function of the bottom tray in the floor area of the safety storage cabinet is to collect leaking substances in the cabinet interior. As a standard feature it cannot be used as additional standing surface.

### 4.8.3 Perforated sheet insert

The bottom tray in the CLASSIC line pure can be used as an additional standing surface in accordance with the regulations if an optionally available perforated sheet is inserted.

#### NOTICE!

##### Perforated sheet insert

The marked collection volume for the largest individual container that can be placed in the cabinet is only available if the perforated sheet that is also integrated is used as a standing surface.

## 5 Transport

### 5.1 Transport measures

The safety storage cabinet is packaged for transport and is protected by transport restraints.

The packaging, transport restraints and safety transport skids protect the cabinet from damage during transport.

#### WARNING!

Risk of crushing from toppling safety storage cabinet if the centre of gravity is unbalanced

The consequences are death or serious injury.

- The safety storage cabinet must be transported by at least two employees.
- Only transport the safety storage cabinet standing and unladen.
- Only drive under the safety storage cabinet using suitable transport means.
- When driving underneath it, only pick up the safety storage cabinet in the centre from the front, rear or from the side.



Fig. 12: Transporting the safety storage cabinet

- 1 Pick up centrally from the front
- 2 Pick up centrally from the side

#### NOTICE!

##### Handling the transport restraints

Incorrect handling can damage the safety transport skids and the safety storage cabinet.

- Transport restraints and safety transport skids should only be removed at the installation location.
- Replace the safety transport skids after transport to the installation location with the enclosed adjustable feet.

## 5.2 Measures for in-house transport

### NOTICE!

#### Tilting the safety storage cabinet

Tilting causes damage to the installed safety technology.

- Secure the doors of the safety storage cabinet using the transport restraints used for shipping in the door gaps.

### NOTICE!

#### Tipping the safety storage cabinet over during transport

Damage to the safety storage cabinet caused by incorrect handling.

- The safety storage cabinet can weigh up to 900 kg.
- Only pick up the safety storage cabinet from the side or back wall.
- Only pick up the safety storage cabinet using special and suitable transport means or lifting equipment.
- Pick up the safety storage cabinet once it is securely lashed and is not at risk of slipping.
- Do not make contact with the adjustable feet during transport.

Please also note the safety information detailed in  *Chapter 5.1 „Transport measures“ on page 79.*

## 6 Install

### 6.1 Attaching the adjustable feet


Fig. 13: Attaching the adjustable feet

- 1 4 x screw SW 19 mm
- 2 Adjustable foot

#### Attaching the adjustable feet

Personal: ☐ Technical specialist employees

1. ➤ Transport the cabinet to its place of use.
2. ➤ Remove packaging.
3. ➤ Remove the adjustable feet from the cabinet roof.
4. ➤ Lift the cabinet and remove the transport skid (A-B).
5. ➤ Screw the adjustable feet completely into the bottom of the cabinet from underneath (C-D).
6. ➤ Position the cabinet and set it down carefully.

## 6.2 Mount the plinth panel


Fig. 14: Mount the plinth panel

1 Three-part plinth panel

### Mount the plinth panel

Personal: ■ Technical specialist employees

1. ➤ Connect the side pieces of the plinth panel to the front cover (A-C).
2. ➤ Push the three-part plinth panel from the front underneath the cabinet (D-F).
3. ➤ Connect the three-part plinth panel using the spring clip to the front adjustable feet.

## 6.3 Align the safety storage cabinet

*The alignment procedure described below is used for precision alignment. Remedy and gross floor unevenness of more than 15 mm on site.*

### 6.3.1 Alignment (with plinth panel)

The safety storage cabinet with plinth panel has long adjustable feet in the corners of the base area.

The adjustable feet are covered and protected by the plinth panel.

#### Alignment with plinth panel

Personal: ■ Technical specialist employees

1. ➤ Pull off the plinth panel.
2. ➤ Lift the cabinet using suitable lifting equipment.
3. ➤ Screw the adjustable feet in or out by hand.

4. ➡ Set the safety storage cabinet back down.
5. ➡ Put the plinth panel back on.

## 6.3.2 Alignment (without plinth panel)

The safety storage cabinet without plinth panel has short adjustable feet in the corners of the base area.

Optionally align the safety storage cabinet from inside or outside via the adjustable feet.

To align from inside, remove the bottom tray.

*For easier accessibility, remove the lowest storage shelf too if necessary in order to remove the bottom tray.*


Fig. 15: Alignment without base

1 Hexagon key

### Outside alignment via the adjustable feet

Personal: ☐ Technical specialist employees ☒

1. ➡ Lift the safety storage cabinet using suitable lifting equipment.
2. ➡ Screw the adjustable feet in or out by hand.
3. ➡ Set the safety storage cabinet down.

### Inside alignment via the adjustable feet

Personal: ☐ Technical specialist employees

1. ➤ Remove the sealing lip around the bottom tray.
2. ➤ Loosen the fastening screws.
3. ➤ Loosen the earthing screw.
4. ➤ Remove the bottom tray.
5. ➤ Screw the adjustable feet in or out using the hexagon key.
6. ➤ Replace the bottom tray.
7. ➤ Tighten the earthing screw.
8. ➤ Tighten the fastening screws.
9. ➤ Position the sealing lip.

### 6.3.3 Check the alignment of the safety storage cabinet


Fig. 16: Checking alignment

Correct alignment of the safety storage cabinet:

- when the doors are closed, the door gaps are of equal width.
- With two doors, the central gap and ceiling gap form an even "T".

## 7 Commissioning

### 7.1 Connection to an exhaust air system

#### NOTICE!

Keep the exhaust air connection socket and supply air opening, as well as the thermocouple on the cabinet roof, constantly accessible and open for maintenance and servicing purposes.

#### NOTICE!

In the event of unplanned use of the cabinet roof as storage area, analyse the function of the safety elements and the weight load and include these in the risk assessment by the responsible operator.

#### NOTICE!

Use of the cabinet roof is not planned in the fire safety concept of the European standard and has also not been tested!

*The installation of industrial ventilation and/or connection to an existing exhaust air system must be carried out by a ventilation company and is not included in the DÜPERTHAL range of products and services.*


Fig. 17: Connection to the exhaust air system

- 1 Exhaust air connection socket
- 2 Supply air opening
- 3 Thermocouple

## Connection to the exhaust air system:

Personal: ■ Technical specialist employees

1. ➤ Connect the exhaust air line to the exhaust air connection socket.
2. ➤ Connect the pipeline with a collar to the exhaust air connection socket.
3. ➤ After installing the safety storage cabinet, check the connection to the exhaust air system with smoke tubes.

## 7.2 Earth the safety storage cabinet


Fig. 18: Earthing connection

### Earthing connection

Personal: ■ Technical specialist employees

1. ➤ Place the screwdriver under the equipotential bonding saddle (A).
2. ➤ Bend up the equipotential bonding saddle.

## 8 Operation

### 8.1 Open the safety storage cabinet

Fig. 19: Door operation

#### Door operation with door handle

- Personal:
- Technical specialist employees
  - Laboratory and warehouse employees

- ➔ Open the safety storage cabinet by pulling on the door handle of the wing door.  
 ⇒ The wing door remains open in any position.

#### NOTICE!

Keep the doors of the safety storage cabinet closed and re-close them manually after each work process. Do not hold the doors of the safety storage cabinet open with a wooden wedge or similar objects. This blocks the safety technology!

*In the event of the safety storage cabinet not being aligned properly, the open wing doors will close themselves or open fully.*

## 8.2 Changing the height of the storage shelves


Fig. 20: Changing the height of the storage shelf

- 1 Fixing screw / Earthing screw
- 2 Storage shelf
- 3 Adjustable support
- 4 Support strip

### Changing the height of the storage shelf

Personal: ☒ Technical specialist employees

1. ➤ Loosen fixing screws / earthing screws.
2. ➤ Remove the storage shelf.
3. ➤ Move the adjustable support into the support strips.
4. ➤ Insert the storage shelf.
5. ➤ Tighten fixing screws / earthing screws.

## 8.3 Fitting additional storage shelves

Installation of additional storage shelves: ➤ Chapter 8.2 „Changing the height of the storage shelves“ on page 88

Personal: ☒ Technical specialist employees

## 8.4 Checking the bottom tray

Personal: ■ Laboratory and warehouse employees

➔ Carry out a visual inspection for extraneous substances on working days.

*Storing objects of any kind in the bottom tray reduces the marked collection volume defined for the safety storage cabinet.*

### NOTICE!

Store hazardous substances so that a visual inspection of the bottom tray for extraneous substances is possible on working days.

### NOTICE!

Remove any leaked liquids from the bottom tray straight away and dispose of them correctly.

### NOTICE!

Following dismantling for cleaning purposes, re-connect the bottom tray to the potential connection.

## 8.5 Inspect the safety labels on the safety storage cabinet to ensure they are complete

*Only clean cabinet stickers with conventional glass cleaners.  
Other cleaners may damage them.*

The following must be attached to the front of the safety storage cabinet and must be clearly visible:

- The instruction "Close the door"
- Prohibitory sign or plaque bearing the words "Unauthorised use prohibited"
- Fire resistance in minutes (e.g. Type 90)
- Name or trademark of the manufacturer
- Serial number and year of manufacture
- Information on the largest individual container volume that can be stored
- Information on the shelves' maximum load capacity

In addition, the following signs must be attached to the front of the safety storage cabinet and must be clearly visible:

Tab. 2: Prohibition sign

| Sign                                                                              | Meaning                                      | Standard             |
|-----------------------------------------------------------------------------------|----------------------------------------------|----------------------|
|  | P02: Fire, open flame and smoking prohibited | DIN EN ISO 7010:2012 |

Tab. 3: Warning sign

| Sign                                                                              | Meaning                              | Standard             |
|-----------------------------------------------------------------------------------|--------------------------------------|----------------------|
|  | W021: Warning – flammable materials! | DIN EN ISO 7010:2012 |

## 8.6 Operating the safety storage cabinet without industrial ventilation

The safety storage cabinet can be operated in accordance with national regulations without industrial ventilation. Labelling through notice!

*If an explosive atmosphere is to be expected, measures must be taken in compliance with national regulations and an explosion protection document must be drawn up.*

### NOTICE!

- Determine the hazardous area in accordance with the hazardous substance regulations under all possible operating conditions and label it clearly, visibly and permanently.
- It is prohibited to handle naked flames or to smoke in areas at risk of explosion.
- Do not use any spark-generating tools in areas at risk of explosion.
- Ensure that no electrostatic charges can form in areas at risk of explosion and that no equipment is operated that has surface temperatures above the ignition temperatures of the stored, flammable liquids.
- Only operate electrical equipment in areas at risk of explosion if they comply with the requirements of European standards, such as DIN EN VDE.

Also note the following points from national regulations:

- Safety storage cabinets without industrial ventilation should, in case of fire, protect the stored goods from impermissible heating and ignition of possibly occurring explosive mixtures.
- There must be no ignition sources inside the safety storage cabinet. If this cannot be fully prevented, depending on the risk assessment, measures must be taken to avoid ignition sources. As a minimum requirement, the measures taken must be in accordance with those designated for Zone 2 areas as defined in the applicable national technical regulations, such as TRGS 727 in Germany.
- Safety storage cabinets without industrial ventilation must be earthed with equipotential bonding.

Clearly label the Zone 2 potentially explosive area according to national accident prevention regulations with the following warning symbol:

The size of all symbols and notices should be appropriate for the size of the safety storage cabinet.

### 9 Opening the safety storage cabinet after a fire

Only open the safety storage cabinet after 24 hours have elapsed, and only with extreme caution.

#### **WARNING!**

##### **Explosive vapour-air mixture**

The consequences are death or serious injury.

- Before opening the safety storage cabinet, remove all ignition sources within a 10 metre radius.
- Open the safety storage cabinet only using spark-free tools.

*Do not use safety storage cabinets that have been damaged by fire.*

## 10 Maintenance

Check the safety storage cabinet for any externally visible damage or defects.

Always perform checks:

- After installation.
- Before commissioning.
- After changes.
- After maintenance.

The safety storage cabinet should also be inspected periodically at the following intervals.

### 10.1 Inspect daily

| Interval | Maintenance work                                                                                                                                                                                                           | Personnel                          |
|----------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|
| Daily    | Bottom tray and pull-out trays <ul style="list-style-type: none"> <li>■ Check in accordance with regulations governing water legislation</li> <li>■ Collect and properly dispose of leaked liquids immediately.</li> </ul> | Laboratory and warehouse employees |

### 10.2 Inspect monthly

| Interval | Maintenance work                                                                                                                                                                                                                              | Personnel                      |
|----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------|
| Monthly  | Closing of the doors <ul style="list-style-type: none"> <li>■ Open the door and inspect the closure.</li> </ul>                                                                                                                               | Technical specialist employees |
|          | Ventilation <ul style="list-style-type: none"> <li>■ Check the effectiveness of the ventilation with a woollen thread or with a smoke tube in the cabinet in front of the left-hand exhaust air duct at the ventilation grilles.</li> </ul>   | Technical specialist employees |
|          | Seals <ul style="list-style-type: none"> <li>■ Check the sealing strips are seated properly in the carcass frame and on the end faces of the doors.</li> <li>■ If visible damage is found, replace the sealing strips immediately.</li> </ul> | Technical specialist employees |

### 10.3 Inspect annually

| Interval | Maintenance work                                                                                                      | Personnel                     |
|----------|-----------------------------------------------------------------------------------------------------------------------|-------------------------------|
| Annually | Safety storage cabinet <ul style="list-style-type: none"> <li>■ Check of the entire safety storage cabinet</li> </ul> | DÜPERTHAL service technicians |

*If faults should occur, assist the technical customer service by providing the cabinet model, production and key number, along with a description of the fault.*

## 11 Faults

| Fault description                     | Cause                                                                  | Remedy                                                                                                                                                                                                                              | Personnel                      |
|---------------------------------------|------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------|
| Doors do not close.                   | Safety storage cabinet is not aligned correctly.                       | Install the safety storage cabinet so it is level.<br>↳ <i>Chapter 6.3.3 „Check the alignment of the safety storage cabinet“ on page 84</i>                                                                                         | Technical specialist employees |
| No extractor present.                 | Venting cut-off flaps closed, as closing mechanism has been triggered. | Replace the locking mechanism.                                                                                                                                                                                                      | DÜPERTHAL service technicians  |
| Doors do not move easily.             | Moving parts, such as hinges, are dirty or corroded.                   | <ul style="list-style-type: none"> <li>■ Remove rust.</li> <li>■ Lubricate parts.</li> <li>■ Remove aggressive substances from the safety storage cabinet.</li> <li>■ Notify the technical customer service.</li> </ul>             | Technical specialist employees |
| Doors open again after being closed.  | Safety storage cabinet is not aligned correctly.                       | <ul style="list-style-type: none"> <li>■ Unscrew the front adjustable feet slightly.</li> <li>■ Set up the safety storage cabinet so it is level.<br/>↳ <i>Chapter 6.3 „Align the safety storage cabinet“ on page 82</i></li> </ul> | Technical specialist employees |
| Doors close again after being opened. | Safety storage cabinet is not aligned correctly.                       | <ul style="list-style-type: none"> <li>■ Unscrew the rear adjustable feet slightly.</li> <li>■ Set up the safety storage cabinet so it is level.<br/>↳ <i>Chapter 6.3 „Align the safety storage cabinet“ on page 82</i></li> </ul>  | Technical specialist employees |

## 12 Spare parts and accessories

- Storage shelves
- Bottom tray
- PP insert
- Anti-slip mat made from rubber
- Door handle
- Perforated sheet insert
- Plinth panels
- Venting connection socket
- Ventilators
- Exhaust air monitoring units

### **13 Disposal**

The safety storage cabinet can be completely dismantled.

Recycle the individual fractions, e.g. metal, insulation panels, door closer, etc. separately.

Comply with national and local disposal regulations.

To save resources, do not place parts of the safety storage cabinet or the whole cabinet in bulky or domestic waste.

## 14 Certificates

### Declaration of Conformity

In accordance with Machinery Directive 2006/42/EC, Annex II A

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hereby declare that the following machine:

Machine designation: Safety cabinet for passive storage of flammable liquids  
Machine model: CLASSIC pure  
Machine size: M, XL  
Fire resistance class: type 90

Key data:

| CLASSIC pure |                                                      |
|--------------|------------------------------------------------------|
| Machine size | Approx. Dimensions<br>(Width x Height x Depth in mm) |
| M            | 594 x 2045 x 612                                     |
| XL           | 1194 x 2045 x 612                                    |

complies with all relevant requirements of Machinery Directive 2006/42/EC.

Institution responsible for review of QS system according to annex X:

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
Additionally, the machine complies with the following harmonised and national standards and specifications:


Transposed harmonised standards: DIN EN ISO 12100:2011

Transposed national standards and technical specifications:  
DIN EN 14470-1:2004  
DIN EN 16121:2017  
DIN EN 16122:2012

Authorised person for compilation of technical documents:

  
(signee)  
Frank Backhaus / CE-authorised person

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TÜV SÜD Management Service GmbH