ORIGINAL DÜPERTHAL

SUPPLY line

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DÜPERTHAL®

Typ 90

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The operating instructions are available to download in other languages at info@dueperthal.com. The German instruction manual is the original version.

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EN SUPPLY line Safety storage cabinets ... 57
1 General information ... 61
2 Safety ... 63
3 Technical specifications ... 69
4 Structure and function ... 73
5 Transport ... 81
6 Installation and commissioning ... 83
7 Operation ... 91
8 Opening the safety storage cabinet after a fire ... 95
9 Maintenance ... 97
10 Faults ... 99
11 Spare parts and accessories ... 101
12 Disposal ... 103
13 Certificates ... 105

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© 11/2019 DÜPERTHAL 71315, 1, en_GB

Table of contents

1 General information ... 61
1.1 Notes for reading ... 61
1.2 Type plate ... 61
2 Safety ... 63
2.1 Function of safety notices ... 63
2.2 Correct use ... 63
2.3 Misuse ... 63
2.4 Obligations of the operator ... 65
2.5 Demands on employees ... 65
2.6 Stored goods ... 65
2.7 Hazardous areas and their labelling ... 66
3 Technical specifications ... 69
3.1 General data ... 69
3.2 Dimensions and equipment ... 70
3.3 Pressure drop during ventilation ... 72
4 Structure and function ... 73
4.1 Construction ... 73
4.2 Earthing options ... 73
4.3 Exhaust air connection and feed opening ... 75
4.4 Doors ... 75
4.5 Safety technology ... 76
4.6 Safety concept for explosion protection ... 76
4.7 Interior fittings ... 78
4.8 Exhaust air monitor ... 79
4.9 Pipe penetration (optional) ... 79
4.10 Extra load adapter (optional) ... 80
5 Transport ... 81
6 Installation and commissioning ... 83
6.1 Requirements on the installation location ... 83
6.2 Attaching the adjustable feet ... 84
6.3 Align the safety storage cabinet ... 85
6.4 Check the alignment of the safety storage cabinet ... 85
6.5 Mount the plinth panel ... 86
6.6 Venting the safety storage cabinet ... 87
6.7 Commissioning the exhaust air monitor ... 89
6.8 Earth the safety storage cabinet ... 90
7 Operation ... 91
7.1 Open the safety storage cabinet ... 91
7.2 Changing the height of the storage shelves ... 92
7.3 Changing storage containers ... 93
7.4 Checking and cleaning the bottom tray ... 93
8 Opening the safety storage cabinet after a fire ... 95
9 Maintenance ... 97
10 Faults ... 99

11 Spare parts and accessories 101
12 Disposal 103
13 Certificates 105

1 General information

1.1 Notes for reading

The following symbols designate specific types of information.

Tab. 1: Explanation of symbol

Symbol	Type of information
	Information for easier and more effective working
	Procedural step
⇨	Result of a procedural step
	Link to another part of the document

1.2 Type plate

The type plate is attached to the outside of the safety storage cabinet door.


Fig. 1: Type plate

1 Model
2 Serial number and year of manufacture
3 Maximum volume of individual containers
4 Maximum load per storage shelf

2 Safety

2.1 Function of safety notices

Safety notices warn against physical or material damage and provide information on how such damage can be avoided.

The following signal words identify the degree of danger and the extent of the risk.

WARNING!

The signal word *'WARNING'* refers to a potential hazard which could result in death or serious injury.

CAUTION!

The signal word *'CAUTION'* refers to an potential hazard which could result in slight or minor injury.

NOTICE!

The signal word *'NOTE'* indicates a situation that can lead to damage to the safety storage cabinet.

2.2 Correct use

Observe the safety instructions in these operating instructions to reduce health risks and avoid dangerous situations.

Any use that is not correct use as define in these operating instructions involves a risk of accidents and a lack of fire protection.

SUPPLY line safety storage cabinets are type tested and classified as *Type 90* in compliance with *'EN14470-1'* with a fire resistance of 90 minutes .

The safety storage cabinet is to be used for passive storage and for active storage of flammable liquids in working spaces.

Passive storage is defined as exclusively storage in a safety storage cabinet in tightly sealed containers without work such as filling, mixing or transfer.

Active storage is defined as follows:

- Storage in open containers or containers that are not tightly closed.
- Filling or emptying non-stationary containers at filling and emptying stations.

2.3 Misuse

Any use that goes beyond the specified correct use is considered to be misuse.

DÜPERTHAL accepts no liability for damage arising from misuse.

In addition, the following safety notices must be observed:

⚠ WARNING!

Storage of living organisms in the safety storage cabinet

Living organisms can come into contact with the stored hazardous substances.

This may result in death or serious injury.

– Use the safety storage cabinet exclusively for storage of flammable liquids.

⚠ WARNING!

Storage of food in the safety storage cabinet

Food can come into contact with the stored hazardous substances.

This may result in death or serious injury.

– Use the safety storage cabinet exclusively for storage of flammable liquids.

⚠ WARNING!

Storage of hazardous substances together

Risk of uncontrolled chemical reactions.

This may result in death or serious injury.

– Only store substances and preparations that are permitted to be stored together in the safety storage cabinet.

⚠ WARNING!

Objects on the cabinet roof

In case of fire, objects on the cabinet roof can impair the function of the safety technology.

The consequences may be death or serious injury.

– Do not store any objects on the cabinet roof.

⚠ WARNING!

Spilled liquids

Inhalation of hazardous vapours from spilled liquids can lead to life-threatening injuries to the respiratory system.

Spilled liquids can lead to painful skin reactions.

– Collect and properly dispose of spilled liquids immediately in accordance with accident prevention regulations.

NOTICE!

Alteration and modifications

Do not alter or modify the safety storage cabinet.

This can lead to a lack of fire protection.

– If alteration or modification of the safety storage cabinet is required, contact DÜPERTHAL.

2.4 Obligations of the operator

The operator is obliged to comply with applicable legal regulations. This includes:

- Issuing operating instructions.
- Carrying out risk assessments.
- Creating explosion protection documents.
- Specifying activities by designated employees.

2.5 Demands on employees

WARNING!

Employees who do not meet these requirements

This may result in death or serious injury.

– Designate employees who meet the requirements to carry out activities.

These operating instructions set out the following employee activities:

- Specialist technical employees
- DÜPERTHAL service technicians

Only people who have been trained by the operator in use of the safety storage cabinet and handling of the stored goods are approved as specialist technical employees.

DÜPERTHAL service technicians

DÜPERTHAL employees are specifically trained by DÜPERTHAL to carry out their activities.

2.6 Stored goods

Storage, handling and use of the stored goods must comply with the applicable national standards and regulations, e.g. *'TRGS 510'* in Germany.

Transfer, filling or emptying of non-stationary containers must comply with the applicable national standards and regulations, e.g. TRGS 509 in Germany.

2.7 Hazardous areas and their labelling

The following must be attached to the front of the safety storage cabinet and must be clearly visible:

- The instruction "Close the door"
- Fire resistance in minutes (e.g. *'type 90'*)
- Name or trademark of the manufacturer
- Serial number and year of manufacture
- Information on the largest individual container volume that can be stored
- Information on the shelves' maximum load capacity

Furthermore, the following signs must be attached to the front of the safety storage cabinet and must be clearly visible:

Tab. 2: Prohibited action sign

Symbol	Meaning	Standard
	P003: No naked flames; fire, open ignition source and smoking prohibited	DIN EN ISO 7010:2012

Tab. 3: Warning sign

Symbol	Meaning	Standard
	W021: Warning: Flammable materials	DIN EN ISO 7010:2012
	W012: Warning: Electricity	DIN EN ISO 7010:2012
	D-W021: Warning: Explosive atmosphere	DIN4844-2:2013/ DIN EN ISO 7010:2012

Symbol	Meaning	Standard
	Warning signs: Electrostatic discharge	ANSI/ESD S8.1:2007

Tab. 4: Mandatory action signs

Symbol	Meaning	Standard
	M002: Read operator's manual	DIN EN ISO 7010:2012
	M004: Wear eye protection	DIN EN ISO 7010:2012
	M009: Wear hand protection	DIN EN ISO 7010:2012

3 Technical specifications

3.1 General data


Fig. 2: General diagram of safety storage cabinet Type 90

1 Exhaust air connection
2 Type plate
3 Exhaust air monitor
4 Exhaust air measuring sensor
5 Bottom tray including standing surface (grating)

3.2 Dimensions and equipment


Fig. 3: SUPPLY line dimensions

H Cabinet height
H1 Adjustable height
B Cabinet width
BS Standing surface width
LB Clear width
T Cabinet depth
T1 Cabinet depth with open doors
TS Standing surface depth
L Maximum volume of largest individual container
M Load-bearing capacity of standing surface (uniformly distributed)
V Maximum collection volume

Tab. 5: SUPPLY line LL - XXL dimensions

	LL (60)	LL (200)	XXL (60)	XXL (200)	XXL (60/200)
H (mm)	2045	2045	2045	2045	2045
B (mm)	894	894	1654	1654	1654
T (mm)	747	747	747	747	747
H1 (mm)	1522	1177	1522	1177	1522/1177
BS (mm)	765	765	702	702	702
TS (mm)	588	588	588	588	588
T1 (mm)	1123	1123	1483	1483	1483
LB (mm)	674	674	2x 685	2x 685	2x 685
L (l)	60	200	60	200	60 / 200
V (l)	66	220	66	220	66 / 220
M (kg)	250	250	250	250	250
Empty weight (kg)	460	480	800	840	820
Max. payload* (kg)	360	360	2x 360	2x 360	2x 360

* When using an extra load adapter, the max. payload is reduced.

3.3 Pressure drop during ventilation

Industrial ventilation of the safety storage cabinet results in a pressure drop at the exhaust air connection, as shown in the following diagram.


Fig. 4: Average pressure drop from SUPPLY line

1 Pressure drop
2 Volumetric flow rate
3 Averaged pressure drop from all cabinet sizes.

Tab. 6: Volumetric flow rate Q and pressure drop ΔP

Model size	Q [m³/h]	ΔP [Pa]
LL	30	30
XXL	60 (corresponds to 30 per cabinet unit)	60

4 Structure and function

4.1 Construction

- Cabinet carcass and doors in multi-layer construction
- External panelling: Discharge-capable, powder-coated sheet steel
- Wall construction: Multi-layer design
- Interior surfaces: Discharge-capable, powder-coated sheet steel
- Elements for transporting exhaust air: Conductive, powder-coated sheet steel and polypropylene (PP)
- Safety technology elements for closure in case of fire: Brass, spring steel (1.410)

4.2 Earthing options


Fig. 5: Earthing options

1 Equipotential bonding saddle on the cabinet carcass

Earthing the safety storage cabinet prevents ignition hazards.

The interior fittings are conductively connected to one another by an equipotential bonding saddle or equipotential bonding screw on the outside of the cabinet carcass.

Correct earthing is stipulated in the applicable national standards and regulations, e.g. *'TRGS 727'* in Germany.

⚠ WARNING!

Risk of fire due to buildup of static electricity

This may lead to death or serious injury.

- Take precautions to prevent buildup of static electricity when transferring hazardous substances.
- Connect conductive parts on the safety storage cabinet to earth.
- Check floating contact points that are created by clamp connections on non-stationary parts for secure seating before every transfer operation and after all maintenance work.
- Avoid introducing electrical charge into the safety storage cabinet by unearthed persons.

4.3 Exhaust air connection and feed opening

Passive storage: The safety storage cabinets can be connected to a technical exhaust air system which ducts outside at a danger-free location. For this purpose, the exhaust air connection and the feed openings are located on the cabinet roof of the safety storage cabinet.

In normal operation, industrial ventilation of safety storage cabinets prevents the occurrence of a potentially explosive atmosphere inside the cabinet.

Exhaust air connection NW 110 mm with reducer adapter NW 75 mm for adaptation to a ventilation system is possible.

The layout of the ventilation ducts in the cabinet means that ventilation takes place directly above the bottom tray and is effective on every cabinet level.

Active storage: Safety storage cabinets for active storage must be connected to a technical exhaust air system and monitored with the accompanying exhaust air monitor.

If the technical ventilation system fails, the exhaust air monitoring unit signals the failure with an acoustic and visual alarm.

The exhaust air monitoring unit has capability for signal transmission to an external control centre via two floating contacts (24 V/2 A).


Fig. 6: Exhaust air system (view from above)

1 Exhaust air connection
2 Air supply opening
3 Thermocouple

4.4 Doors

4.4.1 Door options

SUPPLY line safety storage cabinets have the following door options:

Door operation without arrest system

- The wing door can be opened by pulling on the door handle and remains open in any position.

Door operation with arrest system

- Pulling on the door handle must open the door up until the stop point. The door then locks automatically in the arrest system installed in-house.

Single-handed operation without arrest system

- The wing doors can be opened together by pulling on the door handle and remain open in any position.

4.4.2 Locking cylinder

The door can be locked with the integrated locking cylinder. The key numbers are embossed on the locking cylinder and on the keys supplied, e.g. A007. Locks be subsequently adjusted to the operator's requirements.

4.5 Safety technology

4.5.1 Door closure in case of fire

At an ambient temperature of approx. 50°C, the open doors are closed by the safety technology.

4.5.2 Closure of venting cut-off flaps in case of fire

The cabinet is also equipped with an inspection window for visual inspection of the supply and exhaust air shutters. Above the suspended ceiling, inspection cut-outs identify the ventilation ducts.


Fig. 7: Ventilation ducts open

1 Green inspection cut-outs in the open ventilation ducts

At a temperature of 70°C, the open venting cut-off flaps are closed by the safety technology. The inspection cut-outs are red.

4.6 Safety concept for explosion protection

When all measures of the protection concept are observed, it is assured that the lower explosion protection limit is exceeded only briefly.

When the protection concept is observed, an ex-zone 1 inside the safety storage cabinet and the exhaust air line becomes an ex-zone 2.


Fig. 8: Classification of the potentially explosive areas when the protection concept is observed

1 Exhaust air: 30 m^3/h
2 Air supply
3 Zone 2

Compliance with the protection concept:

- All inspection units that are included in the scope of delivery must be connected and made ready for operation.
- Ensure that the inspection units from the scope of delivery are in 24 h operation.
- Ensure that a volumetric flow rate of at least 30 m^3/h (2.2 m/s) is present when connection is made to the existing exhaust air connection (2x30 m^3/h for cabinet size XXL).
- After commissioning, maintenance work or change of containers, ensure that the volume in the cabinet interior is replaced at least 16 times an hour.

WARNING!

Fire hazard due to unforeseen fault in exhaust air system

This can lead to death or serious injury.

– Do not carry out any filling or decanting operations until the fault has been rectified.
– Immediately end any filling or decanting operations already begun.

4.7 Interior fittings

4.7.1 Bottom tray with grating


Fig. 9: SUPPLY line

1 Bottom tray with grating

The SUPPLY line type safety storage cabinet contains a bottom tray with grating.

The grating serves as a standing surface for the storage containers. The function of the bottom tray in the floor area of the safety storage cabinet is to collect leaking substances in the cabinet interior. The bottom trays are permanently installed in-house and subsequent modification is not possible.

! NOTICE!

Damage due to improper storage

– Do not use the bottom tray as a standing surface without the grating.

4.7.2 Storage shelves (optional)


Fig. 10: Safety storage cabinet with storage shelf

1 Earthing screw
2 Fixing screw
3 Support strip
4 Adjustable support
5 Storage shelf

The SUPPLY line safety storage cabinet contains storage shelves that are distributed uniformly over the cabinet's interior height.

The heights of the storage shelves can be changed.

The highest shelf must not be more than 1.75 m above the floor.

4.8 Exhaust air monitor

The SUPPLY line safety storage cabinet type is equipped with an exhaust air monitor. This may or may not have a ventilator.

4.9 Pipe penetration (optional)

Tested penetrations are optionally available for the safety storage cabinet with Type 90 classification. The penetrations are attached to the safety storage cabinets from the outside and can be provided with holes. Refer to the separate instructions for the pipe penetrations.

When used correctly, it has no negative effect whatsoever on fire resistance. Unused, open holes in the pipe penetrations must be sealed.

! NOTICE!

Later attachment of pipe penetrations

Incorrect attachment in the wrong location can cause damage to the safety storage cabinet.

- They should only be attached on the approved surfaces (see additional instructions for pipe penetrations).
- Holes only based on defined penetration profiles (see separate mounting instructions for each cabinet type).

! NOTICE!

Routing and using pipe penetrations

The routing and use of pipes, cables and hoses are the responsibility of the operator.

Incorrect handling can damage the safety technology in the safety storage cabinet and cause it to fail.

- Carry out a separate risk assessment of the overall setup.

4.10 Extra load adapter (optional)

A tested extra load adapter is available as an option for safety storage cabinets with Type 90 classification. The extra load adapter should be mounted on the cabinet roof. The extra load adapter must be used for supporting loads on the safety storage cabinet. When used correctly, it has no negative effect whatsoever on fire resistance. Refer to the separate instructions for the extra load adapter.

⚠ WARNING!

Objects on the cabinet roof

In case of fire, objects on the cabinet roof can impair the function of the safety technology.

This may result in death or serious injury.

- Do not store any objects on the cabinet roof.
- Using the extra load adapter

The load-bearing capacity of the cabinet is reduced by the extra load (see additional instructions for extra load adapter).

5 Transport

The safety storage cabinet is packaged for transport and is protected against damage by transport restraints. The transport restraints should be refitted before any transport.

⚠ WARNING!

Crush hazard due to safety storage cabinet tipping over

If the safety storage cabinet tips over when not transported with due caution, this can cause potentially fatal crushing.

- Wear personal protective equipment (PPE).
- Only transport with two people.
- Only transport the safety storage cabinet upright and unladen.
- When driving underneath it, only pick up the safety storage cabinet in the centre from the front or from the side.
- Only drive under the safety storage cabinet using suitable transport equipment.


Fig. 11: Transporting the safety storage cabinet

1 Pick up centrally from the front
2 Pick up centrally from the side

! NOTICE!

Handling the transport restraints

Incorrect handling can damage the safety transport skids and the safety storage cabinet.

- Transport restraints and safety transport skids should only be removed at the installation location.
- Replace the safety transport skids after transport to the installation location with the enclosed adjustable feet.

! NOTICE!

Tipping the safety storage cabinet over during transport

Damage to the safety storage cabinet caused by incorrect handling.

- Only pick up the safety storage cabinet from the side or back wall.
- Only pick up the safety storage cabinet using special and suitable transport or lifting equipment.
- Pick up the safety storage cabinet once it is securely lashed and is not at risk of slipping.
- Do not damage the adjustable feet during transport.

6 Installation and commissioning

Install the safety storage cabinet so that the annual maintenance activities can be carried out without restriction.

6.1 Requirements on the installation location

The safety storage cabinet is approved for installation in a building.


Fig. 12: Requirements on the installation location

Consider the following in relation to the installation location:

- The surface must be able to bear the weight of the safety storage cabinet when fully loaded.
- The surface must be horizontal in order to guarantee problem-free functioning of the safety storage cabinet.
- The load-bearing capacity and stability of the surface must be assured both in normal situations and in the event of a fire.
- Do not install the safety storage cabinet near sources of heat.
- Protect the safety storage cabinet against moisture.
 - At a relative humidity of **>**70 **%** use in closed and heated buildings is permissible for a few weeks each year.
- The operating temperature must be between -5°C and **+**40°C.

6.2 Attaching the adjustable feet


Fig. 13: Attaching the adjustable feet

1 Transport equipment for the safety storage cabinet
2 Safety storage cabinet
3 Safety transport skids
4 Attach the safety transport skids (4x screw, size 19 mm)
5 Adjustable foot

Personal:

- Technical specialist employees

1. Transport the cabinet to its place of use.
2. Remove packaging.
3. Remove the adjustable feet from the cabinet roof.
4. Raise the cabinet and loosen the screws for the safety transport skids. (A)
 ⇨ The safety transport skids can be removed (B).
5. Screw the adjustable feet completely into the bottom of the cabinet from underneath (C-D).
6. Position the cabinet and set it down carefully.

6.3 Align the safety storage cabinet

The alignment procedure described below is used for precision alignment. Remedy any major floor unevenness of more than 15 mm on site.

Long adjustable feet are fitted in the corners of the base as standard. These are used to align the safety storage cabinet.

Aligning with adjustable feet

Personal:

- Technical specialist employees

1. Lift the cabinet using suitable lifting equipment.
2. Screw the adjustable feet in or out by hand.
3. Set the safety storage cabinet back down.

Aligning without adjustable feet

Personal:

- Technical specialist employees

Werkzeug:

- Suitable tool

Optionally, the safety storage cabinet can be supplied without adjustable feet.

The safety storage cabinet is supplied without alignment elements. An alignment may be necessary in individual cases.

1. Raise the safety storage cabinet slightly.
2. Place steel or stainless steel spacers underneath the safety storage cabinet.
3. Set the safety storage cabinet down carefully.

6.4 Check the alignment of the safety storage cabinet

If the safety storage cabinet is not aligned properly, the open wing doors will automatically close themselves or open fully, ↳ Chapter 6.3 'Align the safety storage cabinet' on page 85.


Fig. 14: Checking alignment

Correct alignment of the safety storage cabinet:

- when the doors are closed, the door gaps are of equal width.
- With two doors, the central gap and ceiling gap form an even “T”.

6.5 Mount the plinth panel

The adjustable feet are covered and protected by the plinth panel.


Fig. 15: Mounting the plinth panel

1 Three-part plinth panel

Mounting the plinth panel

Personal:

- Technical specialist employees

1. Connect the side pieces of the plinth panel to the front cover (A-C).
2. Push the three-part plinth panel from the front underneath the cabinet (D-F).
3. Connect the three-part plinth panel using the spring clip to the front adjustable feet.

6.6 Venting the safety storage cabinet

6.6.1 Connection to an exhaust air system

WARNING!

Insufficient safety storage cabinet air circulation

A lack of or insufficient air exchange can lead to formation of an explosive atmosphere in the safety storage cabinet.

This may lead to death or serious injury.

– In a safety storage cabinet with ventilation system, an air exchange equivalent to at least 30 m^3/h (for size LL) or 60 m^3/h (for size XXL) must take place with the doors closed.

Installation of industrial ventilation and connection to an existing exhaust air system must be carried out by a qualified company and is not a service provided by DÜPERTHAL.

Connection to the exhaust air system:

Personal:

- Technical specialist employees

1. Connect the exhaust air line to the exhaust air connection socket and secure with a collar.
2. After installing the safety storage cabinet, check the connection to the exhaust air system with smoke tubes.

The power of the exhaust air system can be determined using the technical specifications, Chapter 3.3 'Pressure drop during ventilation' on page 72.

6.6.2 Operating the safety storage cabinet without industrial ventilation

Without industrial ventilation, the safety storage cabinet must only be used for passive storage, *Chapter 2.2 'Correct use' on page 63*.

WARNING!

Safety storage cabinets without industrial ventilation

Risk of fire and explosion due to ignition of explosive mixtures in the safety storage cabinet.

This may result in death or serious injury.

- Do not use any ignition sources in the safety storage cabinet.
- Earth safety storage cabinets using equipotential bonding.

If an explosive atmosphere is to be expected, measures must be taken in compliance with applicable national standards and regulations, e.g. 'TRGS 722' in Germany, and an explosion protection document must be drawn up.

If an explosive atmosphere is to be expected, measures must be taken in compliance with applicable national standards and regulations, e.g. 'TRGS 722' in Germany or 'NFPA 30' in the USA , and an explosion protection document must be drawn up.

WARNING!

Potentially explosive area on safety storage cabinet

This may result in death or serious injury.

- Determine the ex-zone in compliance with the applicable national standards and regulations, e.g. *'TRGS 722'* in Germany, and mark it clearly and permanently.
- Naked flames and smoking in potentially explosive areas are prohibited.
- Do not use any tools that cause mechanically generated sparks.
- Avoid electrostatic charges.
- Avoid equipment with surface temperatures above the ignition temperatures of the flammable liquids stored.
- Only operate electrical equipment in potentially explosive areas if it complies with the requirements of the applicable national standards and regulations, e.g. *'TRGS 722'* in Germany.

Clearly mark the Zone 2 potentially explosive area with the following warning sign in compliance with the applicable national standards and regulations, e.g. the European ATEX product directive *'2014/34/EU'*:

The size of all symbols and notices should be appropriate for the size of the safety storage cabinet.

6.7 Commissioning the exhaust air monitor

Safety storage cabinets for active storage must be connected to a technical exhaust air system and monitored with the accompanying exhaust air monitor. For more information, follow the operating instructions for the individual components.

CAUTION!

Exhaust air monitor on the cabinet roof

In case of fire, objects on the cabinet roof may impair the functioning of the safety storage cabinets.

This may lead to death or serious injury.

- Do not store any objects on the cabinet roof.
- Using the extra load adapter

The load-bearing capacity of the cabinet is reduced by the extra load (see additional instructions for extra load adapter).

6.8 Earth the safety storage cabinet


Fig. 16: Earthing connection

Earthing connection

Personal:

- Technical specialist employees

1. Bend the equipotential bonding saddle upwards.
2. Connect the earthing cable (not included in scope of delivery).

7 Operation

7.1 Open the safety storage cabinet


Fig. 17: Door operation

⚠ WARNING!

Blocked doors

Doors that are held open by objects impair the function of the safety technology.

This may result in death or serious injuries as a result of inadequate fire protection.

– Close the doors after every work process.

Personal:

- Technical specialist employees

→ Open the safety storage cabinet by pulling on the door handle of the wing door.

⇨ The wing door remains open in any position.

7.2 Changing the height of the storage shelves


Fig. 18: Open safety storage cabinet with storage shelf

1 Earthing screw
2 Fixing screw
3 Support strip
4 Adjustable support
5 Storage shelf

Changing the height of the storage shelf

Personal:

- Technical specialist employees

1. Remove the earthing screw.
2. Remove the fixing screws.
3. Remove the storage shelf.
4. Move the adjustable support into the support strips.
5. Insert the storage shelf.
6. Tighten the fastening screws.
7. Tighten the earthing screw.

7.3 Changing storage containers

Personal:
- Technical specialist employees

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1. Open the doors on the safety storage cabinet.
2. Close the storage container tightly.
3. Remove the earth terminal from the storage container.
4. Take the storage container down from the grating using suitable transport equipment.
5. Place a new storage container in the middle of the gratings using suitable transport equipment.
6. Connect the earth terminal to the storage container.
7. Close the doors of the safety storage cabinet.

7.4 Checking and cleaning the bottom tray

Personal:
- Technical specialist employees

→ Perform a daily visual inspection for extraneous substances.

! NOTICE!

Store hazardous substances so that a visual inspection of the bottom tray for extraneous substances is possible on working days.

! NOTICE!

Remove any leaked liquids from the bottom tray straight away and dispose of them correctly.

8 Opening the safety storage cabinet after a fire

After a fire, the safety storage cabinet must not be opened for at least 24 hours, and only with the utmost caution and by specialist employees.

WARNING!

Explosive vapour-air mixture

This may result in death or serious injury.

- Before opening the safety storage cabinet, remove all ignition sources within a 10 metre radius.
- Only open the safety storage cabinet with tools that do not cause any mechanically generated sparks.

WARNING!

Damaged safety storage cabinet due to fire or extinguishing agents

This may result in death or serious injury.

- Do not use safety storage cabinets that have been damaged by fire or extinguishing agents.

9 Maintenance

Check the safety storage cabinet for any externally visible damage or defects.

Always perform checks:

- After installation.
- Before commissioning.
- After changes.
- After maintenance.

The safety storage cabinet should also be inspected periodically at the following intervals.

Interval	Maintenance work	Personnel
Daily	Bottom tray ■ Check in accordance with regulations governing water legislation ■ Collect and properly dispose of leaked liquids immediately.	Laboratory and warehouse employees
	Exhaust air ■ Visual inspection of the exhaust air monitoring unit with reference to the error indication.	Laboratory and warehouse employees
	Earthing ■ Check floating contact points that are created by clamp connections on non-stationary parts such as containers for secure seating.	Laboratory and warehouse employees

Interval	Maintenance work	Personnel
Monthly	Closing of the doors ■ Open the door and inspect the closure.	Technical specialist employees
	Ventilation ■ Check the effectiveness of the ventilation with a woollen thread or with a smoke tube in the cabinet in front of the exhaust air duct at the ventilation grilles. ■ Remove contamination at the exhaust air opening.	Technical specialist employees
	Seals ■ Check the sealing strips are seated properly in the carcass frame and on the end faces of the doors. ■ If visible damage is found, replace the sealing strips immediately.	Technical specialist employees
	Labelling ■ Inspect the safety labels on the safety storage cabinet to ensure they are complete.	Technical specialist employees

Interval	Maintenance work	Personnel
Annually	Safety storage cabinet ■ Check of the entire safety storage cabinet	DÜPERTHAL service technicians

If faults occur, assist the technical customer service by providing the cabinet model, production and key number, along with a description of the fault.

Safety-relevant safety equipment must be tested annually by a qualified person in compliance with BetrSichV and the maintenance interval stipulated by the manufacturer as set out in TRBS 1203.

10 Faults

Fault description	Cause	Remedy	Personnel
Doors do not close.	Safety storage cabinet is not aligned correctly.	Install the safety storage cabinet so it is horizontal. ↳ *Chapter 6.4 'Check the alignment of the safety storage cabinet' on page 85*	Technical specialist employees
	Doors are held open by objects.	Do not wedge or hold doors open with any objects.	Technical specialist employees
	Safety storage cabinet is not correctly filled.	Make sure that containers in the safety storage cabinet are uniformly distributed.	Technical specialist employees
No extractor present.	Venting cut-off flaps closed, as closing mechanism has been triggered.	Replace the locking mechanism.	DÜPERTHAL service technicians
Doors do not move easily.	Moving parts, such as hinges, are dirty or corroded.	■ Remove rust. ■ Lubricate parts. ■ Remove corrosive substances from the safety storage cabinet. ■ Notify technical customer service.	Technical specialist employees
Doors open again after being closed.	Safety storage cabinet is not aligned correctly.	■ Unscrew the front adjustable feet slightly. ■ Align the safety storage cabinet so it is horizontal. ↳ *Chapter 6.3 'Align the safety storage cabinet' on page 85*	Technical specialist employees
Doors close again after being opened.	Safety storage cabinet is not aligned correctly.	■ Unscrew the rear adjustable feet slightly. ■ Align the safety storage cabinet so it is horizontal. ↳ *Chapter 6.3 'Align the safety storage cabinet' on page 85*	Technical specialist employees

11 Spare parts and accessories

Only original parts from DÜPERTHAL are to be used for the safety storage cabinets.

- Storage shelves
- Pull-out shelves
- Bottom tray
- PP insert
- Anti-slip mat made from rubber
- Door handle
- Perforated sheet insert
- Plinth panels
- Venting connection socket
- Ventilators
- Exhaust air monitoring units

12 Disposal

CAUTION!

Dismantling the safety storage cabinet

Risk of injury due to improper dismantling of the safety storage cabinet.

– Ensure that the safety storage cabinet is only dismantled by specialist technical employees.

The safety storage cabinet can be completely dismantled by specialist technical employees.

Recycle the individual material components separately.

Comply with national and local disposal regulations.

To save resources, do not place parts of the safety storage cabinet or the whole cabinet in bulky or domestic waste.

13 Certificates

CE Declaration of Conformity

in accordance with Machinery Directive 2006/42/EC, Annex II A

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hereby declare that the following machine

Machine designation:	Safety cabinet for active and passive storage of flammable liquids
Machine model:	SUPPLY
Machine size:	LL, XXL
Machine type:	49-20BB7u-xyz

Key:

Model size vs. machine type		LL	XXL
Width	BB	09	17
Door design	u	5, 7	7
Colour	x	0, 4, x	
Definition of fittings and features	yz	11, 12, 31, 32, yz	11, 12, 31, 32, 33, 34 und yz

If –xyz is mentioned in the machine number, the order number will define the colour, fittings and special marks.

complies with all relevant requirements of Machinery Directive 2006/42/EC.

Institution responsible for review of QS system according:

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Additionally, the machine complies with the following harmonized and national standards and specifications:

Transposed harmonized standards: DIN EN ISO 12100:2011

Transposed national standards and technical specifications:
DIN EN 14470-1:2004
DIN EN 16121:2017
DIN EN 16122:2012

Authorized person for compilation of technical documents:

(signee)
Frank Backhaus / CE-authorized person

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DÜPERTHAL Sicherheitstechnik GmbH & Co. KG
TÜV SÜD Management Service GmbH