









	Model type
	S
H (mm)	805
B (mm)	601
T (mm)	594
T1 (mm)	1170
Doors	1
Tare weight (kg)	130
Inside cabinet volume (m³)	0,15

* Standard base height (see page 4)

Outer dimensions

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: BENCH standard Type 90, safety cabinet with wing door		
	drawing no: 000106981	page: 1 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 01	date: 2020-07-14	



	Model type
	S
H1 (mm)	554
BS (mm)	480
TS (mm)	424
M (kg)	60
Max. additional load (kg)	60
L (l)	12
V (l)	13,2
Number of storage areas	1

H1: Max. setting height
 BS: Width storage area
 TS: Depth storage area
 M: Load bearing capacity per storage area (uniformly distributed)
 L: Max. volume of the individual container
 V: Max. collection volume

Inner dimensions

This drawing is our property!
 Furnish this information to other
 persons and reproduce or copy is
 not permitted without our
 written permission.

title: BENCH standard
 Type 90, safety cabinet with wing door

drawing no: 000106981


page: 2 of 4

revision: 01

date: 2020-07-14

Düperthal Sicherheitstechnik
 Frankenstrasse 3
 D-63791 Karlstein

 **DÜPERTHAL**




NW: 75
 AB: Exhaust air connection
 ZU: Supply air


Ventilation connection

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: BENCH standard Type 90, safety cabinet with wing door		
	drawing no: 000106981	page: 3 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 01	date: 2020-07-14	





External floor



Classical base



ToeKick-base



	Variable base versions		
	External floor	Classical base	ToeKick-base
HS (mm)	3	20 bis 135	60 bis 135

Base versions

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: BENCH standard Type 90, safety cabinet with wing door		
	drawing no: 000106981	page: 4 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 01	date: 2020-07-14	