





	Model type
	S
H1 (mm)	270
BS (mm)	403
TS (mm)	410
M1 (kg)	20
M2 (kg)	15
Max. additional load (kg)	35
L (l)	2,5
V (l)	2,9
Number of storage areas	2
Tare weight (kg)	145

H1: Max. setting height
 BS: Width storage area
 TS: Depth storage area
 M1: Load bearing capacity per pull-out tray (uniformly distributed)
 M2: Load bearing capacity per additional pull-out tray (uniformly distributed)
 L: Max. volume of the individual container
 V: Max. collection volume

This drawing is our property!
 Furnish this information to other
 persons and reproduce or copy is
 not permitted without our
 written permission.

Düperthal Sicherheitstechnik
 Frankenstrasse 3
 D-63791 Karlstein


title: Additional pull-out tray
 BENCH pro

drawing no.: 000106985

page: 1 of 1

revision: 01

date: 2020-07-15



DÜPERTHAL