






	Model type	
	M	XL
H (mm)	2045	2045
B (mm)	594	1194
T (mm)	612	612
T1 (mm)	1100	1140
Doors	1	2
Tare weight (kg)	275	455
Inside cabinet volume (m³)	0,43	0,95

Outer dimensions

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: CLASSIC line pure Type 90, Safety cabinet with wing door(s)		
	drawing no: 000108001	page: 1 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 00	date: 2020-10-06	




	Model type	
	M	XL
BS (mm)	481	1081
TS (mm)	444	444
M (kg)	75	75
Max. additional load (kg)	360	360
L (l)	10	30
V (l)	11	33
LB (mm)	374	974
Number of storage shelves	4	4

BS: Width storage shelf
TS: Depth storage shelf
M: Load bearing capacity per storage shelf (uniformly distributed)
L: Max. volume of the individual container
V: Max. collection volume
LB: Clear width

Max. setting heights (mm)			
		Model type	
		M	XL
Storage shelves	H1	398	398
	H2	398	398
	H3	398	398
	H4	390	390
Storage shelves adjustable in 16mm pitch			

Inner dimensions


This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: CLASSIC line pure Type 90, Safety cabinet with wing door(s)		
	drawing no.: 000108001	page: 2 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 00	date: 2020-10-06	



NW: 110, with adapter 75
 AB: Exhaust air connection
 ZU: Supply air
 TH: Thermocouple

Ventilation connection

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: CLASSIC line pure Type 90, Safety cabinet with wing door(s)	
	drawing no: 000108001	page: 3 of 4
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 00	date: 2020-10-06
		



● Loaded footprint

Geometry of the adjustable feet

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: CLASSIC line pure Type 90, Safety cabinet with wing door(s)		
	drawing no.: 000108001	page: 4 of 4	
	revision: 00	date: 2020-10-06	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein			