







	Model type			
	LS	M	L	XL
H (mm)	1950	2045	2045	2045
B (mm)	950	594	894	1194
T (mm)	420	612	612	612
T1 (mm)	860	1112	1018	1168
Doors	2	1	2	2
Tare weight (kg)	80	69	91	113
Inside cabinet volume (m³)	0,7	0,7	1,0	1,3

Outer dimensions

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: CHEM standard / CHEM tox Safety cabinet with wing door(s)		
	drawing no: D00104471	page: 1 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 02	date: 2020-02-25	




	Model type			
	LS	M	L	XL
BS (mm)	940	584	884	1184
TS (mm)	370	560	560	560
M (kg)	50	50	50	50
Max. additional load (kg)	200	200	200	200
L (l)	18	20	30	40
V (l)	20	26	34	46
Number of storage shelves	3	3	3	3

BS: Width storage shelf
TS: Depth storage shelf
M: Load bearing capacity per storage shelf (uniformly distributed)
L: Max. volume of the individual container
V: Max. collection volume

Max. container height (mm)			
		Model type	
		LS	M, L, XL
Storage shelves	H1	575	600
	H2	575	600
	H3	508	552
Storage shelves adjustable in 25mm pitch			

Inner dimensions


This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: CHEM standard / CHEM tox Safety cabinet with wing door(s)		
	drawing no.: 000104471	page: 2 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 02	date: 2020-02-25	



NW: 110, with adapter 75
 AB: Exhaust air connection

Ventilation connection

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: CHEM standard / CHEM tox Safety cabinet with wing door(s)		
	drawing no.: 000104471	page: 3 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 02	date: 2020-02-25	 DÜPERTHAL



● Loaded footprint

Geometry of the adjustable feet

This drawing is our property!
Furnish this information to other
persons and reproduce or copy is
not permitted without our
written permission.

Düperthal Sicherheitstechnik
Frankenstrasse 3
D-63791 Karlstein


title: CHEM standard / CHEM tox
Safety cabinet with wing door(s)

drawing no.: 000104471

page: 4 of 4

revision: 02

date: 2020-02-25



 **DÜPERTHAL**