





	Model type
	HSC XL
H (mm)	2045
B (mm)	1194
T (mm)	612
T1 (mm)	1208
Doors	4
Tare weight (kg)	355
Max. total weight (kg)	805
Inside cabinet volume (m³) ACID	0,37
Inside cabinet volume (m³) CLASSIC standard S	0,28
Inside cabinet volume (m³) Storage compartment	0,16

- 1: ACID
2: CLASSIC standard S
3: Storage compartment

Outer dimensions

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: HSC XL Hazardous Substance Centre with wingdoors		
	drawing no.: 000109917	page: 1 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 01	date: 2021-01-18	



	Model type		
	ACID	CLASSIC standard S	Storage compartment
BS (mm)	474	481	555
TS (mm)	404	444	561
M (kg)	40	75	10
Max. additional load (kg)	240	190	20
L (l)	17	10	-
V (l)	19	11	-
LB (mm)	553	374	544
Number of storage areas	6	2	2


Max. setting heights (mm)			
	Model type		
	ACID	CLASSIC standard S	Storage compartment
H1	310	478	286
H2	310	454*	286
H3	294	-	-
H4	310	-	-
H5	310	-	-
H6	294	-	-

* storage shelf is adjustable in 16mm pitch

- 1: ACID
 2: CLASSIC standard S
 3: Storage compartment
 BS: Width storage area
 TS: Depth storage area
 M: Load bearing capacity per storage area (uniformly distributed)
 L: Max. volume of the individual container
 V: Max. collection volume per pull-out tray or bottom tray
 LB: Clear Width

Inner dimensions


This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: HSC XL Hazardous Substance Centre with wingdoors		
	drawing no.: 000109917	page: 2 of 4	
	revision: 01	date: 2021-01-18	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein			



NW: 75
 AB 1: Exhaust air connection ACID
 AB 2: Exhaust air connection CLASSIC standard S
 AB 3: Exhaust air connection storage compartment

Ventilation connection

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: HSC XL Hazardous Substance Centre with wingdoors		
	drawing no: 000109917	page: 3 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 01	date: 2021-01-18	



● Loaded footprint

Geometry of the adjustable feet

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: HSC XL Hazardous Substance Centre with wingdoors		
	drawing no.: 000109917	page: 4 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 01	date: 2021-01-18	