








	Model type			
	M	L	XL	XXL
H (mm)	2030	2030	2030	2030
B (mm)	595	895	1195	1395
T (mm)	596	596	596	596
T1 (mm)	-	936	1236	1336
T2 (mm)	1168	1168	1168	1268
Doors	1	2	2	2
Tare weight (kg)	300	400	495	575
Inside cabinet volume (m³)	0,42	0,67	0,92	1,07

Outer dimensions

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: ECO Type G30, safety cabinet with wing door(s)		
	drawing no.: 000110285	page: 1 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 01	date: 2021-02-10	



	Model type			
	M	L	XL	XXL
H1	1884	1884	1884	1878
BS (mm)	500	800	1100	1294
TS (mm)	367	367	367	361
L (l)	220	220	220	220
LB (mm)	428	705	1005	1199
AH1 (mm)	200	200	200	200
AH1 (mm) Adjustment range	30 - 230	30 - 230	30 - 230	30 - 230
AH2 (mm)	1600	1573	1573	1573
AH2 (mm) Adjustment range	450 - 1600	450 - 1600	450 - 1600	450 - 1600
FH1 (mm)	272	272	272	266
FH1 (mm) Adjustment range	25 - 320	25 - 320	25 - 320	25 - 314
FH2 (mm)	1016	1016	1016	1016
FH2 (mm) Adjustment range	300 - 1300	300 - 1300	300 - 1300	300 - 1300
Max. additional load (kg)	200	300	400	400
Number of retaining belts	2	3	4	4

BS: Width storage area

TS: Depth storage area

L: Max. total volume of the compressed gas cylinder
(according to EN 14470-2)

LB: Clear width



AH: Fitting holder - adjustable

FH: Cylinder holder - adjustable

EK: Roll-in flap

Inner dimensions


This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: ECO Type G30, safety cabinet with wing door(s)		
	drawing no: 000110285	page: 2 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 01	date: 2021-02-10	
			



SA: Flushing gas connection

Flushing gas connection

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: ECO Type G30, safety cabinet with wing door(s)	
	drawing no.: 000110285	page: 3 of 4
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 01	date: 2021-02-10
	 DÜPERTHAL	



	Model type			
	M	L	XL	XXL
A1 (mm)	106	106	106	109
BFD1 (mm)	-	131	281	381

NW: 110, with adapter 75
 AB: Exhaust air connection
 ZU: Supply air connection
 TH: Thermocouple
 BFD: Drilling area ceiling

Ventilation connection / Drilling area ceiling

This drawing is our property!
 Furnish this information to other
 persons and reproduce or copy is
 not permitted without our
 written permission.

Düperthal Sicherheitstechnik
 Frankenstrasse 3
 D-63791 Karlstein

title: ECO
 Type G30, safety cabinet with wing door(s)

drawing no: 000110285

page: 4 of 4

revision: 01

date: 2021-02-10

DÜPERTHAL