








|                                            | Model type                  |                             |
|--------------------------------------------|-----------------------------|-----------------------------|
|                                            | HSC UTS pro LD              | HSC UTS pro XL              |
| 1 / 2                                      | Storage cabinet / UTS pro S | Storage cabinet / UTS pro M |
| H (mm)                                     | 631                         | 631                         |
| B (mm)                                     | 1100                        | 1400                        |
| T (mm)                                     | 594                         | 594                         |
| T1 (mm)                                    | 1068                        | 1068                        |
| T2 (mm)                                    | 1170                        | 1013                        |
| Doors                                      | 2                           | 3                           |
| Tare weight (kg)                           | 180                         | 225                         |
| Max. total weight (kg)                     | 245                         | 290                         |
| Inside cabinet volume (m³) Storage cabinet | 0,11                        | 0,11                        |
| Inside cabinet volume (m³) UTS             | 0,11                        | 0,17                        |

- 1: Storage cabinet  
 2: UTS pro S / UTS pro M  
 \* Standard base height (see page 4)

## Outer dimensions

|                                                                                                                                                   |                                                                                                        |                  |  |
|---------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|------------------|--|
| This drawing is our property!<br>Furnish this information to other persons and reproduce or copy is not permitted without our written permission. | title: HSC UTS<br>Hazardous Substance Center with wing doors                                           |                  |  |
|                                                                                                                                                   | drawing no.: 000111162                                                                                 | page: 1 of 4     |  |
|                                                                                                                                                   | revision: 01                                                                                           | date: 2021-03-24 |  |
| Düperthal Sicherheitstechnik<br>Frankenstrasse 3<br>D-63791 Karlstein                                                                             |  <b>DÜPERTHAL</b> |                  |  |




|                           | Model type      |           |                 |           |
|---------------------------|-----------------|-----------|-----------------|-----------|
|                           | HSC UTS pro LD  |           | HSC UTS pro XL  |           |
|                           | Storage cabinet | UTS pro S | Storage cabinet | UTS pro M |
| H1 (mm)                   | 462             | 340       | 462             | 340       |
| BS (mm)                   | 348             | 403       | 348             | 653       |
| TS (mm)                   | 502             | 418       | 502             | 418       |
| M (kg)                    | 20              | 45        | 20              | 45        |
| Max. additional load (kg) | 20              | 45        | 20              | 45        |
| L (l)                     | 10              | 10        | 10              | 16        |
| V (l)                     | 11              | 11        | 11              | 18        |
| Number of storage areas   | 1               | 1         | 1               | 1         |

- 1: Storage cabinet  
 2: UTS pro S / UTS pro M  
 H1: Max. setting height  
 BS: Width storage area  
 TS: Depth storage area  
 M: Load bearing capacity per storage area (uniformly distributed)  
 L: Max. volume of the individual container  
 V: Max. collection volume

#### Inner dimensions

|                                                                                                                                                            |                                                              |                  |                                                                                       |
|------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------|------------------|---------------------------------------------------------------------------------------|
| This drawing is our property!<br>Furnish this information to other<br>persons and reproduce or copy is<br>not permitted without our<br>written permission. | title: HSC UTS<br>Hazardous Substance Center with wing doors |                  |  |
|                                                                                                                                                            | drawing no.: 000111162                                       | page: 2 of 4     |                                                                                       |
| Düperthal Sicherheitstechnik<br>Frankenstrasse 3<br>D-63791 Karlstein                                                                                      | revision: 01                                                 | date: 2021-03-24 |  |
|                                                                                                                                                            |                                                              |                  |                                                                                       |




|         | Model type     |                |
|---------|----------------|----------------|
|         | HSC UTS pro LD | HSC UTS pro XL |
| A1 (mm) | 148            | 292            |

NW: 75  
AB 1: Exhaust air connection storage cabinet  
AB 2: Exhaust air connection UTS  
ZU: Supply air UTS


#### Ventilation connection

|                                                                                                                                                            |                                                              |                  |                  |
|------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------|------------------|------------------|
| This drawing is our property!<br>Furnish this information to other<br>persons and reproduce or copy is<br>not permitted without our<br>written permission. | title: HSC UTS<br>Hazardous Substance Center with wing doors |                  |                  |
|                                                                                                                                                            | drawing no.: 000111162                                       | page: 3 of 4     |                  |
| Düperthal Sicherheitstechnik<br>Frankenstrasse 3<br>D-63791 Karlstein                                                                                      | revision: 01                                                 | date: 2021-03-24 | <b>DÜPERTHAL</b> |
|                                                                                                                                                            |                                                              |                  |                  |





External floor



Classical base



ToeKick-base



|         | Variable base versions |                |              |
|---------|------------------------|----------------|--------------|
|         | External floor         | Classical base | ToeKick-base |
| HS (mm) | 3                      | 20 bis 135     | 60 bis 135   |

Base versions

|                                                                                                                                                            |                                                                                       |                  |                                                                                       |
|------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|------------------|---------------------------------------------------------------------------------------|
| This drawing is our property!<br>Furnish this information to other<br>persons and reproduce or copy is<br>not permitted without our<br>written permission. | title: HSC UTS<br>Hazardous Substance Center with wing doors                          |                  |                                                                                       |
|                                                                                                                                                            | drawing no: D00111162                                                                 | page: 4 of 4     |  |
|                                                                                                                                                            | revision: 01                                                                          | date: 2021-03-24 |                                                                                       |
| Düperthal Sicherheitstechnik<br>Frankenstrasse 3<br>D-63791 Karlstein                                                                                      |  |                  |                                                                                       |