

Freigegebene Desinfektionsmittel für Liebherr Labor- und Medikamentenkühler der Serie Mediline mit Kunststoff-Innenbehälter

Bezeichnung Desinfektionsmittel für Liebherr Labor- /Medikamentenkühler mit Kunststoff-Innenbehälter:	Hersteller:
Dismozon pur 1 % Lösung	Bode Chemie
Suma Quicksan	Diversey
Incidin Extra N	Ecolab Healthcare
Acrylan	Antiseptica chem.-pharm. Produkte
Buraton 10 F	Schülke und Mayr
Frankocid N	Franken Chemie
Apesin DSR 50	Tana
Nüscosept Spray	Dr. Nüsken
Nüscotan Spezial	Dr. Nüsken
Melsept SF	B. Braun Melsungen
Kohrsolin	Bode Chemie
Neoquat S	Dr. Weigert
Indicin Rapid	Ecolab Healthcare
Bacillocid Spezial	Bode Chemie
Neoform K Spray	Dr. Weigert
Apesin Desinf. Spray	Tana
Nüscosept 100	Dr. Nüsken
Antisept T	Fink Tec
Apesin AP 100 0,50% Lösung	Tana
Perform Pulver 2% Lösung	Schülke und Mayr

Approved disinfectants for Liebherr laboratory and medical coolers of MediLine with plastics inner liner

Name of approved disinfectant Liebherr laboratory / medical cooler with plastics inner liner:	Manufacturer:
Dismozon pur 1 % Lösung	Bode Chemie
Suma Quicksan	Diversey
Incidin Extra N	Ecolab Healthcare
Acrylan	Antiseptica chem.-pharm. Produkte
Buraton 10 F	Schülke und Mayr
Frankocid N	Franken Chemie
Apesin DSR 50	Tana
Nüscosept Spray	Dr. Nüsken
Nüscotan Spezial	Dr. Nüsken
Melsept SF	B. Braun Melsungen
Kohrsolin	Bode Chemie
Neoquat S	Dr. Weigert
Indicin Rapid	Ecolab Healthcare
Bacillocid Spezial	Bode Chemie
Neoform K Spray	Dr. Weigert
Apesin Desinf. Spray	Tana
Nüscosept 100	Dr. Nüsken
Antisept T	Fink Tec
Apesin AP 100 0,50% Lösung	Tana
Perform Pulver 2% Lösung	Schülke und Mayr