








	Model type	
	M	XL
H (mm)	1955	1955
B (mm)	595	1195
T (mm)	601	601
T1 (mm)	1142	1176
Doors	1	2
Tare weight (kg)	155	255
Inside cabinet volume (m³)	0,5	1,05

Outer dimensions

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: ACID pure Safety cabinet with wing door(s)	
	drawing no.: D00112157	page: 1 of 3
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 00	date: 2021-06-23
	 DÜPERTHAL	



	Model type	
	M	XL
BS (mm)	529	1129
TS (mm)	474	474
M (kg)	60	60
Max. additional load (kg)	240	240
L (l)	12	25
V (l)	14	28
LB (mm)	504	1073
Number of storage area	4	4


BS: Width storage area
TS: Depth storage area
M: Load bearing capacity per storage area (uniformly distributed)
L: Max. volume of the individual container
V: Max. collection volume
LB: Clear width

Max. setting heights (mm)	
	Model type M, XL
H1	375
H2	479
H3	479
H4	459
Storage shelves adjustable in 128mm pitch	

Inner dimensions

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: ACID pure Safety cabinet with wing door(s)		
	drawing no.: D00112157	page: 2 of 3	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 00	date: 2021-06-23	

	Model type	
	M	XL
A1 (mm)	300	600



NW: with adapter 100
 AB: Exhaust air connection
 ZU: Supply air connection

Ventilation connection

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: ACID pure Safety cabinet with wing door(s)		
	drawing no.: 000112157	page: 3 of 3	
	revision: 00	date: 2021-06-23	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein			