



Article Number: 4.271.40

Product Range

Safety equipment

Product Category

capacity: approx. 2 x 15 litres
outside dimensions: 610 x 496 x 845 mm (W x D x H)
basin (inside): 2x 298 x 406 x 140 mm(W x D x H)
Basin, lid and frame made from
galvanised sheet steel, red coated.
For cleaning of extremely contaminated components.
1. chamber: pre-cleaning,
2. chamber: final cleaning
3/4"- drain plug in basin bottom for
selective emptying.
Testing: FM-approved

Retail category	Cleaning tanks
Outer width	406 mm
Outer depth	610 mm
Outer height	222 mm
Weight (product)	30,5 kg