



Article Number: 4.272.60

Product Range

Safety equipment

Product Category

capacity: approx. 2 x 34 litres
Outside dimensions: 895 x 496 x 845 mm (WxDxH)
basin (inside): 2x438 x 406 x 222 mm (WxDxH)
Basin, lid and frame made from
galvanised sheet steel, red coated.
For cleaning of extremely contaminated components.
1. chamber: pre-cleaning,
2. chamber: final purification.
3/4"- drain plug in bottom of each chamber
for selective emptying.
Testing: FM-approved

Retail category	Cleaning tanks
Outer width	406 mm
Outer depth	895 mm
Outer height	305 mm
Weight (product)	38,5 kg