



Article Number: 9.31.845

Product Range
Product Category

Safety equipment

With two bottles:
1 x 200 ml of pH-neutral phosphate buffer solution (4.9%) and
1 x 500 ml of sterile sodium chloride solution (0.9%).
With integrated pictogram and mirror.
Incl. eye emergency plan.
Especially suited for workplaces where acids, lye and foreign bodies can get into the eye.
Dimensions wall station (H x W x D):
290 x 230 x 80 mm.
Shelf life: 3 years.
Rinsing time: approx. 7 minutes.
According to DIN EN 15154-4:2009.

Retail category	Emergency products
Retail product filter	Emergency showers
Weight (product)	1,6 kg