

Article Number: 4.270.03

**Product Range** Safety equipment

**Product Category**

Galvanised parts basket for dipping tank  
model 4.271.10 and 4.273.11

Dim. (WxDxH): approx. 406 x 330 x 102 mm.

Offers effective cleaning of multiple parts  
for every cleaning section.

Galvanised parts basket for dipping tank  
model 4.271.10 and 4.273.11

Dim. (WxDxH): approx. 406 x 330 x 102 mm.

Offers effective cleaning of multiple parts  
for every cleaning section.

**HW-category**

Cleaning tanks

**Outer width approx.**

330 mm

**Outer depth approx.**

406 mm

**Outer height approx.**

102 mm