

Article Number: 4.270.04

Product Range Safety equipment

Product Category

Galvanised strainer for dipping tank

model 4.271.10 and 4.273.11

Dim. (WxDxH): approx. 584 x 349 x 19 mm.

Gathers coarse dirt caused by sedimentation.

Offers easy cleaning during solvent exchanging.

Galvanised strainer for dipping tank

model 4.271.10 and 4.273.11

Dim. (WxDxH): approx. 584 x 349 x 19 mm.

Gathers coarse dirt caused by sedimentation.

Offers easy cleaning during solvent exchanging.

HW-category

Cleaning tanks

Outer width approx.

349 mm

Outer depth approx.

584 mm

Outer height approx.

19 mm