

Article Number: 4.270.05

Product Range Safety equipment

Product Category

Galvanised parts basket for dipping tank
model 4.271.20 and 4.273.22

Dim. (WxDxH): approx. 711 x 330 x 102 mm.

Offers effective cleaning of multiple parts
every cleaning section.

Galvanised parts basket for dipping tank
model 4.271.20 and 4.273.22

Dim. (WxDxH): approx. 711 x 330 x 102 mm.

Offers effective cleaning of multiple parts
every cleaning section.

HW-category

Cleaning tanks

Outer width approx.

330 mm

Outer depth approx.

711 mm

Outer height approx.

102 mm