

Article Number: 4.272.20

Product Range Safety equipment

Product Category

capacity: 97 litres

outside dimensions: 895 x 496 x 845 (WxDxH) mm

basin (inside): 899 x 406 x 305 (WxDxH) mm

Basin, lid and frame made from galvanised sheet steel, red coated.

3/4"- drain plug in basin bottom for selective emptying.

Slow, independent closing of lid via pneumatic cylinder.

testing:

capacity: 97 litres

outside dimensions: 895 x 496 x 845 (WxDxH) mm

basin (inside): 899 x 406 x 305 (WxDxH) mm

Basin, lid and frame made from galvanised sheet steel, red coated.

3/4"- drain plug in basin bottom for selective emptying.

Slow, independent closing of lid via pneumatic cylinder.

testing:

HW-category	Cleaning tanks
Outer width approx.	406 mm
Outer depth approx.	895 mm
Outer height approx.	305 mm