

Article Number: 4.272.60

Product Range Safety equipment

Product Category

capacity: approx. 2 x 34 litres

Outside dimensions: 895 x 496 x 845 mm (WxDxH)

basin (inside): 2x438 x 406 x 222 mm (WxDxH)

Basin, lid and frame made from galvanised sheet steel, red coated.

For cleaning of extremely contaminated components.

1. chamber: pre-cleaning,

2. chamber: final purification.

3/4"- drain plug in bottom of each chamber

for selective emptying.

Testing: FM-approved

capacity: approx. 2 x 34 litres

Outside dimensions: 895 x 496 x 845 mm (WxDxH)

basin (inside): 2x438 x 406 x 222 mm (WxDxH)

Basin, lid and frame made from galvanised sheet steel, red coated.

For cleaning of extremely contaminated components.

1. chamber: pre-cleaning,

2. chamber: final purification.

3/4"- drain plug in bottom of each chamber

for selective emptying.

Testing: FM-approved

HW-category

Cleaning tanks

Outer width approx.

406 mm

Outer depth approx.

895 mm

Outer height approx.

305 mm