

Article Number: 9.31.405

**Product Range** Safety equipment

**Product Category**

Sterile sodium chloride solution (0.9%).

Equipped with an ergonomically formed eye cup and dust lid.

The bottle can be used individually or to refill eye rinsing stations.

Shelf life: 3 years.

Rinsing time: approx. 5 minutes.

According to DIN EN 15154-4:2009.

Sterile sodium chloride solution (0.9%).

Equipped with an ergonomically formed eye cup and dust lid.

The bottle can be used individually or to refill eye rinsing stations.

Shelf life: 3 years.

Rinsing time: approx. 5 minutes.

According to DIN EN 15154-4:2009.

**HW-category**

Emergency products